[LBIS Human IL-7 ELISA Kit]

Cat # 637-50441

Please, read this instruction carefully before use.

This kit is manufactured by FUJIFILM Wako Pure Chemical Corporation. Use only the current version of Instruction Manual enclosed with the kit!

Intended use

Interleukin 7 (IL-7) is a glycoprotein cytokine with a molecular weight of approximately 25000 Da consisting of 152 amino acid residues and is an essential cytokine for the adaptive immune system. It is mainly produced by non-hematopoietic cells, including thymocytes, bone marrow and lymphoid tissue stromal cells, keratinocytes, and intestinal epithelial cells.

IL-7 is considered a hematopoietic cytokine that contributes to host defense by its critical and nonredundant role in inducing differentiation and homeostasis of immune cells, including NK cells, B cells, and T cells. Therefore, it is a promising candidate for an immunotherapeutic agent against chronic infectious diseases. It is also associated with rheumatoid arthritis, chronic colitis, and asthma. In addition, the substantial immune reconstitution effect has been demonstrated. Its function in promoting antiviral activity and as a booster factor for anticancer drugs and immune checkpoint inhibitors has been revealed; thus, the high-sensitivity IL-7 measurement is expected to be helpful in immunotherapy.

This kit is a sandwich ELISA system for quantitative measurement of human IL-7 and can precisely measure IL-7 in human serum (plasma) and with high sensitivity.

Research use only. Not for use in diagnostics procedures.

Assay principle

In FUJIFILM Wako Pure Chemical Corporation's LBIS Human IL-7 ELISA Kit, standards or samples are incubated in polyclonal antibody-coated wells to capture IL-7. After 1 hour incubation and washing, biotinylated anti-IL-7 antibody is added and incubated further for 1 hour to bind with captured IL-7. After washing, HRP (horse radish peroxidase)-conjugated streptavidin is added, and incubated for 30 minutes. After washing, bound HRP-conjugated streptavidin is reacted with a chromogen (TMB) for 20 minutes, and reaction is stopped by addition of acidic solution, and absorbance of yellow product is measured spectrophotometrically at 450 nm. The absorbance is nearly proportional to IL-7 concentration. The standard curve is prepared by plotting absorbance against standard IL-7 concentrations. IL-7 concentrations in unknown samples are determined using this standard curve.

Reagents supplied

Components	State	Amount
(A) Anti-IL-7 antibody coated plate	Use after washing.	96 wells/1 plate
(B) Standard human IL-7	Freeze-dried. Use after reconstitution.	1 vial
(C) Buffer solution	Ready for use.	60 mL/1 bottle
(D) Biotinylated anti-IL-7 antibody	Freeze-dried. Use after reconstitution.	1 vial
(E) HRP-conjugated streptavidin	Concentrated. Use after dilution.	150 µL/1 vial
(F) Chromogen (TMB)	Ready for use.	12 mL/1 bottle
(H) Stop solution Be careful!	Ready for use.	12 mL/1 bottle
(I) Wash stock solution (10×)	Concentrated. Use after dilution.	100 mL/1 bottle
Plate seal	—	4 sheets
Instruction Manual	—	1 copy

Storage and expiration

When the complete kit is stored at 2 °C - 8 °C (Do not freeze), the kit is stable until the expiration date shown on the label on the container. Opened reagents should be used as soon as possible to avoid less than optimal assay performance caused by storage environment.

Equipments or supplies required but not supplied Use as a check box

- Deionized water (or Distilled water)
- Test tubes for preparation of standard solution series.
- Glassware for dilution of Wash stock solution (10×) (a graduated cylinder, a bottle)
- Pipettes (disposable tip type). One should be able to deliver 10 µL - 100 µL precisely, and another for 100 µL - 1000 µL.
- Syringe-type repeating dispenser like Eppendorf multipipette plus which can dispense 100 µL.
- Paper towel to remove washing buffer remaining in wells.
- A vortex-type mixer.
- A shaker for 96 well-plate (600 rpm -1200 rpm)
- An automatic washer for 96 well-plate (if available), or a wash bottle with a jet nozzle.
- A 96 well-plate reader (450 nm ± 10 nm, 620 nm: 600 nm - 650 nm)
- Software for data analysis.

Preparation of reagents

- ◆ Bring all reagents of the kit to room temperature (20 °C - 25 °C) before use.
- ◆ Prepare reagent solutions in appropriate volume for your assay. Do not store the diluted reagents.

【Concentrated reagents】

[(B) Standard human IL-7 (Freeze-dried)]


Reconstitute the (B) standard with volume* of deionized water (or distilled water) described in a separate sheet to prepare the Original standard solution (2360 pg/mL). And then prepare dilution series of standard solutions as shown below using the (C) Buffer solution brought to room temperature.

*Because the volume of deionized water (or distilled water) to be added to the freeze-dried standard differs depending on the lot, see the specified volume in a separate sheet.

Below is an example of preparing each standard solution.

Volume of standard solution	(C) Buffer solution	Concentration
Original standard solution 50 µL	450 µL	236 pg/mL
236 pg/mL solution 200 µL	350 µL	85.8 pg/mL
85.8 pg/mL solution 200 µL	350 µL	31.2 pg/mL
31.2 pg/mL solution 200 µL	350 µL	11.3 pg/mL
11.3 pg/mL solution 200 µL	350 µL	4.13 pg/mL
4.13 pg/mL solution 200 µL	350 µL	1.50 pg/mL
Blank	350 µL	0

} Standard Curve

[(D) Biotinylated anti-IL-7 antibody]

Add 150 µL of deionized water (or distilled water) of room temperature to (D) Biotinylated anti-IL-7 antibody (Freeze-dried) (Original (D) solution). Prepare working solution by dilution of original (D) solution with the (C) Buffer solution to **1:100**. 10 mL of the diluted solution is enough for 96 wells.

[(E) HRP-conjugated streptavidin]

Prepare working solution by dilution of (E) with the (C) Buffer solution to **1:100**. 10 mL of the diluted solution is enough for 96 wells.

[(I) Wash stock solution (10×)]

Dilute 1 volume of the concentrated (I) Wash stock solution (10×) to **10 volume** with deionized water (or distilled water) to prepare working solution. Example: 100 mL of concentrated (I) Washing buffer (10×) and 900 mL of deionized water (or distilled water).

【Storage and stability】

[(A) Anti-IL-7 antibody coated plate]

If seal is not removed, put the strip back in a plastic bag with zip lock originally used for well-plate container and store at 2 °C - 8 °C. The strip will be stable until expiration date.

[(B) Standard human IL-7 (Freeze-dried)]

The reconstituted standard solution (original standard solution) should be used within 2 weeks if stored in a refrigerator. Dispose the remaining diluted standard solutions after use.

[(C) Buffer solution] & [(F) Chromogen (TMB)]

Use only volume you need for your assay. Remaining reagents should be stored at 2 °C - 8 °C closing the cap tightly. It maintains stability until expiration date. Once opened, we recommend using as soon as possible to avoid influence by environmental condition.

[(D) Biotinylated anti-IL-7 antibody]

The reconstituted (D) solution should be used within 2 weeks if stored in a refrigerator. Dispose the remaining diluted working solution.

[(E) HRP-conjugated streptavidin]

Remaining working solution (already diluted) should be disposed. The rest of the undiluted solution (unused) : close the cap tightly and store at 2 °C - 8 °C. It is stable until expiration date. Once opened, we recommend using as soon as possible to avoid influence by environmental condition.

[(H) Stop solution]

Close the cap tightly and store at 2 °C - 8 °C. It maintains stability until expiration date.

[(I) Wash stock solution (10×)]

The rest of undiluted solution (unused): close the cap tightly and store at 2 °C -8 °C, it is stable until expiration date. Dispose any remaining diluted buffer.

Preparation of samples

- This kit is intended to measure IL-7 in human serum or plasma.
- We recommend to use EDTA-2Na, EDTA-2K or heparin as anticoagulant.
- Samples should be immediately assayed or stored below –35 °C until assay. Before starting assay, shake thawed samples sufficiently. Do not repeat freeze-and-thaw cycles.
- Centrifuge the sample to remove turbidity or insoluble matters where necessary before use for the assay.
- If presence of interfering substance is suspected, examine by dilution test at more than 2 points.
- Dilution of a sample should be made in a PP or PE test tube using buffer solution prior to adding them to wells.

Assay procedure


Remove the cover sheet of the antibody-coated plate after bringing up to room temperature.

- (1) Wash the (A) anti-IL-7 antibody coated plate by filling the wells with washing buffer and discard 4 times (*①), then strike the plate upside-down onto folded several sheets of paper towel to remove residual buffer in the wells (*⑤).
- (2) Pipette **100 µL** of standard solution to the wells designated for standards.
- (3) Pipette **50 µL** of (C) Buffer solution to the sample wells, and then add **50 µL** of sample to the wells.
- (4) Shake the plate gently on a plate shaker (*②).
- (5) Stick a plate seal (*③) on the plate and incubate for 1 hour at 20 °C - 25 °C.
- (6) Discard the reaction mixture and rinse wells as step (1).

- (7) Pipette **100 µL** of biotinylated anti-IL-7 antibody to all wells, and shake as step (4).
 - (8) Stick a plate seal (*③) on the plate and incubate the plate for 1 hour at 20 °C - 25 °C.
 - (9) Discard the reaction mixture and rinse wells as step (1).
 - (10) Pipette **100 µL** of HRP-conjugated streptavidin to all wells, and shake as step (4).
 - (11) Stick a plate seal (*③) on the plate and incubate the plate for 30 minutes at 20 °C - 25 °C.
 - (12) Discard the reaction mixture and rinse wells as step (1).
 - (13) Pipette **100 µL** of (F) Chromogen (TMB) to wells, and shake as step (4).
 - (14) Stick a plate seal (*③) on the plate and incubate the plate for 20 minutes at 20 °C - 25 °C.
 - (15) Add **100 µL** of the (H) stop solution to all wells and shake as step (4).
 - (16) Measure the absorbance of each well at 450 nm (reference wavelength, 620 nm (*④)) immediately using a plate reader.
- *Refer to the page 7 - 8 for notes of *①, *②, *③, *④, and *⑤.

Calculations

- (1) Create a standard curve by reducing the data using computer software capable of generating 3rd order regression curve, or 4 or 5 parameters. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the human IL-7 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.
- (2) Using the standard curve, read the IL-7 concentration of samples at its absorbance, and multiply the assay value by dilution factor if the sample has been diluted. The standard operation method is 5-fold dilution.

●Precision of assay (Within assay variation)
2 serum samples, 5 replicates assay

n / ID	Sample1 (pg/mL)	Sample2 (pg/mL)
1	2.00	30.4
2	2.13	31.6
3	2.21	31.5
4	2.14	32.1
5	2.14	31.3
mean	2.13	31.4
SD	0.075	0.64
CV(%)	3.5	2.0

Mean CV was within 15 %.

●Reproducibility (Between assay variation)
3 serum samples, 4 days, 3 replicates assay

Day /ID	Sample1 (pg/mL)	Sample2 (pg/mL)	Sample3 (pg/mL)
0	124	7.72	4.08
1	126	7.93	3.99
2	129	7.86	4.18
3	121	7.64	3.88
mean	125	7.79	4.03
SD	3.3	0.13	0.13
CV(%)	2.6	1.7	3.2

Mean CV was within 15 %.

- Recovery test (Standard human IL-7 was added in 4 concentrations to serum or plasma sample and assayed.) The recoveries were 92.0 % - 105 %



Serum sample

Added (pg/mL)	Found (pg/mL)	Recovered (pg/mL)	Recovery (%)
0	8.37	–	–
12.5	19.9	11.5	92.0
25.0	31.9	23.5	94.0
50.0	57.5	49.1	98.2
100	113	105	105

Plasma sample (EDTA)

Added (pg/mL)	Found (pg/mL)	Recovered (pg/mL)	Recovery (%)
0	8.33	–	–
12.5	19.9	11.6	92.8
25.0	32.4	24.1	96.4
50.0	57.7	49.4	98.8
100	109	101	101

- Dilution test (Serum sample and plasma sample were serially diluted by 3 steps.) The dilution curves showed linearity with $R^2=0.998 - 1$

Trouble shooting

- Low absorbance in all wells
Possible explanations:
 - 1)The standard or samples might not be added.
 - 2)Reagents necessary for coloration such as biotinylated anti-IL-7 antibody, HRP-conjugated streptavidin, or chromogen (TMB) might not be added.
 - 3)Wrong reagents related to coloration might have been added. Wrong dilution of biotinylated anti-IL-7 antibody or HRP-conjugated streptavidin.
 - 4)Contamination of enzyme inhibitor(s).
 - 5)Influence of the temperature under which the kits had been stored.
 - 6)Excessive hard washing of the well plate.
 - 7)Addition of chromogenic substrate reagent soon after taking out from a refrigerator might cause poor coloration owing to low temperature.
- Blank OD was higher than that of the lowest standard concentration (1.50 pg/mL).
Possible explanations:

Improper or inadequate washing. (Change washing frequency from 4 times to 5 - 8 times at the constant stroke after the reaction with HRP-conjugated streptavidin.)
- High coefficient of variation (CV)
Possible explanation:
 - 1)Improper or inadequate washing.
 - 2)Improper mixing of standard or samples.
 - 3)Pipetting at irregular intervals.
- Q-1: Can I divide the plate to use it for the other testing?
A-1: Yes, cut off the clear seal on the plate with cutter along strip. Put the residual plate, which is still the seal on, in a refrigerator soon.
- Q-2: I found 96 well-plate is empty when I opened the box.

A-2: As this kit is dried type, not preservation stabilizer is added.

Technical tips / Precautions

- ELISA can be affected by assay environment. Ensure that room temperature at places for assay operation and incubation is strictly controlled at 20 °C to 25 °C. Avoid performing assay under the airstream velocity over 0.4 m/sec. and the humidity less than 30 %, seal the well plate with a plate seal and place the well plate in an incubator or a styrofoam box in each step of incubation.
- To avoid denaturation of the coated antibody, do not let the plate go dry.
- In order to avoid dryness of wells, contamination of foreign substances and evaporation of dispensed reagents, never forget to cover the well plate with a plate seal supplied, during incubation.
- Be careful to avoid any contamination of assay samples and reagents. We recommend the use of disposal pipette tips, and 1 tip for 1 well.
- As the antibody-coated plate is module type of 8 wells × 12 strips, each strip can be separated by cutting the cover sheet with a knife and used independently.
- The chromogen (TMB) should be almost clear pale yellow before use. It turns blue during reaction, and gives yellowish color after addition of stop solution. Greenish color means incomplete mixing.
- Time the reaction from the pipetting of the reagent to the first well.
- Prepare a standard curve for each assay.
- For professional use only, beginners are advised to use this kit under the guidance of experienced person. In manual operation, proficiency in pipetting technique is recommended.
- Wear protective gloves, clothing, eye, and face protection.
- Avoid contact of skin and mucous membranes with kit reagents or specimens. If any reagents come in contact with eyes, skin, or mucous membranes, wash with copious amounts of water and contact a physician.
- Do not drink, eat, or smoke in a place where this kit is used.
- The materials must not be pipetted by mouth.
- Do not use reagents with different lot numbers together.
- The reagents are prepared to give accurate results only when used in combination within the same box. Therefore, do not combine the reagents from kits with different lot numbers. Even if the lot number is the same, it is best not to mix the reagents with those that have been preserved for some period.
- Handle the sample with proper care, being aware that the sample may have an infection risk. This kit contains animal-derived ingredients.
- Residual samples and used tips should be sterilized before disposal.
- Dispose consumable materials and unused contents in accordance with applicable regional/national regulatory requirements.

Summary of assay procedure □ : Use as a check box

*First, read this instruction manual carefully and start your assay after confirmation of details.

- Bring the well-plate and all reagents back to **20 °C - 25 °C for 2 hours.**
- Concentrated washing buffer must be diluted to **10 times** by deionized water (or distilled water) that returned to 20 °C - 25 °C.
- Standard solution dilution example
 Reconstitute the (B) standard with volume* of deionized water (or distilled water) described in a separate sheet to prepare the Original standard solution (2360 pg/mL). And then prepare dilution series of standard solutions as shown below using the (C) Buffer solution brought to room temperature.
 *Because the volume of deionized water (or distilled water) to be added to the freeze-dried standard differs depending on the lot, see the specified volume in a separate sheet.

Preparation of diluted IL-7 standard solutions:

Conc.(pg/mL)	236	85.8	31.2	11.3	4.13	1.50	0
Std. sol.(μL)	50	200*	200*	200*	200*	200*	—
Buffer(μL)	450	350	350	350	350	350	350

*One rank higher standard.

□ Anti-IL-7 antibody coated plate	
□ ↓Washing 4 times (*①) , (*⑤)	
□ Samples (diluted sample : Buffer solution 50 μL + Sample 50 μL), or Standards	100 μL
□ ↓Shaking (*②), Incubation for 1 hour at 20 °C - 25 °C. (Standing (*③)) Add deionized water (or distilled water) of room temperature to (D) Biotinylated anti-IL-7 antibody (Freeze-dried) (Original (D) solution). Prepare working solution by dilution of original (D) solution to 100x with the (C) Buffer solution returned to 20 °C - 25 °C. (This should be prepared during incubation.)	
□ ↓Washing 4 times (*①) , (*⑤)	
□ Biotinylated anti-IL-7 antibody	100 μL
□ ↓Shaking (*②), Incubation for 1 hour at 20 °C - 25 °C. (Standing (*③)) Dilute HRP-conjugated streptavidin to 100x with the (C) Buffer solution returned to 20 °C - 25 °C. (This should be prepared during incubation.)	
□ ↓Washing 4 times(*①) , (*⑤)	
□ HRP-conjugated streptavidin	100 μL
□ ↓Shaking (*②), Incubation for 30 minutes at 20 °C - 25 °C. (Standing (*③))	
□ ↓Washing 4 times(*①) , (*⑤)	
□ Chromogen (TMB) After dispense, the color turns to blue depending on the concentration.	100 μL
□ ↓Shaking (*②), Incubation for 20 minutes at 20 °C - 25 °C. (Standing (*③))	
□ Stop solution After dispense, the color turns to yellow depending on the concentration.	100 μL
□ ↓Shaking (*②) Immediately shake.	
□ Measurement of absorbance (450 nm, Ref 620 nm (*④)) immediately. Ref. wave cancels the dirt in the back of plate.	

*①After dispensing wash buffer to wells, lightly shake the plate on your palm for 10 sec and remove the buffer.
 Guideline of washing volume: 300 μL/well for an automatic washer and for a pipette if the washing buffer is added by pipette. In case of washing by using 8 channel pipette, sometimes the back ground tends to

LBIS Human IL-7 ELISA Kit

be high. If so, change washing frequency from 4 times to 5 - 8 times at the constant stroke after the reaction with HRP-conjugated streptavidin.

Standard of plate-washing pressure: 5 mL/min - 25 mL/min. (Adjust it depending on the nozzle's diameter.)


*② Guideline of shaking: 600 rpm - 1,200 rpm for 10 seconds×3 times.

*③ Seal the plate during the reaction after shaking. Peel off the protective paper from the seal and stick the seal on the plate. Do not reuse the plate seal used once.

*④ 600 nm - 650 nm can be used as reference wavelength.

*⑤ After removal of wash buffer, immediately dispense the next reagent.

Assay worksheet

LBIS Human IL-7 ELISA Kit

[Storage condition] Store the kit at 2 °C - 8 °C (Do not freeze)

[Term of validity] Expiration date is indicated on the container

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