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LabAssay™ ALP

(Alkaline Phosphatase activity assay with *p*-Nitrophenylphosphate as a substrate)

For Research Use Only Not for use in diagnostic procedure

[Introduction]

Alkaline Phosphatase (ALP) is distributed in a variety of tissues such as liver, bone, and small intestine in animals. The change of the enzyme activity in tissues is an important hallmark for physiological phenomena as osteogenesis and so on. This kit is for Alkaline Phosphatase assay in a simultaneous multi-sample assay format with a microplate using *p*-Nitrophenylphosphate as a substrate.

[Kit contents]

	Kit contents	State	Amount
(1)	Substrate Tablet (<i>p</i> -Nitrophenylphosphate Disodium 6.7 mmol/L, after dissolving)	tablet	10 tablets
(2)	Buffer Solution (2.0 mmol/L MgCl ₂ , 0.1 mol/L Carbonate Buffer, pH 9.8)	liquid	50 mL/1 bottle
(3)	Stop Solution (0.2 mol/L Sodium Hydroxide Solution)	liquid	50 mL/1 bottle
(4)	Standard Solution (0.5 mmol/L <i>p</i> -Nitrophenol Solution)	liquid	5 mL/1 bottle

※Assay in microplate format;500 tests.


※Substrate Tablet contains one spare tablet.

[Storage and expiration]

When the complete kit is stored at 2 °C ~ 10 °C, the kit is stable until the expiration date shown on the label on the container.

[Assay principle]

p-Nitrophenylphosphate is hydrolyzed into *p*-Nitrophenol and Phosphoric Acid in the carbonate buffer (pH 9.8) in the presence of Alkaline Phosphatase in sample. Released *p*-Nitrophenol showing yellow color is optically measured at 405 nm wavelength as the enzyme activity.

[Materials and apparatuses prepared]

- 96 well microplate (transparent type) ·Micropipette ·Pipette ·Test tube
- Plate mixer * ·Incubator maintaining at 37 °C *
- Microplate reader with 405 nm wavelength filter
- (* if the microplate reader is not equipped.)

[Preparation of reagents]

- ① Working Assay Solution :

Dissolve a Substrate Tablet into 5 mL of Buffer Solution provided.

* It is easy to dissolve the tablet crashed into pieces.

* Can be stored for two weeks at 2 °C~10 °C in the dark after reconstitution.

② Stop Solution : Ready-to-use

③ Dilution series of Standard Solution : Dilute Standard Solution seriously with distilled water, 0.5, 0.25, 0.125, 0.0625 mmol/L.

Conc.(mmol/L)	0.5	0.25	0.125	0.0625	0
Standard solution(μL)	Original std. sol.: 60	Original std. sol.:60	60*	60*	0
Distilled water(μL)	0	60	60	60	60

* One rank higher standard solution


[Procedure]

Perform the assay in the wells according to the following table scheme.

	Test	Standard	Blank
Working Assay Solution	100 μL	100 μL	100 μL
Sample	Sample 20 μL	Dilution series 20 μL	Distilled water 20 μL
Shake for 1 minute by plate mixer and incubate for 15 minutes at 37 °C			
Stop Solution	80 μL	80 μL	80 μL
Shake for 1 minute by plate mixer and measure the absorbance at 405 nm of wavelength using microplate reader.			

[Standard curve]

The example of standard curve is as follows according to the procedure of standard.

[Unit Definition]

According to the above conditions, one unit of the enzyme activity is defined as release of 1 nmol *p*-Nitrophenol per minute at pH 9.8, 37 °C.

$$\text{Activity (units/}\mu\text{L)} = \frac{C}{15} \times a$$

C : Concentration of *p*-Nitrophenol released by sample, which is calculated as $A_{\text{test}} - A_{\text{blank}}$
(mmol/L=nmol/ μ L)

15 : Reaction Time (min.)

a : Dilutin factor of sample

The relative activity of sample can be determined as follows.

$$\text{units/}\mu\text{g protein or DNA} = \frac{\text{Activity (units/}\mu\text{L)}}{\text{Protein or DNA concentration (}\mu\text{g/}\mu\text{L)}}$$

[Features]

Dynamic assay range : > 0.06 mmol/L

Standard assay range : 0~0.5 mmol/L

(The activity at concentration more than 0.5 mmol/L can be determined by extrapolation due to its linearity.)

Reproducibility : C.V. < 10 %

[Notes]

- Serum sample should be used immediately after collecting blood. Do not use EDTA as an anticoagulant in preparation of blood sample, because of reduction of the activity by the reagent.
- Avoid unnecessary contamination of phosphate in the assay, because of reduction of the activity.
- Bilirubin and materials generated by hemolysis may not significantly affect the assay.
- In dilution of sample, use distilled water or an appropriate buffer, but not phosphate buffer.
- Keep reagents provided indicated condition, before their expiration.
- Avoid exposure of the reagents to direct sunlight.
- Be careful contamination during pipette solutions.
- Be careful of handling of hazardous solutions such as phenols, sodium hydroxide.
- This kit is for biochemical research, but not for diagnostic use.

[Reference]

Yamamoto, M., Takahashi, Y., Tabata, Y.: *Biomaterials*. 24(24), 4375(2003).

LabAssay™ ALP

[Storage condition] Store the kit at 2 °C ~ 10 °C

[Manufacture #] LABALP-M1

[Cat #] 633-51021

This kit is
manufactured by

FUJIFILM Wako Shibayagi Corporation

1062-1 Ishihara, Shibukawa, Gunma, Japan 377-0007

TEL.+81-279(25)0279, FAX.+81-279(23)0313

<E-mail>wksb-info@fujifilm.com <URL><http://www.shibayagi.co.jp>

distributed by

FUJIFILM Wako Pure Chemical Corporation