LabAssay™ Triglyceride (GPO•DAOS method)

For Research Use Only. Not for use in diagnostic procedures.

(Introduction)

Lipids in serum consist of triglycerides, cholesterol, phospholipids, free fatty acids and slight amounts of fat-soluble components such as fat-soluble vitamins and carotenes. Triglycerides, as major components of very low density lipoprotein (VLDL) and chylomicrons, play an important role in metabolism as energy sources and transporters of dietary fat.

LabAssay™ Triglyceride is based on an enzymatic method using *N*-ethyl-*N*-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxyaniline sodium salt (DAOS) as a blue pigment. This kit is used for the quantitative determination of triglycerides in mouse serum. It is a simultaneous multi-sample assay format using a microplate, but measurements can also be made using a test tube.

(Kit contents)


	Kit contents	State	Amount
(1)	Buffer Solution 50 mmol/L Good's (PIPES) buffer, pH6.5	liquid	105 mL/1 bottle
(2)	Chromogen Substrate (when reconstituted) Lipoprotein lipase 99 units/mL Adenosine 5'-Triphosphate Disodium Salt Trihydrate (ATP) 1.7 mmol/L Glycerolkinase 30 units/mL Glycerol-3-phosphate Oxidase (GPO) 1.0 unit/mL Peroxidase (POD) 5.5units/mL <i>N</i> -Ethyl- <i>N</i> -(2-hydroxy-3- sulfopropyl)-3,5-dimethoxyaniline sodium salt (DAOS) 0.54 mmol/L 4-Aminoantipyrene 0.11 mmol/L Ascorbate Oxidase 4.5 units/mL	freeze dry	For 105 mL/1 bottle
(3)	Standard Solution Glycerol 31.2 mg/dL (Corresponding to 300 mg/dL Triolein) Contains Sodium Azide 0.05 %	liquid	4 mL/1 bottle

※Assay in a microplate;350 tests. Assay in a test tube; 35 tests.

(Assay principle)

Triglycerides in a sample are hydrolyzed to glycerol and free fatty acids in a reaction catalyzed by lipoprotein lipase (LPL). Glycerol is converted to glycerol-3-phosphate by glycerolkinase (GK) in the presence of ATP.

Glycerol-3-phosphate so formed is oxidized by glycerol-3-phosphate oxidase (GPO) in a reaction that produces hydrogen peroxide. The hydrogen peroxide produced causes DAOS and 4-Aminoantipyrene to undergo a quantitative oxidative condensation catalyzed by peroxidase (POD), producing a blue pigment. The amount of triglycerides contained in the sample is determined by measuring the absorbance of the blue color.

[Materials and apparatuses to be prepared]

- 96wells microplate (transparent type) ·Micropipette ·Incubator maintaining at 37 °C *
- Plate mixer*
- Microplate reader with 600 nm wavelength filter (* if the microplate reader is not equipped.)

(For test tube method)

- Test tube ·Pipette ·Incubator maintaining at 37 °C
- Spectrophotometer or colorimeter with 600 nm wavelength filter

[Preparation of reagents to be used]

① Chromogen Reagent :

Prepare Chromogen Reagent by Dissolving 1 vial of Chromogen Substrate (for 105 mL) to 105 mL of Buffer Solution. After reconstitution, the solution should be stored at 2 °C~10 °C and used within 10 days.

② Standard solution (Microplate method)

Standard solution is prepared by dilution of the provided Standard

No.	Triglyceride Standard	Distilled or deionized water	Sample volume	Triglyceride Concentration
1	100 μL	200 μL	2 μL	100 mg/dL
2	200 μL	100 μL	2 μL	200 mg/dL
3	undiluted	—	2 μL	300 mg/dL
4	undiluted	—	4 μL	596 mg/dL *1
5	undiluted	—	6 μL	888 mg/dL *1

* 1 The test sample volume is usually 2 μL, but 4 or 6 μL are taken in this case.

The triglyceride concentration must be corrected accordingly as indicated in the table above.

[Procedure]

(1) Assay in a microplate

Perform the assay in the wells according to the following table scheme.

	Test	Standard	Blank
Chromogen reagent	300 μL	300 μL	300 μL
Sample	Serum 2 μL	Standard solution 2 μL	—
	Mix well and incubate at 37 °C for 5 min. Measure the absorbance at 600 nm *2 of the test sample and standard solution with the blank solution as the control.		

* 2 In two wavelength assay, measure using Main wavelength 600 nm/Sub wavelength 700 nm.


(2) Assay in a test tube

Perform the assay in a test tube according to the following table scheme.

	Test	Standard	Blank
Chromogen reagent	3.0 mL	3.0 mL	3.0 mL
Sample	Serum 20 µL	Standard solution 20 µL	–
	Mix well and incubate at 37 °C for 5 min. Measure the absorbance of the test sample and standard solution with the blank solution as the control. Colorimeter with 600 nm. Spectrophotometer : 600 nm ²		

* 2 In two wavelength assay, measure using Main wavelength 600 nm/Sub wavelength 700 nm

[Standard curve]

[Performance]

(1) Sensitivity [assay in a test tube]

- The absorbance is below 0.10, when measuring purified water as a sample.
- The absorbance is 0.09 to 0.25, when measuring 300 mg/dL triglycerides as a sample.

(2) Specificity

The triglycerides concentration is less than ± 12 %, when measuring the known concentration of control serum as a sample.

[Usage Notes]

(1) Sample

- Samples should be used immediately after collecting.
- Hemolysis causes an increase of the triglycerides value.
- Ascorbic acid causes a slightly negative effect on the assay.
- Bilirubin may not significantly affect the assay.

(2) Notes on the assay

- Do not use the reagents described above in any procedures other than those described herein. Performance cannot be guaranteed if the reagents are used in other procedures or for other purposes.
- Operate the measurement apparatuses according to operator’s manuals under appropriate conditions.

LabAssay™ Triglyceride (GPO·DAOS method)

Consult the apparatus manufacturer for details.

- Keep the provided reagents under the indicated conditions before the expiration.
- Do not use reagents, which were frozen in error. Such reagents may give false results.
- After opening the reagents, it is recommended to use them immediately. The opened vials should be capped and kept under the specified conditions.
- Do not use the containers and other materials in the package for any purpose other than those described herein.
- The vial is sealed at reduced pressure. Slowly remove the stopper in order not to blow the powder in the vial.
- Do not use the reagents at other reaction temperatures and reaction times than described herein.
- The reaction induced by Chromogen Reagent is completed in about 5 min at 37 °C. The reaction is also completed within 15 min at room temperature of over 20 °C.
- When the incubation is continued at 37 °C, the absorbance gradually decreases. Therefore, the incubation should be done within 20 min.
- Color level of the reaction solution should show very little change within 1 hour at room temperature.
- This product calculates the concentration of triolein, which exists most within an organism. The concentration of glycerol is 31.2 mg/ dL (corresponding to 300 mg/dL triolein).
- This kit is for research use only. Not for diagnostic use.

(3) Safety precautions

- If a reagent comes into contact with the mouth, eyes, or skin, immediately wash with a lot of water. Consult a physician if necessary.
- A pipette with a safety pipette filler should be used.
- Be careful not to cut yourself with the aluminum cap when removing it from the vial.

(4) Waste

- The waste should be processed appropriately according to regulations.
- All the devices including reagents and vials that come into contact with specimens should be considered potentially infectious.
- Triglyceride Standard Solution contains 0.05 % sodium azide. Sodium azide may react with lead or copper plumbing to form explosive compounds. Even though the reagent contains minute quantity of sodium azide, drains should be flushed well with a large amount of water, when discarding the reagents.

[Reference]

1. Spayd, R. W., Bruschi, B. et al. : *Clin Chem.*, 24, 1343 (1978) .

LabAssay™ Triglyceride (GPO·DAOS method)

[Storage condition] Store the kit at 2 °C~10 °C

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