LabAssay™ Phospholipid (Choline Oxidase•DAOS method)

For Research Use Only. Not for use in diagnostic procedures.

[Introduction]

Phospholipids are a class of lipids, and a major component of all biological membranes, along with glycolipids, cholesterol and proteins. LabAssay™ Phospholipid is based on an enzymatic methods using *N*-ethyl-*N*-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxyaniline sodium salt (DAOS) as a blue pigment. This kit is used for the quantitative determination of phospholipids in mouse serum. It is a simultaneous multi-sample assay format using a microplate, but measurements can also be made using a test tube.

[Kit contents]

	Kit contents	State	Amount
(1)	Buffer Solution ( 50 mmol/L Good’ s buffer, pH7.5)	liquid	50 mL/3 bottles
(2)	Chromogen Substrate After reconstitution Phospholipase D 0.47 units/mL Choline Oxidase 2.0 units/mL Peroxidase 4.2 units/mL <i>N</i> -Ethyl- <i>N</i> -(2-hydroxy -3-sulfopropyl)-3,5-dimethoxyaniline sodium salt (DAOS) 0.77 mmol/L 4-Aminoantipyrine 0.24 units/mL Ascorbate Oxidase 3.9 units/mL	freeze dry	For 50 mL/3 bottles
(3)	Standard Solution Choline Chloride 54 mg/dL (Corresponding to 300 mg/dL Phospholipid)	liquid	5 mL/1 bottles


※Assay in a microplate;500 tests. Assay in a test tube; 50 tests.



[Storage and expiration]

When the complete kit is stored at 2 °C ~ 10 °C, the kit is stable until the expiration date shown on the label on the container.

[Assay principle]

Phospholipids (lecithin, sphingomyelin and lysolecithin) in a sample are hydrolyzed to choline in a reaction catalyzed by phospholipase D. Choline so formed is oxidized by choline oxidase in a reaction that produces hydrogen peroxide. The hydrogen peroxide produced causes DAOS and 4-Aminoantipyrine to undergo a quantitative oxidative condensation catalyzed by peroxidase(POD), producing a blue pigment. The amount of phospholipids contained in the sample is determined by measuring the absorbance of the blue color.

### [Materials and apparatuses to be prepared]

- 96wells microplate (transparent type) ·Micropipette ·Incubator maintaining at 37 °C \*
- Plate mixer \*
- Microplate reader with 600 nm wavelength filter (\* if the microplate reader is not equipped.)

### (For test tube method)

- Test tube ·Pipette ·Incubator maintaining at 37 °C
- Spectrophotometer or colorimeter with 600 nm wavelength filter

### [Preparation of reagents to be used]

#### ①Chromogen Substrate :

Dissolve 1 vial of Chromogen Substrate with 50 mL of Buffer Solution. After reconstitution, the solution should be stored at 2 °C~10 °C and used within one week.

#### ②Standard solution (Microplate method)

Standard solution is prepared by dilution of the provided Standard

No.	Phospholipid Standard Solution	Distilled or deionized water	Sample volume	Phospholipid Concentration
1	50 µL	150 µL	2 µL	75.0 mg/dL
2	100 µL	100 µL	2 µL	150 mg/dL
3	Undiluted	—	2 µL	300 mg/dL
4	Undiluted	—	4 µL	596.1 mg/dL * <sup>1</sup>

\* 1 The test sample volume is usually 2 µL, but 4 µL is taken in this case. The phospholipids concentration must be corrected accordingly as indicated in the table above.

### [Procedure]

#### (1) Assay in a microplate

Perform the assay in the wells according to the following table scheme.

	Test	Standard	Blank
Chromogen reagent	300 µL	300 µL	300 µL
Sample	Serum 2 µL	Standard solution 2 µL	—
	Mix well and incubate at 37 °C for 5 min. Measure the absorbance at 600 nm * <sup>2</sup> of the test sample and standard solution with the blank solution as the control.		

\* 2 In two wavelength assay, measure using Main wavelength 600 nm/Sub wavelength 700 nm.

## (2) Assay in a test tube


Perform the assay in a test tube according to the following table scheme.

	Test	Standard	Blank
Chromogen reagent	3.0 mL	3.0 mL	3.0 mL
Sample	Serum 20 $\mu$ L	Standard solution 20 $\mu$ L	—
	Mix well and incubate at 37 °C for 5 min. Measure the absorbance of the test sample and standard solution with the blank solution as the control. Colorimeter with 600 nm. Spectrophotometer : 600 nm <sup>*2</sup>		

\* 2 In two wavelength assay, measure using Main wavelength 600 nm/Sub wavelength 700 nm.

**[Standard curve]**

[assay in a microplate]

**[Performance]**

## (1) Sensitivity [assay in a test tube]

- The absorbance is below 0.14, when measuring purified water as a sample.
- The absorbance is 0.12 to 0.58, when measuring 300 mg/dL phospholipid as a sample.

## (2) Specificity

The phospholipid concentration is less than  $\pm 15\%$ , when measuring the known concentration of a control serum as a sample.

**[Usage Notes]**

## (1) Sample

- Samples should be used immediately after collecting.
- Ascorbic acid, bilirubin and hemolysis may not significantly affect the assay.
- This method measures lecithin, sphingomyelin, lysolecithin but not cephalin.

## (2) Notes on the assay

- Do not use the reagents described above in any procedures other than those described herein. Performance cannot be guaranteed if the reagents are used in other procedures or for other purposes.

#### LabAssay™ Phospholipid (Choline Oxidase·DAOS method)

- Operate the measurement apparatuses according to operator's manuals under appropriate conditions. Consult the apparatus manufacturer for details.
- Keep the provided reagents under the indicated conditions before the expiration.
- The reaction completes in about 3 min. It is enough to incubate at 37 °C for 5 min. When the incubation is continued, the absorbance gradually decreases. Therefore, the incubation should be done within 15 min.
- Color level of the reaction solution should show very little change within 1 hour at room temperature.
- This kit is for research use only. Not for diagnostic use.

#### (3) Safety precautions

- If a reagent comes into contact with the mouth, eyes, or skin, immediately wash with a lot of water. Consult a physician if necessary.
- A pipette should be used with a safety pipette filler.

#### (4) Waste

- The waste should be processed appropriately according to the law.
- All the devices including reagents and vials that come into contact with specimens should be considered potentially infectious.

#### [Reference]

Takayama, M., Itoh, S., Nagasaki, T. and Tanimizu, I. : *Clin. Chim. Acta* .,79, 93 (1977) .

#### LabAssay™ Phospholipid (Choline Oxidase·DAOS method)

[Storage condition]     Store the kit at   2 °C~10 °C

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