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LabAssay™ Glucose (Mutarotase-GOD method)

For Research Use Only Not for use in diagnostic procedure

[Introduction]

α-D-Glucose and β-D-glucose in solution maintain equilibrium in a constant ratio. Glucose oxidase reacts only with β-D-glucose and does not react with α-D-Glucose. Therefore, α-D-Glucose is converted into β-D-glucose using mutarotase.

LabAssay™ Glucose is the reagent kit for assay of glucose based on an enzymatic method with a combination of mutarotase and glucose oxidase.

The kit is used for the quantitative determination of glucose in mouse serum, plasma or urine. It is a simultaneous multi-sample assay format using a microplate, but measurements can also be made using a test tube.

[Kit contents]

	Kit contents	State	Amount
(1)	Buffer Solution Phosphate buffer, pH 7.1 60 mmol/L Phenol 5.3 mmol/L	liquid	150 mL/1 bottle
(2)	Chromogen Substrate After reconstitution Mutarotase 0.13 units/mL Glucose oxidase(GOD) 9.0 units/mL Peroxidase(POD) 0.65 units/mL 4-Aminoantipyrine 0.50 mmol/L Ascorbate oxidase(AOD) 2.7 units/mL	freeze dry	For 150 mL/1 bottle
(3)	Standard Solution I (200 mg/dL Glucose)	liquid	3 mL/1 bottle
(4)	Standard Solution II (500 mg/dL Glucose)	liquid	3 mL/1 bottle


※Assay in a microplate;500 tests. Assay in a test tube; 50 tests.

[Storage and expiration]

When the complete kit is stored at 2 °C ~ 10 °C, the kit is stable until the expiration date shown on the label on the container.

[Assay principle]

When a sample is mixed with the Choromogen Reagent, the α-form of glucose in the sample is converted to β-form by mutarotase. β-D-Glucose is oxidized and yields hydrogen peroxide by glucose oxidase(GOD). In the presence of peroxidase(POD), the formed hydrogen peroxide yields a red pigment by quantitative oxidation condensation with phenol and 4-aminoantipyrine. The glucose concentration is obtained by measuring absorbance of the red pigment.

[Materials and apparatuses to be prepared]

- 96wells microplate (transparent type) · Micropipette · Incubator maintaining at 37 °C *
- Plate mixer *
- Microplate reader with 505 nm wavelength filter (* if the microplate reader is not equipped.)

(For test tube method)

- Test tube · Pipette · Incubator maintaining at 37 °C
- Spectrophotometer or colorimeter with 505 nm wavelength filter

[Preparation of reagents to be used]

① Chromogen reagent :

Dissolve 1 vial of Chromogen Substrate with 150 mL of Buffer. After reconstitution, the solution should be stored at 2 °C ~ 10 °C and used within 1 month.

② Standard solution (Microplate method)

Standard solution is prepared by dilution of the provided Standard Solution.

No.	Standard Solution I	Standard Solution II	Distilled or deionized water	Sample volume	Concentration
1	50 µL	—	150 µL	2 µL	50 mg/dL
2	100 µL	—	100 µL	2 µL	100 mg/dL
3	undiluted	—	—	2 µL	200 mg/dL
4	—	150 µL	100 µL	2 µL	300 mg/dL
5	—	200 µL	50 µL	2 µL	400 mg/dL
6	—	undiluted	—	2 µL	500 mg/dL

[Procedure]

(1) Assay in a microplate

Perform the assay in the wells according to the following table scheme.

	Test	Standard	Blank
Chromogen reagent	300 µL	300 µL	300 µL
Sample	2 µL	2 µL	— *1
	Mix well and incubate at 37 °C for 5 min. Measure the absorbance at 505 nm*2 of the test sample and standard solution with the blank solution as the control.		

*1 The addition of water is omitted because the absorbance difference in the presence of the addition of 2 µL of water does not exhibit any practical use.

*2 Use 505 nm (Main) and 600 nm (Sub) for two wavelengths.

(2) Assay in a test tube

Perform the assay in a test tube according to the following table scheme.


	Test	Standard	Blank
Chromogen reagent	3.0 mL	3.0 mL	3.0 mL
Sample	20 µL	20 µL	— *1
	Mix well and incubate at 37 °C for 5 min. Measure the absorbance of the test sample and standard solution with the blank solution as the control. Colorimeter with 505 nm. Spectrophotometer : 505 nm*2		

*1 The addition of water is omitted because the absorbance difference in the presence of the addition of 20 µL of water does not exhibit any practical use.

LabAssay™ Glucose (Mutarotase-GOD method) (#LABGLUC-M1)

* 2 Use 505 nm (Main) and 600 nm (Sub) for two wavelengths.

[Standard curve]

[Performance]

(1) Sensitivity [assay in a test tube]

- The absorbance is below 0.03, when measuring purified water as a sample.
- The absorbance is 0.40 to 0.55, when measuring 200 mg/dL glucose as a sample.

(2) Specificity

The glucose concentration is less than $\pm 12\%$, when measuring the known concentration of a control serum or a urine as a sample.

[Usage Notes]

(1) Sample

- Samples should be used immediately after collecting.
- The glucose in whole blood decreases due to glycolysis in erythrocytes when a sample is kept after collecting blood. Therefore, separation of blood cells should be done as quickly as possible.
- Ascorbic acid and bilirubin cause a slightly negative effect on the assay.
- Hemolysis gives a slightly positive effect on the assay.
- Anticoagulants such as heparin, citrate, oxalate and EDTA and sodium fluoride as a glycolytic inhibitor, should not affect the assay when they are used in their usual amounts.

(2) Notes on the assay

- Do not use the reagents described above in any procedures other than those described herein. Performance cannot be guaranteed if the reagents are used in other procedures or for other purposes.
- Operate the measurement apparatuses according to operator's manuals under appropriate conditions. Consult the apparatus manufacturer for details.
- Keep the provided reagents under the indicated conditions before the expiration.
- Do not use reagents, which were frozen in error. Such reagents may give false results.
- After opening the reagents, it is recommended to use them immediately. The opened vials should be capped and kept under the specified conditions.
- Do not use the containers and other materials in the package for any purposes other than those described

LabAssay™ Glucose (Mutarotase-GOD method) (#LABGLUC-M1)

herein.

- Do not use the reagents at other reaction temperatures and reaction times than described herein.
- The reaction completes in about 4 min at 37 °C . When the incubation is continued, the absorbance gradually decreases. Therefore, the incubation should be done within 20 min.
- Glucose can be also measured by a reaction at room temperature (15 °C or more) for 15 min.
- Color level of the reaction solution should show very little change within 1 hour at room temperature.
- This kit is for research use only. Not for diagnostic use.

(3) Safety precautions

- If a reagent comes into contact with the mouth, eyes, or skin, immediately wash with a lot of water. Consult a physician if necessary.
- A pipette with a safety pipette filler should be used.
- Be careful not to cut yourself with the aluminum cap when removing it from the vial.

(4) Waste

- The waste should be processed appropriately according to regulations.
Buffer contains 500 mg/L phenol.
- All the devices including reagents and vials that come into contact with specimens should be considered potentially infectious.

[Reference]

Miwa, I., Okuda, J., Maeda, K. and Okuda, G. : *Clin. Chim. Acta* ., 37, 538 (1972) .

LabAssay™ Glucose (Mutarotase-GOD method)

[Storage condition] Store the kit at 2 °C~10 °C

[Manufacture #] LABGLUC-M1

[Cat #] 638-50971

This kit is
manufactured by

FUJIFILM Wako Shibayagi Corporation

1062-1 Ishihara, Shibukawa, Gunma, Japan 377-0007

TEL.+81-279(25)0279, FAX.+81-279(23)0313

<E-mail>wksb-info@fujifilm.com

<URL><http://www.shibayagi.co.jp>

distributed by

FUJIFILM Wako Pure Chemical Corporation