FUJIFILM**Wako**Code No. 294-84101 (2 purifications)
290-84103 (10 purifications)**For Genetic Research**
MagCapture™ Exosome Isolation Kit PS Ver.2**[Introduction]**

Extracellular Vesicles (EVs) represented by exosomes contain a protein, mRNA, microRNA, DNA, and lipids on the surface or inside. After being secreted from various cells, EVs stably exist in body fluids such as blood, urine, saliva, cerebrospinal fluid, and breast milk. Thus, it is attracting attention as a messenger for cell-to-cell communication and as a biomarker for diseases. As a method for isolating exosomes, an ultracentrifugation method, an affinity method using an antibody against a surface antigen, a precipitation method using a polymer reagent etc. are generally used. However, in all respects such as recovery efficiency, purity, and operability, there was no satisfying way. MagCapture™ Exosome Isolation Kit PS Ver.2 adopts a novel affinity purification method using magnetic beads and phosphatidylserine (PS)-binding protein (PS affinity method). The PS affinity method can easily isolate exosomes and other EVs with PS on the membrane surface from the cell culture medium and body fluids with high purity and high efficiency. Using the 10,000×g supernatant as a sample, PS affinity can get more highly purified exosomes. Besides, PS affinity enables the elution of intact EVs from magnetic beads with a chelating agent because PS affinity captures EVs in a metal ion-dependent manner. Therefore, purified EVs can be used for various applications such as electron microscopic analysis, nanoparticle tracking analysis, uptake assay, and analysis of constituent molecules (proteins, lipids, nucleic acids, etc.).

FUJIFILM Wako Pure Chemical Corporation Products (the Products) are intended for research use only.
If you use for commercial purpose of the Products, please contact us (ffwk-labchem-tec@fujifilm.com).

[Feature]

- Adopted new affinity method (PS affinity method).
- Can obtain EVs with higher yield and higher purity than ultracentrifugation.
- Enables the elution of intact EVs, and can be used for various applications.
- Easy operation with magnetic beads.
- No ultracentrifugation is required.

[Kit contents]**This kit includes 6 components.****• 2 purifications**

- | | |
|------------------------------------|-----------------------------|
| (1) Biotin Capture Magnetic Beads | 120 μ L \times 1 tube |
| (2) Biotin-labeled Exosome Capture | 20 μ L \times 1 tube |

— 1/24 —

- | | |
|--|-----------------------------|
| (3) Exosome Immobilizing/Washing Buffer (10 \times) | 5 mL \times 1 bottle |
| (4) Exosome Binding Enhancer (500 \times) | 300 μ L \times 1 tube |
| (5) Exosome Elution Buffer (10 \times) | 300 μ L \times 1 tube |
| (6) Reaction Tubes | 4 tubes |
- 10 purifications**
- | | |
|--|------------------------------|
| (1) Biotin Capture Magnetic Beads | 600 μ L \times 1 tube |
| (2) Biotin-labeled Exosome Capture | 100 μ L \times 1 tube |
| (3) Exosome Immobilizing/Washing Buffer (10 \times) | 25 mL \times 1 bottle |
| (4) Exosome Binding Enhancer (500 \times) | 1500 μ L \times 1 tube |
| (5) Exosome Elution Buffer (10 \times) | 1500 μ L \times 1 tube |
| (6) Reaction Tubes | 22 tubes |

[Storage]

Store at 2-10°C

[Additional required materials]**Reagents :**

- 1) Distilled Water
- 2) TBS (as necessary)

Equipment :

- 1) Microcentrifuge (Max > 10,000 \times g)
- 2) Vortex mixer
- 3) Tabletop centrifuge
- 4) Magnetic stand
- 5) Rotator or tube mixer
- 6) Micro pipette
- 7) Pipette tip
- 8) Microcentrifuge tube (1.5 mL)
- 9) Centrifuge tube (15 mL) (as necessary)
- 10) Centrifuge tube (50 mL) (as necessary)
- 11) Ultrafiltration unit (Vivaspin6, M.W. cut off : 100 K, Code No. VS0641 or Vivaspin20, M.W. cut off : 100 K, Code No. VS2041) (as necessary)

[Precaution for use]**1. Equipment**

Please use sterile or DNase and RNase free microcentrifuge tubes and pipette tips. We recommend the use of gloves and a mask to avoid contamination of DNase and RNase.

2. Reagents


If you analyze the nucleic acid (RNA, DNA) after isolation of exosomes and other EVs, please handle reagents carefully and keep them as sterile as possible. Additionally, if you analyze the mass spectrometry of proteins after isolation of exosomes and other EVs, please be careful of contamination from other proteins.

3. Handling of biohazardous waste

Please treat human serum, human plasma, and human tissue samples as infectious samples. Biohazardous wastes, which include waste liquid and equipment such as microcentrifuge tubes, pipette tips, and gloves, must be disposed of according to the guidelines of the institution.

— 2/24 —

【Outline of procedure】

【Procedure】

1. Preparation of Sample

This is the section to prepare samples. When exosomes and other large EVs are needed, prepare **1,200 × g supernatant** as a sample. Additionally, when highly purified exosomes are needed, prepare **10,000 × g supernatant** as a sample^{*1}. This sample preparation protocol describes the protocol when using cell culture supernatant, serum/plasma sample and urine sample. When using other body fluid samples, refer to the serum/plasma sample protocol and consider the appropriate pretreatment protocol.

※1 To remove Large EVs such as apoptotic bodies and microvesicles from a sample, use a centrifugal filter unit (Millipore, Ultrafree-MC, GV 0.22 μm sterile, Code No : UFC30GV0S) after centrifuging the supernatant at 10,000 × g for 30 minutes. The filtrate can be used as a test sample.

In the case of cell culture medium

- 1) Culture the cells under the appropriate condition^{*2}.
- 2) Harvest cell culture medium.
- 3) Centrifuge the cell culture medium for 5 minutes at 300 × g at 4°C to remove cells and debris.
- 4) Transfer the supernatant from step 3) into a new tube.
- 5) Centrifuge 4) at 1,200 × g for 20 minutes at 4°C to remove the cell debris.
- 6) Transfer the supernatant from step 5) into a new tube. ... **1,200 × g supernatant**
- 7) Centrifuge 6) at 10,000 × g for 30 minutes at 4°C to remove the large EVs.
- 8) Transfer the supernatant from step 7) into a new tube. ... **10,000 × g supernatant**

(Optional : Concentration of cell culture supernatant by ultrafiltration)

If the volume of **1,200 × g supernatant** or **10,000 × g supernatant** is over 1 mL, concentrating the sample to 1 mL or less is available with ultrafiltration unit. The maximum sample volume

as ultrafiltration is 50 mL (50-fold concentration). However, although EVs can be obtained efficiently by concentrating the sample, it may cause physical damage or inactivation of biological activity. Please consider the appropriate processing method according to the purpose of use.

Precaution for concentration by ultrafiltration

Since EVs may be adsorbed to the centrifugal ultrafiltration unit, and the recovery amount may be reduced, we recommend adding EV-Save™ Extracellular Vesicle Blocking Reagent (Code No : 058-09261) (sold separately) to the sample for preventing loss of EVs. However, this product contains a polymer component, please be careful when using it as a sample for proteomics analysis^{*3}.

※2 Please examine appropriate culture conditions depending on cell lines. Since the amount of EVs in the medium is small, please increase the culture scale as much as possible. (Example : 20 mL or more is better) Besides, to ensure that isolated EVs originate from your cells of interest, please use EV-depleted FBS prepared by ultracentrifugation (Example : at 110,000 × g for 18 hours) or commercialized products.

※3 When performing proteomics analysis, remove the polymer by appropriate preprocessing such as acetone precipitation.

In the case of Serum and Heparin Plasma, Saliva

Please centrifuge serum and heparin plasma used in this section at 1,200 × g beforehand or prepare the processed sample.

When exosomes and Large EVs are needed, use 1,200 × g sup. fraction as a sample for purification. If a floating matter is observed in 1,200 × g centrifuged samples, centrifuge again at 1,200 × g and use as a sample for exosomes and Large EVs purification.

1) Centrifuge a sample at 10,000 × g for 30 minutes at 4°C to remove the large EVs^{*4}.

2) Transfer the supernatant from step 1) into a new tube. ... **(10,000 × g supernatant)**

※4 When Large EVs are needed, use ppt of Large EVs obtained by centrifugation at 10,000 × g as a sample after suspending it with 500 μL to 1 mL of TBS.

In the case of EDTA Plasma and Citrated Plasma

Please centrifuge EDTA plasma and citrated plasma used in this section at 1,200 × g beforehand or prepare the processed sample.

In the case of EDTA and citrated plasma, the anticoagulant contained in the sample inhibits the binding of EVs to Exosome Capture Beads. Therefore, adjust the heparin sodium solution just before the affinity reaction, add it to the centrifuged sample, and proceed to the purification step. Details are described in the section on **EDTA Plasma and Citrated Plasma** in 4. Affinity Reaction.

1) Centrifuge a sample at 10,000 × g for 30 minutes at 4°C to remove the large EVs^{*4}.

2) Transfer the supernatant from step 1) into a new tube. ... **(10,000 × g supernatant)**

In the case of Urine


Please add Exosome Binding Enhancer (500×) to urine used in this section before centrifugation.

Addition of Exosome Binding Enhancer (500×) after centrifugation may result in the formation of a white precipitate. This precipitate may inhibit the elution of extracellular vesicles from Biotin Capture Magnetic Beads. Therefore, after adding Exosome Binding Enhancer (500×), perform centrifugation and proceed to the affinity reaction.


- 1) Add 1/500 volume of Exosome Binding Enhancer (500×) to the sample. Mix by vortexing.
- 2) Centrifuge 1) at 1,200×g for 20 minutes at 4°C.
- 3) Transfer the supernatant from step 2) into a new tube....
1,200×g supernatant of a sample including Exosome Binding Enhancer (500×)
- 4) Centrifuge 3) at 10,000×g for 30 minutes at 4°C.
- 5) Transfer the supernatant from step 4) into a new tube....
10,000×g supernatant of A sample including Exosome Binding Enhancer (500×)
- 6) Proceed **4. Affinity Reaction “In the case of cell culture medium, serum, heparin plasma and urine”**.

[Flowchart of Sample Preparation (Section 1)]


Cell Culture Medium

Serum, Heparin Plasma, Urine, and Saliva

EDTA Plasma and Citrated Plasma

- *1 When **exosomes and large EVs** are needed, use 1,200 × g sup. fraction as sample.
- *2 When **exosomes** are needed, use 10,000 × g sup. fraction as sample.
- *3 When **Large EVs** are needed, use ppt of Large EVs obtained by centrifugation at 10,000 × g as sample after suspending it with TBS.
- *4 The concentration step is an option when using a large volume (~ 50 mL) of cell culture supernatant as a sample for purification.

2. Buffer Preparation

This section is the process of preparing the buffers used in each process. The volume of various buffers required for one sample is 5.0 mL (preparation volume : 5.5 mL) for Exosome Immobilizing/Washing Buffer (1×) and 100 μL (preparation volume : 150 μL) for Exosome Elution Buffer (1×). When performing the affinity reaction overnight, prepare the 3 mL of Exosome Immobilizing/Washing Buffer (1×) used for washing just before use.

Buffer	Required amount	Preparation volume
Exosome Immobilizing/Washing Buffer (1×)	5.0 mL 2.0 mL → for the section 3 3.0 mL → for the section 5	5.5 mL
Exosome Elution Buffer (1×)	100 μL	150 μL

- 1) Add 0.55 mL of Exosome Immobilizing/Washing Buffer (10×) and 4.95 mL of distilled water to a 50 mL centrifuge tube, and add 11 μL of Exosome Binding Enhancer (500×) to prepare the Exosome Immobilizing/Washing buffer (1×)^{*5}.
 - 2) Add 15 μL of Exosome Elution Buffer (10×) and 135 μL of distilled water to a 1.5 mL microcentrifuge tube (not included in the kit) to prepare the Exosome Elution buffer (1×)^{*6, *7}.
 - 3) Mix both 1) and 2) well and store at room temperature until use.
- *5 Make sure to add 1/500 volume of Exosome Binding Enhancer (500×) to Exosome Immobilizing/Washing Buffer (1×). When not added, the recovery efficiency of exosome is significantly reduced.
 - *6 By adding EV-Save™ Extracellular Vesicle Blocking Reagent (Code No : 058-09261) (sold separately) to the Exosome Elution Buffer used in this process, the adsorption of EVs to the tube can be suppressed and the recovery efficiency can be increased. However, since this product contains a polymer component, please carefully consider its use for proteomics analysis before use.
 - *7 When using plasma (heparin, EDTA, and citric acid) as a sample, recovery efficiency may be improved by using 2× Exosome Elution Buffer. In that case, mix 30 μL of Exosome Elution Buffer (10×) and 120 μL of distilled water.

3. Immobilization of Exosome Capture


This is the section to immobilize Biotin-labeled Exosome Capture to Biotin Capture Magnetic Beads. Make sure to use the 1.5 mL Reaction Tube (included) in the downstream experiment using magnetic beads^{*8}.

- 1) Transfer 60 μ L of Biotin Capture Magnetic Beads, which has been well stirred with a vortex mixer, into a new 1.5 mL Reaction Tube.
- 2) Add 500 μ L of Exosome Capture Immobilizing/Washing Buffer (1 \times) into the 1.5 mL Reaction Tube from step 1) and suspend it by vortexing.
- 3) Spin down 2), then places it on the magnetic stand for 1 minute to separate the beads from the solution, and remove the supernatant.
- 4) Add 500 μ L of Exosome Capture Immobilizing/Washing Buffer (1 \times) and 10 μ L of Biotin-labeled Exosome Capture into the 1.5 mL Reaction Tube from step 3) and mix by vortexing^{*9}.
- 5) Mix for 10 minutes at 2-10 $^{\circ}$ C with rotator or tube mixer.
- 6) Spin down 5), then places it on the magnetic stand for 1 minute to separate the beads from the solution, and remove the supernatant.
- 7) Add 500 μ L of Exosome Capture Immobilizing/Washing Buffer (1 \times) into the 1.5 mL Reaction Tube from step 6) and mix by vortexing.
- 8) Spin down 7), then places it on the magnetic stand for 1 minute to separate the beads from the solution, and remove the supernatant.
- 9) Repeat steps7)-8) again. **Exosome Capture-immobilized beads**

*8 The 1.5 mL Reaction Tube included in the kit can suppress the decrease in recovery rate due to bead adsorption because it uses a material that magnetic beads do not easily adsorb to the tube.

*9 Depending on the sample, recovery efficiency may be improved by reducing the amount of Biotin-labeled Exosome Capture. When recovery efficiency was low using the standard protocol, reduce the amount of Biotin-labeled Exosome Capture (Example : 10 μ L \rightarrow 2-5 μ L).

[Flowchart of immobilization of Exosome Capture (Section 3)]

4. Affinity Reaction

This is the affinity reaction section to react the Exosome Capture-immobilized beads with exosomes and other EVs in a sample. The basic protocol use the 1.5 mL Reaction Tube (max volume : 1 mL). When exosomes and other large EVs are needed, use **1,200 \times g supernatant** as a sample. Additionally, when highly purified exosomes are needed, use **10,000 \times g supernatant** as a sample.

Note : The affinity reaction protocol for EDTA plasma and citrated plasma are different.

In the case of cell culture medium, serum, heparin plasma, urine, and saliva

- 1) Transfer a sample (max volume : 1 mL)^{*10} to a new 1.5 mL microcentrifuge tube (not included in the kit) and add 1/500 volume of Exosome Binding Enhancer (500 \times) to the sample. Mix by vortexing. **A sample including Exosome Binding Enhancer (500 \times)**
- 2) Spin down 1), then transfer it into the 1.5 mL Reaction Tube containing **Exosome Capture-immobilized beads**. Mix by vortexing.
- 3) Mix for 1 hours \sim at RT or 4 $^{\circ}$ C with rotator^{*11}.
- 4) Spin down 3), then places it on the magnetic stand for 1 minute to separate the beads from the solution, and remove the supernatant^{*12}. **EVs-binding beads**

In the case of EDTA/citrated plasma

- 1) Dissolve 10,000 U of sodium heparin (Code No : 085-00134) with distilled water to make 1,000 U/mL heparin sodium solution^{*13}.
- 2) Transfer a sample (max volume : 1 mL)^{*10} to a new 1.5 mL microcentrifuge tube (not included in the kit) and add 1/200 volume of heparin sodium solution to the sample. (final conc. : 5 U/mL)
- 3) Add 1/100 volume of Exosome Binding Enhancer (500 \times) to 2) and mix by vortexing. **A sample including Exosome Binding Enhancer (500 \times)**
- 4) Spin down 3), then transfer it into the 1.5 mL Reaction Tube containing **Exosome Capture-immobilized beads**. Mix by vortexing.
- 5) Mix for 1 hours \sim at RT or 4 $^{\circ}$ C with rotator^{*11}.
- 6) Spin down 5), then places it on the magnetic stand for 1 minute to separate the beads from the solution, and remove the supernatant^{*12}. **EVs-binding beads**

*10 When the sample volume is less than 0.5 mL, add TBS to scale up the sample volume to 0.5 mL or more to obtain the better mixture of the Exosome Capture-immobilized beads with a sample. (Example : 100-200 μ L \rightarrow 500 μ L)

*11 Appropriate reaction time varies depending on the sample and sample volume. For obtaining a sufficient recovery amount, extend the reaction time as appropriate. Regarding the temperature, select either RT or 4 $^{\circ}$ C depending on the purpose of use of the purified exosome.

*12 When further recovery of remaining EVs is needed, please store the removed supernatant in another tube.

*13 The maximum usage per kit (for 10 purification) is 100 μ L. After preparation, please store at 4 $^{\circ}$ C.

5. Washing of EVs-binding beads

This is the washing section of EVs-binding beads.

Note : Make sure to add Exosome Binding Enhancer (500×) to Exosome Capture Immobilizing/Washing Buffer (1×).

- 1) Add 1 mL of Exosome Capture Immobilizing/Washing Buffer (1×) into the 1.5 mL Reaction Tube containing **EVs-binding beads** and mix by vortexing.

In the case of EDTA/citrated plasma

After adding 1/200 volume of heparin sodium solution to 1 mL of Exosome Capture Immobilizing/Washing Buffer (1×) only at the first washing (final conc. : 5 U/mL), suspend with a vortex mixer.

- 2) Spin down 1), then places it on the magnetic stand for 1 minute to separate the beads from the solution, and remove the supernatant.
- 3) Add 1 mL of Exosome Capture Immobilizing/Washing Buffer (1×) and mix by vortexing.
- 4) Spin down 3), then places it on the magnetic stand for 1 minute to separate the beads from the solution, and remove the supernatant.
- 5) Repeat steps 3)-4) one more time.
- 6) Spin down again the 1.5 mL Reaction Tube from step 6), then places it on the magnetic stand for 1 minute to separate the beads from the solution, and remove the supernatant completely^{*14}. **Washed EVs-binding beads**

※14 Please remove Exosome Capture Immobilizing/Washing Buffer (1×) in the tube completely, because the remaining buffer causes a reduction of elution efficiency of captured exosomes and other EVs.

6. Elution of Extracellular Vesicles

This section is the elution process of EVs from Washed EVs-binding beads.

- 1) Add 50 μL of Exosome Elution Buffer (1×) into the 1.5 mL Reaction Tube containing **Washed EVs-binding beads** and mix by vortexing, then spin down^{*15, *16, *17}.
- 2) Places 1) on the magnetic stand for 1 minute to separate the beads from the solution, and transfer the supernatant into a new 1.5 mL microcentrifuge tube (not included in the kit).
- 3) Add 50 μL of Exosome Elution Buffer (1×) into the 1.5 mL Reaction Tube from step 2) again, and mix by vortexing, then spin down^{*15, *16, *17}.
- 4) Places 3) on the magnetic stand for 1 minute to separate the beads from the solution, and transfer the supernatant into the 1.5 mL microcentrifuge tube from step 2) to pool together^{*18, *19}. (Total : 100 μL)

※15 If the magnetic beads aggregate and can not be uniformly suspended with only a vortex mixer, suspend the beads evenly while applying tapping or pipetting operation. Please avoid excessive suspension using a tube mixer.


※16 By adding 25 μL of Exosome Elution Buffer (1×), you can obtain a more concentrated sample.

※17 When using plasma (heparin, EDTA, and citric acid) as a sample, recovery efficiency may be improved by using 2 × Exosome Elution Buffer.

※18 Exosome Capture-immobilized beads are reusable after eluting EVs (up to 4 times). Therefore, when you need to recover more remaining EVs from samples, please repeat the section 4-6 to increase the recovered amount. After recovering the eluate, please add the stored supernatant of section 4 (refer to ※12) into the 1.5 mL Reaction Tube containing Exosome Capture-immobilized beads from step 4) of section 6 and repeat from the section 4-3) to the section 6-4) as necessary. When using EDTA/Citrated plasma, repeat from the section 4-5) to the section 6-4). This kit contains the required reagents necessary for reuse up to 4 times.

※19 There is a possibility that the pooled eluate have a small amount of magnetic beads. When EVs are analyzed by Nanoparticle Tracking Analysis (NTA) or electron microscopic analysis, filter pooled eluate with 0.45 μm pore size filter unit before analysis as necessary. (Example : Ultrafree MC 0.45 μm, Code No. UFC30HV25, Millipore)

[Flowchart of Affinity Purification (Section4-6)]

[Optional]**1. RNA extraction from Purified EVs**

When extracting RNA from eluted EVs, add 100 μ L of Exosome Elution Buffer (1 \times) to pooled eluate to prepare 200 μ L of sample in total. Then, extract RNA from EVs with microRNA Extractor[®] SP Kit (Code No : 295-71701) in accordance with the instruction manual.

2. Recycling of Exosome Capture-immobilized beads

Used Exosome Capture-immobilized beads can be recycled for "A. Repeated extraction of the remaining EVs from the same sample (refer to ※12)" and "B. Purification of EVs from the same lot of culture supernatant sample and body fluid sample." Recycling is up to 4 times. When purifying EVs from the same lot of sample, prepare samples first and proceed to section 4-1) of the affinity reaction.

When storing used Exosome Capture-immobilized beads, add 1 \times TBS containing 0.05 w/v% of sodium azide to the Reaction Tube containing the used Exosome Capture-immobilized beads. After suspending it, store at 4°C. We recommend using Exosome Capture immobilized-beads as soon as possible.

[Related Products]

Code No.	Description	Size
297-79201	PS Capture [™] Exosome ELISA Kit (Anti Mouse IgG POD)	96 reactions
298-80601	PS Capture [™] Exosome ELISA Kit (Streptavidin HRP)	96 reactions
299-77603	MagCapture [™] Exosome Isolation Kit	2 purifications
293-77601	PS	10 purifications
290-80301	PS Capture [™] Exosome Isolation Resin Kit	1 Kit (0.5 mL Slurry)
297-79701	PS Capture [™] Exosome Flow Cytometry Kit	300 reactions
290-83601	CD63-Capture Human Exosome ELISA Kit (Streptavidin HRP)	96 reactions
296-83701	CD9-Capture Human Exosome ELISA Kit (Streptavidin HRP)	96 reactions
292-83801	CD81-Capture Human Exosome ELISA Kit (Streptavidin HRP)	96 reactions
014-27763	Anti CD9, Monoclonal Antibody (1K)	100 μ L
019-28173	Anti CD9, Rat Monoclonal Antibody (77B)	100 μ L
017-28211	Anti CD9, Rat Monoclonal Antibody (77B), Biotin-conjugated	50 μ L
019-27953	Anti CD9, Rat Monoclonal Antibody (30B), Biotin Conjugated	100 μ L
012-27063	Anti CD63, Monoclonal Antibody (3-13)	100 μ L
014-27643	Anti CD63, Monoclonal Antibody (3-13), Fluorescein Conjugated	100 tests

Code No.	Description	Size
017-27753	Anti CD63, Monoclonal Antibody (3-13), Red Fluorochrome (635) Conjugated	100 tests
019-27713	Anti CD63, Monoclonal Antibody (3-13), Biotin Conjugated	100 μ L
011-27773	Anti CD81, Monoclonal Antibody (17B1)	100 μ L
010-28223	Anti CD81, Rat Monoclonal Antibody (9B)	100 μ L
011-28111	Anti CD81, Rat Monoclonal Antibody (9B), Biotin-conjugated	50 μ L
052-09301	Exosomes, from COLO201 cells, purified	50 μ L
058-09261	EV-Save [™] Extracellular Vesicle Blocking Reagent	1 mL
290-35591	Magnet Stand	1 each
295-71701	microRNA Extractor [®] SP Kit	50 reactions
085-00134	Heparin Sodium	10,000 U
317-90175	10 \times TBS Buffer (pH 7.4)	500 mL

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
 Telephone : + 81-6-6203-3741
 Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
 Richmond, VA 23237
 U.S.A.
 Telephone : + 1-804-271-7677
 Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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