**Anti SARS-CoV-2 S-RBD IgG ELISA Kit Wako****[Introduction]**

This ELISA kit quantifies IgG antibody against SARS-CoV-2 S protein receptor-binding domain (RBD) in serum or plasma. The concentration of the antibody in samples can be calculated by measuring the absorbance. Assay time is about three hours. This kit contains all the reagents used in the assay. This product is for research use only.

**[Kit Performance]**

Standard curve range	2-250 U/mL
Target	IgG antibody specific for SARS-CoV-2 RBD
Sample	Serum and plasma
Sample volume	5 $\mu$ L
Assay time	About three hours
Detection method	Coloring method

**[Kit Content]**

Component	Condition	Volume or quantity
SARS-CoV-2 S Protein-coated Plate (SARS-CoV-2 S protein RBD recombinant antigen coated)	Ready-to-use	1 plate (96 wells)
Peroxidase (POD)-labeled Antibody Solution (POD-labeled anti human IgG monoclonal antibody (mouse))	Ready-to-use	6 mL $\times$ 1 tube
Concentrated Wash Solution (10 $\times$ )	Use after dilution	100 mL $\times$ 1 tube
Coloring Reagent (TMB Solution)	Ready-to-use	11 mL $\times$ 1 tube
Stop Solution (1 N HCl)	Ready-to-use	11 mL $\times$ 1 tube
Sample Dilution Buffer (0 U/mL)	Ready-to-use	100 mL $\times$ 2 tubes
Anti SARS-CoV-2 IgG Standard-1 (62.5 U/mL)	Ready-to-use	500 $\mu$ L $\times$ 1 tube
Anti SARS-CoV-2 IgG Standard-2 (125 U/mL)	Ready-to-use	500 $\mu$ L $\times$ 1 tube
Anti SARS-CoV-2 IgG Standard-3 (187.5 U/mL)	Ready-to-use	500 $\mu$ L $\times$ 1 tube
Anti SARS-CoV-2 IgG Standard-4 (250 U/mL)	Ready-to-use	500 $\mu$ L $\times$ 1 tube
Plate Seal	Ready-to-use	3 sheets
Instruction Manual		1 copy

**[Measurement Principle]**

SARS-CoV-2 protein RBD recombinant antigen has been immobilized on the plate. Sample is added to each well, and anti SARS-CoV-2 antibody in the sample binds to the antigen on the plate. Then POD-labeled anti human IgG monoclonal antibody (mouse) is added, forming "Antigen-Anti SARS-CoV-2 IgG-POD labeled anti human IgG antibody" complex on the plate. The amount of POD bounded on the plate is proportional to the amount of anti SARS-CoV-2 antibody in the sample. Therefore by measuring of POD activity in the well, the concentration of anti SARS-CoV-2 antibody in the sample can be quantified.

**[Equipment and materials required but not supplied]**

- Purified water (distilled water)
- Tubes for diluting samples
- Glass utensils for diluting Wash Solution (graduated cylinder and beaker)
- Vortex mixer
- Aspirator
- Pipettes with disposable tips
- Multichannel pipette
- Paper towels (for removal of liquid)
- Microplate shaker
- 96-well plate washer (preferable if available)
- 96-well plate reader (Measurement wavelength 450 nm/620 nm)
- Data processing software

**[Reagent Preparation]**

1. Wash Solution  
Dilute Concentrated Wash Solution ( $\times$ 10) 10 times with purified water.
2. Standard  
Use Sample Dilution Buffer and Anti SARS-CoV-2 IgG Standard (1-4) as they are.
3. Sample Dilution  
Dilute 5  $\mu$ L of sample with 1,000  $\mu$ L of Sample Dilution Buffer.

\* Use other reagents as they are.

**[Stability and Storage of Reagents]**

Store the kit at 2-10°C. The kit is stable until the expiration date under the condition (The expiration date is printed on the label). The kit should not be used beyond the expiration date. Please use opened reagents as soon as possible since they may be affected by storage conditions.


**[Assay Procedure]**

Prepare all the reagents as described in [Reagent Preparation] before starting the assay.

- |   |   |
|---|---|
| Reference Standard/Sample<br>50 $\mu$ L     | (1) Add 50 $\mu$ L of diluted samples or standard to each well of the antigen-coated plate (standard are assayed in duplicate). |
| ↓ RT, 1 hour                                |   |
| Wash<br>300 $\mu$ L $\times$ 3 times        | (2) Cover with a plate seal and shake it on a microplate shaker. Incubate an hour at room temperature.                          |
| ↓   |   |
| POD-labeled Antibody Solution<br>50 $\mu$ L | (3) Discard the samples and standard using an aspirator.  |
| ↓ RT, 1 hour                                |   |
| Wash<br>300 $\mu$ L $\times$ 5 times        | (4) Add 300 $\mu$ L of Wash Buffer and discard it.  |
| ↓   |   |
| Coloring (TMB solution)<br>100 $\mu$ L      | (5) Repeat wash as in (4) 3 times.  |
| ↓ RT, 10 min                                | (6) Invert the plate and gently blot it against clean paper towels.   |
| Stop Solution (1 N HCl)<br>100 $\mu$ L      | (7) Add 50 $\mu$ L of POD-labeled Antibody Solution to each well.   |
| ↓   | (8) Cover with a new plate seal and shake it on a microplate shaker. Incubate for an hour at room temperature.                  |
| Absorbance measurement<br>450 nm/620 nm     | (9) Aspirate each well.   |
- (10) Add 300  $\mu$ L of Wash Buffer and discard it.  
 (11) Repeat wash as in (10) 5 times.  
 (12) Add 100  $\mu$ L of TMB Solution to each well.  
 (13) Cover with a new plate seal and shake it on a microplate shaker. Incubate for 10 minutes at room temperature (the color of the wells turn blue).  
 (14) Add 100  $\mu$ L of Stop Solution (1 N HCl) to each well (the color of the wells turn yellow).  
 (15) Within 30 minutes after the reaction stops, measure the absorbance (450 nm/620 nm) using 96-well plate reader.

#### [Calculation of Results]

- After the measurement according to [Assay Procedure], prepare a standard curve by plotting the absorbance of each standard on Y-axis against the concentration on X-axis.
- Determine the anti SARS-CoV-2 antibody concentration of the samples from the standard curve and multiply the value by the dilution factor.
- Human IgG recombinant antibody against SARS-CoV-2 S protein RBD is used as standard.

#### [Note]

In case non-specific reactants (e.g. heterophile antibody) exist in samples, correct measurement results may not be obtained.

#### [Precautions]

- (1) If any reagents touch eyes, skin or mucous membranes, wash eyes/skin thoroughly with copious amounts of water and contact physician, when necessary.
- (2) Handle the samples carefully with the attention of virus infection.
- (3) When the samples or waste liquid scatter, wipe them off using disinfectant such as sodium hypochlorite (effective chlorine concentration : 1,000 ppm) or ethanol.
- (4) Stop Solution (1 N HCl) contains hydrochloric acid. Handle it carefully not to spill.
- (5) Do not use the reagents if accidentally frozen.
- (6) Avoid using products from different lot together.
- (7) Use the containers and accessories in this product only for the designated purpose.
- (8) ELISA can be affected by the laboratory environment. Avoid the environment where there is current of air including the current of air from air conditioner and low humidity.
- (9) Dispose properly according to the rules and regulations of your site.
- (10) Treat reagents and containers that had contact with the samples as the one having infectious risk. Incinerate them or dispose of them according to the regulations.

〈Example of disposal method〉

Autoclave processing :

Sterilize with high pressure steam for 20 minutes or more at more than 121°C.

Hypochlorous acid processing :

Soak the wastes in aqueous solution that contains 1,500-2,000 ppm or more of free chlorine for more than 60 minutes.

#### [Package]

Code No.	Size
290-84201	96 reactions

#### [References]

- (1) The Japanese Association for Infectious Diseases Significant Scientific Evidences about COVID-19 (24/07/2020) ([https://www.kansensho.or.jp/uploads/files/topics/2019ncov/covid19\\_sse\\_0724.pdf](https://www.kansensho.or.jp/uploads/files/topics/2019ncov/covid19_sse_0724.pdf))
- (2) The Japanese Association for Infectious Diseases 新型コロナウイルス感染症に対する検査の考え方－遺伝子診断、抗体、抗原検査の特長と使い分け－ ([https://www.kansensho.or.jp/uploads/files/topics/2019ncov/covid19\\_kensaguide\\_0526.pdf](https://www.kansensho.or.jp/uploads/files/topics/2019ncov/covid19_kensaguide_0526.pdf))

### FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
 Telephone : + 81-6-6203-3741  
 Facsimile : + 81-6-6201-5964  
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road  
 Richmond, VA 23237  
 U.S.A.  
 Telephone : + 1-804-271-7677  
 Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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