Oxytocin ELISA Kit Wako

[1. Introduction]

This item is an ELISA kit designed for the quantitative detection and measurement of oxytocin in biological fluids such as saliva, urine, serum and plasma. The kit can be measured by a simple pretreatment, and does not require purification using a C18 column and organic solvent, which was conventionally performed as a pretreatment method. For Research Use Only. Not for use in diagnostic procedures.

Oxytocin is called as a love or happy hormone consisting of 9 amino acids. It also has stress relief, learning ability improvement, anxiolytic, and anti-fear effects, and is expected to have therapeutic effects on mental illness (autism, depression, etc.).

[2. Performance]

Standard curve range	4.00-1,024 pg/mL
Measuring object	Oxytocin
Sample	Human saliva, urine, serum, plasma Mouse serum, plasma Rat serum, plasma
Sample volume	50 μ L (n = 1, minimum volume) 200 μ L (n = 2, recommended volume)
Assay time	Approx. 2.5 hours
Detection method	Luminescent detection *

* Requires luminescence plate reader for measurement.

[3. Material supplied]

Components	State	Volume or quantity
Antibody-coated Plate	Use after washing	1 plate 96 wells (8 \times 12)
Oxytocin Standard	Use after dilution	100 μ L/ 1 bottle
Buffer	Ready-to-use	60 mL/ 1 bottle
Biotin-conjugated Oxytocin	Use after reconstitution and dilution	1 bottle
Peroxidase-conjugated Streptavidin Solution	Concentrated Use after dilution	100 μ L/ 1 bottle
Luminescent Reagent 1	Ready-to-use	6 mL/1 bottle
Luminescent Reagent 2	Ready-to-use	6 mL/1 bottle
Wash Solution (10 \times)	Concentrated Use after dilution	100 mL/ 1 bottle
Sample Buffer 1	Ready-to-use	5 mL/1 bottle

Sample Buffer 2	Ready-to-use	12 mL/ 1 bottle
Plate Seal	—	4 sheets
Instruction manual	—	1 book

[4. Measurement principle]

The microplate is coated with anti-oxytocin antibody. In each well, the standard solution or sample and biotin-conjugated oxytocin are incubated to proceed with the antigen-antibody reaction. Furthermore, peroxidase-conjugated streptavidin is added to proceed with the biotin-streptavidin binding. Finally, peroxidase activity in each well is measured to determine oxytocin in the sample.

[5. Equipment or materials required but not supplied]

- ☐ Purified water (distilled water)
- ☐ Test tube for dilution of a standard solution and samples
- ☐ Glass utensils for dilution of Wash Buffer (e.g. Graduated cylinders and beakers)
- ☐ Pipettes with disposable tips (one capable of pipetting 10 μ L of liquid accurately and one capable of pipetting 200 to 500 μ L)
- ☐ Dispenser capable of dispensing
- ☐ Water-absorbable material such as paper towel (to remove liquid remaining on a plate after washing)
- ☐ Mixer (Vortex type)
- ☐ Shaker for 96-well plate (approximately 600 to 800 rpm)
- ☐ Automatic washer for 96-well plate (If available) or washing bottle
- ☐ Microplate reader for luminescence detection
- ☐ Software for data analysis

[6. Reagent preparation]

Bring the reagents to room temperature about 2 hours prior to use.

6.1. Standard Solution

Standard solutions are prepared using Standard solution (102.4 ng/mL) and Buffer as shown below.

Concentration of the standard solution (pg/mL)	Volume of the standard solution	Volume of Buffer
1024	Standard solution (102.4 ng/mL) : 30 μ L	2970 μ L
512	Standard solution (102.4 pg/mL) : 150 μ L	150 μ L
256	Standard solution (512 pg/mL) : 150 μ L	150 μ L
64.0	Standard solution (256 pg/mL) : 50 μ L	150 μ L
16.0	Standard solution (64 pg/mL) : 50 μ L	150 μ L
4.00	Standard solution (16 pg/mL) : 50 μ L	150 μ L
0 (B0)	—	150 μ L

6.2. Biotin-conjugated Oxytocin

Reconstitute Biotin-conjugated Oxytocin with 100 μ L of distilled water and completely dissolved. Dilute the solution 100-fold with Buffer.

6.3. Peroxidase-conjugated Streptavidin Solution

Peroxidase-conjugated Streptavidin Solution is diluted 100-fold with Buffer.

6.4. Luminescent Reagent 1/2

Mix Luminescent Reagent 1 and Luminescent Reagent 2 in 1:1 (vol./vol.) and incubate for 30 minutes before use. Keep in the dark.

6.5 Wash Solution (10×)

Wash Solution (10×) is diluted 10-fold with purified water (distilled water).

(e.g.) When solution for 96-well reaction is prepared.

Add 100 mL of Wash Solution (10×) to 900 mL of purified water (distilled water) to prepare 1,000 mL of Wash Solution (1×).

Other reagents are ready-to-use.

7. Stability and storage method of each reagent

7.1. Antibody-coated Plate

Unused antibody-immobilized strips (kept refrigerated and sealed) should be returned into the zip-seal bag provided in the kit and stored at 2-10°C. These are stable until the expiration date.

7.2. Oxytocin Standard

Store at 2-10°C. It is stable until the expiration date. Discard the remaining diluted standard solution after use.

7.3. Buffer

If a part of Buffer is used, transfer the volume slightly greater than needed to another container, immediately close cap lid tightly without bringing the remaining Buffer to room temperature, and store at 2-10°C. It is stable until the expiration date.

7.4. Biotin-conjugated Oxytocin

Store at 2-10°C. It is stable until the expiration date.

Reconstituted standard solution should be stored at 2-10°C and should be used up within 2 weeks. Discard the remaining diluted solution after use.

7.5. Peroxidase-conjugated Streptavidin Solution

If the kit is divided for multiple assays, prepare these solutions by dilution just after taking out of the refrigerator, immediately close the cap tightly without bringing the remained stock solution to room temperature, and store at 2-10°C. It is stable until the expiration date. Discard the remaining diluted solution after use.

7.6. Luminescent Reagent 1, Luminescent Reagent 2

If the kit is divided for multiple assays, prepare these solutions just after taking out of the refrigerator, immediately close the cap tightly, and store at 2-10°C. It is stable until the expiration date. Discard the remaining mixed Luminescent Reagent.

7.7. Wash Buffer (10×)

Store the Wash Buffer (10×) with the cap tightly closed at 2-10°C if applicable.

It is stable until the expiration date. Discard the remaining diluted washing solution.

7.8. Sample Buffer 1, Sample Buffer 2

Store with the cap tightly closed at 2-10°C. It is stable until the expiration date.

8. Sample preparation

8-1. Add Sample Buffer 1 to sample in the table below. Shake it and incubate 5 minutes. Shake it again and incubate 5 minutes.

8-2. Add Sample Buffer 2 which gel is mixed evenly to sample in the table below. Shake it and incubate 5 minutes. Shake it again and incubate 5 minutes.

8-3. After shake it, centrifuge (5000-6000 g, 10 minutes, 4°C) and use the supernatant as a sample for measurement. Don't suck the gel in the buffer when collecting it.

	Urine			Saliva/Serum/plasma		
Sample volume	50 μ L	100 μ L	200 μ L	50 μ L	100 μ L	200 μ L
Sample Buffer 1	5 μ L	10 μ L	20 μ L	5 μ L	10 μ L	20 μ L
Sample Buffer 2	27.5 μ L	55 μ L	110 μ L	27.5 μ L	75 μ L	150 μ L
Dilution factor	1.375			1.475		

(Precaution for sample)

- Add 100 KIU/mL aprotinin and 1/750 amount (v/v) Proclin950 to sample at the time of the sample collection.
- It is recommended to store the sample with 100 KIU/mL aprotinin and 1/750 amount (v/v) Proclin950 at -80°C or lower for long-term storage. Avoid repeated freezing and thawing. Thaw the frozen sample just before assay and thoroughly agitate. Prepare the sample before use.
- If presence of interfering substances is suspected for the sample, dilute the concerned one at multiple dilution ratios to check the dilution linearity. Dilute the prepared sample with the Buffer.
- Samples with turbidity and insoluble matter should be used for measurement after removal by centrifugation.
- When diluting a sample, dilute it with physiological saline or Buffer in this kit in advance using a test tube (PP, PE), and dispense it into the measurement well after the preparation process.
- For assay values below the detection limit, try concentrating the sample using the following method.
Note, however, that the amount of sample required will be greater than that used for a normal assay.
For human saliva, human serum and human urine samples, we confirmed concentrating samples up to 3-fold.

[Reference]

Sample Concentration Method

Example of concentration protocol (3-fold concentration, duplicate assay)

- (1) Collect samples of 300 μ L (3 times the amount of sample required for a normal duplicate assay).
- (2) Add 30 μ L of Sample Buffer 1. Stir, then allow to stand for 5 minutes. Stir again, and allow to stand for further 5 minutes.

- (3) Add 225 μL (for saliva/serum samples) or 165 μL (for urine samples) of Sample Buffer 2. Stir, then allow to stand for 5 minutes. Stir again, and allow to stand for further 5 minutes.
- (4) Stir, then centrifuge the sample (5000-6000 g, 10 min, 4°C) and separate 360 μL of the supernatant.
- (5) Concentrate and dry the separated supernatant in accordance with the following conditions.
Equipment : miVac Duo LV (Cat. No. DPP-10000-G00, SP Scientific)
Concentration conditions : 30°C, 3 hours
Note : Store dried samples at 4°C.
- (6) Add 120 μL of the Buffer supplied with the kit to the dried sample (3-fold concentrated by volume).

[9. Assay Procedure]

1. Wash 4 times the plate with Wash Buffer (1×). Invert the plate and gently blot it against clean paper towels.
2. Add 50 μL of biotin-conjugated oxytocin solution to each well.
3. Agitate the plate on a microplate shaker.
4. Add 50 μL of standard solution to each well.
5. Add 50 μL of samples after preparation (shown as 8. Sample preparation) to each well.
6. Agitate the plate on a microplate shaker.
7. Cover with a Plate Seal. Incubate for 2 hour at room temperature (20-25°C).
8. Discard the solution and wash 4 times with Wash Buffer (1×). Invert the plate and gently blot it against clean paper towels.
9. Add 100 μL of diluted peroxidase-conjugated streptavidin solution to each well.
10. Agitate the plate on a microplate shaker.
11. Cover with a new Plate Seal. Incubate for 30 minutes room temperature (20-25°C).
12. Repeat the wash as in step 8.
13. Add 100 μL of mixed luminescent reagent 1 and 2 (1 : 1) to each well.
14. Shake the plate for 1 minute using the microplate shaker.
15. After agitation, determine the luminescent intensity (RLU) using the 96-well microplate reader (for luminescence measurement). It is recommended to perform the measurement 10 to 20 minutes after addition of the luminescent reagent.

(Calculation)

Create a standard curve by plotting the concentrations (pg/mL) of the standard solutions along the X axis and the RLU (B/B₀ (%)) along the Y axis. Read the concentration (pg/mL) corresponding to the RLU (B/B₀ (%)) of the diluted sample. Multiply the concentration read with the sample dilution factor to obtain a measured value.

*Use of a tertiary polynomial equation with 4 or 5 parameters is recommended for the arithmetic processing using computer software.

*The oxytocin concentration using this kit is determined with reference to the WHO International Standard (Oxytocin 4th

International Standard ; NIBSC code : 76/575).


Examples

Oxytocin (pg/mL)	RLU	B/B ₀ (%)
1024	17269	16.4%
512	35560	33.7%
256	56127	53.2%
64.0	84849	80.4%
16.0	93648	88.8%
4.00	100273	95.1%
0 (B ₀)	105488	100%

	RLU	B/B ₀ (%)	Measured (pg/mL)
Sample 1	53147	50.4%	283
Sample 2	79755	75.6%	85.9

$$B/B_0(\%) = (\text{RLU of standard or sample} / \text{RLU of 0 pg/mL (B}_0\text{)}) \times 100$$

[10. Standard curve]

[11. Precautions]

1. This kit shall be used by an operator who has completed training for ELISA method or under appropriate supervision.
2. This kit shall be used by an operator who has a stably reproducible pipetting skill if the measurement is manually operated.
3. Wear protective gloves, clothing, eye, and face protection.
4. Avoid skin contact with any reagent. In case where any reagent in this kit comes in contact with eyes, mouth, wound, or skin by mistake, provide emergency care, for instance, rinsing immediately with plenty of water, and seek medical advice if necessary.
5. Do not drink, eat, or smoke in a place where this kit is used.
6. Handle the sample with proper care, being aware that the sample may have an infection risk. This kit contains animal-derived ingredients.
7. Soak any used samples or consumables in a solution containing 1% formalin and 2% glutaraldehyde or $\geq 0.1\%$ sodium hypochlorite solution for at least 1 hour, or autoclave them before disposition. Dispose used consumables and unused reagents in accordance with rules of the facility and regional regulations.

8. Do not use reagents with different lot numbers together. When allowing the plate to stand in each step, always affix a plate sealer to protect the wells from drying, contamination with foreign matters, uneven temperature, and evaporation of a dispensed reagent.
9. ELISA can be affected by assay environment. Ensure that room temperature at places for assay operation and incubation is strictly controlled at 20 to 25°C (on a bench or in an incubator). Avoid performing assay under air flow (including air flow from an air-conditioning equipment) or in a low-humidity environment.

[12. Assay Procedure Summary]

Be sure to read the instruction manual to check the sample conditions, measurement conditions, and measurement method before performing the measurement operation.

- ☐ Bring the plate and reagents to room temperature (20-25°C) about 2 hours before use.
- ☐ Dilution of Wash Solution (10×) : Dilute Wash Solution (10×) 10-folds with purified water.
- ☐ Dilution of the standard solution : prepare the standard solutions by mixing with the Buffer as shown in a table below.

Example	Concentration(pg/mL)	1024	512	256	64.0	16.0	4.00	0
	Standard (μL)	Original: 30	150*	150*	50*	50*	50*	-
	Buffer (μL)	2970	150	150	150	150	150	150

*Standard solution at a 1-step higher concentration

- ☐ Preparation of biotin-conjugated oxytocin solution
Reconstitute Biotin-conjugated Oxytocin with 100 μL of distilled water and completely dissolved. Dilute the solution 100-fold with Buffer.
- ☐ Do sample preparation (shown in 8. Sample preparation)
- ☐ **Antibody-coated Plate**
- ☐ ↓ Washing 4 times (*1)
Diluted biotin-conjugated oxytocin solution 50 μL/well
↓ Agitate (*2)
- ☐ **Standard solution or sample after preparation 50 μL/well**
- ☐ ↓ Agitate (*2) and shake at room temperature (20-25°C) for 2 hours (500 rpm)
- ☐ Preparation of peroxidase-conjugated streptavidin solution (dilute 100-folds with Buffer)
- ☐ ↓ Washing 4 times (*1)
- ☐ **Peroxidase-conjugated streptavidin solution 100 μL/well**
- ☐ ↓ Agitate (*2) and incubate at room temperature (20-25°C) for 30 minutes (*3)
- ☐ Preparation of the luminescent reagent (mix Luminescent Reagent 1 and Luminescent Reagent 2 in 1 : 1 (vol./vol.))
- ☐ ↓ Washing 4 times (*1)
- ☐ **Luminescent reagent 100 μL/well**
- ☐ ↓ Agitate at room temperature (20-25°C) for 1 minute.
- ☐ **Measurement of luminescence intensity (to be measured between 10 and 20 minutes)**

- (*1) In each washing operation, dispense the wash solution (1×) into wells, gently agitate the filled plate on the palm for about 10 seconds, and then empty the wells. After washing the wells 4 times consecutively, reverse the plate, and tap it against paper towel to remove the washing solution completely. After removal of the washing solution, immediately dispense the next solution with care not to dry the wells. Use of a pipet set at the liquid volume of 300 μL may be appropriate for dispensing the washing solution into each well.
- (*2) Three repeats of agitation at 600 to 800 rpm for 10 seconds may be appropriate.
- (*3) After agitation, cover with a plate sealer. Remove the liner from the plate sealer, and apply its adhesive side to the plate for affixation. Do not re-use any plate sealer.

[Storage] Store at 2-10°C

[Expiration date] Indicated on the label

[Package] For 96 assays

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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