EV-Save™ Extracellular Vesicle Blocking Reagent for *in vivo*

EV-Save™ Extracellular Vesicle Blocking Reagent for *in vivo* is a polymer reagent to prevent the absorption of extracellular vesicles (EVs) derived from cell culture supernatants and/or body fluids to labware such as tubes and pipette tips. It contributes to reduce the loss due to the absorption of EVs during experiments and storage. The combination of this and MagCapture™ Exosome Isolation Kit PS Ver.2 (Code No. 290-84103) is recommended because of reducing the loss of EVs effectively*1. Furthermore this has protective effect against freeze-thaw of EVs. Adding this can suppress the damage from freeze-thaw to EVs.

*1 It is possible to get the same effect to against EVs using by MagCapture™ Exosome Isolation Kit PS (Code No. 293-77601).

[Procedure]

⟨For isolation of EVs from cell culture supernatants using MagCapture™ Exosome Isolation Kit PS Ver.2 (V2.0 kit)⟩

1. Add 1/100 volume of this product to cell culture supernatants and mix it by inverting or tapping.
2. Proceed to the point just before the step of "6. Elution of Extracellular Vesicles" described in the instruction of V2.0 kit.
3. Add 1/100 volume of this product to "Exosome Elution Buffer" accompanied with V2.0 kit and mix it by inverting or tapping.
4. Elute EVs according to the instruction of V2.0 kit.
5. Store the eluted EVs*2 or proceed to the next experiment.

*2 Although the eluted sample that is added EV-Save™ can be stored at -20°C or less, avoid repeated freeze-thaw cycles.

⟨For isolation of EVs from serum or plasma using MagCapture™ Exosome Isolation Kit PS Ver.2 (V2.0 kit)⟩

1. Proceed to the point just before the step of "6. Elution of Extracellular Vesicles" described in the instruction of V2.0 kit.
2. Add 1/100 volume of this product to "Exosome Elution Buffer" accompanied with V2.0 kit and mix it by inverting or tapping.
3. Elute EVs according to the instruction of V2.0 kit.
4. Store the eluted EVs*2 or proceed to the next experiment.

⟨For isolation of EVs from cell culture supernatants using other methods⟩

1. Add 1/100 volume of this product to cell culture supernatants and mix it by inverting or tapping.
2. After isolation of EVs, add 1/100 volume of this product to isolated EVs and mix it by inverting or tapping.

3. Store the eluted EVs*2 or proceed to the next experiment.

⟨Option : For remove EDTA in EVs that are isolated by MagCapture™ Exosome Isolation Kit PS Ver.2 (V2.0 kit)⟩

1. Elute EVs according to the instruction of V2.0 kit.
2. Add 1/100 volume of this product to equilibration buffer that is described the instruction of PD SpinTrap™ G-25 (cytiva, Code No. 28-9180-04).
3. Carry out the gel filtration according to the instruction of PD SpinTrap™ G-25.
4. Store the obtained EVs*2 or proceed to the next experiment.

[Attention]

- This product can't reduce the loss by ultrafiltration. When the ultrafiltration is necessary, EV-Save™ Extracellular Vesicle Blocking reagent (Code No. 058-09261) is useable.
- When performing protein quantification, a BCA assay is recommended.
- Regarding EV staining with ExoSperkler and PKH67, an increase of background noise have not been confirmed. However, it may increase depending on the sample, we recommend the comparison with the control to accurately determine the uptake of the fluorescent dye.

[Storage]

-20°C

[Package]

1 mL

[Related Products]

Code No.	Product Name	Package
299-77603	MagCapture™ Exosome Isolation Kit PS	2 tests
293-77601		10 tests
294-84101	MagCapture™ Exosome Isolation Kit PS Ver.2	2 tests
290-84103		10tests
290-80301	PS Capture™ Exosome Isolation Resin Kit	1 kit (0.5 mL Slurry)
297-79201	PS Capture™ Exosome ELISA Kit (Anti Mouse IgG POD)	96 tests
298-80601	PS Capture™ Exosome ELISA Kit (Streptavidin HRP)	96 tests
296-83701	CD9-Capture Human Exosome ELISA Kit (Streptavidin HRP)	96 tests
290-83601	CD63-Capture Human Exosome ELISA Kit (Streptavidin HRP)	96 tests
292-83801	CD81-Capture Human Exosome ELISA Kit (Streptavidin HRP)	96 tests
297-79701	PS Capture™ Exosome Flow Cytometry Kit	300 tests
014-27763	Anti CD9, Monoclonal Antibody (1K)	100 μL

Code No.	Product Name	Package
019-28173	Anti CD9, Rat Monoclonal Antibody (77B)	100 μ L
017-28211	Anti CD9, Rat Monoclonal Antibody (77B), Biotin-conjugated	50 μ L
019-27953	AntiCD9, Rat Monoclonal Antibody (30B), Biotin Conjugated	100 μ L
012-27063	Anti CD63, Monoclonal Antibody (3-13)	100 μ L
014-27643	Anti CD63, Monoclonal Antibody (3-13), Fluorescein Conjugated	100 tests
017-27753	Anti CD63, Monoclonal Antibody (3-13), Red Fluorochrome (635) Conjugated	100 tests
019-27713	Anti CD63, Monoclonal Antibody (3-13), Biotin Conjugated	100 μ L
011-27773	Anti CD81, Monoclonal Antibody (17B1)	100 μ L
010-28223	Anti CD81, Rat Monoclonal Antibody (9B)	100 μ L
011-28111	Anti CD81, Rat Monoclonal Antibody (9B), Biotin Conjugated	50 μ L
052-09301	Exosomes, from COLO201 cells, purified	50 μ L
053-09451	EV-Up TM EV Production Basal Medium for MSC, AF	95 mL
298-84001	EV-Up TM MSC EV Production Supplement, AF	for 100 mL
295-71701	microRNA Extractor [®] SP Kit	50 tests
058-09261	EV-Save TM Extracellular Vesicle Blocking Reagent	1 mL

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
 Telephone : + 81-6-6203-3741
 Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
 Richmond, VA 23237
 U.S.A.
 Telephone : + 1-804-271-7677
 Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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