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[LBIS Human Insulin ELISA Kit (Luminescent type)]

Cat # 634-52911 (Manufacture # AKH-INSL)

Please, read this instruction carefully before use.

This kit is manufactured by FUJIFILM Wako Shibayagi Corporation. Use only the current version of Instruction Manual enclosed with the kit! For the detailed assay procedure, refer to [Key points for ELISA by movie] on our website.

1. Intended use

LBIS Human Insulin ELISA Kit (Luminescent type) is a sandwich ELISA system for quantitative measurement of human insulin. This is intended for research use only. Not for use in diagnostics procedure.

2. Storage and expiration

When the complete kit is stored at 2 °C -8 °C (Do not freeze), the kit is stable until the expiration date shown on the label on the container. Opened reagents should be used as soon as possible to avoid less than optimal assay performance caused by storage environment.

3. Introduction

Insulin is a peptide hormone secreted from B cells of islet of Langerhans in the pancreas with a molecular weight of about 5800 and pI 5.4. It is consisted of 2 chains, A and B. It has 3 disulfide bonds formed between A6 and A11, A7 and B7, and A20 and B19. Insulin exists as a dimer molecule in acidic to neutral solution without Zn ion, and as a hexamer including two Zn ions in neutral solution if Zn ions are present. Main targets of insulin are liver, muscle, and adipose tissue. Insulin actions in these targets are as follows. In the liver, it promotes glycogenesis, protein synthesis, fatty acid synthesis, carbohydrate utilization, and inhibition of gluconeogenesis. In the muscle, it promotes membrane permeability for carbohydrates, amino acids and K ion, glycogenesis, protein synthesis, while inhibits protein degradation. In the adipose tissue, it promotes membrane permeability for glucose and fatty acid synthesis.

A precursor of insulin, called proinsulin with a single polypeptide chain, is first synthesized in the cell, then sulfide bonds are formed, and finally by enzymatic cutting at two sites, active insulin and c-peptide (connecting peptide) are formed. Potency of an insulin preparation was originally determined by bioassay. However, whole body bioassay inevitably shows poor precision owing to individual variation.

WHO issued 1st International Standard for human insulin in 1986 which has the potency of 26 IU/mg (0.038 mg/IU). In the same year, 1st International Standard of bovine insulin, the potency of which is 25.7 IU/mg, and Porcine insulin 1st International Standard, 26 IU/mg, were provided. Before these standards, in 1974, 1st International Reference Preparation of human insulin for immunoassay was provided as 3 IU/ampoule. Based on the above data, if the biological activity of insulin per molecule is the same among various animal species, potencies of animal insulin might be calculated from their molecular weights. But, so far, we do not have experimental proof about this. As the molecular weights of insulin of various animals are nearly the same, and the differences are within 1 %, there may be no critical fault if we think that the general potency of insulin is 26 IU/mg. Rat and mouse have two molecular species of insulin, type 1 and type 2. Amino acid sequences of these molecular species are same between rat and mouse. But as their ratios are different between these two animal species, it is recommended to use standard preparation derived from each animals.

4. Assay principle

In Shibayagi's LBIS Human Insulin ELISA Kit (Luminescent type), HRP-conjugated anti-insulin antibody, and standard or sample are incubated in monoclonal anti-insulin-coated wells to capture insulin bound with HRP-conjugated anti-insulin antibody. After 2 hours' incubation and washing, HRP remaining in wells are reacted with a luminescent reagent. Finally, peroxidase activity in each well is measured to determine insulin in the sample. The Relative Light Units (RLU) is proportional to insulin concentration. The standard curve is prepared by plotting RLU against standard insulin concentrations. Insulin concentrations in unknown samples are determined using this standard curve.

5. Precautions

- For professional use only, beginners are advised to use this kit under the guidance of experienced person. In manual operation, proficiency in pipetting technique is recommended.
- Wear protective gloves, clothing, eye, and face protection.
- Avoid contact of skin and mucous membranes with kit reagents or specimens. If any reagents come in

LBIS Human Insulin ELISA Kit (Luminescent type) (AKH-INSL)

contact with eyes, skin, or mucous membranes, wash with copious amounts of water and contact a physician.

- Do not drink, eat, or smoke in a place where this kit is used.
- Handle the sample with proper care, being aware that the sample may have an infection risk. This kit contains animal-derived ingredients.
- Do not use reagents with different lot numbers together.
- It is recommended to store the sample frozen at -35 °C or lower if long-term storage is intended. Avoid repeated freezing and thawing. Thaw the frozen sample just before assay and thoroughly agitate. Prepare the sample before use.
- Do not use a sample with hemolysis or containing high lipid.
- Centrifuge the sample to remove turbidity or insoluble matters where necessary before use for the assay.
- If presence of interfering substances is suspected for the sample, dilute the concerned one at multiple dilution ratios to check the dilution linearity. Dilute the prepared sample with the Diluent Buffer.
- When allowing the plate to stand in each step, always affix a plate sealer to protect the wells from drying, contamination with foreign matters, uneven temperature, and evaporation of a dispensed reagent.
- In order to avoid dryness of wells, contamination of foreign substances and evaporation of dispensed reagents, never forget to cover the well plate with a plate seal supplied, during incubation.
- The materials must not be pipetted by mouth.
- Residual samples and used tips should be sterilized before disposal.
- Dispose consumable materials and unused contents in accordance with applicable regional/national regulatory requirements.

6. Reagents supplied

Components	State	Amount
(A) Anti-Insulin antibody coated plate	Use after washing	96 wells/1 plate
(B) Standard human insulin solution (200 ng/mL)	Concentrated. Use after dilution	50 µL/1 vial
(C) Buffer solution	Ready to use.	60 mL/1 bottle
(D) HRP-conjugated anti-insulin antibody	Concentrated. Use after dilution.	100 µL/1 vial
(E) Luminescent Reagent 1	Mix (E) Luminescent Reagent 1 and (F) Luminescent Reagent 2 in 1 : 1 (vol./vol.)	6 mL/1 bottle
(F) Luminescent Reagent 2		6 mL/1 bottle
(I) Wash stock solution (10×)	Concentrated. Use after dilution.	100 mL/1 bottle
Plate seal	—	1 sheet
Instruction Manual	—	1 copy

(D) HRP-conjugated anti-insulin antibody: Vials contain more than volumes shown in the list. You can easily take out 100 µL, respectively, from vials.

7. Equipments or supplies required but not supplied ☐ Use as a check box

- ☐ Deionized water (or Distilled water). ☐ Test tubes for preparation of standard solution series.
- ☐ Glassware for dilution of Wash stock solution (10×) (a graduated cylinder, a bottle). ☐ Pipettes (disposable tip type). One should be able to deliver 5 µL precisely, and another for 50 µL-100 µL, 100 µL -1000 µL.
- ☐ Syringe-type repeating dispenser like Eppendorf multipette plus which can dispense 50 µL. ☐ Paper towel to remove washing buffer remaining in wells. ☐ A vortex-type mixer. ☐ A shaker for 96 well-plate (500-1200 rpm). ☐ An automatic washer for 96 well-plate (if available), or a wash bottle with a jet nozzle. ☐ A 96 well-plate reader for luminescence detection ☐ Software for data analysis.

8. Preparation of reagents

- ◆ Bring the reagents to room temperature (20 °C -25 °C) about 2 hours prior to use.
- ◆ Reagents described as “Use after dilution” in **【6. Reagent supplied】** shall be prepared as follows.
- ◆ Prepare reagent solutions in appropriate volume for your assay. Do not store the diluted reagents.

【Concentrated reagents】

[(B) Standard human insulin solution (200 ng/mL)]

Make a serial dilution of master standard (200 ng/mL) solution to prepare each standard solution (0.0381-20 ng/mL). * Unit reduction for µIU/mL is 26 IU/mg. (Refer to 3. Introduction.)

Volume of standard solution	Buffer solution	Concentration (ng/mL)	Concentration (μ IU/mL)*
Original solution: 10 μ L	90 μ L	20.0	520.0
20.0 ng/mL solution: 50 μ L	125 μ L	5.71	148.6
5.71 ng/mL solution: 50 μ L	125 μ L	1.63	42.45
1.63 ng/mL solution: 50 μ L	125 μ L	0.466	12.13
0.466 ng/mL solution: 50 μ L	125 μ L	0.133	3.465
0.133 ng/mL solution: 50 μ L	125 μ L	0.0381	0.9901
0 (Blank)	125 μ L	0	0

[(D)HRP-conjugated anti-insulin antibody]

Prepare working solution by dilution of (D) with the buffer solution (C) to **1:100**.

6 mL of the diluted solution is enough for 96 wells.

[(E) Luminescent Reagent 1]&[(F) Luminescent Reagent 2]

Mix Luminescent Reagent 1 and Luminescent Reagent 2 in **1 : 1 (vol./vol.)** and incubate for 15-30 minutes before use. **Keep in the dark.** Dispose any unused mixed luminescent reagent.

6 mL of the diluted solution is enough for 96 wells.

[(I) Wash stock solution (10 $\times$)]

Dilute 1 volume of the concentrated Wash stock solution (10 $\times$) to 10 volume with deionized water (or distilled water) to prepare working solution. Example: 100 mL of concentrated washing buffer (10 $\times$) and 900 mL of deionized water (or distilled water).

【Storage and stability】

[(A)Anti-Insulin coated plate]

If seal is not removed, put the strip back in a plastic bag with zip-seal originally used for well-plate container and store at 2 °C -8 °C. The strip will be stable until expiration date.

[(B) Standard human insulin solution (200 ng/mL)]

Standard solutions prepared above should be used as soon as possible, and should not be stored.

[(C)Buffer solution]

If not opened, store at 2 °C -8 °C. It maintains stability until expiration date. Once opened, we recommend using as soon as possible to avoid influence by environmental condition.

[(D)HRP-conjugated anti-insulin antibody]

If the kit is divided for multiple assays, prepare these solutions just after taking out of the refrigerator, immediately close the cap tightly, and store at 2 °C,-8 °C. It is stable until the expiration date. Unused working solution (already diluted) should be disposed.

[(E) Luminescent reagent 1, (F) Luminescent reagent 2]

If the kit is divided for multiple assays, prepare these solutions just after taking out of the refrigerator, immediately close the cap tightly, and store at 2 °C,-8 °C. It is stable until the expiration date. Discard the remaining mixed luminescent reagent.

[(I) Wash stock solution (10 $\times$)]

The rest of undiluted buffer: if stored tightly closed at 2 °C -8 °C, it is stable until expiration date.

Dispose any unused diluted buffer.

9. Technical tips

- Be careful to avoid any contamination of assay samples and reagents. We recommend the use of disposal pipette tips, and 1 tip for 1 well.
- The reagents are prepared to give accurate results only when used in combination within the same box. Therefore, do not combine the reagents from kits with different lot numbers. Even if the lot number is the same, it is best not to mix the reagents with those that have been preserved for some period.
- Optimally, the reagent solutions of the kit should be used immediately after reconstitution. Otherwise, store them in a dark place at 2 °C -8 °C.
- Time the reaction from the pipetting of the reagent to the first well.
- Dilution of the assay sample must be carried out using the buffer solution provided in the kit.
- To avoid denaturation of the coated antibody, do not let the plate go dry.
- As the antibody-coated plate is module type of 8 wells $\times$ 12 strips, each strip can be separated by cutting the cover sheet with a knife and used independently.
- When ELISA has to be done under the airstream velocity over 0.4 m/sec. and the humidity less than 30 %,

seal the well plate with a plate seal and place the well plate in an incubator or a styrofoam box in each step of incubation. For more details, watch our web movie [Assay circumstance].

10. Preparation of samples

This kit is intended to measure insulin in human serum, plasma (preferably obtained with EDTA or heparin), culture medium and tissue/cell extracts. The necessary sample volume for the standard procedure is 5 μ L. Samples should be immediately assayed or stored below -35°C for several days. Defrosted samples should be mixed thoroughly for best results.

Hemolytic and hyperlipemic serum samples are not suitable.

* To avoid influence of blood (high lipid or hemolysis, etc.), if your original samples have heavy chyle or hemolysis, do not use them for assay.

If presence of interfering substance is suspected, examine by dilution test at more than 2 points. Dilution of a sample should be made in a test tube using buffer solution prior to adding them to wells. Turbid samples or those containing insoluble materials should be centrifuged before testing to remove any particulate matter.

Storage and stability

Insulin in samples will be inactivated if stored at 2°C - 8°C . If it is necessary to store sample in refrigerator (2°C - 8°C), add aprotinin at final concentration of 100-500 KIU/mL. (KIU: kallikrein inhibitor unit).

If you have to store assay samples for a longer period, snap-freeze samples and keep them below -35°C . Avoid repeated freeze-thaw cycles.

11. Assay procedure


Remove the cover sheet of the anti-Insulin-coated plate after bringing up to room temperature.

- (1) Wash the anti-Insulin-coated plate (A) by filling the wells with 300 μ L of washing buffer and discard 4 times(*①), then strike the plate upside-down onto several sheets of paper towel to remove residual buffer in the wells.
- (2) Pipette 50 μ L of HRP-conjugated anti-insulin antibody to all wells. Shake the plate gently on a plate shaker(*②).
- (3) Pipette 5 μ L of standard solution to the wells designated for standards.
- (4) Pipette 5 μ L of sample to the designated sample wells..
- (5) Shake the plate gently on a plate shaker(*②).
- (6) Stick a plate seal(*③) on the plate and incubate for 2 hours at room temperature (20°C - 25°C).
- (7) Discard the reaction mixture. Rinse wells by filling the wells with 300 μ L of washing buffer and discard 4 times(*①), then strike the plate upside-down onto several sheets of paper towel to remove residual buffer in the wells.
- (8) Pipette 50 μ L of mixed luminescent reagent to all wells, and shake for 1 minute using a plate shaker. (*②)
- (9) After agitation, determine the luminescent intensity using a 96-well microplate reader (for luminescence measurement). It is recommended to perform the measurement 10 to 20 minutes after addition of the luminescent reagent.

*Refer to the page 14-15 for notes of *①, *② and *③.

12. Calculations

- (1) Prepare a standard curve for each assay. Prepare a standard curve using two-way logarithmic section paper by plotting RLU(Y-axis) against insulin concentration (ng/mL) on X-axis.
 - (2) Using the standard curve, read the insulin concentration of a sample at its RLU, and multiply the assay value by dilution factor if the sample has been diluted. Though the assay range is wide enough, in case the RLU of some samples is higher than that of the highest standard, please repeat the assay after proper dilution of samples with the buffer solution.
- We recommend the use of 3rd order regression curve for log-log plot, or 4 or 5 parameters method for log-normal plot in computer calculation.
 - Physiological or pathological situation of animals should be judged comprehensively taking other examination results into consideration.

13. Performance characteristics

- Assay range; The assay range of the kit is 0.0381 ng/mL ~ 20 ng/mL.

NIBSC/WHO International Standard Insulin, Human (code: 83/500) (4th) approximate value (μIU/mL)
 = 26 x 「LBIS Human Insulin ELISA Kit (Luminescent type)」 value (ng/mL)
 1 IU = 0.03846 mg, 26.0 IU/mg

- Precision of assay (Within assay variation)

3 samples, 5 replicated assay

n / ID	Sample1 (ng/mL)	Sample2 (pg/mL)	Sample3 (pg/mL)
1	5.84	0.980	0.207
2	5.92	0.943	0.210
3	5.79	0.962	0.207
4	5.91	0.940	0.201
5	5.74	0.973	0.201
mean	5.84	0.960	0.205
SD	0.078	0.018	0.004
CV(%)	1.3	1.9	2.0

Mean CV was within 15 %.

- Reproducibility (Between assay variation)

3 samples, 4 days, assayed in triplicate

Day / ID	Sample1 (ng/mL)	Sample2 (ng/mL)	Sample3 (ng/mL)
Day 0	5.52	0.472	0.742
Day 1	5.72	0.469	0.723
Day 2	5.55	0.467	0.785
Day 3	5.55	0.440	0.772
mean	5.58	0.462	0.756
SD	0.090	0.015	0.028
CV(%)	1.6	3.2	3.7

Mean CV was within 15 %

- Recovery test

Standard insulin was added in 6 concentrations to 3 samples and were assayed.

The recoveries were 91.4 ~ 103 %, 91.4 % ~ 103 %

Serum

Added (ng/mL)	Found (ng/mL)	Recovered (ng/mL)	Recovery (%)
—	0.409	—	—
0.143	0.548	0.139	97.2
0.358	0.746	0.337	94.1
0.896	1.24	0.830	92.6
2.24	2.57	2.16	96.0
5.60	5.64	5.23	93.4
14.0	13.7	13.3	95.0


Plasma (EDTA)

Added (ng/mL)	Found (ng/mL)	Recovered (ng/mL)	Recovery (%)
—	0.484	—	—
0.143	0.632	0.148	103
0.358	0.834	0.350	97.8
0.896	1.33	0.843	94.1
2.24	2.62	2.13	95.1
5.60	5.87	5.39	96.3
14.0	13.3	12.8	91.4

● Dilution test

Two serum samples were serially diluted by 6 steps.

The dilution curves showed linearity with $R^2 = 0.9998 \sim 0.9999$.

14. Trouble shooting

● Low RLU in all wells

Possible explanations:

- 1) The standard or samples might not be added.
- 2) Wrong reagents related to luminescence might have been added. Wrong dilution of HRP-conjugated anti-insulin antibody.
- 3) Contamination of peroxidase enzyme inhibitor(s).
- 4) Influence of the temperature under which the kits had been stored.
- 5) Excessive hard washing of the well plate.
- 6) Addition of luminescent reagent soon after taking out from a refrigerator might cause poor luminescence owing to low temperature.

● Intense luminescence in all wells including blank

Possible explanations:

- 1) Improper or inadequate washing. (Change washing frequency from 4 times to 5-6 times at the constant stroke after the reaction with HRP-conjugated anti-insulin / sample.)
- 2) Too high incubation temperature. Adjust the temperature to 20 °C-25 °C.

● High coefficient of variation (CV)

Possible explanation:

- 1) Improper or inadequate washing.
- 2) Improper mixing of standard or samples.
- 3) Pipetting at irregular intervals.

● Q-1: Can I divide the plate to use it for the other testing?

A-1: Yes, cut off the clear seal on the plate with cutter along strip. Put the residual plate, which is still the seal on, in a refrigerator soon

● Q-2: I found there contains liquid in 96 well-plate when I opened the box. What is it?

A-2: When we manufacture 96 well-plate, we insert preservation stabilizer in wells.

For detailed FAQs and explanations, refer to "Trouble shooting and Important Points in Shibayagi's ELISA kits" on our website.

Summary of assay procedure ☐ : Use as a check box


* First, read this instruction manual carefully and start your assay after confirmation of details.

For more details, watch our web movie [ELISA by MOVIE] on our website.

☐ Bring the well-plate and all reagents back to room temperature at 20 °C -25 °C for 2 hours.

☐ Wash stock solution (10×) must be diluted to 10 times by deionized water (or distilled water that returned to 20 °C -25 °C).

☐ Standard human insulin solution dilution example:

☐ (D)HRP-conjugated anti-insulin antibody : Dilute to **100 times** by using buffer solution(C) and use.

☐ Anti-Insulin-coated plate		
☐ ↓Washing 4 times(* ①)		* ④
☐ HRP-conjugated anti-insulin antibody	50 µL	* ⑤
☐ ↓Shaking(* ②)		
☐ Samples/Standards	5 µL	* ⑤
☐ ↓Shaking(* ②), Incubation for 2 hours at 20 °C -25 °C. (Standing(* ③))		* ⑥
☐ Mix Luminescent Reagent 1 and Luminescent Reagent 2 in 1 : 1 (vol./vol.) and incubate for 15-30 minutes before use.	Keep in the dark.	
☐ ↓Washing 4 times(* ①)		* ④
☐ Mixed luminescent reagent	50 µL	* ⑤
☐ ↓Shaking for 1 minute(* ②)		
☐ Measurement of luminescence intensity.		
☐ It is recommended to perform the measurement 10 to 20 minutes after addition of the luminescent reagent.		

* ① After dispensing wash buffer to wells, lightly shake the plate on your palm for 10 seconds and remove the buffer. Guideline of washing volume: 300 µL/well for an automatic washer and for a pipette if the washing buffer is added by pipette. In case of washing by using 8 channel pipette, sometimes the back ground tends to be high. If so, change washing frequency from 4 times to 5-6 times at the constant stroke after the reaction with HRP-conjugated anti-insulin antibody/sample.
Standard of plate-washing pressure: 5-25 mL/min. (Adjust it depending on the nozzle's diameter.) Refer to our web movie [Washing of microplate].

* ② Guideline of shaking: 500-1200 rpm for 10 seconds × 3 times.

* ③ Seal the plate during the reaction after shaking. Peel off the protective paper from the seal and stick the seal on the plate. Do not reuse the plate seal used once.

* ④ After removal of wash buffer, immediately dispense the next reagent.


* ⑤ Refer to our web movie [Handling of pipetting].

* ⑥ Refer to our web movie [Assay circumstance].

Worksheet example

	Strip 1&2	Strip 3&4	Strip 5&6	Strip 7&8	Strip 9&10	Strip 11&12
A	20.0 ng/mL	Sample 2	Sample 10	Sample 18	Sample 26	Sample 34
B	5.71 ng/mL	Sample 3	Sample 11	Sample 19	Sample 27	Sample 35
C	1.63 ng/mL	Sample 4	Sample 12	Sample 20	Sample 28	Sample 36
D	0.466 ng/mL	Sample 5	Sample 13	Sample 21	Sample 29	Sample 37
E	0.133 ng/mL	Sample 6	Sample 14	Sample 22	Sample 30	Sample 38
F	0.0381 ng/mL	Sample 7	Sample 15	Sample 23	Sample 31	Sample 39
G	0(Blank)	Sample 8	Sample 16	Sample 24	Sample 32	Sample 40
H	Sample 1	Sample 9	Sample 17	Sample 25	Sample 33	Sample 41

Assay worksheet

LBIS Human Insulin ELISA Kit (Luminescent type)

[Storage condition] Store the kit at 2 °C -8 °C (Do not freeze).

[Term of validity] Expiration date is indicated on the container.

[Manufacture #] AKH-INSL

[Cat #] 634-52911

This kit is
manufactured by

FUJIFILM Wako Shibayagi Corporation

1062-1 Ishihara, Shibukawa, Gunma, Japan 377-0007

TEL.+81-279(25)0279, FAX.+81-279(23)0313

<E-mail>wksb-info@fujifilm.com

<URL><http://www.shibayagi.co.jp>

distributed by

FUJIFILM Wako Pure Chemical Corporation