Human s-IgA ELISA Kit Wako

«For Research Use Only»

[1. Introduction]

Secretory IgA (s-IgA) is a type of immunoglobulin that is most abundant in exocrine secretions. It plays a major role in mucosal immunity and functions as the front line of immune response in the gastrointestinal and respiratory tracts. s-IgA levels in saliva are useful as a biomarker in various mental stress and psychological studies.

This ELISA kit is designed for specific and simple measurement of s-IgA in human saliva.

[2. Kit Performance]

Calibration curve range	0.082-20 $\mu\text{g}/\text{mL}$
Sample	Human saliva
Sample amount required	5 μL
Assay time	About 2.5 hours
Intra-assay CV (%)	CV < 6%
Inter-assay CV (%)	CV < 12%

[3. Materials supplied]

Components	State	Amount
Antibody-coated Plate	Use after washing	96 wells (8×12)/plate
s-IgA Standard	Lyophilized product, use after reconstitution	1 vial
Anti-Human s-IgA Antibody, HRP-conjugated	Ready-to-use	12 mL×1 vial
TMB Solution	Ready-to-use	12 mL×1 vial
Stop Solution	Ready-to-use	12 mL×1 vial
Buffer (2×)	Use after preparation	25 mL×1 vial
Wash Solution (20×)	Use after preparation	50 mL×1 vial
Plate Seal	Ready-to-use	3 sheets
Instruction Manual	—	1 copy

[4. Assay Principle]

Mouse anti-human s-IgA antibody is immobilized in each well of the assay plate. The standard solutions or samples are added to these wells and allowed to react with the antibody. Next, Anti-Human s-IgA Antibody, HRP-conjugated is added to the wells. Finally, the concentrations of s-IgA in the samples can be determined by measuring the peroxidase activity in the wells.

[5. Supplies and Equipment]

- Micropipettes and disposable tips (5 to 1,000 μL) (8- or 12-channel pipettes are recommended.)
- Microplate reader which can read extinction 3.0 at 450 nm

- Microplate shaker or conventional shaker
- Test tubes for dilution of standard solutions and samples (glass or polypropylene)
- Microplate washer (For manual operation, a continuous dispenser, needle dispenser, and an aspirator or vacuum pump are recommended.)
- Graduated cylinder (1,000 mL)
- Distilled or deionized water

[6. Reagent preparation]

The reagents in the kit must be brought to room temperature (20-30°C) before use (about 2 hours).

6.1 Preparation of Standard Solution

To prepare a 20 $\mu\text{g}/\text{mL}$ (0.5 mL) standard solution, add 0.5 mL of the Buffer to the vial containing the standard, allow to stand for about 10 minutes to dissolve, and mix thoroughly. Take 0.1 mL of the standard stock solution and dilute with 0.2 mL of the Buffer to make a 6.67 $\mu\text{g}/\text{mL}$ standard solution. Repeat the same dilution procedure to prepare standard solutions of 2.22, 0.74, 0.25, and 0.082 $\mu\text{g}/\text{mL}$. For 0 $\mu\text{g}/\text{mL}$ standard solution, use the Buffer.

Concentration ($\mu\text{g}/\text{mL}$)	Volume of standard solution	Buffer
6.67	Standard stock solution : 0.1 mL	0.2 mL
2.22	6.67 $\mu\text{g}/\text{mL}$ standard solution : 0.1 mL	0.2 mL
0.74	2.22 $\mu\text{g}/\text{mL}$ standard solution : 0.1 mL	0.2 mL
0.25	0.74 $\mu\text{g}/\text{mL}$ /standard solution : 0.1 mL	0.2 mL
0.082	0.25 $\mu\text{g}/\text{mL}$ standard solution : 0.1 mL	0.2 mL
0 (Blank)	—	0.2 mL

6.2 Buffer (2×)

Dilute necessary amount of the Buffer (2×) 2-fold with distilled water.

6.3 Wash Solution (20×)

Dilute 50 mL (the whole amount) of the Wash Solution (20×) with 950 mL of distilled water.

If diluted wash solution is stored in a refrigerator, bring it to room temperature before use.

- Other reagents are ready to use.

[7. Sample Preparation]

Frozen saliva samples should be thawed at room temperature and stirred well. Dilute saliva samples 80-fold with the Buffer prepared in 6.2, "Buffer (2×)."

1. Prepare 1 test tube for each saliva sample and dispense 395 μL of the Buffer into each test tube.
2. Add 5 μL of saliva sample to each tube and mix well.


[8. Assay Procedure]

1. Bring the kit components to room temperature (20-30°C). Prepare Standard Solutions, Buffer (2×), Wash Solution

(20×), and diluted samples as described above.

2. Add 350 μL of the Wash Solution into each well and aspirate the solution. Alternatively, invert the plate, discard the buffer, and then remove the remaining solution by blotting on a paper towel. Repeat this washing procedure further twice (total 3 times).
3. Add 100 μL of the Buffer into each well, then add 50 μL of standard solutions (0, 0.082, 0.25, 0.74, 2.22, 6.67, and 20 $\mu\text{g}/\text{mL}$) or diluted samples.
4. Agitate the plate on a microplate shaker for 1 minute (500 rpm). Cover with the plate seal, and incubate for 1 hour at room temperature with shaking (approximately 100 rpm).
5. After the reaction is complete, discard the solution from each well and wash as in step 2, 6 times in total.
6. Add 100 μL of Anti-Human s-IgA Antibody, HRP-conjugated into each well.
7. Cover with the plate seal and incubate for 1 hour at room temperature with shaking (approximately 100 rpm).
8. Discard the solution from each well and wash as in step 2, 6 times in total.
9. Add 100 μL of TMB Solution into each well and allow to stand in a dark place for 30 minutes at room temperature.
10. Add 100 μL of the Stop Solution to each well.
11. Measure the absorbance at 450 nm (sub-wavelength : 620 nm) using a microplate reader.
12. Using commercially available software, create a standard curve from the measured values of each s-IgA standard solution using a 5 (or 4)-parameter regression equation and determine s-IgA concentrations of the diluted samples. Multiply the obtained concentrations by 80 to calculate the s-IgA concentrations of the original saliva samples. When using semi-log graph paper, plot the concentration of the standard solution on the x axis (log scale) and the absorbance of the standard solution at each concentration on the y axis (linear scale) to create a standard curve. Determine the s-IgA concentration that corresponds to the absorbance of each sample using the standard curve.

[9. Typical Standard Curve]

[10. Precautions]

1. Saliva samples should be stored frozen at -30°C or below.

Avoid repeated freezing and thawing of samples.

2. In general, reagents should be freshly prepared. In particular, the standard solutions should be used immediately after preparation. If the kit is used in portions, the prepared standard solutions should be divided into small portions and stored frozen at -30°C or below.
3. Precipitates may form in the Wash Solution (20×) during storage. The precipitates will dissolve with dilution, and the quality will not be affected.
4. When diluting standard solution, be sure to use a new tip for each dilution step. Dispensing to each well should be performed accurately, as it affects the accuracy of the assay. When adding samples into wells, use a new tip for each sample and take care to avoid cross-contamination.
5. If the concentration of a diluted sample exceeds 20 $\mu\text{g}/\text{mL}$, dilute the 80-fold diluted sample further with the Buffer and use for the assay.
6. For incubation at room temperature, always use a microplate shaker to shake the assay plate (except for the chromogenic reaction). Shaking should be done slowly (about 100 rpm) to prevent splashing of the reaction solution on the plate seal.
7. All measurements should be performed in duplicate.
8. Absorbance should be taken immediately after the reaction is stopped.
9. The coloration level of the substrate may be slightly affected by the reaction temperature, time, degree of shaking of the assay plate, etc. Therefore, the standard curve must be prepared for each assay.
10. Avoid exposure to strong light during storage or use of each reagent.
11. Do not combine kits from different lots.
12. TMB Solution should be brought to room temperature in dark place before use.
13. Some of the reagents contain biological materials of human origin and should be handled accordingly.

[11. Measurement Examples]

11.1 Spiked Recovery Test

	Spiked s-IgA ($\mu\text{g}/\text{mL}$)	Assay value ($\mu\text{g}/\text{mL}$)	Recovery (%)
Saliva sample 1	0	64.8	—
	40	108.8	103.8
	160	220.8	98.2
	400	404.8	87.1
Saliva sample 2	0	46.4	—
	40	96.0	111.1
	160	196.0	95.0
	400	432.0	96.8
Saliva sample 3	0	120.8	—
	40	189.6	117.9
	160	274.4	97.7
	400	476.0	91.4

	Spiked s-IgA ($\mu\text{g/mL}$)	Assay value ($\mu\text{g/mL}$)	Recovery (%)
Saliva sample 4	0	112.0	—
	40	157.6	103.7
	160	244.8	90.0
	400	441.6	86.3
Saliva sample 5	0	105.6	—
	40	152.0	104.4
	160	250.4	94.3
	400	600.8	118.8

11.2 Dilution Linearity Test

	Dilution	Measured value ($\mu\text{g/mL}$)
Saliva sample 1	1/1	363.5
	1/2	201.7
	1/4	107.7
	1/8	53.6
Saliva sample 2	1/1	237.4
	1/2	123.9
	1/4	65.3
	1/8	32.5
Saliva sample 3	1/1	310.0
	1/2	167.3
	1/4	89.5
	1/8	42.7
Saliva sample 4	1/1	293.3
	1/2	145.3
	1/4	74.7
	1/8	38.3
Saliva sample 5	1/1	188.2
	1/2	97.5
	1/4	47.4
	1/8	26.8

11.3 Cross-reactivity

Related antibody	Cross-reactivity (%)
Serum IgG (Human)	< 0.1
Serum IgA (Human)	0.3
Serum IgM (Human)	0.0
Serum IgE (Human)	0.0

[12. Assay Procedure Summary]

Ensure that you read the instruction manual and check the sample conditions, assay conditions, and assay procedure before beginning the assay.

<input type="checkbox"/>	Thoroughly equilibrate the plates and reagents to room temperature (20°C to 30°C) before use (approx. 2 hours)
<input type="checkbox"/>	Dilution of Buffer (2×) and Wash Solution (20×) : Dilute the Buffer (2×) 2-fold and Wash Solution (20×) 20-fold with distilled water.

<input type="checkbox"/>	Preparation of Standard Solutions (example) : To prepare a 20 $\mu\text{g/mL}$ Standard Stock Solution, add 0.5 mL of the Buffer to the vial of s-IgA standard, allow to stand for about 10 minutes to dissolve thoroughly, and mix. Dilute with the Buffer as below to prepare standard solutions of 6.67, 2.22, 0.74, 0.25, and 0.082 $\mu\text{g/mL}$.																					
Dilution example	<table border="1"> <thead> <tr> <th>Concentration ($\mu\text{g/mL}$)</th> <th>6.67</th> <th>2.22</th> <th>0.74</th> <th>0.25</th> <th>0.082</th> <th>0</th> </tr> </thead> <tbody> <tr> <td>Standard Solution (mL) Original</td> <td>0.1</td> <td>0.1*</td> <td>0.1*</td> <td>0.1*</td> <td>0.1*</td> <td>—</td> </tr> <tr> <td>Buffer (mL)</td> <td>0.2</td> <td>0.2</td> <td>0.2</td> <td>0.2</td> <td>0.2</td> <td>0.2</td> </tr> </tbody> </table>	Concentration ($\mu\text{g/mL}$)	6.67	2.22	0.74	0.25	0.082	0	Standard Solution (mL) Original	0.1	0.1*	0.1*	0.1*	0.1*	—	Buffer (mL)	0.2	0.2	0.2	0.2	0.2	0.2
Concentration ($\mu\text{g/mL}$)	6.67	2.22	0.74	0.25	0.082	0																
Standard Solution (mL) Original	0.1	0.1*	0.1*	0.1*	0.1*	—																
Buffer (mL)	0.2	0.2	0.2	0.2	0.2	0.2																
	* : Standard solution at 1 step higher concentration																					

<input type="checkbox"/>	Assay plate with immobilized antibody
<input type="checkbox"/>	↓ Wash 3 times.
<input type="checkbox"/>	Buffer 100 μL
<input type="checkbox"/>	Samples or Standard Solutions 50 μL
<input type="checkbox"/>	↓ Agitate the plate in on a microplate shaker for 1 minute (approximately 500 rpm).
<input type="checkbox"/>	↓ Cover with the plate seal and incubate for 1 hour at room temperature with shaking (approximately 100 rpm).
<input type="checkbox"/>	↓ Wash 6 times.
<input type="checkbox"/>	Anti-Human s-IgA Antibody, HRP-conjugated 100 μL
<input type="checkbox"/>	↓ Cover with the plate seal and incubate for 1 hour at room temperature with shaking (approximately 100 rpm).
<input type="checkbox"/>	↓ Wash 6 times.
<input type="checkbox"/>	TMB Solution 100 μL
<input type="checkbox"/>	↓ Allowed to stand for 30 minutes at room temperature in dark.
<input type="checkbox"/>	Stop Solution 100 μL
<input type="checkbox"/>	Absorbance measurement (primary wavelength at 450 nm, secondary wavelength at 620 nm)

[Storage]

Store at 2-10°C

[Expiration date]

Indicated on the label

[Package]

For 96 assays

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
 Telephone : + 81-6-6203-3741
 Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
 Richmond, VA 23237
 U.S.A.
 Telephone : + 1-804-271-7677
 Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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