Oxytocin ELISA Kit Wako

[1. Introduction]

This item is an ELISA kit designed for the quantitative detection and measurement of oxytocin in biological fluids such as saliva, urine, serum and plasma. The kit can be measured by a simple pretreatment, and does not require purification using a C18 column and organic solvent, which was conventionally performed as a pretreatment method. For Research Use Only. Not for use in diagnostic procedures. Oxytocin is called as a love or happy hormone consisting of 9 amino acids. It also has stress relief, learning ability improvement, anxiolytic, and anti-fear effects, and is expected to have therapeutic effects on mental illness (autism, depression, etc.).

[2. Assay Principle]

The microplate is coated with anti-oxytocin antibody. In each well, the standard solution or sample and biotin-conjugated oxytocin are incubated to proceed with the antigen-antibody reaction. Furthermore, peroxidase-conjugated streptavidin is added to proceed with the biotin-streptavidin binding. Finally, peroxidase activity in each well is measured to determine oxytocin in the sample.

[3. Kit Performance]

Calibration curve range	4.00 ~ 12,500 pg/mL
Analyte	Oxytocin
Sample	Human saliva/urine/serum/plasma Mouse serum/plasma Rat serum/plasma
Sample amount required	50 μ L (n = 1, minimum volume) 200 μ L (n = 2, recommended volume) ※ 1
Assay time	Approx. 2.5 hours
Detection method	Luminescent detection ※ 2

※ 1 A sample volume of 200 μ L is recommended for duplicate measurements. The volume of 200 μ L is necessary to separate the supernatant without disturbing the "gel-containing pellet" generated after adding Pretreatment Solution 2 and centrifugation. (In theory, if the gel is not disturbed, sample volumes less than 200 μ L can be used.)

※ 2 Requires luminescence plate reader for measurement.

[4. Precautions for Use]

1. It is important to bring the kit reagents to room temperature (20-25°C) before use (approximately 2 hours).
2. Prepare the reagents to be used in the next step in advance, before starting each step.
3. Use this kit after completing training in ELISA methods or

— 1/20 —

under the guidance of an instructor.

4. For manual measurements, individuals who have demonstrated consistent reproducibility in pipette operation should perform the assay.
5. Wear gloves, glasses, and protective clothing during preparation and when using this kit.
6. Prevent reagents coming into contact with the skin. In the event the kit's reagents accidentally come into contact with eyes, mouth, wounds, skin, etc., immediately wash thoroughly with running water and seek medical attention from a physician if necessary.
7. Do not eat, drink, or smoke in the area where this kit is being used.
8. Handle samples with particular care as they may pose a risk of infection. This kit contains animal-derived components.
9. Used specimens and consumables should be soaked in 1% formaldehyde, 2% glutaraldehyde or sodium hypochlorite solution (>0.1%) for at least one hour, or autoclaved and disposed of. Used consumables and unused chemicals should be disposed of in accordance with the regulations of the facility and local laws.
10. Do not mix reagents from different lots.
11. It is important to cover the plate with a plate seal to prevent drying of the wells, contamination, temperature fluctuations, and evaporation of dispensed reagents during incubation at each step.
12. ELISA methods can be affected by the assay environment. Strictly adhere to the following: Room temperature should be 20-25°C (on the bench or in the incubator) during operation and incubation. Do not perform the assay in environments with wind speeds (including air conditioning) exceeding 0.4 m/sec and humidity levels below 30%.

[5. Materials supplied]

5-1. Kit Components

Item	State	Amount
Antibody-coated Plate	Use after washing	96 wells (8 × 12) / 1 plate
Oxytocin Standard	Use after preparation	100 μ L / 1 vial
Buffer	Ready-to-use	60 mL / 1 vial
Biotin-conjugated Oxytocin	Lyophilized / Use after reconstitution	1 vial
Peroxidase-conjugated Streptavidin Solution	Use after dilution	100 μ L / 1 vial
Luminescent Reagent 1	Ready-to-use	6 mL / 1 vial
Luminescent Reagent 2	Ready-to-use	6 mL / 1 vial
Wash Solution (10 ×)	Use after preparation	100 mL / 1 vial
Sample Buffer 1	Ready-to-use	5 mL / 1 vial
Sample Buffer 2	Ready-to-use	12 mL / 1 vial
Plate Seal	Ready-to-use	4 sheets
Manual	—	1 copy

5-2. Handling of Unused Reagents

When using the kit in divided portions, refer to the following

— 2/20 —

for the stability and storage of unused reagents. If unopened, they remain stable through the expiration date.

- (1) **Antibody-coated Plate**
The plate can be used in divided portions. Unused antibody-coated strips should be returned into the zip-lock bag provided in the kit and stored at 2-10°C.
- (2) **Oxytocin Standard**
Store unused portion at 2-10°C.
Note : Do not store diluted standard solutions.
- (3) **Buffer**
When using a portion of the solutions, transfer a slightly larger amount than necessary to separate containers. The remaining unused solutions should be tightly capped and stored at 2-10°C immediately; do not bring them to room temperature.
- (4) **Biotin-conjugated Oxytocin**
Store unused portion at 2-10°C.
After reconstitution, store at 2-10°C and use within 2 weeks.
Do not store Biotin-conjugated Oxytocin solution diluted 100-fold with Buffer.
- (5) **Peroxidase-conjugated Streptavidin Solution**
When using the kit in divided portions, take out the solution from the refrigerator when ready to prepare dilutions. The remaining undiluted solution should be tightly capped and stored at 2-10°C immediately ; do not bring it to room temperature.
Discard any remaining diluted solutions.
- (6) **Luminescent Reagent 1/Luminescent Reagent 2**
When using the kit in several divided portions, take out the solution from the refrigerator, prepare, immediately close the lid tightly, and store the remaining stock solution at 2-10°C (do not return to room temperature).
- (7) **Wash Solution (10×)**
When storing the Wash Solution (10×), close the lid tightly, and store at 2-10°C.
Discard the unused portions of diluted Wash Solutions.
- (8) **Sample Buffer 1/Sample Buffer 2**
When storing the Sample Buffer 1/Sample Buffer 2, close the lid tightly, and store at 2-10°C.

[6. Apparatus, equipment, and materials required]

- Purified water (distilled water)
- Tubes for dilution of standard solutions and samples
- Glassware for dilution of Wash Solution (graduated cylinders, beakers)
- Micropipettes (for 10 µL disposable tips and 200 to 500 µL disposable tips)
- Continuous dispensing pipette, capable of dispensing 100 µL continuously
- Paper towels or other absorbent material (to remove liquid left on the plate after washing)
- Vortex-type mixer
- Microplate shaker (range approx. 600 to 800 rpm)
- Automatic washer for 96-well plate (if available) or washing

bottle

- Microplate reader for 96-well plates (for absorbance measurement)
- Software for data processing

[7. Sample Preparation]

7-1. Recommendations and precautions before sample preparation

- **Sample collection**
Add 100 KIU/mL aprotinin and 1/750 amount (v/v) Proclin950 to sample.
- **Sample storage**
 - It is recommended to store the sample with 100 KIU/mL aprotinin and 1/750 amount (v/v) Proclin950 at -80°C or lower for long-term storage. Avoid repeated freezing and thawing.
 - Frozen samples should be thawed just before the assay and thoroughly mixed.
- **Sample dilution**
 - Please dilute the sample before the following preparation process.
 - Use the Buffer provided in the kit or saline as the dilution solvent.
 - Use test tubes (PP, PE) to dilute the samples.
 - Dilute the samples just before use.
- **Other issues**
 - When turbidity or insoluble materials are present in the samples, remove them by centrifugation or other methods before the assay.
 - If it is suspected that the sample contains an interfering substance, perform a dilution linearity check with two or more different dilutions of the same sample.

7-2. Sample Preparation

- (1) Add Sample Buffer 1 to sample in the table below. Shake it and incubate 5 minutes. Shake it again and incubate 5 minutes.
- (2) Add Sample Buffer 2 which gel is mixed evenly to sample in the table below. Shake it and incubate 5 minutes. Shake it again and incubate 5 minutes.
- (3) After shake it, centrifuge (5,000-6,000 g, 10 minutes, 4°C) and use the supernatant as a sample for measurement. Don't suck the gel in the buffer when collecting it.

	saliva/urine/serum/plasma		
Sample volume	50 µL	100 µL	200 µL
Sample Buffer 1	5 µL	10 µL	20 µL
Sample Buffer 2	25 µL	50 µL	100 µL
Dilution factor	1.35		

[8. Reagent Preparation]

Reagents in the kit must be brought to room temperature (20-25°C) before use (about 2 hours).

8-1. Preparation of Standard Solution

Before starting the procedure, Prepare a oxytocin standard stock solution (100 ng/mL) with the Buffer supplied with the

kit pre-equilibrated to room temperature.

Examples :

Concentration of the standard solution (pg/mL)	Volume of the standard solution	Volume of Buffer
12,500	Standard stock solution : 30 μ L	210 μ L
2,500	12,500 pg/mL solution : 50 μ L	200 μ L
500	2,500 pg/mL solution : 50 μ L	200 μ L
100	500 pg/mL solution : 50 μ L	200 μ L
20.0	100 pg/mL solution : 50 μ L	200 μ L
4.00	20.0 pg/mL solution : 50 μ L	200 μ L
0 (B ₀)	–	200 μ L

8-2. Biotin-conjugated Oxytocin

Before starting the procedure, dissolve the lyophilized product in 100 μ L of purified water and dilute the required amount 100-fold with Buffer.

The stock solution after dissolution can be stored at 2-10°C for up to 2 weeks.

Discard any remaining diluted solutions.

8-3. Peroxidase-conjugated Streptavidin Solution

During the 2 hours competitive reaction, dilute 100-fold with room temperature Buffer.

8-4. Luminescent Reagent 1/Luminescent Reagent 2

※ extemporaneously prepared

Mix Luminescent Reagent 1 and Luminescent Reagent 2 in a 1 : 1 ratio (vol./vol.) and stored protected from light 30 minutes before use*.

(e.g.) When solution for 96-well reaction is prepared.

Mix 6 mL of Luminescent Reagent 1 and 6 mL of Luminescent Reagent 2

* This means during the reaction of Peroxidase-conjugated Streptavidin Solution.

8-5. Wash Solution (10 \times)

Wash Solution (10 \times) is diluted 10-fold with purified water (distilled water) .

(e.g.) When solution for 96-well reaction is prepared.

Add 100 mL of Wash Solution (10 \times) to 900 mL of purified water (distilled water) to prepare 1,000 mL of Wash Solution (1 \times).

Other reagents are ready-to-use.

[9. Assay Procedure]

• Prepare the reagents to be dispensed in the next step before starting each step.

• For dispensing reagents in steps (2), (9) and (13), continuous dispensing using a multichannel pipette is recommended.

(1) Fill each well of the Antibody-coated Plate with the diluted Wash Solution and remove. Repeat four times.

Then, turn the plate upside down on a paper towel and gently tap to remove any remaining liquid in the wells.

Note : Steps (2) through (5) below should be completed within 10-15 minutes.

(2) Dispense 50 μ L of Biotin-conjugated Oxytocin solution to each well.

(3) Agitate the plate on a microplate shaker. Lightly agitation (for about 2 to 3 seconds once) until the liquid level in the well is even.

(4) Dispense 50 μ L of standard solution into the wells designated for standards.

(5) Dispense 50 μ L of pretreated-samples (ref : 7. Sample Preparation) into the wells designated for samples.

(6) Agitate the plate using a microplate shaker.

Note : Shaking conditions : 600-800 rpm for 10 seconds. Repeat 3 times.

After each cycle of shaking, stop and then restart again.

(7) Cover with a Plate Seal.

Incubate at room temperature (20-25°C) for 2 hours with shaking 500 rpm by microplate shaker.

Note : Peel the protective paper off the plate seal and seal with the adhesive side facing the plate.

Do not reuse a plate seal that has already been used.

(8) After the reaction is complete, remove the reaction mixture and fill each well with Wash Solution and wash 4 times.

Then, turn the plate upside down on a paper towel and gently tap to remove any remaining liquid in the wells.

(9) Dispense 100 μ L of diluted Peroxidase-conjugated Streptavidin Solution to each well.

(10) Agitate the plate on a microplate shaker.

Note : Shaking conditions : 600-800 rpm for 10 seconds. Repeat 3 times.

After each cycle of shaking, stop and then restart again.

(11) Cover with a Plate Seal.

Incubate at room temperature (20-25°C) for 30 minutes without shaking.

Note : Peel the protective paper off the plate seal and seal with the adhesive side facing the plate.

Do not reuse a plate seal that has already been used.

(12) After the reaction is complete, remove the reaction mixture and fill each well with Wash Solution and wash 4 times.

Then, turn the plate upside down on a paper towel and gently tap to remove any remaining liquid in the wells.

(13) Dispense 100 μ L of Luminescent Reagent prepared at room temperature into each well.

(14) Shake the plate for 1 minute using the microplate shaker.

(15) After agitation, determine the luminescent intensity (RLU) using the 96-well microplate reader (for luminescence measurement).

It is recommended to perform the measurement 10 to 20 minutes after addition of the Luminescent Reagent.

[Supplementary Notes]

Washing Procedure

For each wash, dispense the Wash Solution into the wells, then agitate gently by hand for 10 seconds and discard.

After 4 consecutive washes, place the plate upside down on a paper towel and tap to completely remove the Wash Solution.

After removing the Wash Solution, immediately dispense the next solution to prevent the wells from drying out. The recommended volume of Wash Solution to be dispensed is 300 μL /well.

If foam remains during washing, when manually washing using a wash bottle, allow the Wash Solution to overflow during the final cycle of washing to prevent the formation of foam.

[10. Calculation Method]

- (1) Create a standard curve by plotting the concentration of the standard solution (pg/mL) on the X-axis against the B/B₀ (%) value on the Y-axis.
- (2) Read the concentration (pg/mL) corresponding to the B/B₀ (%) of the diluted sample.
- (3) Multiply the read concentration by the sample dilution factor to obtain the measured concentration.

[Note]


- Luminescence intensity varies depending on the manufacturer and model of the plate reader.
- For calculations using computer software, we recommend the use of a cubic polynomial, and 4 or 5 parameters.
- The oxytocin concentration using this kit is determined with reference to the WHO International Standard (Oxytocin 4th International Standard ; NIBSC code : 76/575).

Examples :

Oxytocin (pg/mL)	RLU ($\times 100$)	B/B ₀ (%)
12,500	12,764	18.8
2,500	38,585	56.9
500	52,102	76.8
100	59,618	87.9
20.0	63,615	93.8
4.00	65,255	96.2
0 (B ₀)	67,820	100.0

$$B/B_0(\%) = (\text{RLU of standard or sample} / \text{RLU of } 0 \text{ pg/mL (B}_0\text{)}) \times 100$$

[11. Typical standard Curve]

[12. Assay Procedure Summary]

Ensure that you read the instruction manual and check the sample conditions, assay conditions, and assay procedure before beginning the assay.

- Plates and reagents must be brought to room temperature (20 to 25°C) before use (about 2 hours).
- Dilution of Wash Solution (10 \times) : Dilute 10-fold with purified water pre-equilibrated to room temperature.
- Preparation of standard solutions (example) : Prepare the standard solution (100 ng/mL) with the Buffer supplied with the kit pre-equilibrated to room temperature. The following is an example.

Dilution example	Concentration (pg/mL)	12,500	2,500	500	100	20.0	4.00	0
	Standard (μL)	stock solution	30	50*	50*	50*	50*	50*
Buffer (μL)		210	200	200	200	200	200	200

* : Standard solution at 1 step higher concentration

- *Preparation of Biotin-conjugated Oxytocin : Dissolve the lyophilized product in 100 μL of purified water and dilute the required amount 100-fold with the Buffer.
- *Prepare samples according to [7. Sample Preparation].

- Antibody-coated Plate**
- ↓ Wash 4 times (*1).
- Diluted biotin-conjugated oxytocin solution 50 μL /well**
- ↓ Agitate using a microplate shaker
*Lightly agitation (for about 2 to 3 seconds once)
- Standard solution or sample after preparation 50 μL /well**
- ↓ Agitate using a microplate shaker (*2), room temperature (20-25°C) (*3)
- ↓ cover with a plate seal.
- ↓ Agitate using a microplate shaker for 2 hours (500 rpm)
*Preparation of peroxidase-conjugated streptavidin solution (dilute 100-folds with Buffer equilibrated to room temperature)
- ↓ Washing 4 times (*1).
- Peroxidase-conjugated streptavidin solution 100 μL /well**
- ↓ Agitate using a microplate shaker (*2)
- ↓ incubate at room temperature (20-25°C) for 30 minutes (*3)
*Preparation of the luminescent reagent (mix Luminescent Reagent 1 and Luminescent Reagent 2 in 1 : 1 (vol./vol.))
- ↓ Washing 4 times (*1).
- Luminescent reagent 100 μL /well**
- ↓ Agitate at room temperature (20-25°C) for 1 minute
- Measurement of luminescence intensity (to be measured between 10 and 20 minutes)**

- (*1) For each wash, dispense the Wash Solution into the wells, then agitate gently by hand for 10 seconds and discard. After 4 consecutive washes, place the plate upside down on a paper towel and tap to completely remove the Wash Solution. After removing the Wash Solution, immediately dispense the next solution to prevent the wells from drying out. The recommended volume of Wash Solution to be dispensed is 300 μL /well. If foam remains during washing,

when manually washing using a wash bottle, allow the Wash Solution to overflow during the final cycle of washing to prevent the formation of foam.

- (*2) Shaking should be done at 600-800 rpm for 10 seconds. Repeat 3 times.
- (*3) After shaking is complete, cover with a plate seal and incubate without shaking. Peel the protective paper off the plate seal and seal with the adhesive side facing the plate. Do not reuse a plate seal that has already been used.

[Storage]

Store at 2°C to 10°C

[Expiration date]

Indicated on the label

[Package size]

For 96 assays

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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