## LBIS™ Mouse/Rat Insulin ELISA Kit (Luminescent)

Please, read this instruction carefully before use.

### 1. Intended use

LBIS™ Mouse/Rat Insulin ELISA Kit (Luminescent) is a sandwich ELISA system for quantitative measurement of mouse/rat insulin. This is intended for research use only. Not for use in diagnostics procedure.

### 2. Introduction

Insulin is a peptide hormone secreted from B cells of islet of Langerhans in the pancreas with a molecular weight of about 5800 and pI 5.4. It is consisted of 2 chains, A and B. It has 3 disulfide bonds formed between A6 and A11, A7 and B7, and A20 and B19. Insulin exists as a dimer molecule in acidic to neutral solution without Zn ion, and as a hexamer including two Zn ions in neutral solution if Zn ions are present. Main targets of insulin are liver, muscle, and adipose tissue. Insulin actions in these targets are as follows. In the liver, it promotes glycogenesis, protein synthesis, fatty acid synthesis, carbohydrate utilization, and inhibition of gluconeogenesis. In the muscle, it promotes membrane permeability for carbohydrates, amino acids and K ion, glycogenesis, protein synthesis, while inhibits protein degradation. In the adipose tissue, it promotes membrane permeability for glucose and fatty acid synthesis.

A precursor of insulin, called proinsulin with a single polypeptide chain, is first synthesized in the cell, then sulfide bonds are formed, and finally by enzymatic cutting at two sites, active insulin and c-peptide (connecting peptide) are formed. Potency of an insulin preparation was originally determined by bioassay. However, whole body bioassay inevitably shows poor precision owing to individual variation.

WHO issued 1<sup>st</sup> International Standard for human insulin in 1986 which has the potency of 26 IU/mg (0.038 mg/IU). In the same year, 1<sup>st</sup> International Standard of bovine insulin, the potency of which is 25.7 IU/mg, and Porcine insulin 1<sup>st</sup> International Standard, 26 IU/mg, were provided. Before these standards, in 1974, 1st International Reference Preparation of human insulin for immunoassay was provided as 3 IU/ampoule. Based on the above data, if the biological activity of insulin per molecule is the same among various animal species, potencies of animal insulin might be calculated from their molecular weights. But, so far, we do not have experimental proof about this. As the molecular weights of insulin of various animals are nearly the same, and the differences are within 1%, there may be no critical fault if we think that the general potency of insulin is 26 IU/mg. Rat and mouse have two molecular species of insulin, type 1 and type 2. Amino acid sequences of these molecular species are same between rat and mouse. But as their ratios are different between these two animal species, it is recommended to use standard preparation derived from each animals.

### 3. Assay Principle

In FUJIFILM Wako Pure Chemical Corporation's LBIS™ Mouse/Rat Insulin ELISA Kit (Luminescent), Peroxidase-conjugated Antibody Solution, and standard or sample are incubated in monoclonal anti-insulin-coated wells to capture insulin bound with Peroxidase-conjugated Antibody Solution. After 2 hours' incubation and washing, HRP remaining in wells are reacted with a luminescent reagent. Finally, peroxidase activity in each well is measured to determine insulin in the sample. The Relative Light Units (RLU) is proportional to insulin concentration. The standard curve is prepared by plotting RLU against standard insulin concentrations. Insulin concentrations in unknown samples are determined using this standard curve.

### 4. Performance Characteristics

- Assay Range ; The assay range of the kit is 0.0391 ng/mL - 20.0 ng/mL.

1 IU=0.038 mg, 26 IU/mg

• Precision of assay (Within assay variation)  
3 samples, 5 replicated assay  
Mouse\* Serum

n/ID	Sample1 (ng/mL)	Sample2 (ng/mL)
1	5.23	0.776
2	5.27	0.768
3	5.36	0.844
4	5.38	0.782
5	5.35	0.809
mean	5.31	0.795
SD	0.065	0.031

• Reproducibility (Between assay variation)  
3 samples, 4 days, assayed in triplicate  
Mouse\* Serum

Day/ID	Sample1 (ng/mL)	Sample2 (ng/mL)	Sample3 (ng/mL)
Day 0	2.48	0.664	0.175
Day 1	2.52	0.642	0.171
Day 2	2.55	0.618	0.160
Day 3	2.42	0.621	0.163
mean	2.49	0.636	0.167
SD	0.056	0.021	0.0069
CV (%)	2.2	3.3	4.1

• Recovery Test

Standard insulin was added in 6 concentrations to 3 samples and were assayed.

The recoveries were 93.3% - 97.7%, 93.8% - 97.8%

Mouse\* Serum

Added (ng/mL)	Found (ng/mL)	Recovered (ng/mL)	Recovery (%)
-	0.294	-	-
0.256	0.544	0.250	97.7
0.640	0.891	0.597	93.3
1.60	1.85	1.56	97.5
4.00	4.05	3.76	94.0



Mouse\* Plasma (EDTA)

Added (ng/mL)	Found (ng/mL)	Recovered (ng/mL)	Recovery (%)
-	0.346	-	-
0.256	0.592	0.246	96.1
0.640	0.949	0.603	94.2
1.60	1.85	1.50	93.8
4.00	4.26	3.91	97.8

• Dilution Test

Two serum samples were serially diluted by 3 steps.

The dilution curves showed linearity with  $R^2 = 0.9932 - 0.9994$ ,  $R^2 = 0.9912 - 0.9998$ .

\* : Mouse : Male CD-1 (ICR), Rat : Male Sprague Dawley

5. Precautions

- For professional use only, beginners are advised to use this kit under the guidance of experienced person. In manual operation, proficiency in pipetting technique is recommended.
- Wear protective gloves, clothing, eye, and face protection.
- Avoid contact of skin and mucous membranes with kit reagents or specimens. If any reagents come in contact with eyes, skin, or mucous membranes, wash with copious amounts of water and contact a physician.

- Do not drink, eat, or smoke in a place where this kit is used.
- Handle the sample with proper care, being aware that the sample may have an infection risk. This kit contains animal-derived ingredients.
- Do not use reagents with different lot numbers together.
- It is recommended to store the sample frozen at  $-35^{\circ}\text{C}$  or lower if long-term storage is intended. Avoid repeated freezing and thawing. Thaw the frozen sample just before assay and thoroughly agitate. Prepare the sample before use.
- Do not use a sample with hemolysis or containing high lipid.
- Centrifuge the sample to remove turbidity or insoluble matters where necessary before use for the assay.
- If presence of interfering substances is suspected for the sample, dilute the concerned one at multiple dilution ratios to check the dilution linearity. Dilute the prepared sample with the Diluent Buffer.
- When allowing the plate to stand in each step, always affix a plate sealer to protect the wells from drying, contamination with foreign matters, uneven temperature, and evaporation of a dispensed reagent.
- In order to avoid dryness of wells, contamination of foreign substances and evaporation of dispensed reagents, never forget to cover the well plate with a plate seal supplied, during incubation.
- The materials must not be pipetted by mouth.
- Residual samples and used tips should be sterilized before disposal.
- Dispose consumable materials and unused contents in accordance with applicable regional/national regulatory requirements.

## 6. Reagents supplied

Components	Use Status	Amount
(A) Antibody-coated Plate	Use after washing.	96 wells/1 plate
(B) Mouse/Rat Insulin Standard	Concentrated. Use after dilution.	50 $\mu\text{L}$ /1 bottle
(C) Buffer	Ready to use.	60 mL/1 bottle
(D) Peroxidase-conjugated Antibody Solution	Concentrated. Use after dilution.	100 $\mu\text{L}$ /1 bottle
(E) Luminescent Reagent 1	Mix (E) Luminescent Reagent 1 and (F) Luminescent Reagent 2 in 1 : 1 (vol./vol.)	6 mL/1 bottle
(F) Luminescent Reagent 2		6 mL/1 bottle
(I) Wash Solution (10 $\times$ )	Concentrated. Use after dilution.	100 mL/1 bottle
(J) Plate Seal	–	1 sheet

(D) Peroxidase-conjugated Antibody Solution : Vials contain more than volumes shown in the list. You can easily take out 100  $\mu\text{L}$ , respectively, from vials.

### **【Storage and Stability】**

#### [(A) Antibody-coated Plate]

If seal is not removed, put the strip back in a plastic bag with zip-seal originally used for well-plate container and store at  $2^{\circ}\text{C} - 10^{\circ}\text{C}$ .

#### [(B) Mouse/Rat Insulin Standard (200 ng/mL)]

Standard solutions prepared above should be used as soon as possible, and should not be stored.

#### [(C) Buffer]

If not opened, store at  $2^{\circ}\text{C} - 10^{\circ}\text{C}$ . Once opened, we recommend using as soon as possible to avoid influence by environmental condition.

#### [(D) Peroxidase-conjugated Antibody Solution]

If the kit is divided for multiple assays, prepare these solutions just after taking out of the refrigerator, immediately close the cap tightly, and store at  $2^{\circ}\text{C} - 10^{\circ}\text{C}$ . Unused working solution (already diluted) should be disposed.

#### [(E) Luminescent Reagent 1, (F) Luminescent Reagent 2]

If the kit is divided for multiple assays, prepare these solutions just after taking out of the refrigerator, immediately close the cap tightly, and store at  $2^{\circ}\text{C} - 10^{\circ}\text{C}$ . Discard the remaining mixed luminescent reagent.

#### [(I) Wash Solution (10 $\times$ )]

The rest of undiluted buffer : stored tightly closed at  $2^{\circ}\text{C} - 10^{\circ}\text{C}$ . Dispose any unused diluted buffer.

## 7. Equipments or supplies required but not supplied Use as a check box

Deionized water (or Distilled water).  Test tubes for preparation of standard solution series.

Glassware for dilution of Wash Solution (10 $\times$ ) (a graduated cylinder, a bottle).  Pipettes (disposable tip type). One should be able to deliver 5  $\mu\text{L}$  precisely, and another for 50  $\mu\text{L}$  - 100  $\mu\text{L}$ , 100  $\mu\text{L}$  - 1000  $\mu\text{L}$ .

Syringe-type repeating dispenser like Eppendorf multipette plus which can dispense 50  $\mu\text{L}$ .  Paper towel to remove washing buffer remaining in wells.  A vortex-type mixer.  A shaker for 96 well-plate (500 rpm - 1200 rpm).  An automatic washer for 96 well-plate (if available), or a wash bottle with a jet nozzle.  A 96 well-plate reader for luminescence detection  Software for data analysis.

## 8. Preparation of Samples

This kit is intended to measure insulin in mouse/rat serum, plasma (preferably obtained with EDTA or heparin), culture medium and tissue/cell extracts. The necessary sample volume for the standard procedure is 5  $\mu$ L. Samples should be immediately assayed or stored below -35°C for several days. Defrosted samples should be mixed thoroughly for best results. Hemolytic and hyperlipemic serum samples are not suitable.

\*To avoid influence of blood (high lipid or hemolysis, etc.), if your original samples have heavy chyle or hemolysis, do not use them for assay.

If presence of interfering substance is suspected, examine by dilution test at more than 2 points. Dilution of a sample should be made in a test tube using Buffer prior to adding them to wells. Turbid samples or those containing insoluble materials should be centrifuged before testing to remove any particulate matter.

## Storage and Stability

Insulin in samples will be inactivated if stored at 2°C - 10°C. If it is necessary to store sample in refrigerator (2°C - 10°C), add aprotinin at final concentration of 100 KIU/mL - 500 KIU/mL (KIU : kallikrein inhibitor unit).

If you have to store assay samples for a longer period, snap-freeze samples and keep them below -35°C. Avoid repeated freeze-thaw cycles.

## 9. Preparation of Reagents

◆ Bring the reagents to room temperature (20°C - 25°C) about 2 hours prior to use.

◆ Reagents described as "Use after dilution" in **[6. Reagents supplied]** shall be prepared as follows.

◆ Prepare reagent solutions in appropriate volume for your assay. Do not store the diluted reagents.

### [Concentrated reagents]

#### [(B) Mouse/Rat Insulin Standard (200 ng/mL)]

Make a serial dilution of master standard (200 ng/mL) solution to prepare each standard solution (0.0391 - 20.0 ng/mL).

\*Unit reduction for  $\mu$ IU/mL is 26 IU/mg. (Refer to 2. Introduction.)

Volume of standard solution	Buffer	Concentration (ng/mL)	Concentration ( $\mu$ IU/mL)*
Original solution : 10 $\mu$ L	90 $\mu$ L	20.0	520
20.0 ng/mL solution : 50 $\mu$ L	75 $\mu$ L	8.00	208
8.00 ng/mL solution : 50 $\mu$ L	75 $\mu$ L	3.20	83.2
3.20 ng/mL solution : 50 $\mu$ L	75 $\mu$ L	1.28	33.3
1.28 ng/mL solution : 50 $\mu$ L	110 $\mu$ L	0.400	10.4
0.400 ng/mL solution : 50 $\mu$ L	110 $\mu$ L	0.125	3.25
0.125 ng/mL solution : 50 $\mu$ L	110 $\mu$ L	0.0391	1.02
0 (Blank)	110 $\mu$ L	0	0

#### [(D) Peroxidase-conjugated Antibody Solution]

Prepare working solution by dilution of (D) with the Buffer (C) to **1 : 100**.

6 mL of the diluted solution is enough for 96 wells.

#### [(E) Luminescent Reagent 1] & [(F) Luminescent Reagent 2]

Mix Luminescent Reagent 1 and Luminescent Reagent 2 in **1 : 1 (vol./vol.)** and incubate for 15 - 30 minutes before use.

**Keep in the dark.** Dispose any unused mixed luminescent reagent.

6 mL of the diluted solution is enough for 96 wells.

#### [(I) Wash Solution (10 $\times$ )]

Dilute 1 volume of the Wash Solution (10 $\times$ ) to 10 volume with deionized water (or distilled water) to prepare working solution. Example : 100 mL of concentrated washing buffer (10 $\times$ ) and 900 mL of deionized water (or distilled water).

## 10. Assay Procedure

Remove the cover sheet of the anti-Insulin-coated plate after bringing up to room temperature.

- (1) Wash the (A) Antibody-coated Plate by filling the wells with 300  $\mu$ L of washing buffer and discard 4 times (\*①), then strike the plate upside-down onto several sheets of paper towel to remove residual buffer in the wells.
- (2) Pipette 50  $\mu$ L of (D) Peroxidase-conjugated Antibody Solution to all wells. Shake the plate gently on a plate shaker (\*②).
- (3) Pipette 5  $\mu$ L of standard solution to the wells designated for standards.
- (4) Pipette 5  $\mu$ L of sample to the designated sample wells.
- (5) Shake the plate gently on a plate shaker (\*②).
- (6) Stick a Plate Seal (\*③) on the plate and incubate for 2 hours at room temperature (20°C - 25°C).


- (7) Discard the reaction mixture. Rinse wells by filling the wells with 300  $\mu$ L of washing buffer and discard 4 times (\*①), then strike the plate upside-down onto several sheets of paper towel to remove residual buffer in the wells.
  - (8) Pipette 50  $\mu$ L of mixed luminescent reagent to all wells, and shake for 1 minute using a plate shaker (\*②).
  - (9) After agitation, determine the luminescent intensity using a 96-well microplate reader (for luminescence measurement). It is recommended to perform the measurement 10 to 20 minutes after addition of the luminescent reagent.
- \*Refer to 14. Summary of Assay Procedure for notes of \*①, \*② and \*③.

## 11. Technical Tips

- Be careful to avoid any contamination of assay samples and reagents. We recommend the use of disposal pipette tips, and 1 tip for 1 well.
- The reagents are prepared to give accurate results only when used in combination within the same box. Therefore, do not combine the reagents from kits with different lot numbers. Even if the lot number is the same, it is best not to mix the reagents with those that have been preserved for some period.
- Optimally, the reagent solutions of the kit should be used immediately after reconstitution. Otherwise, store them in a dark place at 2°C - 10°C.
- Time the reaction from the pipetting of the reagent to the first well.
- Dilution of the assay sample must be carried out using the Buffer provided in the kit.
- To avoid denaturation of the coated antibody, do not let the plate go dry.
- As the Antibody-coated Plate is module type of 8 wells  $\times$  12 strips, each strip can be separated by cutting the cover sheet with a knife and used independently.
- When ELISA has to be done under the airstream velocity over 0.4 m/sec. and the humidity less than 30%, seal the well plate with a plate seal and place the well plate in an incubator or a styrofoam box in each step of incubation.

## 12. Calculations

- (1) Prepare a standard curve for each assay. Prepare a standard curve using two-way logarithmic section paper by plotting RLU (Y-axis) against insulin concentration (ng/mL) on X-axis.
  - (2) Using the standard curve, read the insulin concentration of a sample at its RLU, and multiply the assay value by dilution factor if the sample has been diluted. Though the assay range is wide enough, in case the RLU of some samples is higher than that of the highest standard, please repeat the assay after proper dilution of samples with the Buffer.
- We recommend the use of 3rd order regression curve for log-log plot, or 4 or 5 parameters method for log-normal plot in computer calculation.

## 13. Trouble Shooting

- Low RLU in all wells  
Possible explanations :
  - 1) The standard or samples might not be added.
  - 2) Wrong reagents related to luminescence might have been added. Wrong dilution of Peroxidase-conjugated Antibody Solution.
  - 3) Contamination of peroxidase enzyme inhibitor(s).

- 4) Influence of the temperature under which the kits had been stored.
- 5) Excessive hard washing of the well plate.
- 6) Addition of luminescent reagent soon after taking out from a refrigerator might cause poor luminescence owing to low temperature.

• Intense luminescence in all wells including blank

Possible explanations :

- 1) Improper or inadequate washing. (Change washing frequency from 4 times to 5 - 6 times at the constant stroke after the reaction with Peroxidase-conjugated Antibody Solution/sample.)
- 2) Too high incubation temperature. Adjust the temperature to 20°C - 25°C.

• High coefficient of variation (CV)

Possible explanation :

- 1) Improper or inadequate washing.
- 2) Improper mixing of standard or samples.
- 3) Pipetting at irregular intervals.

• Q-1 : Can I divide the plate to use it for the other testing?

A-1 : Yes, cut off the clear seal on the plate with cutter along strip. Put the residual plate, which is still the seal on, in a refrigerator soon


• Q-2 : I found there contains liquid in 96 well-plate when I opened the box. What is it?

A-2 : When we manufacture 96 well-plate, we insert preservation stabilizer in wells.

#### 14. Summary of Assay Procedure : Use as a check box

\* First, read this instruction manual carefully and start your assay after confirmation of details.

- Bring the well-plate and all reagents back to room temperature at 20°C - 25°C for 2 hours.
- Wash Solution (10×) must be diluted to 10 times by deionized water (or distilled water) that returned to 20°C - 25°C.
- Standard mouse/rat insulin solution dilution example :

(D) Peroxidase-conjugated Antibody Solution : Dilute to **100 times** by using (C) Buffer and use.

- Antibody-coated Plate
- ↓ Washing 4 times (\*①) \*④
- Peroxidase-conjugated Antibody Solution 50 µL
- ↓ Shaking (\*②)
- Samples/Standards 5 µL
- ↓ Shaking (\*②), Incubation for 2 hours at 20°C - 25°C. (Standing (\*③))
- Mix (E) Luminescent Reagent 1 and (F) Luminescent Reagent 2 in 1 : 1 (vol./vol.) and incubate for 15-30 minutes before use. Keep in the dark.
- ↓ Washing 4 times (\*①) \*④
- Mixed luminescent reagent 50 µL
- ↓ Shaking for 1 minute (\*②)
- Measurement of luminescence intensity.
- It is recommended to perform the measurement 10 to 20 minutes after addition of the luminescent reagent.

\*① After dispensing wash buffer to wells, lightly shake the plate on your palm for 10 seconds and remove the buffer. Guideline of washing volume : 300  $\mu$ L/well for an automatic washer and for a pipette if the washing buffer is added by pipette. In case of washing by using 8 channel pipette, sometimes the back ground tends to be high. If so, change washing frequency from 4 times to 5 - 6 times at the constant stroke after the reaction with Peroxidase-conjugated Antibody Solution/sample.

Standard of plate-washing pressure : 5 - 25 mL/min. (Adjust it depending on the nozzle's diameter.)

\*② Guideline of shaking : 500 rpm - 1200 rpm 3 times for 10 seconds.


\*③ Seal the plate during the reaction after shaking. Peel off the protective paper from the seal and stick the seal on the plate. Do not reuse the plate seal used once.

\*④ After removal of wash buffer, immediately dispense the next reagent.

#### Worksheet Example

	Strip 1&2	Strip 3&4	Strip 5&6	Strip 7&8	Strip 9&10	Strip 11&12
A	20.0 ng/mL	Sample 1	Sample 9	Sample 17	Sample 25	Sample 33
B	8.00 ng/mL	Sample 2	Sample 10	Sample 18	Sample 26	Sample 34
C	3.20 ng/mL	Sample 3	Sample 11	Sample 19	Sample 27	Sample 35
D	1.28 ng/mL	Sample 4	Sample 12	Sample 20	Sample 28	Sample 36
E	0.400 ng/mL	Sample 5	Sample 13	Sample 21	Sample 29	Sample 37
F	0.125 ng/mL	Sample 6	Sample 14	Sample 22	Sample 30	Sample 38
G	0.0391 ng/mL	Sample 7	Sample 15	Sample 23	Sample 31	Sample 39
H	0 (Blank)	Sample 8	Sample 16	Sample 24	Sample 32	Sample 40

#### Assay Worksheet

#### 15. Storage and Expiration

The complete kit is stored at 2°C - 10°C (Do not freeze). Opened reagents should be used as soon as possible to avoid less than optimal assay performance caused by storage environment.

LBIS™ Mouse/Rat Insulin ELISA Kit (Luminescent)

[Storage] Store at 2°C - 10°C (Do not freeze).

[Expiration date] Indicated on the label.

[Package] For 96 tests

[Cat #] 296-89301

---

## **FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://fwk.fujifilm.co.jp>

### **FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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