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Read this operation manual before using this product.

Research Use Only
For Amino Acid Automated Analysis

Code No. 014-23841


Amino Acid Analysis Reagent (for LC/MS) (APDSTAG®)

【Precaution】

This reagent is for laboratory use only.
Measurement results obtained using reagent cannot be used for diagnostic procedures.

【Measurement principle】

When reaction reagent solution is added to a sample under the alkaline conditions and they are warmed, derivatized compounds of 3-Aminopyridyl-*N*-hydroxysuccinimidyl Carbamate bound to amino group are produced. This solution is separated by liquid chromatography to detect the area value of amino acid derivatized compound per mass-to-charge ratio. Then concentration of each amino acid in the sample is calculated as a ratio of the area value obtained from the sample to that of the standard solution.

【Operation procedures】

(1) Preparation of reagent

Reaction reagent solution

: Add 5 mL of Acetonitrile (not supplied with this product) to 100 mg of this product and mix them to prepare reaction reagent solution. After preparation, store at a temperature of 2 to 10°C and use within 2 weeks.

(2) Preparation of material not supplied with this product

Amino Acids Mixture Standard Solution

: It should be obtained separately. Use it following the operation manual attached to the product.

(3) Measurement operation

< Pretreatment for reagent >

Centrifuge the reaction reagent solution and use the supernatant.

If the supernatant is turbid before use, centrifuge it and use the supernatant.

< Standard operation >

- [1] Add 10 μ L of the sample and 5 μ L of reaction reagent solution to 185 μ L of the APDSTAG® Wako Borate Buffer (not supplied with this product) and mix thoroughly. Then, warm the solution at a temperature of 55 to 60°C for 5 to 15 minutes.
 - [2] Transfer the solution obtained in the step [1] to a column of liquid chromatography. Perform elution in the gradient mode with APDSTAG® Wako Eluent (not supplied with this product) and Acetonitrile (for LC/MS) (not supplied with this product).
 - [3] Feed the substances eluted in the step [2] into mass analyzer and allow area values to be recorded automatically.
- 2) Operate Amino Acids Mixture Standard Solution (not supplied with this product) in the same procedure as described above in 1). Calculate the concentration of each amino acid in the sample as a ratio to the area value of each amino acid obtained from the standard solution. If internal standard method is used, calculate the concentration of each amino acid using a ratio of the area value of internal standard substance to that of the measurement target component, not using area values.

3) Calculation method

Concentration of each amino acid in the sample (μ mol/L) = concentration of amino acid in Amino Acids Mixture Standard Solution (μ mol/L) \times (area value of amino acid in the sample) / (area value of amino acid in Amino Acids Mixture Standard Solution)
(With internal standard method)

Concentration of each amino acid in the sample (μ mol/L) = concentration of amino acid in Amino Acids Mixture Standard Solution (μ mol/L) \times (ratio of area value of amino acid in the sample)^{*1} / (ratio of area value of amino acid in Amino Acids Mixture Standard Solution)^{*2}

* 1 : Area value of amino acid in the sample / area value of internal standard substance

* 2 : Area value of Amino Acids Mixture Standard Solution / area value of internal standard substance

(Note) When these operations are applied to an analysis device, follow the instructions for each device.

【Precautions for use/handling】

< Precaution for handling (hazard prevention) >

When reagent is accidentally splashed into mouth or eyes, or attaches to the skin, flush immediately with plenty of water and seek the medical attention of a doctor if necessary.

< Precautions for use >

- (1) This product easily absorbs moisture. Put in an aluminum bag and store at a temperature of 2 to 10°C. Do not use after expiration date.
- (2) Use as soon as possible after opened. Close the cap and store under the designated conditions if stored.
- (3) Do not use the bottle and attachment of this product for any other purpose.
- (4) APDSTAG® Wako Eluent, Acetonitrile, APDSTAG® Wako Borate Buffer and Amino Acids Mixture Standard Solution are not supplied with this product. Use a separate product.

< Precautions for disposal >

Follow the Waste Disposal and Public Cleaning Act and the standard of water discharge for proper disposal.

【Storage and expiration date】

Storage : Inactive gas included. Store at a temperature of 2 to 10°C in a light-proof container.

Expiration date : Indicated on the label.

【Package】

Code No.	Product Name	Package Size
014-23841	Amino Acid Analysis Reagent (for LC/MS) (APDSTAG®)	100 mg

【References】

- (1) Kazutaka, S., Takashi, O., Akihisa, Y., Kazuo, H. and Hiroshi, M. : *Rapid Commun. Mass. Spectrom.*, **23** (10), 1483-1492 (2009).
- (2) Rafael, Venta : *Clinical Chemistry*, **47** (3), 575-683 (2001).

[Related Products]

Code No.	Product Name	Package Size
010-23061	APDSTAG® Wako Eluent	1 L
018-19853	Acetonitrile (for LC/MS)	3 L
019-23151	APDSTAG® Wako Borate Buffer	1 L
015-27891	Amino Acid Mixture Standard Solution, Type AN [CRM]	1 mL × 5A
011-27871	Amino Acid Mixture Standard Solution, Type B [CRM]	1 mL × 5A
018-27881	Amino Acid Mixture Standard Solution, Type H [CRM]	1 mL × 5A
016-28161	Amino Acids Mixture Standard Solution, Type AN (High Range)	1 mL × 5A
010-28164	Amino Acids Mixture Standard Solution, Type AN (High Range)	5 mL
012-28141	Amino Acids Mixture Standard Solution, Type B (High Range)	1 mL × 5A
016-28144	Amino Acids Mixture Standard Solution, Type B (High Range)	5 mL
019-28151	Amino Acids Mixture Standard Solution, Type H (High Range)	1 mL × 5A
013-28154	Amino Acids Mixture Standard Solution, Type H (High Range)	5 mL

Manufactured by FUJIFILM Healthcare Manufacturing Corporation

2613-2, Ogohara, Komono-Cho, Mie-Gun, Mie 510-1222, Japan

Distributed by FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2407KA4