Code No. 012-25841 (2 mL, Net volume 1 mL)
018-25843 (10 mL, Net volume 5 mL)

For immunochemistry

Anti PA tag Antibody Beads

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PA tag is consisted of 21 amino acids (GVAMPGAEDDVV) and has been used as the tag for detection and/or purification of PA tag fusion recombinant proteins.

Anti PA tag Antibody Beads is the slurry of agarose beads immobilized Anti PA tag rat monoclonal antibody. This product is used for the purification of PA tag fusion proteins by immunoprecipitation method and able to purify PA tag fusion proteins with neutral pH condition by the competitive elution of PA tag Peptide. In addition, This product is easily regenerated by using PA tag Washing Solution (Code No. 169-27261) on neutral pH condition.

[Formulation]

1 × PBS aqueous solution with 0.05w/v% sodium azide, pH 7.2.

[Immobilized antibody]

Anti PA tag, Rat Monoclonal Antibody (NZ-1)

[Immobilized antibody subclass]

IgG_{2a}

[Antibody bonding amount]

Indicated on the label.

[Antigen binding capacity]

0.015-1.0 mg protein/1 mL beads

The binding capacity depends on the target protein and the experiment.

[Usage]

10mL conditioned culture medium/50-100 μL beads

200-1,000 μL cell lysate/20 μL beads

Please examine appropriate usage according to the target protein and cell lines etc.

[Additional required reagents and materials]

Centrifuge, Rotator, Vortex mixer, Tube mixer, Centrifuge tube, 1.5 mL microcentrifuge tube, Cell lysis buffer (see below), PBS(-), TBS, Protease inhibitors, Phosphatase inhibitors

Recommended cell lysis buffer

●RIPA Buffer (Code No. 182-02451)

50 mmol/L Tris-HCl, pH 8.0, 150 mmol/L Sodium Chloride, 0.5w/v% Sodium Deoxycholate, 0.1w/v% Sodium Dodecyl Sulfate, 1.0w/v% NP-40 substitute.

●Cell Lysis Buffer M (Code No. 038-21141)

20 mmol/L Tris-HCl, pH 7.4, 200 mmol/L Sodium Chloride, 2.5 mmol/L Magnesium Chloride, 0.05w/v% NP-40 substitute.

In the case of adding a protease inhibitor etc. to a cell lysis buffer, please prepare it according to its instruction manual.

[Protocol]

The immunoprecipitation of PA tag fusion proteins and the regeneration of Anti PA tag Antibody Beads

1. The preparation of samples

1-A In the case of adherent cells

- 1) Culture the cell lines expressing PA tag fusion proteins ($1 \times 10^6 \sim 1 \times 10^7$ cells).^{※1}
- 2) Remove the medium from culture dishes and wash cultured cells twice with PBS(-).
- 3) Detach the cells from dishes using a cell scraper or trypsin-EDTA solution and transfer the cell suspension into a new 1.5 mL centrifugation tube. In the case of trypsin-EDTA solution, add the medium including serum to the cell suspension immediately to avoid the damage due to trypsin treatment.
- 4) Centrifuge the cell suspension at $200 \times g$ for 5 minutes at 4 °C.
- 5) Remove the supernatant.
- 6) After adding 1 mL of PBS(-) and suspending the cell pellet by pipetting, centrifuge the cell suspension at $200 \times g$ for 5 minutes at 4 °C.
- 7) Remove the supernatant.
- 8) Repeat 6)-7) steps.
- 9) Add 1 mL of cell lysis solution to the cell pellet and suspend by pipetting or vortex mixer.
- 10) Incubate on ice for 10 minutes.
- 11) Centrifuge the cell lysate at $20,000 \times g$ for 10~20 minutes at 4 °C.^{※2}
- 12) Transfer the supernatant into a new 1.5 mL microcentrifuge tube. ... **Cell lysate**

1-B In the case of suspension cells

- 1) Culture the cell lines expressing PA tag fusion proteins ($1 \times 10^6 \sim 1 \times 10^7$ cells).^{※1}
- 2) Centrifuge the culture solution at $200 \times g$ for 5 minutes at 4 °C.
- 3) Remove the supernatant.
- 4) After adding 1 mL of PBS(-) and suspending the cell pellet by pipetting, transfer the cell suspension into a new 1.5 mL microcentrifuge tube.
- 5) Centrifuge the cell suspension at $200 \times g$ for 5 minutes at 4 °C.
- 6) Remove the supernatant.
- 7) After adding 1 mL of PBS(-) and suspending the cell pellet, centrifuge the cell suspension at $200 \times g$ for 5 minutes at 4 °C.
- 8) Remove the supernatant.
- 9) Add 1 mL of cell lysis solution to the cell pellet and suspend by pipetting or vortex mixer.
- 10) Incubate on ice for 10 minutes.
- 11) Centrifuge the cell lysate at $20,000 \times g$ for 10~20 minutes at 4 °C.^{※2}
- 12) Transfer the supernatant into a new 1.5 mL microcentrifuge tube. ... **Cell lysate**

1-C In the case of conditioned culture medium

- 1) Culture the cell lines expressing PA tag fusion proteins ($1 \times 10^6 \sim 1 \times 10^7$ cells).^{※1}
- 2) Recover the cell culture medium into a new centrifuge tube and centrifuge at $2,000 \times g$ for 10 minutes at 4 °C for cell and cell debris removal.^{※3}

- 3) Transfer the supernatant into a new centrifuge tube and adjust its pH to 7.5 with 1 mol/L Tris. ... **Conditioned culture medium**

※1 Please examine some appropriate culture conditions according to cell lines.
 ※2 If a white floating layer appeared in supernatant after centrifugation, please filtrate it using a filter (pore size : 0.45 μ m).
 ※3 Cells etc. are removable by using a filter (pore size : 0.22 μ m etc.).

2. Prewashing of Antibody Beads

- 1) Suspend Anti PA tag Antibody Beads adequately by vortex mixer.
- 2) Transfer 40-200 μ L (net beads volume 20-100 μ L) of Anti PA tag Antibody Beads into a new 1.5 mL microcentrifuge tube.※4
- 3) Centrifuge at 5,000-8,000 $\times$ g for 30 seconds at 4 $^{\circ}$ C.
- 4) Incubate for 1-2 minutes on ice.
- 5) Remove the supernatant carefully not to aspirate beads.
- 6) Add 0.5-1.0 mL of chilled TBS or chilled PBS(-) and suspend by vortex mixer.※5
- 7) Centrifuge at 5,000-8,000 $\times$ g for 30 seconds at 4 $^{\circ}$ C and remove the supernatant carefully not to aspirate beads.
- 8) Repeat 6)-7) step. ... **Prewashed beads**

※4 Please examine an appropriate beads quantity and centrifugation tube for your purpose.
 ※5 Please choose an appropriate washing buffer for your purpose.

3. Immunoprecipitation

3-A In the case of cell lysate

- 1) Add 200-1,000 μ L of **Cell lysate** into **Prewashed beads**.※6
- 2) Incubate with rotation for more than 2 hours at 4 $^{\circ}$ C.
- 3) Centrifuge at 5,000-8,000 $\times$ g for 30 seconds at 4 $^{\circ}$ C.
- 4) Remove the supernatant carefully not to aspirate beads.
- 5) Add 0.5-1.0 mL of chilled TBS or chilled PBS(-) and suspend by vortex mixer.
- 6) Centrifuge at 5,000-8,000 $\times$ g for 30 seconds at 4 $^{\circ}$ C and remove the supernatant carefully not to aspirate beads.
- 7) Repeat 5)-6) step. ... **Antigen binding beads**

※6 Fill up to 1 mL by cell lysis buffer as needed.

3-B In the case of conditioned culture medium

- 1) Add 10 mL of **Conditioned culture medium** into **Prewashed beads**.※7
- 2) Incubate with rotation for more than 2 hours at 4 $^{\circ}$ C.
- 3) Centrifuge at 5,000-8,000 $\times$ g for 1 minute at 4 $^{\circ}$ C.
- 4) Remove the supernatant carefully not to aspirate beads.
- 5) Transfer the beads into a new 1.5 mL microcentrifuge tube using 1.0 mL of chilled TBS or chilled PBS(-).※8
- 6) Centrifuge at 5,000-8,000 $\times$ g for 30 seconds at 4 $^{\circ}$ C and remove the supernatant carefully not to aspirate beads.
- 7) Add 0.5-1.0 mL of chilled TBS or chilled PBS(-) and suspend by vortex mixer.
- 8) Centrifuge at 5,000-8,000 $\times$ g for 30 seconds at 4 $^{\circ}$ C and remove the supernatant carefully not to aspirate beads.
- 9) Repeat 7)-8) step 1-3 times. ... **Antigen binding beads**

※7 The use of 50-100 μ L beads per 10 mL conditioned culture medium is recommended.
 ※8 The yield of PA tag fusion proteins may increase by reuse of a conditioned culture medium removed at step 4).

4. The elution of PA tag fusion proteins

Please choose an appropriate elution method depending on protein character and downstream experiment.

4-A The elution by PA tag Peptide (non-denaturing condition)

- 1) Prepare the PA tag Peptide stock solution (final concentration 5mg/mL). Dissolve 5 mg of PA tag Peptide (Code No. 167-25501) in 1 mL of TBS (10 mmol/L Tris-HCl, pH 7.5, 150 mmol/L Sodium Chloride).※9
- 2) Prepare the PA tag Peptide working solution (final concentration 100-200 μ g/mL) by using the PA tag Peptide stock solution.
- 3) Add above PA tag Peptide working solution of 1-5 times the amount of **Antigen binding beads** into the 1.5 mL microcentrifuge tube containing **Antigen binding beads**.
- 4) Incubate with rotation for 30 minutes at 4 $^{\circ}$ C.※10
- 5) Centrifuge at 5,000-8,000 $\times$ g for 30 seconds at 4 $^{\circ}$ C and transfer the supernatant into a new 1.5 mL microcentrifuge tube carefully not to aspirate beads.
- 6) Store the collected supernatant on ice.※11

4-B The elution by 0.1 mol/L Glycine-HCl, pH 3.0 (acidic condition)

- 1) Add 0.1 mol/L Glycine-HCl, pH 3.0 of 1-5 times the amount of **Antigen binding beads** into the 1.5 mL microcentrifuge tube containing **Antigen binding beads**.
- 2) Incubate with rotation for 5 minutes at room temperature.
- 3) Centrifuge at 5,000-8,000 $\times$ g for 30 seconds at room temperature and transfer the supernatant into a new 1.5 mL microcentrifuge tube carefully not to aspirate beads.
- 4) Neutralize by neutralization buffer. (Example : 1 mol/L Tris-HCl, pH 9.0)
- 5) Store the neutralized supernatants on ice.

4-C The elution by SDS sample buffer (denaturing condition)

In the case of using 2 $\times$ SDS sample buffer including reducing reagents such as 2-mercaptoethanol or DTT, heavy chain (50kDa) and light chain (25kDa) of Anti PA tag Rat Monoclonal Antibody is denatured.

- 1) Add 40-200 μ L of 2 $\times$ SDS sample buffer (0.125 mol/L Tris-HCl, pH 6.8, 4w/v% SDS, 20w/v % Glycerol, 10% 2-mercaptoethanol or 100 mmol/L DTT, 0.01% Bromophenol Blue) into a 1.5 mL microcentrifuge tube containing **Antigen binding beads** and suspend by vortex mixer.
- 2) Boil for 5 minutes at 95 $^{\circ}$ C.
- 3) Centrifuge at 5,000-8,000 $\times$ g for 30 seconds at room temperature and transfer the supernatant into a new 1.5 mL microcentrifuge tube.
- 4) Store the collected supernatants on ice.

- ※9 Please store the used PA tag Peptide stock solution at -20°C , and use it as promptly as possible.
- ※10 In order to increase the elution efficiency of PA tag fusion proteins, please incubate with rotation for more than 5 minutes.
- ※11 Please remove the peptide from the solution by dialysis or gel filtration as needed.

5. The regeneration of the used Anti PA tag Antibody Beads by PA tag Washing Solution (Code No. 169-27261).

- 1) Add 10 times the volume of PA tag Washing Solution (Code No. 169-27261) into the used Anti PA tag Antibody Beads.
- 2) Incubate with rotation for 10 minutes at room temperature.
- 3) Centrifuge at $5,000-8,000 \times g$ for 30 seconds at room temperature and remove the supernatant carefully not to aspirate beads.
- 4) Repeat 1)-3) steps 3 times.
- 5) Add 10 times the volume of TBS or PBS(-) and suspend by vortex mixer.
- 6) Centrifuge at $5,000-8,000 \times g$ for 30 seconds at room temperature and remove the supernatant carefully not to aspirate beads.
- 7) Repeat 5)-6) step 3 times.
- 8) Add 2 times the volume of TBS or PBS(-).
- 9) Store $2-10^{\circ}\text{C}$.

[Storage]
 $2-10^{\circ}\text{C}$

[Reference]
 Fujii, Y. *et al.*, "PA tag : A versatile protein tagging system using a super high affinity antibody against a dodecapeptide derived from human podoplanin.", *Protein Exp. Purif.*, **95**, 240-247 (2014).

[Related products for Anti PA tag Antibody Beads]

Code No.	Product Name	Package
169-27261	PA tag Washing Solution	50 mL
167-25501	PA tag Peptide	5 mg
163-25503		25 mg
182-02451	RIPA Buffer	100 mL
188-02453		250 mL
038-21141	Cell Lysis Buffer M	100 mL
034-21143		250 mL

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
 Telephone : + 81-6-6203-3741
 Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
 Richmond, VA 23237
 U.S.A.
 Telephone : + 1-804-271-7677
 Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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