

FUJIFILM

Wako

Code No. 012-29241 ( 2,000,000 units)  
018-29243 (20,000,000 units)

**Achromopeptidase™, Crude,  
Lytic Enzyme**  
(TBL-1)

**アクロモペプチダーゼ™, 粗製品, 溶菌酵素**  
(TBL-1)

Bacteriolytic enzymes are hydrolases that act on bacterial cell wall and cause cell-lysis. The most investigated bacteriolytic enzyme has been lysozyme. Lysozyme is now produced on an industrial scale and has found many applications. Lysozyme, however, is bactericidal only against some Gram-positive organisms (e. g., *Micrococci*, *Bacilli*, *Sarcinae*, *Pediococci*, *Kurthiae*) and not against pathological organisms such as *Staphylococci*, *Streptococci* and *Clostridia*. Many new bacteriolytic enzymes have been described and several originate from microorganisms. The soil bacterium first isolated in 1966 has been found to produce a highly active bacteriolytic enzyme with a broad specificity. Like lysozyme, the enzyme causes lysis of *Micrococci*, *Bacilli* and *Sarcinate* but unlike lysozyme will also cause lysis of *Staphylococci*, *Streptococci*, *Clostridia* and *Leuconostoc* and, in addition, is active against some Gram-negative organisms. FUJIFILM Wako has investigated the properties of the bacteriolytic enzyme and has succeeded in producing it on a large scale. This bacteriolytic enzyme is now available as Achromopeptidase (TBL-1) .

Source : *Lysobacter enzymogenes*

Activity : indicated on the label

Molecular weight : about 20,000

Stability at various pH values : Most stable at pH 6.0 and virtually stable at  
pH 5.5 ~ 10

Appearance : Pale yellowish brown ~ brown, powder

Heat stability : incubation in 50 mmol/L phosphate buffer at pH 6.0 for 10 minutes

Temperature	Activity
Less than 40°C	Stable
50°C	95% remaining
60°C	75% remaining
70°C	Lost

Optimal pH : 7.5 ~ 8.5

Optimal ionic strength : 10 mmol/L sodium chloride (in 10 mmol/L Tris-HCl  
buffer solution (pH 8.0))

### Bacteriolytic Spectrum

Microorganism	Degree of lysis*	
	TBL-1 (300 U/mL)**	Egg lysozyme (1,000 U/mL)
<i>Micrococcus caseolyticus</i>	+++	-
<i>Micrococcus luteus</i>	+++	-
<i>Pediococcus acidilactici</i>	+++	+
<i>Staphylococcus aureus</i>	+++	±
<i>Microbacterium arborescens</i>	++	+
<i>Bacillus alvei</i>	+++	±
<i>Clostridium acetobutylicum</i>	+++	-
<i>Brevibacterium leucinophagum</i>	+	±
<i>Lactobacillus sakei</i>	+++	+
<i>Enterococcus faecalis</i>	+++	±
<i>Leuconostoc dextranicum</i>	++	±
<i>Achromobacter liquidum</i>	+	-
<i>Nocardioides simplex</i>	+	+
<i>Beijerinckia indica</i>	+	±
<i>Mycobacterium diernhoferi</i>	±	±
<i>Kurthia zopfii</i>	++	+

\*determined by decrease of the optical density at 600 nm in each microorganism suspension under the standard assay conditions incubated for indicated times.

+++ ; ≥ 80% decrease within 10 min and ++ ; ≥ 80%, + ; 80-40, ± ; 40-10%, - ; 10-0% decrease within 60 min.

\*\*The activity value of the conventional product (014-09661).

### Assay Method

#### 1) Reagents

A : 0.01 mol/L Tris-HCl + 0.01 mol/L NaCl buffer solution (pH 8.0)

B : Substrate

Suspend 31-32 mg of lyophilized *Micrococcus lysodeikticus* IFO 3333 in Reagent A to bring the volume up to 30 mL. Transfer 5 mL of the solution and dilute with Reagent A to bring the volume up to 25 mL. Place it on ice for 1 hour.

C : Enzyme solution

Dissolve 10 mg of TBL-1 with Reagent A to bring the volume up to 20 mL. Place it on ice for 30 minutes with stirring occasionally.

### Procedure :

Reagent	Test	Blank
B	3.0 mL	3.0 mL
	Pre-incubate in 10 mm absorption cell (37°C 5 min.) [Trans] * 10mm 吸収セルに入れ、37°Cで5分間予備加温する。	
C	0.05 mL	-
A	-	0.05 mL

Immediately, record the absorbance of 600 nm at 37°C for 1 ~ 6 minutes after addition of Reagent C and Reagent A with deionized water as a reference.

[Trans] \* 直ちに、37°Cで波長600nmにおける吸光度を水を対照液として試液C及び試液A添加後1~6分間自記記録する。

### Enzymatic activity

One unit of enzyme is defined as that which when present in 1 mL of reaction solution causes turbidity to decrease at the rate of 0.001 per minute under the above conditions.

Calculation formula for determining activity.

$$\text{Units/mg} = \frac{\Delta E - \Delta E_0}{0.001} \times 3.05 \times \frac{1}{S} \times \frac{20}{0.05} \times f$$

$\Delta E$ : The change of absorbance of Test for 1 minute.  
 $\Delta E_0$ : The change of absorbance of Blank for 1 minute.  
S: Amount of TBL-1 (mg)  
f: 0.9 (correction factor)

**Trans**\*

### 活性単位の計算式

$$\text{Units/mg} = \frac{\Delta E - \Delta E_0}{0.001} \times 3.05 \times \frac{1}{S} \times \frac{20}{0.05} \times f$$

$\Delta E$ : 本試験の1分間当りの吸光度変化  
 $\Delta E_0$ : 空試験の1分間当りの吸光度変化  
S: 量りとした本品の質量 (mg)  
f: 0.9 (補正係数)

### Cautions in Use

- (1) The optimal pH value for TBL-1 is between 7.5 and 8.5. It should be used in this range.
- (2) TBL-1 is greatly influenced by ionic strength and its activity is markedly decreased with increasing ionic strength. In 10 mmol/L Tris-HCl buffer solution (pH 8.0) containing 10 mmol/L sodium chloride the enzyme is optimally active.
- (3) Aqueous solutions of TBL-1 at concentrations up to 0.4% (w/v) remain transparent. At higher concentrations, aqueous solutions of TBL-1 become turbid due to impurities. The turbidity can be cleared by centrifugation. If large amounts of insoluble substances are noted, do not use the enzyme as it may have deteriorated.
- (4) TBL-1 is not purified to homogeneity and therefore may contain enzymes other than TBL-1 and this should be taken into consideration when this product is used. Purified TBL-1 (Achromopeptidase<sup>®</sup>, Purified, Lytic Enzyme (TBL-1), 019-28531) is available for use when required.

### 【Storage】

2 ~ 10°C

### 【Package】

Code No.	Packaging
012-29241	2,000,000 units
018-29243	20,000,000 units

**[References]**

- 1) Isono, M., Takahashi, K. and Yamazaki, Y., : Patent Notice, Sho 46-42953.
- 2) Nakamura, K., Okazawa, Y., Soejima, M. and Masaki, T. : *Agri. Biol. Chem.*, **37**, 2667 (1973).
- 3) Horinouchi, S., Uozumi, T., Beppu, T. and Arima, K. : *ibid.*, **41**, 2487 (1977).
- 4) Kotani, S. : *Protein, Nucleic Acid, Enzyme*, **13**, 1136 (1968).
- 5) Kotani, S. : *Protein, Nucleic Acid, Enzyme*, **14**, 38 (1969).

\* : **Trans** is the Japanese translation.

---

**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : +81-6-6203-3741  
Facsimile : +81-6-6201-5964  
<http://fwvk.fujifilm.co.jp>

Group Companies



Distributors

