Code No. 029-19161 (1 mL × 10)

Brain Tissue Freezing Medium

Brain Tissue Freezing Medium is appropriate for cryopreservation of rat or mouse embryonic brain tissue. Brain tissues may simply be cryopreserved by thawing the product and placing the tissue in the media. When culturing the neurons derived from cryopreserved tissues, we recommend that you use Neuron Dissociation Solutions (Code No. 291-78001, 297-78101) and Neuron Culture Medium (Code No. 148-09671).

[Component]

Component	Package
Brain Tissue Freezing Medium	1 mL × 10

[Reagent volume]

Animal	Kinds of tissue	Amount of tissue	Volume
Rat	Cerebral Hemisphere (embryonic day 17)	1 embryo	1 mL
	Cerebral Cortex (embryonic day 17)	2 embryos	
	Hippocampus (embryonic day 19)	2.5 embryos	
	Cerebral Striatum (embryonic day 17)	2 embryos	
Mouse	Cerebral Hemisphere (embryonic day 15)	1 embryo	
	Cerebral Cortex (embryonic day 15)	2 embryos	
	Hippocampus (embryonic day 16)	2.5 embryos	
	Cerebral Striatum (embryonic day 15)	2 embryos	

[Procedures]

< Precautions for use >

Please shorten the time from the dissection to Primary-cooling. The approximate time is within one hour. The cell viability is decreased as time proceeds after the dissection. Then the dead cells are removed when you isolate cells with Neuron Dissociation Solutions. Therefore the isolated cells are decreased. Please perform the procedures on ice.

< Preparation >

1. Fill a cooler (with a lid) with crushed ice.
2. Thaw this product at room temperature or 4 °C and store this product at 4 °C. Mix the medium before use.

< Dissection >

1. Narcotize the animals.
2. Sterilize the whole body (from whiskers to the tail) of narcotized animals with 70% ethanol.

3. Bring the animals in the clean booth with the stainless steel tray laid out by waste.
4. Take the embryos out from the animals and transfer them to the dish filled with L-15 solution. Cool them on ice until the next procedure.
※L-15 solution : Leibovitz's L-15 medium with 0.05 mg/mL gentamicin sulfate.
5. Transfer the embryo to the cover of the dish and pick up the brain tissues from the embryo quickly.
6. Transfer the brain tissues to the new dish filled with L-15 solution and cool them on ice.

< Primary-Cooling >

1. Place the Brain Tissue Freezing Medium and cooler near the clean bench.
2. In the clean bench, cover the brain tissues with the media and close the cap tightly.
3. Confirm that the tissues are immersed in this product. (Check the tissues from below. They should be easily visible.)
4. Put the tubes in the cooler so the tubes are in the crushed ice with only the lids showing.
5. Place the lid on the cooler and let the tubes sit for 90 minutes.

< Secondary-Cooling >

1. After Primary-Cooling, transfer the tubes from the ice to a freezer box.
2. Immediately place the freezer box in a -80 °C freezer. It is best to store directly on the bottom of a rack inside the freezer.
※Finish this procedure as quickly as possible to minimize temperature change.
※The tubes should be handled by the lids. Avoid touching the body of the tubes.
3. Store the tubes for 3 hours.
※Avoid opening and closing the freezer's door often during this time.

< Storage >

1. After Secondary-Cooling, transfer the tubes to a liquid nitrogen storage container from the -80 °C freezer as quickly as possible.
※The tubes should be handled by the lids. Avoid touching the body of the tubes.

< Thawing and Washing >

1. Remove the tubes from the liquid nitrogen storage container.
2. Thaw at room temperature. If the thawing time is shorter than 30 min, neurons may not thaw enough.
3. Remove Brain Tissue Freezing Medium from the tubes about 0.8 mL.
4. Add 1 mL of Washing Solution (Leibovitz's L-15 Medium or HBSS (+)) to the tubes, then transfer the thawed tissues and Washing Solution to a 15 mL centrifuge tube.
5. Add 4 mL of Washing Solution to the centrifuge tube.

- Let stand at room temperature for 5 min, and centrifuge at 300 rpm at room temperature for 1 min.
- Remove the supernatant, and add 5 mL of Washing Solution.
- Let stand at room temperature for 5 min, and centrifuge at 300 rpm at room temperature for 1 min.
- Dissociate and isolate neurons with Neuron Dissociation Solutions and culture them in Neuron Culture Medium.

[Storage]

Store at -80°C .

Avoid repeating freeze-thaw.

[Related Reagents]

Code No.	Product Name
291-78001	Neuron Dissociation Solutions
297-78101	Neuron Dissociation Solutions S
148-09671	Neuron Culture Medium
128-06075	Leibovitz's L-15 Medium with L-Glutamine, Phenol Red and Sodium Pyruvate
082-09365	HBSS(+) with Phenol Red

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