FUJIFILM

Wako

〈For Research Use Only〉 Code No. 031-19781 (5 mL)

for Genetic Research

Co-Agarose

Co-Agarose is designed for affinity purification of 6 × His fusion protein. Cobalt is used as cation for this type of purification because it provides effective binding and selectivity for a wide range of proteins. Iminodiacetic acid (IDA) is used as ligand.

【Features】

- One step purification
- High yield of purified protein
- High selectivity
- Hydrophilic linkage- low non specific binding

【Agarose】

Cross-linked 6% agarose

【Binding Capacity】

20-40 $\mu\text{mol Co}^{2+}$ /mL gel

【Ligand】

Iminodiacetic Acid (IDA)

【Bead size】

40-180 μm

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50% suspension in 20% ethanol

【General purification method】

Linkage condition :

The linking buffer generally used is acetate (50 mmol/L) or phosphate (10-150 mmol/L). The linking pH comes close to neutrality (normally pH 7.0-8.0) but can vary in the following range : pH 5.5-8.5. Generally, to avoid the effects of ionic interchange, 0.15-0.5 mol/L NaCl is added. It is important to avoid agents such as EDTA or citrate in this buffer.

Elution condition :

- a. Addition of a competitive ligand (generally, imidazole) allows the elution of the retained proteins. In general, 0.5 mol/L imidazole is sufficient to desorb a particular protein. It is also appropriate to use concentration gradients of this agent (0-0.5 mol/L). Most of the proteins are desorbed by concentrations around 250 mmol/L. Other agents

that may be used as competitive ligands are histidine and ammonium chloride.

- b. The reduction of pH (with or without gradient) also permits the elution of the protein of interest (normally pH 3.0-4.0).

A more drastic method uses agents such as EDTA or EGTA (0.05 mol/L) which cause the desorption of the protein as well as of the chelating metal.

【Storage】

2 ~ 10°C

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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