Code No. 031-25151 (50 mL)

For Plant Optical Clearing Reagent

ClearSee™**【Introduction】**

Dr. Daisuke Kurihara *et al.*¹⁾ developed the clearing solution, termed ClearSee™, to diminish chlorophyll autofluorescence while maintaining fluorescent protein stability and tissue structures of plant samples. ClearSee™ is applicable to multicolor deep imaging of various plant samples without sectioning.

【Features】

- Procedures are only 3 steps.
- Every fluorescence microscopes can be used.
- ClearSee™-treated samples can be store for long term (at least 6 months).
- The ClearSee™-treated samples can be observed for many times.

【Additional required materials】**Reagents :**

- 1) Ultra pure water
- 2) Paraformaldehyde (Code No.160-16061)
- 3) 2 mol/L Sodium Hydroxide Solution (Code No.194-05631)
- 4) 10 × PBS(-) (Code No. 163-25265)
- 5) 1 × PBS(-) (Code No. 164-25511)

(Equipment)

- 1) Autoclave
- 2) Draft chamber
- 3) Vacuum pump
- 4) Desiccator
- 5) Incubator for 60 °C
- 6) Fluorescence microscopes
- 7) Conical tubes (50 mL)
- 8) Microtubes (1.5 mL)
- 9) Parafilm
- 10) Needle
- 11) Forcep
- 12) Cover glass
- 13) Slide glass
- 14) Vaseline

【Procedures】**· Fixative solution**

1. Add 8 mL of sterilized water in a conical tube.
2. Add 0.4 g of paraformaldehyde (in the draft chamber).
3. After adding 50 μ L of 2 mol/L Sodium Hydroxide Solution, close and seal a tube with parafilm (in the draft chamber).
4. Incubate at 60 °C until clearing with occasionally inversion.
5. After dissolving and cooling, add 1 mL of 10 × PBS to adjust pH (in the draft chamber).
6. Adjust 10 mL by sterilized water (in the draft chamber). (Preparation freshly is better.)

(Fixation)

1. Put plant samples in a microtube and add 1.3 mL of fixative solution (in the draft chamber).


Fig.1 Arabidopsis leaves and seedling in fixative solution

< Notes >

The sample volume should not exceed 20% of the volume of fixative solution. If the sample volume is larger, prepare large centrifuge tube and increase fixative solution. Because the samples are floated, soak the samples into fixative solution as possible.

2. To prevent the escape of samples, seal the microtube with parafilm and open holes by a needle.


Fig.2 Open holes on the parafilm

- Put the samples in the desiccator and turn on the vacuum pump. After closing a desiccator and turn off the vacuum pump, fix for 30 minutes at room temperature.

< Notes >

Vacuum degree depends on plant species and tissue type (eg. -90 kPa for Arabidopsis). Strong vacuum induces the damage of the samples. Adjust vacuum degree to appear slowly bubbles from the samples.

Fig.3 Bubbles from the samples in a microtube

- Slowly open the desiccator and turn on again the vacuum pump. After closing a desiccator and turn off the vacuum pump, fix for 30 minutes at room temperature.

< Notes >

Take care to prevent the damage of samples by open the desiccator. Two times of vacuum treatment enhance the penetration of fixative solution into the samples. When fixative solution penetrates into the samples, the samples are soaked after brief centrifugation (Fig.4). If the samples are floated even after centrifugation, additional vacuum treatment are required.

Fig.4 Arabidopsis leaves and seedling after fixative treatment

(Wash)

- Slowly open the desiccator. After removing the fixative solution, add 1 mL of 1 × PBS and store for 1 minute (in the draft chamber).
- After removing 1 × PBS, add new 1 mL of 1 × PBS and store for 1 minute.

(Clearing)

- After removing 1 × PBS, add 1.3 mL of ClearSee™.
- Seal the microtube with parafilm and open holes by a needle. Put the samples in the desiccator and turn on the vacuum pump. After closing a desiccator and turn off the vacuum pump, store for 1 hour at room temperature.


Fig.5 Arabidopsis leaves and seedling in ClearSee™ solution

- Slowly open the desiccator. After removing the parafilm and close the microtube, store at room temperature in the dark for clearing.

< Notes >

Invert the microtube every 1~2 days to diffuse the chlorophyll. Required incubation times depends on plant species and tissue type. For examples, the incubation times for clearing are 1~2 days for roots, 4~7 days for leaves and seedlings, 2 weeks for pistils, and 4 weeks for mature tissues.


Fig.6 Chlorophyll are observed around the samples (1~2 days incubation)

Fig.7 Invert the tube

Fig.8 Substitute new ClearSee™


Fig.9 Cleared Arabidopsis leaves and seedling (3~5 days incubation)

(Observation)

10. Put vaseline on a slide glass.
11. Transfer the ClearSee™-treated samples into ClearSee™ on a slide glass.
12. Cover with a cover glass and seal with the cover glass with vaseline.
13. Observe the ClearSee™-treated samples under a fluorescence microscope.


Fig.10 The ClearSee™-treated samples are cleared under a bright field. Autofluorescence of chlorophyll are diminished. The fluorescent proteins are detected even in the vascular bundles. Scale bar : 200 μ m.

< Notes >

Take care to prevent the damage of samples by transfer the samples. To prevent the evaporation of ClearSee™ by vaseline because of deposition. The samples can be stored into ClearSee™ in a microtube after observation.

【Points】

(Fixation)

Fixation step is important for clearing without the damage of samples. To investigate the fixation conditions of plant samples for the first time, observe the fluorescent proteins in the samples after fixation.

(Various plant species)

ClearSee™ is applicable to various plant species other than Arabidopsis. Fig.11 shows ClearSee™-treated leaves of torenia, petunia, tobacco, tomato, cucumber, and rice after 6-days incubation. In the case of rice, removal of cuticular wax are required for penetration of solution. Immerse the rice samples into organic solvents such as chloroform for 10~30 seconds and then fix with a fixative solution.


Fig.11 ClearSee™-treated leaves of various plant species (6-days incubation)

(Fluorescent dye staining)

ClearSee™ is compatible with staining by fluorescent dyes. Add Hoechst 33342 or DAPI to a final concentration of 10 μ g/mL for nuclear staining, or Calcofluor White to a final concentration of 1 mg/mL for cell wall staining for 1 hour. After washing with new ClearSee™ solution for 1 hour, observe the samples under a fluorescence microscope.

【Caution】

Waste each solution with compliance to regulations of institutions and governments.

【Storage】

Stored in the dark.

【Package】

50 mL

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Code No.	Description	Size
160-16061	Paraformaldehyde	100 g
194-05631	2 mol/L Sodium Hydroxide Solution	100 mL
163-25265	10 $\times$ PBS(-)	500 mL
164-25511	1 $\times$ PBS(-)	5 L

【Reference】

- 1) Kurihara, D. *et al.* : *Development*, **142**, 4168 (2015).

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
 Telephone : + 81-6-6203-3741
 Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
 Richmond, VA 23237
 U.S.A.
 Telephone : + 1-804-271-7677
 Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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