FUJIFILM

Wako

Code No. 053-09451

## EV-Up™ EV Production Basal Medium for MSC, AF

EV-Up™ EV Production Basal Medium for MSC, AF is a serum-free and animal component-free medium for the production of extracellular vesicles (EV) with high bioactivity while maintaining the viability of mesenchymal stem cells (MSC). This medium is available for MSCs derived from various tissue sources (bone marrow, umbilical cord, adipose tissue etc.) regardless of the type of growth medium.

To ensure the best performance, mix this medium well with EV-Up™ MSC EV Production Supplement, AF (code No. 298-84001)\*.

This product is for laboratory research use only.

\* EV-Up™ MSC EV Production Supplement, AF contains MSC supplement A and B.

### [Procedures]

(Preparation of complete medium)

- Thaw Supplement A and B at 2-10°C or room temperature.  
*NOTE : Do not thaw the supplements at 37°C . Supplement A may develop a white precipitate after thawing. Please use after completely dissolving the precipitate with a vortex mixer etc.*
- Add all of thawed Supplement A and B to EV-Up™ MSC EV Production basal medium and mix gently to prepare complete medium (EV-Up™ MSC EV Production medium).  
EV-Up™ EV Production Basal Medium for MSC, AF should be used immediately.
- Bring the EV-Up™ MSC EV Production medium to room temperature before use.  
*NOTE : Do not use a water bath.*

(Example : Culturing MSC with EV-Up™ MSC EV Production medium)

- Prepare the 80-90% confluent MSCs.
- Aspirate the culture medium, and rinse the MSCs once with D-PBS(-).
- Remove D-PBS(-), and replace with EV-Up™ MSC EV Production medium.

Culture container	EV-Up™ MSC EV Production medium (mL)
6 well-plate	4
T25 flask	10
T75 flask	30
T150 flask	60

- Incubate at 37°C with 5% CO<sub>2</sub> for 3-5 days.  
*NOTE : No need to change medium.*
- After incubating for 3-5 days, collect the culture supernatant into a centrifugal tube.

— 1/4 —

- Centrifuge at 2,000×g at 4°C for 20 minutes.
- Collect the culture supernatant into a new tube after centrifugation. Be careful not to aspirate the residue.

*NOTE : Add the 1/100 amount of EV-Save™ Extracellular Vesicle Blocking Reagent (code No. 058-09261) in the culture supernatant, it can suppress the adsorption of EVs to the tubes and tips.*

- Isolate the EVs from the culture supernatant using MagCapture™ Exosome Isolation Kit PS (code No. 293-77601) etc.

### [Storage]

Store at 2-10°C.

If you wish to use EV-Up™ EV Production Basal Medium for MSC, AF in small quantities, freeze the supplements in small portions and prepare only the required amount on demand.

### [Package]

Code No.	Packaging
053-09451	95 mL

### [Related reagents]

Code No.	Product name
298-84001	EV-Up™ MSC EV Production Supplement, AF
058-09261	EV-Save™ Extracellular Vesicle Blocking Reagent
293-77601	MagCapture™ Exosome Isolation Kit PS

## FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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