Exosomes, from COLO201 cells, purified (CD9)

This product contains highly purified exosomes derived from COLO201 cells (ATCC® CCL-224™, American Type Culture Collection), originating from human colorectal adenocarcinoma cancer tissue, and isolated using a PS affinity method. These exosomes express CD9, CD63, and CD81, which are well-known exosomal surface markers.

They are suitable for use as a positive control in qualitative analyses such as Western blotting, ELISA, and flow cytometry. When measuring CD9 signals on exosomes with the PS Capture™ Exosome ELISA Kit series or the CD9-Capture Human Exosome ELISA Kit, this product can serve as a calibrator. By generating a standard curve with this product, the measured absorbance values can be converted to concentrations, enabling correction of inter-assay variations.

[Storage]

Store at -20°C *1

*1 Up to 10 freeze-thaw cycles at -20°C or -80°C do not affect the quality of this product. However, it is recommended to store the product in working aliquots to maintain stability. For example, prepare and store 10 aliquots of 5 µL each at -20°C, and thaw at room temperature before use.

[Concentration]

10 µg/mL

The concentration indicated is a protein concentration corrected based on the CD9 signal value, which serves to minimize lot-to-lot variation. This product is intended for use as a positive control and calibrator.

[Precaution for use]

This product consists of exosomes derived from COLO201 cells. The source cells have tested negative for human immunodeficiency virus (HIV), hepatitis B virus (HBV), human papillomavirus (HPV), Epstein-Barr virus (EBV), and cytomegalovirus (CMV). Other potential adventitious agents have not been tested; therefore, this product should be handled as potentially infectious material.

[Procedure]

〈Use as a Positive Control for Human CD9, CD63, and CD81 Exosomes〉

Dilute this product to an appropriate concentration and use it as a positive control for the qualitative detection of exosomes in Western blotting, ELISA, and flow cytometry assays.

〈Use as a Calibrator for PS Capture™ Exosome ELISA Kit Series and CD9-Capture Human Exosome ELISA Kit〉

This product can be used as a calibrator with the following ELISA kits to measure CD9 signals on exosomes, employing an anti-CD9 antibody*2 as the detection antibody:

(ELISA Kits)

- (A) 297-79201 PS Capture™ Exosome ELISA Kit (Anti-Mouse IgG POD)
- (B) 298-80601 PS Capture™ Exosome ELISA Kit (Streptavidin HRP)
- (C) 296-83701 CD9-Capture Human Exosome ELISA Kit (Streptavidin HRP)

By generating a standard curve using this product, measured absorbance values in samples can be converted into concentrations, allowing for inter-assay value normalization.

Dilute this product with the provided Reaction Buffer to ensure the ELISA absorbance at 450 nm (with absorbance at 620-650 nm subtracted as reference) falls within an optical density range of 0.1 to 3.0 O.D.

*2 The anti-CD9 antibody used for detection is not included in kits (A) or (B). When using these kits, please provide an appropriate anti-CD9 antibody independently. Our antibody lineup is available at the following link : <https://labchem-wako.fujifilm.com/category/01115.html>

-Usage example-

To prepare a standard curve for the PS Capture™ Exosome ELISA Kit (Streptavidin HRP), perform a two-fold serial dilution series over seven steps, starting from a concentration of 10 ng/mL:

- (1) Prepare a 100-fold dilution (100 ng/mL) by adding 2 µL of this product to 198 µL of Reaction Buffer provided in the kit.
- (2) Perform a two-fold serial dilution series using Reaction Buffer according to the table below, preparing concentrations from 10 ng/mL down to 0.15625 ng/mL in seven dilution steps.

Final Conc. (ng/mL)	Starting Conc. (ng/mL)	Volume to dilute (µL)	Reaction Buffer (µL)
10.0	100	50	450
5.0	10.0	250	250
2.50	5.0	250	250
1.25	2.50	250	250
0.625	1.25	250	250
0.3125	0.625	250	250
0.15625	0.3125	250	250
0.00 (Blank)	—	—	250

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FUJIFILM Wako Pure Chemical Corporation