EV-Save™ Extracellular Vesicle Blocking Reagent for TFF

[Product Information]

Code No.	Product Name	Volume	Storage
059-09671	EV-Save™ Extracellular Vesicle Blocking Reagent for TFF	20 mL	Refrigerated

[Overview]

This product is an additive designed to minimize the loss of extracellular vesicles (EVs) during buffer exchange and concentration performed using tangential flow filtration (TFF). Adding 1% (v/v) of the reagent to the buffer during buffer exchange reduces nonspecific adsorption to the TFF membrane and prevents clogging, thereby improving EV yield.

[Additional Materials Required (☐ : Check Box)]

- ☐ TFF system or peristaltic pump
- ☐ TFF module
- ☐ Tubing for fluid transfer
- ☐ Equilibration buffer
- ☐ Exchange buffer (only for buffer exchange)

[TFF Conditions]

When using Repligen Microkros mPES 100 kDa (C02-E100-05-N), the following are the recommended conditions.

Flow rate : 18 mL/min

Transmembrane pressure (TMP) : 3 psi

[Operating Procedure]

When using the Repligen KR2 system as the TFF device and the Repligen Microkros mPES 100 kDa (C02-E100-05-N) as the TFF module, the following are the recommended operating conditions.

A. EV Concentration Using the TFF System

- (1) Add EV-Save™ Extracellular Vesicle Blocking Reagent for TFF to the sample to be concentrated at a 1 : 100 dilution. Mix thoroughly.
- (2) Add EV-Save™ Extracellular Vesicle Blocking Reagent for TFF to 15 mL of equilibration buffer at a 1 : 100 dilution. Mix thoroughly. The equilibration buffer should have a composition equivalent to that of the sample solvent.
- (3) Set up the TFF apparatus, including the module and tubing.
- (4) Pre-wash the TFF module according to the Repligen instruction manual.
- (5) In Manual Mode, perform equilibration by preparing 15 mL of the equilibration buffer in the feed reservoir and

circulating it at a flow rate of 18 mL/min with a TMP of 3 psi until 10 mL has been discharged.

- (6) Load the sample to be concentrated into the feed reservoir.
- (7) In Auto Mode (Concentration Mode), concentrate the EV sample to the desired volume ratio maintaining a flow rate of 18 mL/min and a TMP of 3 psi.
- (8) After concentration, use Manual Mode to recover any remaining EV sample from the tubing and module.

B. EV Buffer Exchange Using the TFF System

- (1) Add EV-Save™ Extracellular Vesicle Blocking Reagent for TFF to the sample to be concentrated at a 1 : 100 dilution. Mix thoroughly.
- (2) Add EV-Save™ Extracellular Vesicle Blocking Reagent for TFF to 15 mL of equilibration buffer at a 1 : 100 dilution. Mix thoroughly. The equilibration buffer should have a composition equivalent to that of the sample solvent.
- (3) Add EV-Save™ Extracellular Vesicle Blocking Reagent for TFF to the exchange buffer at a 1 : 100 dilution. Mix thoroughly.
- (4) Set up the TFF apparatus, including the module and tubing.
- (5) Pre-wash the TFF module according to the Repligen instruction manual.
- (6) In Manual Mode, perform equilibration by preparing 15 mL of equilibration buffer in the feed reservoir and circulating it at a flow rate of 18 mL/min with a TMP of 3 psi until 10 mL has been discharged.
- (7) Load the sample to be concentrated into the feed reservoir.
- (8) Fill the supply reservoir with the exchange buffer.
- (9) In Auto Mode (Diafiltration Mode), perform buffer exchange at a flow rate of 18 mL/min with a TMP of 3 psi until the desired volume ratio is reached, preferably within 7 diafiltration volume exchanges (equivalent to ≤ 128 -fold buffer exchange).
- (10) After buffer exchange, use Manual Mode to recover any remaining EV sample from the tubing and module.

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