FUJIFILM

Wako

Code No. 095-07891 (20 mL)
097-07895 (500 mL)

For Blotting ImmunoOneStep™ Western

ImmunoOneStep™ Western is a ready-to-use reagent that enables rapid and specific immunodetection by combining blocking, primary antibody incubation and secondary antibody incubation into a single step in Western blotting. For example, compared to protocols using skim milk, this product reduces the time required for blocking and antibody incubation.

[Features]

1. Ready-to-Use
2. Blocking, primary antibody incubation and secondary antibody incubation in Western blotting are completed in a single one-step reaction in as little as 1 hour
3. Protein-free formulation reduces background signals

[Storage]

Store at room temperature.

[Required materials but not supplied]


- Primary antibody
- HRP-conjugated secondary antibody
- Wash buffer : PBS-T (Phosphate Buffered Saline with Tween-20) or TBS-T (Tris Buffered Saline with Tween-20)
- Chemiluminescent or colorimetric substrate
- Orbital or rocking shaker
- Detection system (CCD imaging system or X-ray film)

[Example Protocol]

Preparation

Separate the target protein by SDS-PAGE. Then, transfer the protein onto a PVDF membrane or a nitrocellulose (NC) membrane using a blotting apparatus.

1. After transfer in Western blotting, immerse the membrane in PBS-T or TBS-T for 5 minutes.


2. Dilute the primary and secondary antibody separately in different tubes using this product.*

*For example, if the primary antibody dilution is 1 : 1,000 and the secondary antibody dilution is 1 : 10,000, prepare each diluted antibody solution separately by adding 20 μL of the primary antibody or 2 μL of the secondary antibody

— 1/6 —

to 10 mL of ImmunoOneStep™ Western.

Dilution Example	Antibody	ImmunoOneStep™ Western
Primary Antibody (1 : 1,000)	20 μL	10 mL
Secondary Antibody (1 : 10,000)	2 μL	10 mL

3. Add 10 mL of the diluted primary antibody solution to a container, then add 10 mL of the diluted secondary antibody solution, and gently mix.

4. **Within 10 minutes**, immerse the membrane in the solution prepared in Step 3 and incubate with gentle shaking at room temperature for 1-2 hours.

5. Wash the membrane 3 times with PBS-T or TBS-T while shaking for 5 minutes each time.

6. Perform chemiluminescent or chromogenic detection.

[Note]

1. When diluting the secondary antibody with ImmunoOneStep™ Western, use a dilution of at least 1 : 10,000. Excess secondary antibody increases background.
2. Do not immerse the membrane in ImmunoOneStep™ Western for more than 4 hours. Prolonged incubation may increase background.
3. If high background is observed, consider optimizing the concentrations of the primary and secondary antibodies and/or the incubation time.

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FUJIFILM Wako Pure Chemical Corporation