## Luminol Reaction Reagent Set

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Luminol is a compound that emits a blue-violet luminescence when it reacts with hydrogen peroxide. This luminescence reaction is catalyzed by transition metals such as iron, copper, cobalt, their complexes, and oxidizing enzymes such as peroxidases. The amount of luminescence can be used to measure trace amounts of hydrogen peroxide and metals, as well as for stability testing. The iron complex in hemoglobin contained in blood acts as a catalyst and activates the luminescence reaction, which making it suitable for detecting bloodstains.

The luminol reaction produces chemiluminescence when luminol (5-amino-2,3-dihydro-1,4-phthalazinedione) is oxidized to form 2-aminophthalate dianion. Luminol is the most important compound among chemiluminescent substances, and its luminescence mechanism has been extensively studied.

In an alkaline aqueous solution, luminol emits a blue-violet light in the presence of various substances with peroxidase-like catalytic activity (iron, copper, cobalt, plants, hemoglobin, etc.).


### 【Contents】

Each pair of pouches contains :

- Luminol 0.5 g
- Sodium peroxide 2.5 g

### 【Preparation Method】

1. Prepare 500 mL of distilled water in a glass beaker or wide-mouth polyethylene container.
2. Cut along the black line and open the aluminum pouch containing the luminol and sodium peroxide.
3. Add the two reagents simultaneously to the distilled water in the glass beaker or wide-mouth polyethylene container.
4. Stir the solution with a glass rod if using a glass beaker. If using a polyethylene container, close the lid and mix the solution.

**[Examples of Detection Methods]**

- Add the sample dropwise into a small amount of the luminol reagent mixture in a beaker or test tube.
- Transfer the dissolved luminol reagent mixture into a spray bottle and spray onto the area for detection.

Note : The spray bottle should be made of plastic.

To prevent deterioration of the reagent, a transparent container should be covered with aluminum foil to protect from light, and store in a refrigerator until immediately before use.

**[Storage]**

Store in the dark at 2 ~ 10°C, under inert gas.

**[Precautions for Use]**

- Wear protective goggles, gloves, and a mask.
- Sodium peroxide may generate heat when in contact with a small amount of water. It is recommended to add it to a large volume of water when dissolving. Using chilled distilled water is safer because the heat will be absorbed. Exercise caution while handling it, as it can ignite on impact.
- If the two reagents cannot be added simultaneously, separate the aluminum pouches in the middle and add the sodium peroxide first, followed by luminol, to the distilled water. The sodium peroxide will make the solution alkaline and help to dissolve the luminol.
- The mixed reagent solution is alkaline. If the solution comes into contact with the eyes, immediately rinse with large amounts of water for several minutes and seek medical attention.

## &lt; FAQ &gt;

Q : What is the shelf life of the prepared solution?

A : It depends on the storage conditions, but it is recommended to use the solution within one week.

Since the solution begins to deteriorate immediately after preparation, it is recommended to prepare the solution just before use to ensure that bright luminescence can be observed. When storing, keep the solution refrigerated (2 ~ 10°C) and protected from light.

Q : How should I store the reagents before and after preparation?

A : This product before preparation contains sodium peroxide itself. After preparation, it is an alkaline solution containing 5% or less of sodium peroxide.

Store them in accordance with the appropriate laws and regulations of the facility, region, and country.

Q : What is the disposal method?

A : Residual waste (waste liquid) : Do not pour into the drainage system. Dispose of it in accordance with the appropriate laws and regulations of the facility, region, and country.

Contaminated containers and packaging : Dispose of them in accordance with the appropriate laws and regulations of the facility, region, and country.

---

## FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://flwk.fujifilm.co.jp>

### FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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