MSCulture™ High Growth Supplement**[Product and Storage]**

Product	Cat No.	Volume	Storage	Note
MSCulture™ High Growth Supplement	133-19331	5 mL	-20°C	Supplemental reagent
MSCulture™ High Growth Basal Medium	132-19345	500 mL	2-10°C	Basal Medium

If you wish to use MSCulture™ High Growth Basal Medium in small quantities, freeze the supplements in small portions and prepare only the required amount on demand.

[Note]

This supplement should be mixed with two more products, MSCulture™ High Growth Basal Medium (Cat No. 132-19345, **Optional**) and animal serum (i.e. FBS), for preparation of MSCulture-based MSC complete medium.

[About MSCulture™-based MSC complete medium]

MSCulture™-based MSC complete medium is optimized for mesenchymal stem cells (MSCs), and can be used for high quality culture and well proliferation. This medium can be used for several tissue-derived MSCs, including human bone marrow, human adipose tissue, and human umbilical cord.

After expansion of MSC by MSCulture™-based MSC complete medium, switching to an exosome producing medium, EV-Up™ (basal medium (Cat No. 053-09451) and supplement (Cat No. 298-84001)), facilitates efficient exosome production.

This product is for research use only.

[Required equipment and reagents]

Plastic labware :

- Centrifuge tube (15 mL)
- Culture dishes (depending on culture scale)

Reagents :

- MSCulture™ High Growth Supplement (Cat No. 133-19331)
- MSCulture™ High Growth Basal Medium (Cat No. 132-19345)
- Animal serum (i.e. Wako, FBS, Cat No. 554-04855)
- D-PBS (-) (i.e. Wako, Cat No. 045-29795)
- Trypsin (i.e. Wako, Cat No. 205-20255, 209-16941)

[Protocol]

■ Optimal density of MSCs on each culture plate

Prepare and seed optimal number of MSCs on your chosen culture plate below.

[Table 1]

Plate	Cells / vessel
6-well plate	0.45×10^5 cells
10 cm dish	2.75×10^5 cells
T25 flask	1.25×10^5 cells
T75 flask	3.75×10^5 cells
T225 flask	1.12×10^6 cells

■ Procedure**1. Preparation of MSCulture™-based MSC complete medium**

- Use laminar flow hood for refraining and preventing contamination.

- (1) Completely thaw MSCulture™ High Growth Supplement at room temperature (RT).

- Do not thaw MSCulture™ High Growth Supplement using a water bath.
- If some precipitation is observed after thawing MSCulture™ High Growth Supplement, vortex until precipitates have disappeared.

- (2) Add MSCulture™ High Growth Supplement and animal serum (i.e. 10-15%) to MSCulture™ High Growth Basal Medium.

- Prepared MSCulture™-based MSC complete medium should be used immediately.
- During storage, DO NOT freeze MSCulture™-based MSC complete medium.
- MSCulture™ High Growth Supplement can be re-frozen for long-term storage.

2. Culture MSCs from cryopreserved stock

- (1) Add 10 mL of MSCulture™-based MSC complete medium to a 15 mL centrifuge tube.

- Medium should be at RT before cells are suspended.

- (2) Thaw MSCs using water bath and add to (1) immediately.
- (3) Centrifuge at $200 \times g$ (1,000 rpm) for 3 min, and discard supernatant.
- (4) Add 1 mL of MSCulture™-based MSC complete medium and suspend by pipetting gently.
- (5) Count cell number and seed cells on suitable culture labware (see. Table 1).

- Recommended cell density is over 5×10^3 cells/cm²

- (6) Incubate MSCs up to 80-90% confluency at 37°C, 5% CO₂

- Culture period without passage is less than 96 hrs. If MSCs CANNOT proliferate around 80-90% confluency within 96 hrs, increase cell density to over 5×10^3 cells/cm². It is not necessary to change the medium during culture.

3. Passage and expansion of MSCs

- (1) Wash with D-PBS (-) twice, before MSC is trypsinized.
- (2) Transfer trypsinized MSCs into a 15 mL centrifuge tube.
- (3) Centrifuge at $200 \times g$ (1,000 rpm) for 3 min, and discard supernatant.

- (4) Add 1 mL of MSCulture™-based MSC complete medium and suspend by pipetting gently.
- (5) Count cell number and seed cells on suitable culture labware (see, Table 1).

• Recommended cell density is over 5×10^3 cells/cm²

- (6) Incubate MSCs up to 80-90% confluency at 37°C, 5% CO₂

• Culture period without passage is less than 96 hrs. If MSCs CANNOT proliferate around 80-90% confluency within 96 hrs, increase cell density. It is not necessary to change the medium during culture.

4. Optional : Exosome generation and purification using EV-Up™

- (1) Incubate MSCs up to 80-90% confluency at 37°C, 5% CO₂
- (2) Discard MSCulture™-based MSC complete medium, and wash with D-PBS (–) twice.
- (3) Discard D-PBS (–), and add suitable volume of EV-Up™ (see below).

Culture Labware	EV-Up™ MSC EV Medium volume (mL)
6-well plate	4
10 cm dish	20
T25 flask	10
T75 flask	30
T225 flask	90

- (4) Incubate MSCs for 3-5 days at 37°C, 5% CO₂

• It is not necessary to change the medium during culture.

- (5) Transfer EV-Up™ into a centrifuge tube.
- (6) Centrifuge culture medium at $2,000 \times g$ for 20 min
- (7) Transfer supernatant into a new centrifuge tube.

• Purify EVs using MagCapture™ Exosome Isolation Kit PS Ver.2 (Wako Cat No. 290-84103).

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : +1-804-271-7677
Facsimile : +1-804-271-7791
<http://www.wakousa.com>

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