FUJIFILM

Wako

Code No. 137-18751 (2 mL)

For immunochemistry
**MagCapture™ HP Anti PA tag Antibody
Magnetic Beads**

MagCapture™ HP Anti PA tag Antibody Magnetic Beads are composed of the Anti PA tag, rat monoclonal antibody attached to magnetic beads, which specifically recognize the PA tag which is consisted of 12 amino acids (GVAMPGAEDDVV), and can be used for purification and immunoprecipitation of PA-tagged recombinant proteins.

PA-tagged recombinant proteins can be eluted by competitive elution with PA tag Peptide under neutral pH condition after immunoprecipitation. In addition, used MagCapture™ HP Anti PA tag Antibody Magnetic Beads can also be recycled easily under mild condition at pH 6.0 by PA tag Washing Solution (Code No. 169-27261) after elution of PA-tagged recombinant proteins.

This product is intended for laboratory research use only.

[Formulation]

1 × TBS aqueous solution with 50% glycerol, 0.05w/v% sodium azide, pH 7.4

[Immobilized antibody]

Anti PA tag, Rat Monoclonal Antibody (NZ-1)

[Subclass of Immobilized antibody]

IgG2a

[Binding capacity of antigen]

16 µg protein / 1 mL beads slurry

Note : the binding capacity is for reference. It varies depending on protein of interest.

[Usage]

10 mL of cell culture medium / 50-100 µL beads slurry

200-1,000 µL of cell lysate / 20 µL beads slurry

Note : the above usage of beads slurry is for reference. Please optimize the appropriate usage in accordance with goal, protein of interest, and cell lines etc.

[Additional required materials]

Magnetic stand, Centrifuge, Rotator, Vortex mixer, Tube mixer, Centrifuge tube, Microcentrifuge tube (1.5 mL), Cell lysis buffer (see below), Wash buffer such as PBS (-), TBS etc., Protease inhibitors, Phosphatase inhibitors

Recommended cell lysis buffer

●RIPA Buffer (Code No.182-02451)

50 mmol/L Tris-HCl, pH 8.0, 150 mmol/L Sodium Chloride, 0.5w/v% Sodium Deoxycholate, 0.1w/v% Sodium Dodecyl Sulfate, 1.0w/v% NP-40 substitute

●Cell Lysis Buffer M (Code No. 038-21141)

20 mmol/L Tris-HCl, pH 7.4, 200 mmol/L Sodium Chloride, 2.5 mmol/L Magnesium Chloride, 0.05w/v% NP-40 substitute

When inhibitors are added to a cell lysis buffer, please use it in accordance with manufacturer's instruction.

[Protocol]

The immunoprecipitation method for PA-tagged recombinant proteins and the recycling method for used MagCapture™ HP Anti PA tag Antibody Magnetic Beads

1. Sample preparation : cell lysate and cell culture medium

1-A In the case of adherent cells

- 1) Incubate the PA-tagged-protein-expressing cell lines under appropriate conditions.^{※1}
- 2) Remove the medium from culture dishes and wash cultured cells twice with PBS(-).
- 3) Detach the cells from culture dishes with a cell scraper or trypsin-EDTA solution and transfer the cell suspension into a new 1.5 mL microcentrifuge tube. In the case of trypsin-EDTA solution, add the medium including serum to the cell suspension immediately to avoid the damage due to trypsin treatment.
- 4) Centrifuge 3) at 200 × g for 5 minutes at 4 °C.
- 5) Remove the supernatant.
- 6) Add PBS(-) and suspend the cell pellet, then centrifuge the cell suspension at 200 × g for 5 minutes at 4 °C.
- 7) Remove the supernatant.
- 8) Repeat 6)-7) steps once again.
- 9) Add 1mL of cell lysis buffer to the cell pellet and suspend by pipetting or vortex mixer.
- 10) Incubate on ice for 10 minutes.
- 11) Centrifuge 10) at 20,000 × g for 10-20 minutes at 4 °C.^{※2}
- 12) Transfer the supernatant to a new 1.5 mL microcentrifuge tube.
... **Cell lysate**

1-B In the case of suspension cells

- 1) Incubate the PA-tagged-protein-expressing cell lines under appropriate condition.^{※1}
- 2) Centrifuge the culture medium at 200 × g for 5 minutes at 4 °C.
- 3) Remove the supernatant.
- 4) Add 1 mL of PBS(-) and suspend the cell pellet, then transfer all to a new 1.5 mL microcentrifuge tube.
- 5) Centrifuge 4) at 200 × g for 5 minutes at 4 °C.
- 6) Remove the supernatant.
- 7) Add 1 mL of PBS(-) and suspend the cell pellet again, then centrifuge the cell suspension at 200 × g for 5 minutes at 4 °C.
- 8) Remove the supernatant.

- 9) Add 1 mL of cell lysis buffer to the cell pellet and suspend by pipetting or vortex mixer.
- 10) Incubate on ice for 10 minutes.
- 11) Centrifuge 10) at $20,000 \times g$ for 10-20 minutes at 4 °C.^{※2}
- 12) Transfer the supernatant to a new 1.5 mL microcentrifuge tube. ... **Cell lysate**

1-C In the case of cell culture medium

- 1) Incubate the PA-tagged-protein-expressing cell lines under appropriate condition.^{※1}
- 2) Collect the cell culture medium to a centrifuge tube, then centrifuge at $2,000 \times g$ for 10 minutes at 4 °C to remove cell.^{※3}
- 3) Transfer the supernatant to a new centrifuge tube.^{※4} ... **Cell culture medium**

- ※1 Please examine appropriate culture conditions depending on cell lines.
- ※2 When a white-floating layer appears in supernatant after centrifugation, please filter it with a centrifugal filter (pore size: 0.45 μ m).
- ※3 Cells etc. are removable by using the filtration filter (pore size: 0.22 μ m etc.).
- ※4 It is possible that adding a 1 mol/L Tris-HCl into collected cell culture medium to adjust its pH to 7.5.

2. Prewashing of Antibody Magnetic Beads

- 1) Suspend MagCapture™ HP Anti PA tag Antibody Magnetic Beads sufficiently by vortex mixer until it reaches a uniform consistency.
- 2) Transfer 50-100 μ L of MagCapture™ HP Anti PA tag Antibody Magnetic Beads into a 1.5 mL centrifuge tube.^{※5}
- 3) Add 1 mL of cell lysis buffer and suspend by vortex mixer.
- 4) Spin-down 3), then places it on the magnetic stand for 30 seconds to 1 minute.
- 5) Remove the supernatant after separation of magnetic beads and supernatant.
- 6) Repeat 3)-5) of prewashing steps twice. ... **Prewashed beads**

- ※5 Please examine appropriate amount of beads slurry depending on the purpose.

3. Immunoprecipitation

3-A In the case of cell lysate

- 1) Add 200-1,000 μ L of **Cell lysate** into **Prewashed beads**.^{※6}
- 2) Incubate with rotator for 2 hours~at 4 °C.
- 3) Spin-down 2) and places it on the magnetic stand until magnetic beads and supernatant separate sufficiently, then remove the supernatant.
- 4) Add 1.0 mL of cold cell lysis buffer and suspend by vortex mixer.
- 5) Spin-down 4), then places it on the magnetic stand for 30 seconds to 1 minute.
- 6) Remove the supernatant after separation of magnetic beads and supernatant.
- 7) Repeat 4)-6) of washing steps twice. ... **Antigen binding beads**

- ※6 Fill up to 1mL with cell lysis buffer as needed.

3-B In the case of cell culture medium

- 1) Add 1mL of **Cell culture medium** into **Prewashed beads**.^{※7}
- 2) Incubate with rotator for 2 hours~at 4 °C.
- 3) Spin-down 2) and places it on the magnetic stand until magnetic beads and supernatant separate sufficiently, then remove the supernatant.
- 4) Add 1.0 mL of cold cell lysis buffer and suspend by vortex mixer.
- 5) Spin-down 4), places it on the magnetic stand for 30 seconds to 1 minute.
- 6) Remove the supernatant after separation of magnetic beads and supernatant.
- 7) Repeat 4)-6) of washing steps twice. ... **Antigen binding beads**

- ※7 Please examine appropriate amount of beads slurry depending on the purpose.

4. Elution of PA-tagged recombinant Proteins

Please select an appropriate elution method in accordance with goal and property of protein.

4-A Elution by PA tag Peptide (non-denaturing condition)

- 1) Prepare a stock solution of PA tag Peptide (final concentration: 5 mg/mL) by dissolving 5 mg of PA tag Peptide (Code No. 167-25501) in 1mL of 1 $\times$ TBS (10 mmol/L Tris-HCl, pH 7.5, 150 mmol/L Sodium Chloride).^{※8}
- 2) Prepare a diluted solution of PA tag Peptide (final concentration : 100-200 μ g/mL) by diluting the stock solution from step 1).
- 3) Add 30-100 μ L of the diluted solution of PA tag Peptide into **Antigen binding beads**.^{※9}
- 4) Incubate with rotator or tube mixer for 30 minutes at 4 °C.^{※10}
- 5) Spin-down 4) and places it on the magnetic stand until magnetic beads and supernatant separate, then transfer the supernatant into a new 1.5 mL microcentrifuge tube. ... **Eluted sample**
- 6) Store the **Eluted sample** on ice. In the case of long storage, store at -20 °C.^{※11}

Point : Repeating 3)-5) of elution steps with PA tag peptide can increase the yield of eluted protein. Please examine the appropriate number of times for elution with PA tag peptide.

4-B Elution by 0.1 mol/L Glycine-HCl, pH 3.0 (acidic condition)

- 1) Add 30-100 μ L of 0.1 mol/L Glycine-HCl, pH 3.0 into **Antigen binding beads**.^{※9}
- 2) Incubate with rotator or tube mixer for 5-10 minutes at room temperature.^{※12}
- 3) Spin-down 2) and places it on the magnetic stand until magnetic beads and supernatant separate, then transfer the supernatant into a new 1.5 mL microcentrifuge tube. ... **Eluted sample**
- 4) Neutralize **Eluted sample** immediately with neutralization buffer. (Example : 1 mol/L Tris-HCl, pH 9.0) ... **Neutralized sample**
- 5) Store the **Neutralized sample** on ice. In the case of long storage, store at -20 °C.

4-C Elution by SDS sample buffer (denaturing condition)

Note : when 2 × SDS sample buffer includes reducing agent such as 2-mercaptoethanol or DTT, heavy chain (50 kDa) and light chain (25kDa) may dissociate from anti PA tag antibody and are eluted into a recovered sample.

- 1) Add 50 μ L of 2 × SDS sample buffer (0.125 mol/L Tris-HCl, pH 6.8, 4w/v% SDS, 20w/v% Glycerol, 10% 2-mercaptoethanol or 100 mmol/L DTT, 0.01% Bromophenol Blue) into **Antigen binding beads** and mix by vortex mixer.
- 2) Boil for 5 minutes at 95 °C.
- 3) Spin-down 2) and places it on the magnetic stand until magnetic beads and supernatant separate, then transfer the supernatant into a new 1.5 mL microcentrifuge tube. . . . **Eluted sample**
- 4) Store the **Eluted sample** on ice. In the case of long storage, store at -20 °C.

- ※8 Please store the used stock solution of PA tag Peptide at -20 °C, and use it as soon as possible.

※9 Please examine appropriate amount depending on the purpose.

※10 In order to increase the efficiency of competitive elution, please be sure to incubate with rotator or tube mixer for 5 minutes~.

※11 Please remove the peptide in the solution by dialysis or gel filtration as needed.

※12 There is a possibility that proteins of interest are denatured if reaction time is too long.

5. Recycling of the used Anti PA tag Antibody Magnetic Beads

The used Anti PA tag Antibody Magnetic Beads can be recycled with PA tag Washing Solution (Code No. 169-27261). Here, the recycling protocol for 50 μ L of beads slurry is described.

Note : denatured magnetic beads cannot be recycled.

- 1) 500 μ L of PA tag Washing Solution into the used Anti PA tag Antibody Magnetic Beads.
- 2) Incubate with rotator for 10 minutes at room temperature.
- 3) Spin-down 2) and places it on the magnetic stand until magnetic beads and supernatant separate, then remove the supernatant.
- 4) Repeat 1)-3) steps three times.
- 5) Add 500 μ L of TBS and suspend sufficiently by vortex mixer.
- 6) Spin-down 5) and places it on the magnetic stand until magnetic beads and supernatant separate, then remove the supernatant.
- 7) Repeat 5)-6) steps three times.
- 8) Add 2 times the amount of TBS or PBS(-).
- 9) Add 50 μ L of storage buffer (1 × TBS aqueous solution with 50% glycerol, 0.05w/v% sodium azide, pH 7.4).
- 10) Store 2-10 °C.

[Storage]

2-10 °C

[References]

1. Fujii, Y. *et al.*, *Protein Expr. Purif.*, **95**, 240-7 (2014)
2. 藤井勇樹ら, 蛋白質科学会アーカイブ, **7**, e075 (2014)
3. Kaneko, M. K. *et al.*, *Biochem. Biophys. Res. Commun.*, **432**, 40-45 (2013)
4. Fujii, Y. *et al.*, *J. Cell Sci.*, **129**, 1512-1522 (2016)

[Related products]

Code No.	Product Name	Package
169-27261	PA tag Washing Solution	50 mL
167-25501	PA tag Peptide	5 mg
163-25503		25 mg
182-02451	RIPA Buffer	100 mL
188-02453		250 mL
038-21141	Cell Lysis Buffer M	100 mL
034-21143		250 mL

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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