FUJIFILM**Wako**

<For Research Use Only> Code No. 141-07983 (10 mL)
149-07984 (100 mL)

for genetic research

Ni-Agarose

Ni-Agarose is designed for affinity purification of 6 × His fusion protein. This is used iminodiacetic acid(IDA) as ligand. Nickel is the most widely used cation for this type of purification because it provides effective binding and selectivity for a wide range of proteins.

〔Features〕 · One step purification
· High yield of purified protein
· High selectivity
· Hydrophilic linkage low non specific binding

〔Agarose〕 Cross linked 6% agarose

〔Binding/Loading Capacity〕 20-40 μmol Me²⁺/mL gel

〔Ligand〕 Iminodiacetic Acid (IDA)

〔Bead size〕 40-180 μm

〔Storage buffer〕 50% suspension solution in 20% ethanol

〔General purification method〕

Linkage condition

The linking buffer is used acetate (50 mmol/L) or phosphate (10-150 mmol/L). Its pH comes close to neutrality (normally pH 7.0-8.0) but can vary in the following range: pH 5.5-8.5. Generally, to avoid the effects of ionic interchange, 0.15-0.5 mol/L NaCl is added.

It is important to avoid agents such as EDTA or citrate in this buffer.

Elution condition

- Addition of a competitive ligand (generally, imidazole) allowing the elution of the retained proteins. In general, 0.5 mol/L imidazole is sufficient to desorb a particular protein. It is also appropriate to use concentration gradients of this agent (0-0.5 mol/L). Most of the proteins are desorbed by concentrations around 250 mmol/L. Other agents that may be used as competitive ligands are histidine and ammonium chloride.
- The reduction of pH (with or without gradient) also permits the elution of the protein of interest (normally pH 3.0-4.0). A more drastic method uses agents such as EDTA or EGTA (0.05 mol/L) which cause the desorption of the protein as well as of the chelating metal.

〔Storage〕 2~10 °C

〔Literatures〕

- Bestida, A., Latorre, M. and Garcia-Junceda, E., *Chem-BioChem.*, **4**, 531 (2003)
- Porath, J., Carlsson, J., Olsson, I. and Belfrage, G., *Nature.*, **258**, 598 (1975)

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://ffwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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