FUJIFILM

Wako

Code No. 166-24231 (1 mg)

## Plasmin Solution, from Human Plasma

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Plasmin exhibits preferential cleavage at the carboxyl side of lysine and arginine residues with higher selectivity than trypsin. This product dissolves the fibrin of blood clots and acts as a proteolytic factor in a variety of other processes including embryonic development, tissue remodeling, tumor invasion and inflammation.

This product is for laboratory research use only ; use in any such application is the responsibility of the user.

**Formulation** : 100 mmol/L sodium phosphate buffer (pH 7.3) containing 1mmol/L 6-aminohexanoic acid and 25% glycerol

**Source** : Human Plasma

**Activity** : Indicated on the label


**Protein Concentration** : Indicated on the label

**CAS No.** : 9001-90-5

## Assay Method

## 1. Principle

*N*-(*p*-Tosyl)-Gly-Pro-Lys *p*-nitroanilide

## 2. Reagents

A. 0.1 mol/L Phosphate Buffer, pH 7.8

B. Substrate Solution

Dissolve 17.8 mg of *N*-(*p*-Tosyl)-Gly-Pro-Lys *p*-nitroanilide acetate in 2 mL of water

C. Enzyme Solution

Add Reagents A to 40  $\mu$ L of this product to bring the final volume up to 10 mL. Store on ice.

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- 1/3 -

## 3. Procedure

| Reagent                                                                                  | Test    | Blank   |
|------------------------------------------------------------------------------------------|---------|---------|
| A                                                                                        | 2.70 mL | 2.70 mL |
| B                                                                                        | 0.15 mL | 0.15 mL |
| Pre-incubate at 25 °C for 5 minutes.                                                     |         |         |
| C                                                                                        | 0.10 mL | —       |
| A                                                                                        | —       | 0.10 mL |
| Measure the absorbance at 405 nm wavelength with Reagent A as the control for 5 minutes. |         |         |

## 4. Unit Definition

One Plasmin unit is the amount of enzyme which causes the production of 1  $\mu$ mol of *p*-Nitroaniline for 1 minute at 25 °C and pH 7.8 with *N*-(*p*-Tosyl)-Gly-Pro-Lys *p*-nitroanilide acetate as the substrate.

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- 2/3 -

**Storage** : Store at  $-20^{\circ}\text{C}$

**Warning** : This product uses blood ingredients which have tested negative for the HBs antibody, the HIV antibody and the HCV antibody. However, precautions should be taken when handling this product to prevent possible infection.

\* : Trans is the Japanese translation.

---

## FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://fwk.fujifilm.co.jp>

### FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

### FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12  
D-41468 Neuss  
Germany  
Telephone : + 49-2131-3111-0  
Facsimile : + 49-2131-311100  
<http://www.wako-chemicals.de>

2306KA2