FUJIFILM

Wako

Code No. 160-28911 (25 mg)
166-28913 (100 mg)
164-28914 (500 mg)
160-28916 (1 g)

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Source : *Pichia pastoris* expressed *Parengyodontium album*
(*Tritirachium album*) proteinase K

Appearance : White crystalline powder

Activity : Indicated on the label.

Measurement of enzyme activity :

1. Principle

Hemoglobin is reacted with proteinase K to hydrolyze. The obtained peptides are colored with a Phenol reagent.

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2. Reagents

- A. 0.3 mol/L Trichloroacetic acid solution
Dissolve 4.9 g of trichloroacetic acid in distilled water for the final volume of 100 mL.
- B. Phosphate buffer (pH 7.5)
(Containing 6 mol/L urea and 5 mmol/L sodium chloride)
Dissolve 0.680 g of potassium dihydrogen phosphate in distilled water for a volume of 50 mL. Add 36.04 g of urea and 0.03 g of sodium chloride, adjust pH to 7.5 with 1 mol/L sodium hydroxide solution, and dilute with distilled water to bring the volume up to 100 mL.
- C. Hemoglobin solution
Dissolve 500 mg of hemoglobin from bovine blood in solution B for the final volume of 25 mL.
- D. 0.02 mol/L Calcium chloride solution
Dissolve 0.55 g of calcium chloride in distilled water for the final volume of 250 mL.
- E. Enzyme solution
Measure 4 mg of Proteinase K, recombinant exactly and dissolve in solution D for the final volume of 200 mL.
- F. Sodium carbonate solution
Dissolve 6 g of sodium carbonate in distilled water for the final volume of 100 mL.
- G. L-Tyrosine standard solution
Measure 10 mg of L-tyrosine exactly and dissolve in 0.2 mol/L hydrochloric acid for the final volume of 100 mL. (100 μ g tyrosine/mL)
- 1) G₁ solution
Pipet 5 mL of this solution and dilute with distilled water

- to bring the volume up to 50 mL. (10 μ g tyrosine/mL)
- 2) G₂ solution
Pipet 15 mL of this solution and dilute with distilled water to bring the volume up to 100 mL. (15 μ g tyrosine/mL)
- 3) G₃ solution
Pipet 5 mL of this solution and dilute with distilled water to bring the volume up to 25 mL. (20 μ g tyrosine/mL)
- 4) G₄ solution
Pipet 5 mL of this solution and dilute with distilled water to bring the volume up to 20 mL. (25 μ g tyrosine/mL)

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3. Procedure

1) Test

Reagent	Test	Blank
C	5 mL	—
D	0.8 mL	0.8 mL
—	Keep at 37°C.	—
E	0.2 mL	0.2 mL
A	—	10 mL
—	Incubate at 37°C for 10 minutes.	Shake.

C	—	5 mL
A	10 mL	—
Incubate at 37°C for 30 minutes.		
Filtration.		
Filtrate	2 mL	2 mL
F	5 mL	5 mL
Phenol reagent(1+2)	1 mL	1 mL
Incubate at 37°C for 20 minutes.		
Measure the absorbance at 660 nm with water as a reference using a 10 mm cell.		
[Trans] *		
吸収セル 10mm を用い、波長 660nm における吸光度を水を対照液として測定する。		
E ₁ : Absorbance of test, E ₀ : Absorbance of blank		

2) Calibration curve

Reagent	Blank	10 μg Tyr/mL	15 μg Tyr/mL	20 μg Tyr/mL	25 μg Tyr/mL
—	Water 2 mL	G ₁ solution 2 mL	G ₂ solution 2 mL	G ₃ solution 2 mL	G ₄ solution 2 mL
F	5 mL	5 mL	5 mL	5 mL	5 mL
Phenol reagent(1+2)	1 mL	1 mL	1 mL	1 mL	1 mL
Incubate at 37°C for 20 minutes.					
Measure the absorbance at 660 nm with water as a reference using a 10 mm cell, and prepare the calibration curve.					
[Trans] *					
吸収セル 10mm を用い、波長 660nm における吸光度を水を対照液として測定し、検量線を作成する。					

4. Unit Definition

One unit is defined as the enzyme amount that generates peptides equivalent to 1 μmol tyrosine as the coloring substance reacted with a Folin reagent at 37°C for one minute under the above conditions.

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A : Activity (units/mg)

B : Amount of tyrosine that was calculated with calibration curve and (E₁ - E₀) (μg/mL)

S : Amount of sample (mg)

181.19 : Molecular weight of tyrosine

0.02 : Sample concentration of Enzyme solution (Reagent E) (mg/mL)

0.2 : Volume of Enzyme solution (Reagent E) (mL)

16 : Total volume of reaction solution (mL)

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FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

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FUJIFILM Wako Pure Chemical Corporation