FUJIFILM

Wako

Code No. 178-00551 (1 L)
174-00553 (250 mL)

Quick-CBB PLUS

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Quick-CBB PLUS stains protein bands on polyacrylamide gel with CBB G-250. This staining method is simple and fast, and does not require a fixing process. Therefore waste of reagents such as methanol and acetic acid is reduced. Additionally, the protein bands are clearly visible in contrast with low background.

【Features】

1. Quick-CBB PLUS does not use organic solvents such as methanol or acetic acid.
2. The protocol is simple and fast.
3. Stains protein bands clearly with low background.

【Procedure】

1. Rinse the SDS-PAGE gel with deionized water for 5 minutes in a staining tray. Repeat 2 more times. Drain off the deionized water from the tray.
2. Add 25 mL of Quick-CBB PLUS to the tray and shake it for 30-60 minutes until the protein bands are visible.
3. Drain off the solution and rinse the gel with distilled water for 60 minutes.

【Usage note】

1. Use a tray with a smooth inner surface.
2. Immerse the gel completely in staining solution while staining.

【Cautions for handling】

1. If the Quick-CBB PLUS appear to cause precipitate due to low temperature, warm it to resolubilize before use.
2. Quick-CBB PLUS does not contain any toxic substances. However, the solution is blue, so disposal in a specific vessel may be necessary.

【Reference】

Blakesley, R., W. and Boezi, A., J. : *Anal. Biochem.*, **82**, 580 (1977)

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : +1-804-271-7677
Facsimile : +1-804-271-7791
<http://www.wakousa.com>

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