FUJIFILM**Wako**

< For research Use Only >

Code No. 180-01891 (50 reactions)

for genetic research

**RNastabilizer, Irreversible RNase
Inactivation Reagent**

[Introduction]

RNastabilizer is a reagent that effectively stabilizes RNA prepared from organs containing a large amount of RNase, such as pancreas, liver and kidney. RNase contamination in the isolated RNA sample can induce experimental failure. A minute small amount of residual RNase may cause degradation of RNA.

RNastabilizer is a unique non-protein reagent that will irreversibly inactivate RNase in the RNA preparation step. RNastabilizer is applicable to commercial RNA isolation kits based on a solid-phase method with silica matrix. RNastabilizer can aid in the preparation of high quality, stable RNA preparation by virtue of inactivating RNase irreversibly.

[Advantages]

- Irreversible inactivation of RNase.
- Inhibition of reactivating RNase, which is inactivated once by a reducing agent such as mercaptoethanol.
- Preparing very stable RNA from organs such as pancreas, liver and kidney, which contain a large amount of RNase.
- Applicable to the commercial RNA isolation kits based on a silica-matrix method, such as Invitex GmbH's and QIAGEN's kits.

[Contents]

3.5 mL (in methanol) × 1

[Storage and Stability]

Store at -20 °C (in a dark place).

Reagent expiration dates are listed on the product label.

RNastabilizer solution is guaranteed for 1 year from the manufacture day under the above-mentioned conditions. If the color of the solution changes to light yellow or yellow, do not use it. If any precipitation of RNastabilizer is seen, redissolve it by vortex before use. This precipitation does not affect the quality

[Reactions]

50 reactions (for Invitek's kit)

[Method]

1. Note before starting
If any precipitation of RNastabilizer is seen, redissolve by vortex, and then place at room temperature.
2. Standard protocol
RNastabilizer is applicable to the washing step in RNA isolation kits based on silica matrix method. Please refer to the below-mentioned protocol, and apply to the applicable RNA isolation kit.
 - (1) Add 0.1 volume of RNastabilizer to the washing solution^(*) for the first of the washing steps.
(*) Generally, the RNA isolation kit uses two types of washing solution in the column washing steps. For example:
Invitek's kit : "Wash Buffer R1" of Invisorb Spin Tissue RNA mini Kit
QIAGEN's Kit : "Buffer RW1" of RNeasy Mini Kit

The mixture of washing solution and RNastabilizer can not be stored. This mixture should be prepared just before use.

- (2) Please refer to the manufacture's manual of your RNA isolation kit until RNA-adsorbing to silica matrix.
- (3) Add the RNastabilizer-containing washing solution (adequate volume as described in the protocol of your kit) to the spin column.
- (4) Incubate for 20 minutes at room temperature. Then, centrifuge at a recommended speed in the manufacture's manual of your RNA isolation kit.
- (5) Wash the column with the second washing solution, which does not contain RNastabilizer, according to the manufacture's manual of your RNA isolation kit.

After RNA-binding to the silica column, be sure to carry out both times of washing steps, which are the 1st washing step with RNastabilizer and the 2nd step without it. The effect of the RNastabilizer does not appear except for the addition of RNastabilizer to the first washing step.

- (6) Carry out the following elution step according to the manufacture's manual of your RNA isolation kit.

[The example of use]

Application protocol of RNastabilizer to Invisorb Spin Tissue RNA Mini Kit (Invitek GmbH : Code No. 10607002) :

- (1) Add 50 μ L of RNastabilizer to 500 μ L of Wash Buffer R1 and mix.


- (2) Add 500 μ L of Wash Buffer R1 containing RNastabilizer to the spin filter, and place at room temperature for 20 min.

The mixture of washing solution and RNastabilizer can not be stored. This mixture should be prepared just before use. The effect of the RNastabilizer does not appear except for the addition of RNastabilizer to the first washing step.

- (3) Centrifuge at $10,000 \times g$ for 30 sec. at room temperature.
(4) After the next washing step with Wash Buffer R2, please refer to the manufacturer's manual of RNA isolation kit.

[Result] Preparation of total RNA from mouse liver and pancreas using RNastabilizer :

The extracted total RNA was run on denaturing agarose gel after incubation at 37 °C over night. Total RNA extracted using RNastabilizer has high purity and superior storage stability by virtue of irreversible RNase inactivation.

Related products

Catalog No.	Description	Package Size
10607002	Invisorb TM Spin Tissue RNA Mini Kit	50 reactions
10607003	Invisorb TM Spin Tissue RNA Mini Kit	250 reactions

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://ffwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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