

FUJIFILM**Wako**Code No. 184-03511 (100 μ L)
180-03513 (100 μ L \times 5)

rBC2LCN-FITC (AiLecS1-FITC), AF

rBC2LCN is a recombinant protein of N-terminal domain of BC2L-C and has been reported as a marker of undifferentiated human pluripotent stem cells (hPSCs), human ES cells and human iPS cells. This product is the rBC2LCN labeled with fluorescein isothiocyanate (FITC).

rBC2LCN specifically recognize a sugar chain which exists on the surface of hPSCs. When you add rBC2LCN-FITC to the medium which is used for hPSCs culture, you can stain the living hPSCs. rBC2LCN-FITC is not toxic to cells, therefore you can culture hPSCs without any effect on the growth and pluripotency in the medium containing rBC2LCN-FITC. Also, this product can be used for flow cytometry.

This product is produced using non-animal ingredients as raw materials.

This product is for laboratory research use only ; use in any such application is the responsibility of the user.

[Formulation]

Phosphate Buffered Saline sterilized with 0.1 μ m filter

[Working Dilution]

Live Cell Imaging 1 : 100 ~ 1 : 1,000
Flow Cytometry 1 : 100 ~ 1 : 1,000

[Protocol for Live Cell Imaging]

- 1) Culture hPSCs in appropriate culture conditions.
- 2) Add 1 ~ 10 μ L of rBC2LCN-FITC per medium 1 mL.
- 3) Incubate the cells over 30 min under appropriate conditions.
- 4) Replace the medium containing rBC2LCN-FITC to a fresh medium, HBSS (+) or D-PBS (-).
- 5) View and analyze the cells on an imaging instrument (excitation 495 nm, emission 520 nm).
* Step 4 may be omitted.

[Protocol for Flow Cytometry]

- 1) Dissociate the hPSCs and wash with FCM* buffer once. Centrifuge and discard the supernatant.
- 2) Resuspend the cells with FCM* buffer to prepare cell suspension of approximately 5×10^6 cells/mL.
- 3) Add 1 ~ 10 μ L of rBC2LCN-FITC per 1 mL of cell suspension.
- 4) Incubate the cells in the dark over 30 min at room temperature.
- 5) Centrifuge the cell suspension at 1,000 rpm for 3 min, and discard the supernatant.
- 6) Wash the cells with 1 mL of FCM* buffer, centrifuge the cell suspension at 1,000 rpm for 3 min and discard the supernatant.

- 7) Suspend the washed cells in appropriate volume of FCM* buffer.
- 8) Analyze the cells by flow cytometry.
* FCM buffer : D-PBS (-) containing 10 mmol/L EDTA and 1% FBS, and so on.

[Attention]

- 1) Staining with the rBC2LCN-FITC persists for a few days even after removing the rBC2LCN-FITC.
- 2) For continuous observations, it might be better to add the same concentration of rBC2LCN-FITC at the daily medium change. Cytotoxicity of the rBC2LCN-FITC has not been observed.
- 3) rBC2LCN-FITC is applicable to not only live cell imaging, but also staining fixed cells.
- 4) The signal of background might become higher when the cells are cultured with the medium containing serum.
- 5) When sorting cells by flow cytometry, please add Y-27632 at a final concentration of 10 μ mol/L to the FCM buffer at the time of cell dispersion.

[Storage]

Store at -20°C, in the dark.

After thawing, store at 2~10°C. Avoid repeating freeze-thaw.

[References]

- 1) Onuma, Y., Tateno, H., Hirabayashi, J., Ito, Y. and Asashima, M. : *Biochem. Biophys. Res. Commun.*, **431**, 524 (2013).
- 2) Tateno, H., Matsushima, A., Hiemori, K., Onuma, Y., Ito, Y., Hasehira, K., Nishimura, K., Ohtaka, M., Takayasu, S., Nakanishi, M., Ikehara, Y., Nakanishi, M., Ohnuma, K., Chan, T., Toyoda, M., Akutsu, H., Umezawa, A., Asashima, M and Hirabayashi, J. : *Stem Cells Transl. Med.*, **2**, 265 (2013).
- 3) Tateno, H., Onuma, Y., Ito, Y., Hiemori, K., Aiki, Y., Shimizu, M., Higuchi, K., Fukuda, M., Warashina, M., Honda, S., Asashima, M. and Hirabayashi, J. : *Sci. Rep.*, **4**, 4069 (2014).

[Package]

Code No.	Package
184-03511	100 μ L
180-03513	100 μ L \times 5

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggenstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>