

Code No. 192-18001
199-18011**SuperSep™ Phos-tag™ (50 μmol/L),
100×100×6.6mm**



SuperSep™ Phos-tag™ are precast polyacrylamide gels containing Phos-tag™ with zinc ion (Zn^{2+}). Phos-tag™ is a functional molecule that binds specifically to the phosphate group^{1), 2), 3)}. This product can trap phosphorylated proteins during SDS-PAGE, and it allows detection of phosphorylated and non-phosphorylated proteins as different bands. The gels have a neutral pH to enhance stability.

This product is designed to be used with XCell SureLock® Mini-Cell (Catalog No. EI0001, Life Technologies Corporation). XCell SureLock® is a registered trademark of Life Technologies Corporation.


■ Function of Phos-tag™

- ◇ Selectivity of binding of a phosphate ion is much higher than that of other anions.
- ◇ Stable complex is formed under physiological conditions (pH 5 to 8).

Basic structure of Phos-tag™

 M^{2+} : manganese or zinc ion**■ Principle of Phos-tag™ SDS-PAGE****■ Example of Phos-tag™ SDS-PAGE****Time course of phosphorylation by using the Tyrosine kinase Abl**

Phosphorylated tyrosine was prepared by GST binding protein of tyrosine kinase Abl and the substrate peptide (Abltide) and separated with conventional SDS-PAGE and Phos-tag™ SDS-PAGE, respectively.



< Note >

- Handle with care. The glass plates are fragile.
- This product is not suitable for measurement of protein molecular weight.
- The gel becomes clouded when being incubated with any buffer including methanol such as acetate-methanol solution for fixing proteins and transfer buffer. It does not affect the quality of the product.
- Handle the gel carefully not to tear it.
- It is not recommended to apply ordinary prestained protein ladders to polyacrylamide gels containing Phos-tag™ such as this product because the ladder patterns is disturbed. A prestained ladder recommended is WIDE-VIEW™ Prestained Protein Size Marker III (Code No. 230-02461), which is designed to reduce disturbance of a ladder pattern in Phos-tag SDS-PAGE.

[Format]

Cassette size	100×100×6.6 (mm)
Max well volume	17 wells : 25 μL/well

[Package]

5 gels

[Storage]

Keep at 2 ~ 10°C. Do not freeze.

[Additional materials required]


- Running Buffer (×10) (Code No. 184-01291)
 - ▶ 0.25 mol/L Tris, 1.92 mol/L Glycine, 1% SDS.
- Sample Buffer (2ME+) (×4) (Code No. 191-13272)
 - ▶ 0.25 mol/L Tris-HCl, pH 6.8, 40 w/v% Glycerol, 8% SDS, 20 v/v% 2-Mercaptoethanol, 0.02 w/v% BPB.
- Electrophoresis apparatus -XCell SureLock® Mini-Cell (Catalog No. EI0001, Life Technologies Corporation).
- Disodium Dihydrogen Ethylenediamine Tetraacetate (EDTA)
- Transfer buffer (×10) : 10×Tris-Glycine Buffer (Code No. 201-18601)
 - ▶ 0.25 mol/L Tris, 1.92 mol/L Glycine.
- Methanol (Code No. 131-01826)
- Blotting papers
- PVDF or nitrocellulose membrane (Code No. 033-22453, 032-22663)

[Procedure]**1. Sample preparation for Phos-tag SDS-PAGE**


- 1) Remove contaminant in protein samples by TCA precipitation, dialysis or desalting*.
 - * Phos-tag SDS-PAGE is vulnerable to contaminant in protein samples, especially chelating reagent, vanadic acid, inorganic salts, surfactants. Cleaning them up is strongly recommended before Phos-tag SDS-PAGE.
- 2) Mix sample with Sample Buffer Solution (2ME+) (×4) and heat for 5 min at 95°C. Cool it to room temperature.

2. Phos-tag SDS-PAGE


- 1) Remove a gel cassette from a pouch and set it into the electrophoresis apparatus.
- 2) Fill the reservoir with 1× running buffer.
- 3) Carefully remove the comb from the gel.
- 4) Load the samples into lanes. Blank lanes may cause distortion of bands. Load the same amount of 1× sample buffer into blank lanes.
- 5) Start electrophoresis, according to the instruction of the electrophoresis apparatus.
- 6) Remove the gel cassette from the electrophoresis apparatus and rinse the gel cassette with water.
- 7) Insert a gel knife into a bottom slit of gel cassette. Rever up the one side gel plate.




- 8) Carefully remove the one side gel plate, while the gel remains on the other plate.



- 9) Insert a gel knife into the slit between the gel and gel plate, and cut pressingly the gel. In the same way, cut pressingly the gel at left and right sides.



- 10) Turn it upside down over a dish and detach the gel with the help of the knife.




— 3/4 —

3. Pre-treatment for Transfer

An additional procedure, elimination of the zinc ion (Zn^{2+}) from the gel using EDTA, is necessary before transfer. This procedure increases transfer efficiency of proteins from a gel to a membrane.

- 1) Prepare 1× transfer buffer with 10 mmol/L EDTA and without EDTA.
 - ▶ 1× transfer buffer with 10 mmol/L EDTA
25 mmol/L Tris, 0.192 mol/L Glycine, 10 v/v% methanol and 10 mmol/L EDTA.
 - ▶ 1× transfer buffer without 10 mmol/L EDTA
25 mmol/L Tris, 0.192 mol/L Glycine and 10 v/v% methanol.
- 2) Soak the gel in 1× transfer buffer with 10 mmol/L EDTA for a minimum of 20 minutes with gentle agitation. Repeat it 3 times with buffer exchanges.

For 20 minutes, 3 times



- 3) Soak the gel in 1× transfer buffer without 10 mmol/L EDTA for 10 minutes with gentle agitation.

For 10 minutes, 1 time



Transfer buffer without EDTA

- 4) Transfer the proteins from the gel to a membrane*.

* A Wet-tank method is strongly recommended for effective protein transfer.

[References]

- 1) Kinoshita, E., Kinoshita-Kikuta, E., Takiyama, K. and Koike, T. : *Mol. Cell. Proteomics*, **5**, 749 (2006).
- 2) Yamada, S., Nakamura, H., Kinoshita, E., Kinoshita-Kikuta, E., Koike, T. and Shiro, Y. : *Anal. Biochem.*, **360**, 160 (2007).
- 3) Kinoshita-Kikuta, E., Aoki, Y., Kinoshita, E. and Koike, T. : *Mol. Cell. Proteomics*, **6**, 356 (2007).

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://ffwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2303KA3

— 4/4 —