<For Research Use Only> Code No. 193-16711 (5 Gels)
190-16721 (5 Gels)

SuperSep™ Phos-tag™ (50 μmol/L), 10%, 13well

SuperSep™ Phos-tag™ (50 μmol/L), 10%, 17well

SuperSep™ Phos-tag™ are precast polyacrylamide gels containing Phos-tag™. Phos-tag™ is a functional molecule that binds specifically to the phosphate group^{1), 2), 3)}. This product can trap phosphorylated proteins during SDS-PAGE, and it allows detection of phosphorylated and non-phosphorylated proteins as different bands. The gels have a neutral pH to obtain sharp bands.

[Notes]

- It is not recommended to apply ordinary prestained protein ladders to polyacrylamide gels containing Phos-tag™ such as this product because the ladder patterns is disturbed. A prestained ladder recommended is WIDE-VIEW™ Prestained Protein Size Marker III (Code No. 230-02461), which is designed to reduce disturbance of a ladder pattern in Phos-tag SDS-PAGE.
- The gel becomes clouded when incubating it with any buffer including methanol such as acetate-methanol solution for fixing proteins and transfer buffer. It does not affect the quality of the product.
- Handle the gel carefully not to tear it.

[Features]

- Separates phosphorylated and non-phosphorylated proteins
- High resolution and sharp bands
- Long-term stability

[Format]

Plate Size	100 × 100 × 3 (mm)
Gel Size	90 × 85 × 1 (mm)
Well Volume	13wells : 30 μL/well 17wells : 25 μL/well

[Package]

5 gels

[Additional materials required]

- Running Buffer Solution (× 10) (Code No. 184-01291)
0.25 mol/L Tris, 1.92 mol/L Glycine, 1% SDS
- Sample Buffer Solution (2ME +) (× 2) (Code No. 196-11022)
0.125 mol/L Tris-HCl, pH 6.8, 20w/v% Glycerol, 4w/v% SDS, 10vol% 2-Mercaptoethanol, 0.002w/v% BPB
- Electrophoresis apparatus-EasySeparator™
(Code No. 058-07681)
This product is dedicated for EasySeparator™.

[Procedure]

- Remove the gel plate and set it into the apparatus.
- Fill the reservoir with running buffer.
- Carefully remove the comb from the gel.
- Load the samples* into the wells.
* Prepare the sample using the following method. Mix sample with Sample Buffer Solution (2ME +) (× 2) and heat for 5 min at 95 °C. Cool it to room temperature.
- Start the electrophoresis at 20 mA/gel for 1 hour until the BPB reaches the bottom of the gel.
- Remove the gel and proceed to the next step.

[Q & A]

Q1. What gel staining method can be used ?	A1. It is possible to use CBB staining, Silver staining, negative staining and fluorescent staining.
Q2. Can I destain the gel after CBB staining ?	A2. You can destain by using the following microwave method. [Processing method] Transfer the stained gel to a container with 100 mL of deionized water. Place a folded Kimwipe on the gel and heat 1min with microwave. Change the deionized water and repeat 3 times.
Q3. Can Western blot be performed ?	A3. It can be used for Western blot. However, transfer efficiency is poor, so it is necessary to remove zinc by EDTA processing. [Processing method] The gel is soaked in a general transfer buffer containing 5 mmol/L EDTA for 10 minutes with gentle agitation. Repeat the same process again three times. Finally, the gel is soaked in a general transfer buffer without EDTA for 10 minutes with gentle agitation. Transfer on a membrane. If transfer efficiency is poor, increase the treatment time or frequency of exchange of the buffer containing EDTA.
Q4. Why are the protein bands distorted ?	A4. EDTA, inorganic salts, surfactants, etc., cause distortion or tailing of bands. Desalinate by precipitation with TCA or dialysis. Blank lanes also cause distortion. Load the same amount of sample buffer in blank lanes.
Q5. Why aren't the phosphorylated and non-phosphorylated of target proteins separated ?	A5. Load control samples such as β-casein and dephosphorylated β-casein and check that the band shift is occurred. When it doesn't occur, it may not be separated at the Phos-tag™ and/or acrylamide concentration of this product.
Q6. Can cell lysate be used ?	A6. It can be used. It is possible that Rf value is small and separation of bands is unclear. In that case, gels with low-concentration Phos-tag™ should be used.

Q7. What is the amount of protein loaded per well ?	A7. Purified proteins : 1~5 μ g (For CBB staining) Tissue and cell lysate : 10~30 μ g (Adjust the amount of sample according to the expression level of the target protein.)
Q8. What proteins size markers can be used ?	A8. It is impossible to measure molecular weight with size markers accurately. It is recommended to use de-phosphorylated sample of the target protein as a negative control.
Q9. Which Mn complex or Zn complex Phos-tag TM acrylamide is formed in this product ?	A9. It forms a Zinc complex.
Q10. How do I verify if the band shift is due to phosphorylation ?	A10. Perform electrophoresis with the same acrylamide concentration gel without Phos-tag TM and confirm that the band shift is not occurred.

〔Storage〕

Store at 2~10℃ in the dark.

〔References〕

- 1) Kinoshita, E., Kinoshita-Kikuta, E., Takiyama, K. and Koike, T. : *Mol. Cell. Proteomics*, **5**, 749 (2006).
- 2) Yamada, S., Nakamura, H., Kinoshita, E., Kinoshita-Kikuta, E., Koike, T. and Shiro, Y. : *Anal. Biochem.*, **360**, 160 (2007).
- 3) Kinoshita-Kikuta, E., Aoki, Y., Kinoshita, E. and Koike, T. : *Mol. Cell. Proteomics*, **6**, 356 (2007).

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://ftwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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