FUJIFILM**Wako**Code No. 198-18841 (50 μ L)
194-18843 (500 μ L)**For Genetic Research
SAFELOOK™ Green Nucleic Acid Stain**

SAFELOOK™ is a non-carcinogenic and less harmful nucleic acid stain developed to replace ethidium bromide that has mutagenic and/or toxic properties.

SAFELOOK™ Green Nucleic Acid Stain has an excitation wavelength at 490 nm and an emission wavelength at 520 nm (bound to DNA) or 635 nm (bound to RNA). It can be used as both precast and post-staining.

Nucleic Acid : dsDNA, ssDNA, RNA
Applicable Gel : Agarose
Excitation wavelength : 490 nm
Emission wavelength : 520 nm (DNA), 635 nm (RNA)
Sensitivity : 0.1-0.3 ng
Light source : UV, LED

[Protocol]**Precast**

1. Prepare 100 mL of molten agarose gel solution.
2. Add 5 μ L of SAFELOOK™ Green Nucleic Acid Stain to 100 mL of the agarose gel solution and mix gently.
3. Cast the gel and allow it to solidify.
4. Load samples with loading-dye into wells of the gel and run the gel.
If fluorescent signal is weak, add 5-10 μ L of SAFELOOK™ Green Nucleic Acid Stain to 200 mL of running buffer.
5. Visualize the gel using UV or LED light.

Post-staining

1. Run gel according to your protocol.
2. Dilute SAFELOOK™ Green Nucleic Acid Stain
1 : 5000 - 1 : 10000 in TAE or TBE buffer to make staining solution.
Do not use a glass container. Absorption of the dye on the glass surface may cause decrease the staining efficiency.
3. Place the gel in the staining solution.
4. Agitate the gel gently and incubate at room temperature for 10-40 minutes, protected from light.
5. Visualize the stained gel using UV or LED light source.
Although destaining step is not required basically, washing the gel with water may be effective to reduce the background.

[Storage]

In the dark at 2-10°C.

[Packaging]

Code No.	Packaging
198-18841	50 μ L
194-18843	500 μ L

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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