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
Source	: Porcine Pancreas
Appearance	: Lyophilized
Specific Activity	: Indicated on the label
Optimal pH	: pH 7.6

Assay method**1. Principle****2. Procedure**

- Set up dual beam spectrophotometer with temperature control as follows :
Cuvette Temperature : 30 °C
Wavelength : 253 nm
Kinetic program capable of measuring absorbance over 60 seconds.
- Resuspend Trypsin with 200 μL of Dilution Buffer (200-15952) (0.1 mg/mL).
Mix well by vortexing, and incubate at 30 °C for 15 minutes, then place on ice.
- Set up Activity Assays as follows :
 - Prepare a Blank Solution by adding 25 μL of Dilution Buffer (200-15952) to 3.15 mL of Reaction Buffer (see below).
 - Blank spectrophotometer by adding 1 mL of Blank Solution to a quartz cuvette, placing in cuvette sample holder, pressing Auto Zero, and monitoring A_{253} for 1 minute. If the ΔA_{253} of the Blank is $> \pm 0.005$, clean cuvette(s) and repeat.
 - Assay Resuspended Trypsin as follows :
25 μL of Resuspended Trypsin + 3.15 mL Reaction Buffer.
Mix by inverting several times.
Immediately place 1 mL into a quartz cuvette and monitor ΔA_{253} /minute over a period of 1 minute.

(Solutions)

Reaction Buffer : Make these two solutions first-
Buffer A (50 mmol/L Tris-HCl, pH 7.6 (25 °C), 1 mmol/L CaCl_2)
100 mmol/L BAEE (N-Benzoyl-L-arginine Ethyl Ester)-this should be prepared in Buffer A. For 10 mL, dissolve 342.8 mg BAEE in 9 mL of Buffer A. Mix well, then bring to 10 mL. This solution should be aliquoted into 200 μL samples and stored at -20 °C .
Just prior to use, add 150 μL of 100 mmol/L BAEE to 31 mL of Buffer A. Mix thoroughly.

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3. Unit definition

One unit is the amount of enzyme, which increases 0.001 of absorbance at 253 nm per minute at 30 °C in pH 7.6 with BAEE (N-Benzoyl-L-arginine Ethyl Ester) as substrate.

(Calculation)

Calculate Activity and Specific Activity as follows, and record the average value :

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〔Protocols for In-Gel Digestion〕

Use a siliconized microcentrifuge tube and Pipette-Tip to prevent trapping any proteins. Use a gel staining kits for Mass Spectrometry, such as Wako's Silver Stain MS Kit (Code No. 299-58901) and Negative Gel Stain MS Kit (Code No. 293-57701).

1. Separate protein samples by electrophoresis.
2. Cut the protein bands from the gel and slice them into a microcentrifuge tube.
3. Destain the gel. (Use Destaining solution which is included in the gel staining kit for Mass Spectrometry)
4. Add 300 μ L Acetonitrile(ACN) to the tube. Shake the gel pieces for 30 minutes with a mixer(Dehydration).
5. Remove the ACN and cover the microcentrifuge tube with Parafilm.
6. Make some hole on the parafilm and vacuum dry for 15 minutes.
7. Reduce the proteins of the gel pieces with 100 μ L of 10 mmol/L Dithiothreitol (DTT) in 100 mmol/L Ammonium Bicarbonate and incubate for 1 hour at 56℃.
8. After cooling to room temperature, replace the DTT solution with the same volume of 50 mmol/L Iodoacetamide in 100 mmol/L Ammonium Bicarbonate and incubate the gel pieces for 45 minutes in the dark with occasional vortexing.
9. Wash the gel pieces for 10 minutes with 100 μ L of 100 mmol/L Ammonium Bicarbonate.
10. Dehydrate the gel pieces for 15 minutes with 300 μ L ACN.
11. Swell the gel pieces for 15 minutes with 100 μ L of 100 mmol/L Ammonium Bicarbonate.
12. Dehydrate again the gel pieces for 15 minutes with 300 μ L ACN.
13. Remove the liquid phase, and vacuum dry the gel pieces for 15 minutes.

14. Swell the gel pieces for 45 minutes with 100 μ L Trypsin solution* in an ice-cold bath.

※ Stock Solution : Dissolve a bottle of Trypsin in 200 μ L Dilution Buffer (200-15952).

※ Reaction Solution : Dilute 10 μ L Stock Solution in 50 mmol/L Tris-HCl, pH8.5 to 10 μ g/mL.

15. Remove the Trypsin solution, and incubate the gel pieces overnight at 37℃ with 10 μ L of 50 mmol/L Tris-HCl, pH8.5.
16. Shake the gel pieces for 20 minutes with 50 μ L of 20 mmol/L Ammonium Bicarbonate and extract the peptides.
17. Shake the gel pieces three times for 20 minutes with 5% Formic acid in 50% ACN and extract the peptides.
18. If necessary, concentrate the peptides with Speed Vac.
19. Desalt and purify the peptides with ZipTip®.
20. If necessary, concentrate the peptides to 2 μ L with weak vacuum.
21. Add the matrix and analyze it by mass spectroscopy.

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* **Trans** is the Japanese translation.

ZipTip® is registered trademark of Millipore.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
 Telephone : + 81-6-6203-3741
 Facsimile : + 81-6-6201-5964
<http://fwvk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
 Richmond, VA 23237
 U.S.A.
 Telephone : + 1-804-271-7677
 Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
 D-41468 Neuss
 Germany
 Telephone : + 49-2131-3111-0
 Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>