FUJIFILM

Wako

Code No. 208-17371 (100 mL)

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This product contains 3, 3', 5, 5'-Tetramethylbenzidine (TMB) and hydrogen peroxide in a mildly acidic buffer.

Unreacted substrate should be colorless or very light yellow in appearance. When this substrate system is reacted with peroxidase, a soluble blue reaction product is obtained. The reaction can be stopped using the stop solution, producing a soluble yellow reaction product.

This product is not used for membrane or immunohistochemical applications that require a precipitating reaction product.

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[Usage]

This product is supplied as a ready-to-use solution.

The product should be allowed to equilibrate to room temperature (25°C) prior to use. For microwell applications, 100 μ L of substrate solution is added to each well. The substrate reaction can be stopped using equal volumes of 1 mol/L HCl or 0.3 mol/L H₂SO₄. After stopping, a soluble yellow product develops which is read in the 450 nm range.

Dilution of the substrate is not recommended. To reduce the intensity of a reaction, it is recommended that the antibody or conjugates be diluted.

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— 1/2 —

[Expiration date] Indicated on the label.

[Storage] Store at 2 ~ 10°C in the dark.

[Package] 100 mL

* : **[Trans]** is the Japanese translation.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-311-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

— 2/2 —

2308KA2