FUJIFILM

Wako

<For Research Use Only>

Code No. 296-62701 (20 reactions)
290-62704 (500 reactions)

for Genetic Research

S. cerevisiae Direct Transformation Kit *Wako*

[Introduction]

The *S. cerevisiae* Direct Transformation Kit *Wako* is specially designed for rapid and easy transformation of the budding yeast, *Saccharomyces cerevisiae*.

The *S. cerevisiae* Direct Transformation Protocol allows for successful transformation simply by mixing a plasmid DNA and the kit reagent with overnight yeast cell cultures. No complicated competent cell preparation is required. The traditional yeast transformation methods require repeated centrifugation or washing steps, thus are problematic for high-throughput transformation. The *S. cerevisiae* Direct Transformation Kit *Wako* is well suited for high-throughput transformation of a large number of yeast strains grown in 96-well plates. You may choose either a Tube protocol or a 96-well plate protocol for your purpose.

[Features]

- Direct transformation : One-step procedure by mixing your plasmid and the kit reagent with yeast culture.
- Simply : No need to prepare competent cells.
- High-throughput : Suitable for transformation of a large number of yeast strains using a 96-well plate.

[Transformation Efficiency]

Table 1. Transformation efficiency of BY4743 strain with pRS316 plasmid (*URA3*, *CEN/ARS*).

	Transformation efficiency (cfu/ μ g)
Tube protocol	$> 5 \times 10^3$
96-well plate protocol	$> 5 \times 10^2$

Transformation efficiency is calculated as the number of colonies per microgram of plasmid DNA

[Kit contents]

	(20 reactions)	(500 reactions)
1) Sc Transformation Reagent	500 μ L	12.5 mL
2) Carrier DNA (5 μ g/ μ L)	50 μ L	1.25 mL

[Storage]

Store at -20°C

※The number of times of freeze-thaw should not exceed 10 times. If you do transformation more than 10 times, dispense the reagents and store at -20°C .

[Protocol]

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I-1. Preparation of Yeast culture

<Required Reagents and Equipments>

- Shaking incubator ($28-30^{\circ}\text{C}$)
- Sterile two position type tube (14 mL)
- YPD medium

[Preparation of YPD medium, 1L]

Dissolve the following in distilled water to 1 L. Add 20 g agar for plates.

Yeast extract	10 g
Peptone	10 g
Glucose	20 g

Autoclave for 20 minutes at 121°C .

<Before Starting>

Streak a YPD plate with your *S. cerevisiae* strain and incubate the plate at $28-30^{\circ}\text{C}$ for 2-3 days.

<Yeast Culture Protocol>

1. Inoculate the colonies of your *S. cerevisiae* strain into 2 mL of YPD medium. Measure its OD at 600 nm. After adjusting the OD₆₀₀ to 0.01 with medium, transfer the 1 mL of its suspension to the 14 mL two position type tube. Cover the tube with the lid loosely.

Note : One mL of culture can be used for approximately 30 transformations (25 μ L/transformation). If you want to prepare over 30 transformations, increase the number of culture tubes but not the volume. The volume of culture significantly affects the aeration and growth conditions.

- We highly recommend to perform yeast culture with a reciprocal or rotary shaking incubator at 150 rpm for 18-21 hours. However, yeast cell growth may vary from strain to strain, the value of OD₆₀₀ should be between 3.5 and 5.0 at the end of culture. Cell concentration plays an important role in transformation efficiency.
- The OD₆₀₀ should be measured by 10-fold to 50-fold dilution of the culture.

2. Incubate at $28-30^{\circ}\text{C}$ for ~ 24 hours with vigorous shaking. The culture should reach a cell density to an OD₆₀₀ of 3.5-5.0.

I-2. Transformation

<Required Reagents and Equipments>

- Incubator or water bath (42°C)
- Incubator ($28-30^{\circ}\text{C}$)
- Sterile microcentrifuge tube (1.5 mL)
- Appropriate selective plates for your strain and plasmid
- Sterile water

<Before Starting>

Thaw completely the Sc Transformation Reagent and Carrier DNA at room temperature.

<Transformation Protocol>

1. Add 1 μg of plasmid DNA¹⁾ and 2 μL of Carrier DNA to 20 μL of Sc Transformation Reagent²⁾ in a microtube, and then add sterile water to give a final volume of 25 μL . Mix by vortexing.

- 1) The quality of DNA used for this method plays an important role in transformation efficiency. We recommend the QIAGEN plasmid Kit, EndoFree Plasmid Kit, or CONCERT™ High Purity Plasmid Prep Purification System for plasmid DNA purification.
 - 2) Slow pipette operation is required because of the viscosity of this reagent.

2. Add 25 μL of yeast culture to the mixture and mix by vortexing.
3. Incubate at 42 °C for 2 hours.
4. Add 150 μL of sterile water to the mixture and spread on an appropriate selection plate using a sterile spreader. If you want to obtain single colonies, dilute the mixture to 20-200 fold with sterile water and spread 100 μL of the diluted mixture on an appropriate selection plate.
5. Incubate the plate at 28-30 °C for 2-3 days.

II. Protocol 2 (96-well plate protocol)

II-1. Preparation of Yeast culture

<Required Reagents and Equipments>

- Incubator (28-30 °C)
- Sterile 96-well microtiter plate (Nunc cat #: 267245)
- Microplate sealing tape
- Sterile 96-pin replicator (GENETIX, UK. Cat #: X5052)
- Square plate (Nissui Pharmaceutical Co., LTD: Nissui Kaku Schale)
- Sterile water
- YPD medium (See Section I for preparation (p.2))

<Before Starting>

Inoculate a YPD plate with your *S. cerevisiae* strain from glycerol stock using a sterile 96-pin replicator. Incubate the plate at 28-30 °C for 2-3 days.

<Yeast culture Protocol>

1. Aliquot 25 μL YPD medium into each well of a 96-well microtiter plate.
2. Pick *S. cerevisiae* strains from a square YPD plate culture using a 96-pin replicator and inoculate YPD medium in a 96-well plate.
3. Seal the plate with a microplate sealing tape and incubate at 28-30 °C for 18-24 hours without shaking.

II-2. Transformation

<Required Reagents and Equipments>

- Incubator (28-37 °C)
- Multi channel pipettor
- Sterile 25 mL reservoir
- Appropriate selective plates for your strain and plasmid
- Sterile water

<Before Starting>

Thaw completely Sc Transformation Reagent and Carrier DNA at room temperature.

<Transformation Protocol>

1. Add 1 μg of plasmid DNA¹⁾ and 2 μL of Carrier DNA to 20 μL of Sc Transformation Reagent²⁾, and then add sterile water to give a final volume of 25 μL per 1 well. In the case for preparing multiple wells, the premix can be prepared as follows.

	1 well	96 wells	X wells
Sc Transformation Reagent ¹⁾	20 μL	2.4 mL	$20 \times 1.25 \times X \mu\text{L} =$
plasmid DNA ²⁾	1 μg	120 μg	$1 \times 1.25 \times X \mu\text{g} =$
Carrier DNA (5 $\mu\text{g}/\mu\text{L}$)	2 μL	240 μL	$2 \times 1.25 \times X \mu\text{L} =$
Sterile water	up to 25 μL	3 mL	$25 \times 1.25 \times X \mu\text{L} =$

2. Suspend precipitated cells in the cultured 96-well microtiter plate using a plate mixer³⁾. Immediately add 25 μL of the DNA mixture (Step 1) to each well of yeast culture using multi-channel pipettor. Seal the plate and mix using a plate mixer.
3. Incubate at 37 °C for 2 hours.
4. Spot 5-10 μL of the mixture on the appropriate selection plate.
5. Incubate the selection plate at 28-30 °C for 2-3 days.

- 1) The quality of DNA used for this method plays an important role in transformation efficiency. We recommend QIAGEN plasmid Kit, EndoFree Plasmid Kit, or CONCERT™ High Purity Plasmid Prep Purification System for plasmid DNA purification.
 - 2) Slow pipette operation is required because of the viscosity of this reagent.
 - 3) If you do not have a plate mixer, a multi-channel pipettor can be used.

[Appendix (Flow Chart of the Experiments)]

I. Tube protocol

2 mL of YPD medium
↓ ← Inoculate yeast colonies
Adjust OD₆₀₀ : 0.01 by medium
↓
1 mL of the suspension
↓ Incubate at 28-30 °C with shaking
Yeast culture

1.5 mL of microcentrifuge tube
↓ { 20 µL of Sc Transformation Reagent
1 µg of plasmid DNA
2 µL of Carrier DNA
Up to 25 µL by sterile water
Mix by vortexing
← 25 µL of yeast culture
Incubate at 42 °C for 2 hr.
↓ ← 150 µL sterile water
Plating
↓ Incubate at 28-30 °C for 2-3 days
Transformed yeast

II. 96-well plate protocol

96 microtiter well
↓ ← 25 µL of YPD
Inoculate yeast strains into each well
Incubate at 28-30 °C for 18-24 hr. without shaking
Suspend precipitated cells using plate mixer
↓ { 20 µL of Sc Transformation Reagent
1 µg of plasmid DNA
2 µL of Carrier DNA
Up to 25 µL per 1 well by sterile water
Mix with plate mixer
Incubate at 37 °C for 2 hr.
↓
Spotting 5-10 µL
↓ Incubate at 28-30 °C for 2-3 days
Transformed yeast

Allied products

Wako catalog No.	Description	Package Size
045-31162	D(+)-Glucose	25 g

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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