Code No. 294-80201(for 100mL)
290-80203(for 1L)

ScreenFect™UP-293

ScreenFect™UP-293 is a transfection reagent kit which is highly efficient for production of antibodies and proteins in Expi293F™ or FreeStyle™ 293-F cells. ScreenFect™ UP-293 Booster, which is added to transfected cells, enhances remarkably protein expression level in cells.

【Kit contents】

	For 100mL	For 1L
ScreenFect™UP-293 Transfection Reagent	0.2mL	2x1mL
ScreenFect™UP-293 Dilution Buffer	6.7mL	67mL
ScreenFect™UP-293 Booster	0.5mL	5mL

【Storage】

Store at 2-10°C. Do not freeze.

【 Additional required reagents 】

	For Expi293™ Expression System	For FreeStyle™ 293 Expression System
Cells	Expi293F™ Cells	FreeStyle™ 293-F Cells
Medium	Expi293™ Expression Medium	FreeStyle™ 293 Expression Medium
Dilution Buffer	-	Opti-MEM™ I Reduced Serum Medium

【 Attention】

- Please be careful of bacteria contamination under experiment.
- After storage of several months without using the reagent, slight precipitation in ScreenFect™ UP-293 Transfection Reagent might occurs. Vortex it well before use.
- Optimization study of amounts of ScreenFect™UP-293 Transfection Reagent, plasmid DNA and cell culture time is required for successful protein expression.


Expi293F™, Freestyle™ and Opti-MEM™ are registered trademarks of Thermo Fisher Scientific.


Transfection protocol for 30 mL culture volume


- ①  One day prior to transfection, prepare culture cells to be appropriate cell numbers for transfection using a 125mL or 250mL of Erlenmeyer shake flask .
Note: Recommended seeding cell density is $>0.8 \times 10^6$ viable cells/mL.

- ②  Prepare cell suspension
Just before transfection, prepare 28mL of cell suspension with the cell densities below in a 125 mL of Erlenmeyer shake flask.
【For Expi293F™ Cells】
2.7x10⁶ viable cells/mL medium (7.5x10⁷ viable cells/28mL of medium)
【For FreeStyle™ 293-F Cells】
1.1x10⁶ viable cells/mL medium (3.0x10⁷ viable cells/28mL of medium)
Note: cell viability must not be < 90 %.

- ③  Mix ScreenFect™UP-293 Transfection Reagent well with a vortex mixer.

- ④  Prepare a diluted transfection reagent
【For Expi293F™ Cells】
Add 60μL of ScreenFect™UP-293 Transfection Reagent into 940μL of ScreenFect™UP-293 Dilution Buffer in a sterile microfuge tube and mix it using pipette strokes.
【For FreeStyle™ 293-F Cells】
Add 30μL of ScreenFect™UP-293 Transfection Reagent into 970μL of Opti-MEM™ I Reduced Serum Medium in a sterile microfuge tube and mix it using pipette strokes.

- ⑤  Prepare a diluted DNA solution
【For Expi293F™ Cells】
Add 30μg of plasmid DNA with a density of 1μg/μL into 970μL of ScreenFect™UP-293 Dilution Buffer in a sterile microfuge tube and mix it using pipette strokes.
【For FreeStyle™ 293-F Cells】
Add 30μg of plasmid DNA with a density of 1μg/μL into 970μL of Opti-MEM™ I Reduced Serum Medium in a sterile microfuge tube and mix it using pipette strokes.

- ⑥  Prepare a lipoplex solution
Mix the diluted transfection reagent and the diluted DNA solution with rapid pipette strokes and spin down it. Incubate it for 5-20minutes at room temperature.
Note: Do not vortex the mixed solution.

- ⑦  Add the 2mL of the lipoplex solution into the cell suspension in the Erlenmeyer Shaker Flask with slowly swirling.
Note: The final volume should be approximately 30mL in the flask.

- ⑧  Incubate the cells at 37°C(8% CO₂), rotating at 125rpm on an orbital shaker for 16hours.

- ⑨  Add 150μL of ScreenFect™UP-293 Booster to the cells.

- ⑩  Continue to incubate the cells at 37°C(8% CO₂), rotating at 125rpm on an orbital shaker .
Collect cells or culture supernatants containing target proteins in 48-120 hours after transfection.
Note: The incubation period required for successful transfection depends on the target protein and cell viability.

Expi293F™, Freestyle™ and Opti-MEM™ are registered trademarks of Thermo Fisher Scientific.

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