Code No. 290-80301 (1 Kit : 0.5 mL Slurry)

For Genetic Research PS Capture™ Exosome Isolation Resin Kit

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This kit can easily isolate high purity extracellular vesicles (EVs) such as exosomes having phosphatidylserine (PS) on the membrane surface by using PS Affinity method, which use a resin immobilizing PS-binding protein.

This kit uses nonmagnetic resin and captures extracellular vesicles metallic ion-dependently and elute extracellular vesicles in nearly intact state under neutral condition using a chelating agent similarly to MagCapture™ Exosome Isolation Kit PS (Code No. 293-77601) that utilizes magnetic beads.

【Feature】

- Highly pure extracellular vesicles are obtained by PS affinity method.
- Extracellular vesicles are obtained in nearly intact state.
- Purification is possible simply by centrifugation without a magnetic stand.
- An optimal protocol can be set by the amount of a sample (1.2 mL to 40 mL).
- Ultrafiltration concentration of a sample is unnecessary.

【Kit Contents】

This kit includes 5 components.

- 1 Kit : 0.5 mL Slurry
- | | |
|--------------------------------|------------------|
| (1) Exosome Capture Resin | 0.5 mL × 1 tube |
| (2) Exosome Basic Buffer (10×) | 50 mL × 1 bottle |
| (3) Binding Enhancer (100×) | 10 mL × 1 bottle |
| (4) Washing Enhancer (100×) | 5 mL × 1 bottle |
| (5) Eluting Enhancer (100×) | 0.5 mL × 1 tube |


【Storage】

Store at 2-10°C

【Necessary Equipment】

- Centrifuge
 - For preparation of samples : A cooling type at 10,000×g centrifugation force
 - For isolation of exosome : With a swing rotor at 2,000 to 3,000×g centrifugation force (applicable to 5 mL centrifuge tube, 15 mL centrifuge tube or 50 mL centrifuge tube)
At 2,000 to 3,000×g centrifugation force (1.5 to 2 mL microtube)
- Tabletop centrifuge
- Vortex mixer
- Rotator

【Outline of Procedure】

【Preparation of sample】

A simple centrifugation shown below is recommended with this kit to obtain each main fraction.

- [1] Use a supernatant obtained by centrifugation at 10,000×g as a sample to isolate Small Extracellular Vesicles (Small EVs) such as exosomes.
- [2] To isolate Large Extracellular Vesicles (Large EVs) such as microvesicles, obtain a supernatant centrifuged at 1,200×g first and, then, centrifuge this at 10,000×g and suspend the obtained precipitate in TBS and then use it as a sample.
- [3] Use the supernatant centrifuged at 1,200×g as the sample to isolate both Small EVs and Large EVs together.


Centrifugation Condition

1. Recover the cell culture medium, centrifuge at 4°C, 300×g, for 5 minutes, and transfer the supernatant to a new tube.
2. Centrifuge the supernatant from step-1 at 4°C, 1,200×g, for 20 minutes, and transfer the supernatant to a new tube (**1,200×g supernatant**).
3. Centrifuge the supernatant from step-2 at 4°C, 10,000×g, for 30 minutes, and transfer the supernatant to a new tube (**10,000×g supernatant**).

Note [1] : If a fetal bovine serum (FBS) is added to a medium for cell culture, use a FBS after elimination of extracellular vesicles by ultracentrifugation (example: 18 hours centrifugation at 110,000×g) or commercially available extracellular vesicle-free FBS.

Note [2] : Use fresh cell culture supernatant as a sample.

(Flowchart of Sample Preparation)

[Isolation Process]

We recommend the following conditions with **a sample volume equivalent to 80% of the container capacity as ideal to avoid an extra space in the binding reaction container** since extracellular vesicles including exosomes are likely to adsorb to containers.

Table-1 Recommended conditions by a container

Binding reaction container	Sample volume	Amount of resin (slurry)	Volume of washing buffer after reaction	Volume of elution buffer
Remarks	80% of the binding reaction container	Sample volume × 1/200 to 1/100	Sample volume/time	Amount of resin (slurry) × 2.5 times/time
1.5 mL tube	1.2 mL	20 μ L (minimum volume)	1.2 mL/time × 3 times	50 μ L (minimum volume)/time × 2 times
2 mL tube	1.6 mL	20 μ L (minimum volume)	1.6 mL/time × 3 times	50 μ L (minimum volume)/time × 2 times
5 mL tube	4.0 mL	20 μ L ~ 40 μ L	4 mL/time × 3 times	50 μ L ~ 100 μ L/time × 2 times
15 mL tube	12 mL	60 μ L ~ 120 μ L	12 mL/time × 3 times	150 μ L ~ 300 μ L/time × 2 times
50 mL tube	40 mL	200 μ L ~ 400 μ L	40 mL/time × 3 times	0.5 mL ~ 1.0 mL/time × 2 times

Preparation of each reagent

· Preparation of washing buffer

Prepare the “washing buffer” of a necessary volume by adding 1/10-volume of Exosome Basic Buffer (10×), 1/100-volume of Binding Enhancer (100×), and 1/100-volume of Washing Enhancer (100×) to purified water.

Note [3] : The washing buffer is also used for washing the resin in the initial process in addition to the volumes in the above table. Prepare the washing buffer of that volume (20 times the amount of resin) as well.

· Preparation of elution buffer

Prepare the “elution buffer” of a necessary volume by adding 1/10-volume of Exosome Basic Buffer (10×) and 1/100-volume of Eluting Enhancer (100×) to purified water.

Operation

In the isolation of extracellular vesicles, perform centrifugation in each process under the following conditions according to the sample volume.

Table-2 Centrifugation conditions by a container

Container	Recommended Sample volume	Centrifugation condition
1.5 mL tube	1.2 mL	3,000 × g, 1 minute, room temperature (※)
2 mL tube	1.6 mL	3,000 × g, 1 minute, room temperature (※)
5 mL tube	4.0 mL	2,000 - 3,000 × g, 5 minutes, room temperature, swing rotor
15 mL tube	12 mL	2,000 - 3,000 × g, 10 minutes, room temperature, swing rotor
50 mL tube	40 mL	2,000 - 3,000 × g, 10 minutes, room temperature, swing rotor

(※) Tabletop centrifuge can be used alternatively.

1. Add the washing buffer of 20 times the volume of a resin slurry in a tube (example : 400 μ L of the washing buffer to 20 μ L of the resin slurry), add Exosome Capture Resin (resin slurry homogenized by vortex immediately before), vortex lightly and centrifuge (see Table-2 above for the centrifugation condition).
2. Remove the supernatant as much as possible (aspirate loosely and, then, remove the supernatant using a pipet of a small volume, taking care not to aspirate the resin as much as possible).
3. Add the sample to the tube, add 1/100-volume of the Binding Enhancer (100 ×), and agitate using a rotator for binding of the resin and extracellular vesicles in the sample (refrigerated, 3 hours to overnight).
4. Centrifuge the sample-resin mixture after binding to precipitate the resin (see Table-2 above for centrifugation condition).
5. Remove the supernatant as much as possible (aspirate loosely and, then, remove the supernatant using a pipet of a small volume, taking care not to aspirate the resin as much as possible).

Note [4] : When further recovery of remaining extracellular vesicles is needed, please store the used supernatant.

6. Add the washing buffer and vortex lightly (see Table-1 above for the volume of the washing buffer to be added), and centrifuge to precipitate the resin (see Table-2 above for centrifugation condition).
7. Remove the supernatant.
8. Repeat steps (6 and 7) once.
9. Add the washing buffer and vortex lightly (see Table-1 above for the volume of the washing buffer to be added), and centrifuge to precipitate the resin (see Table-2 above for centrifugation condition).

Note [5] : When a binding reaction container other than 1.5 mL tube is used, transfer the resin to a new 1.5 mL tube for elution step. Leave the centrifuged supernatant of about the same amount as the resin and suspend them, and transfer it to a new 1.5 mL tube, then centrifuge for 1 min to precipitate the resin.

10. Remove the supernatant as much as possible (aspirate loosely and, then, remove the supernatant using a pipet of a small volume, taking care not to aspirate the resin as much as possible).
11. Add the “elution buffer” to the resin after washing and vortex (see Table-1 above for the volume of the elution buffer to be added), and leave to stand for 10 minutes at room temperature.
12. Vortex again after 10 minutes, centrifuge at 3,000 × g for 1 minute at room temperature (tabletop centrifuge can be used alternatively), and transfer the supernatant to a new 1.5 mL or 2 mL tube, taking care not to aspirate the resin as much as possible.

13. To the remaining resin add the “elution buffer” again and vortex (see Table-1 above for the volume of the elution buffer to be added), and leave to stand for 10 minutes at room temperature.


14. Vortex again after 10 minutes, centrifuge at $3,000 \times g$ for 1 minute at room temperature (tabletop centrifuge can be used alternatively), and transfer the supernatant to combine with the supernatant from step 12, taking care not to aspirate the resin as much as possible.

Note [6] : Exosome Capture Resin is reusable after eluting EVs (up to 2 times). Therefore, when you need to recover more remaining EVs from samples, please repeat from the step-3 in **[Isolation Process]** to increase the recovered amount. After collecting the eluate, add the washing buffer and transfer the Exosome Capture Resin (resin slurry homogenized by vortex immediately before) to the appropriate reaction container. After removing the washing buffer, add the used sample stored in Step-5 and start with Step-3. In this case, it is unnecessary to add 1/100-volume of the Binding Enhancer (100×) again. In addition, this kit contains reagents necessary for reuse twice.

15. As necessary, filter out with a $0.2 \mu\text{m}$ filter to remove the resin (particle size : about $34 \mu\text{m}$) completely, or for bacteria removal (centrifugal filtration filter is preferable).

Recycling of Exosome Capture Resin after eluting of extracellular vesicles
Used Exosome Capture Resin can be recycled for “Repeated extraction of extracellular vesicles from the same sample (refer to note [6])” and “Purification of extracellular vesicles from sample of the same lot”. When you need to purify extracellular vesicles from sample of the same lot, prepare sample first and start with the step-1 in **[Isolation Process]**. Up to the second reuse, extracellular vesicles can be recovered almost equally. In the third to fourth reuse, it is confirmed that the amount of extracellular vesicles in flow-through increases and the recovery performance decays. When you need to store the used Exosome Capture Resin, please add the Washing buffer in the kit or $1 \times \text{TBS}$ containing 0.05 w/v% sodium azide, and suspend it. After suspension, keep it at $2-10^\circ\text{C}$. Please use it as soon as possible.

[Flowchart of Isolation Process]

Binding reaction container	Recommended Sample volume	Amount of resin (slurry)	Volume of washing buffer after reaction	Volume of elution buffer
Remarks	80% of the binding reaction container	Sample volume \times 1/200 to 1/100	Sample volume / time	Amount of resin (slurry) \times 2.5 times / time
1.5 mL tube	1.2 mL	$20 \mu\text{L}$ (minimum)	$1.2 \text{ mL} / \text{time} \times 3 \text{ times}$	$50 \mu\text{L}$ (minimum) / $\text{time} \times 2 \text{ times}$
2 mL tube	1.6 mL	$20 \mu\text{L}$ (minimum)	$1.6 \text{ mL} / \text{time} \times 3 \text{ times}$	$50 \mu\text{L}$ (minimum) / $\text{time} \times 2 \text{ times}$
5 mL tube	4.0 mL	$20 \mu\text{L} \sim 40 \mu\text{L}$	$4 \text{ mL} / \text{time} \times 3 \text{ times}$	$50 \mu\text{L} \sim 100 \mu\text{L} / \text{time} \times 2 \text{ times}$
15 mL tube	12 mL	$60 \mu\text{L} \sim 120 \mu\text{L}$	$12 \text{ mL} / \text{time} \times 3 \text{ times}$	$150 \mu\text{L} \sim 300 \mu\text{L} / \text{time} \times 2 \text{ times}$
50 mL tube	40 mL	$200 \mu\text{L} \sim 400 \mu\text{L}$	$40 \text{ mL} / \text{time} \times 3 \text{ times}$	$0.5 \text{ mL} \sim 1.0 \text{ mL} / \text{time} \times 2 \text{ times}$

Container	Recommended Sample volume	Centrifugation condition
1.5 mL tube	1.2 mL	3,000 × g, 1 minute, room temperature (※)
2 mL tube	1.6 mL	3,000 × g, 1 minute, room temperature (※)
5 mL tube	4.0 mL	2,000 - 3,000 × g, 5 minutes, room temperature, swing rotor
15 mL tube	12 mL	2,000 - 3,000 × g, 10 minutes, room temperature, swing rotor
50 mL tube	40 mL	2,000 - 3,000 × g, 10 minutes, room temperature, swing rotor

(※) Tabletop centrifuge can be used alternatively.

[Related Products]

Code No.	Description	Size
299-77603	MagCapture™ Exosome Isolation Kit PS	2 purifications
293-77601		10 purifications
297-79201	PS Capture™ Exosome ELISA Kit (Anti Mouse IgG POD)	96 reactions
297-79701	PS Capture™ Exosome Flow Cytometry Kit	300 reactions
016-27061	Anti CD63, Monoclonal Antibody (3-13)	20 μL
012-27063	Anti CD63, Monoclonal Antibody (3-13)	100 μL
290-35591	Magnet Stand	1 each
295-71701	microRNA Extractor SP Kit	50 reactions
317-90175	10 × TBS (pH 7.4)	500 mL

[References]

- 1) W. Nakai, *et al.*, *Sci Rep* **6**, 33935 (2016).
- 2) S. Osawa, *et al.*, *Biochem. Biophys. Res. Commun.*, **488** (1), 232-238 (2017).
- 3) S. Nagashima, *et al.*, *J. Virol.*, **91** (22), Published Online (2017).
- 4) T. Yoshida, *et al.*, *Curr. Protoc. Cell Biol.*, Published Online (2017).
- 5) S. Saito, *et al.*, *Sci Rep* **8**, 3997 (2018).

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
 Telephone : + 81-6-6203-3741
 Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
 Richmond, VA 23237
 U.S.A.
 Telephone : + 1-804-271-7677
 Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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