LBIS™ Der f 2 ELISA Kit

Please, read this instruction carefully before use.

1. Intended use

LBIS™ Der f 2 ELISA Kit is a sandwich ELISA system for quantitative measurement of Der f 2 (*Dermatophagoides farina*). This is intended for research use only.

2. Assay principle

In FUJIFILM Wako Pure Chemical Corporation's LBIS™ Der f 2 ELISA Kit, standards or diluted samples are incubated in monoclonal antibody coated wells to capture Der f 2. After 1 hour incubation and washing, Peroxidase-conjugated anti-Der f 2 antibody is added and incubated for 1 hour. After washing, HRP-complex remaining in wells is reacted with a TMB Solution for 20 minutes, and reaction is stopped by addition of acidic solution, and absorbance of yellow product is measured spectrophotometrically at 450 nm. The absorbance is nearly proportional to Der f 2 concentration. The standard curve is prepared by plotting absorbance against standard Der f 2 concentrations. Der f 2 concentrations in unknown samples are determined using this standard curve.

3. Performance characteristics

- Assay range

The assay range of the kit is 0.78 ng/mL - 50 ng/mL.

- Specificity [(D) Peroxidase conjugated Antibody Solution]

The kit uses two monoclonal antibodies specific to Der f 2.

Sample	Cross reaction
Der f 1	less than the detection limit.

- Precision of assay

Within assay variation (2 samples, 5 replicates assay) Mean CV was within 10%.

- Reproducibility

Between assay variation (3 samples, 4 days, duplicate assay) Mean CV was within 10%.

- Recovery test

Der f 2 was added in 3 concentrations to 2 samples and was assayed.

The recoveries were 95.1% - 104%

- Dilution test

2 samples were serially diluted by 3 steps.

The dilution curves showed excellent linearity. ($R^2=0.9959\&0.9996$)

4. Precautions / Technical tips

- For professional use only. Beginners are advised to use this kit under the guidance of experienced person.
- Wear gloves and laboratory coats when handling assay materials.
- Use clean laboratory glassware.
- Do not drink, eat or smoke in the areas where assays are carried out.
- Handle the standard solution and assay samples with great care as ones potentially containing allergenic substances.
- This kit contains components of animal origin. These materials should be handled as potentially infectious.
- Be careful not to allow the reagent solutions of the kit to touch the skin, eyes and mucus membranes. Especially be careful for the stop solution because it is 1 M sulfuric acid. The stop solution and the substrate solution may cause skin/eyes irritation. In case of contact with these wash skin/eyes thoroughly with water and seek medical attention, when necessary.
- Avoid contact with the acidic stop solution and TMB Solution, which contains hydrogen peroxide and tetramethylbenzidine. Wear gloves and eye and clothing protection when handling these reagents.
- The materials must not be pipetted by mouth.
- Residual samples and used tips should be rinsed in 1% formalin, 2% glutal aldehyde, or more than 0.1% sodium hypochlorite solution for more than 1 hour, or be treated by an autoclave before disposal.
- Dispose consumable materials and unused contents in accordance with applicable regional / national regulatory requirements.
- In order to avoid dryness of wells, contamination of foreign substances and evaporation of dispensed reagents, never forget to cover the well plate with a plate seal supplied, during incubation.

- ELISA can be easily affected by your laboratory environment. Room temperature should be at 20°C - 25°C strictly. Avoid airstream velocity over 0.4 m/sec. (including wind from air conditioner), and humidity less than 30%.
- In manual operation, proficiency in pipetting technique is recommended.
- The reagents are prepared to give accurate results only when used in combination within the same box. Therefore, do not combine the reagents from kits with different lot numbers. Even if the lot number is the same, it is best not to mix the reagents with those that have been preserved for some period.
- Be careful to avoid any contamination of assay samples and reagents. We recommend the use of disposal pipette tips, and 1 tip for 1 well.
- Optimally, the reagent solutions of the kit should be used immediately after reconstitution. Otherwise, store them in a dark place at 2°C - 10°C.
- Time the reaction from the pipetting of the reagent to the first well.
- Prepare a standard curve for each assay.
- Dilution of the assay sample must be carried out using the buffer solution provided in the kit.
- The TMB Solution should be almost clear pale blue before use. It turns blue during reaction, and gives yellowish color after addition of stop solution. Greenish color means incomplete mixing.
- To avoid denaturation of the coated antibody, do not let the plate go dry.
- As the antibody-coated plate is module type of 8 wells × 12 strips, each strip can be separated by cutting the cover sheet with a knife and used independently.
- When ELISA has to be done under the airstream velocity over 0.4 m/sec. and the humidity less than 30%, seal the well plate with a plate seal and place the well plate in an incubator or a styrofoam box in each step of incubation.

5. Reagents supplied

Components	Use Status	Amount
(A) Antibody-coated plate	Use after washing	96 wells/1 plate
(B) Der f 2 Standard	Concentrated. Use after dilution	200 μ L/1 vial
(C) Buffer	Ready for use.	60 mL/1 bottle
(D) Peroxidase-conjugated Antibody Solution	Concentrated. Use after dilution.	200 μ L/1 vial
(F) TMB Solution	Ready for use.	12 mL/1 bottle
(G) Buffer solution for extraction	Ready for use.	100 mL/1 bottle
(H) Stop Solution	Ready for use.	12 mL/1 bottle
(I) Wash Solution (10 \times)	Concentrated. Use after dilution.	100 mL/1 bottle
(J) Plate seal	—	3 sheets

[Storage and stability]

[(A) Antibody-coated Plate]

If seal is not removed, put the strip back in a plastic bag with zip-seal originally used for well-plate container and store at 2°C - 10°C. The strip will be stable until expiration date.

[(B) Der f 2 Standard]

Standard solutions prepared above should be used as soon as possible, and should not be stored.

[(C) Buffer] & [(F) TMB Solution] & [(G) Buffer solution for extraction]

If not opened, store at 2°C - 10°C. It maintains stability until expiration date. Once opened, we recommend using them as soon as possible to avoid influence by environmental condition.

[(D) Peroxidase-conjugated Antibody Solution]

Unused working solution (already diluted) should be disposed.

[(H) Stop Solution]

Close the stopper tightly and store at 2°C - 10°C. It maintains stability until expiration date.

[(I) Wash Solution (10 \times)]

The rest of undiluted buffer : if stored tightly closed at 2°C - 10°C, it is stable until expiration date.

Dispose any unused diluted buffer.

6. Equipments or supplies required but not supplied ☐ Use as a check box

- ☐ Deionized water (or Distilled water) ☐ Test tubes for preparation of standard solution series.
- ☐ Glassware for dilution of Wash Solution (10 \times) (a graduated cylinder, a bottle)
- ☐ Pipettes (disposable tip type). One should be able to deliver 10 μ L - 50 μ L precisely, and another for 50 μ L - 500 μ L.
- ☐ Syringe-type repeating dispenser like Eppendorf multipipette plus which can dispense 100 μ L.
- ☐ Paper towel to remove washing buffer remaining in wells. ☐ A vortex-type mixer.
- ☐ A shaker for 96 well-plate (600 rpm - 1200 rpm) ☐ An automatic washer for 96 well-plate (if available), or a wash bottle

with a jet nozzle (refer to our web movie [Washing of microplate]).

☐ A 96 well-plate reader (450 nm \pm 10 nm, 620 nm : 600 nm - 650 nm) ☐ Software for data analysis.

7. Preparation of samples

- Perform pretreatment (extraction) of the collected assay sample using (G) Buffer solution for extraction.

Dilute the collected assay sample 20 folds (W/V) with (G) Buffer solution for extraction, and stir at room temperature for 30 minutes. Then, remove insoluble matters from the diluted assay sample with centrifugal filtration tube (0.45 μ m) or centrifugation, and use the clarified sample as the extracted sample. Dilute the extracted sample with (C) Buffer so that the result can fall within the standard curve, and use 100 μ L of the diluted extracted sample for each well.

* To collect the sample, use commercially available collection filter, etc.

(no filters for sample collection supplied in this kit)

Make sure to dilute samples more than 1 : 10. Dilution should be carried out with the (C) Buffer solution of the kit using small test tubes before assay.

Example of dilution :	Rate	1 : 10	1 : 100	1 : 200	1 : 400	...
Extracted sample :		50 μ L	* 50 μ L	* 250 μ L	* 250 μ L	
(C)Buffer :		450 μ L	450 μ L	250 μ L	250 μ L	

*one rank higher diluted sample

[Handling of samples]

- For extraction of samples, use of (G) Buffer solution for extraction supplied in this kit is recommended.
- Samples should be diluted with (C) Buffer before use. Sample dilution should be carried out with the (C) Buffer using small test tube such as PP, PE or glass, before assay.
- Turbid samples or ones containing insoluble matters should be clarified before use.
- Either (G) Buffer solution for extraction supplied in this kit or the extraction method described here, or both do not guarantee complete extraction for all the samples.
- For samples diluted without using (G) Buffer solution for extraction supplied in this kit, effects of interfering substances should be checked.
- If presence of interfering substance is suspected, linearity of dilutions should be confirmed using at least 2 different dilution levels.

[Storage of samples] [Storage and stability]

- Use extracted samples for measurement as soon as possible. Freeze the extracted samples at -20°C or lower for long-term storage. Repeated freezing and thawing should be avoided. Extracted samples should be immediately assayed or stored below -20°C for several days. If you have to store assay samples for a longer period, snap-freeze samples and keep them below -20°C . Defrosted samples should be mixed thoroughly for best results. Avoid repeated freezing and thawing.

8. Preparation of Reagents

- Bring all reagents of the kit to room temperature (20°C - 25°C) before use.
- Prepare reagent solutions in appropriate volume for your assay. Do not store the diluted reagents.

[Concentrated reagents]

[(B) Der f 2 Standard (500 ng/mL)]

Make a serial dilution of original standard solution to prepare each standard solution. An example is shown below.

Volume of standard solution	(C) Buffer	Concentration (ng/mL)
Original solution : 50 μ L	450 μ L	50
50 ng/mL solution : 250 μ L	250 μ L	25
25 ng/mL solution : 250 μ L	250 μ L	12.5
12.5 ng/mL solution : 250 μ L	250 μ L	6.25
6.25 ng/mL solution : 250 μ L	250 μ L	3.13
3.13 ng/mL solution : 250 μ L	250 μ L	1.56
1.56 ng/mL solution : 250 μ L	250 μ L	0.78
0 (Blank)	250 μ L	0

[(D) Peroxidase conjugated Antibody Solution]

Prepare working solution by dilution of (D) with the (C) Buffer to **1 : 100**.

[(I) Wash Solution (10×)]

Dilute 1 volume of the concentrated (I) Wash Solution (10×) to **10 volume** with deionized water (or distilled water) to prepare working solution. Example : 100 mL of concentrated (I) Wash Solution (10×) and 900 mL of deionized water (or distilled water).

9. Assay procedure


Remove the cover sheet of the Antibody-coated Plate after bringing up to room temperature.

- (1) Wash the Antibody-coated Plate (A) by filling the well with washing buffer and discard 3 times (*①), then strike the plate upside-down onto several sheets of paper towel to remove residual buffer in the wells.
- (2) Pipette 100 μ L of diluted samples to the designated sample wells.
- (3) Pipette 100 μ L of standard solution to the wells designated for standards.
- (4) Shake the plate gently on a plate shaker (*②)
- (5) Stick a Plate Seal (*③) on the plate and incubate for 1 hour at 20°C - 25°C.
- (6) Discard the reaction mixture and rinse wells as step (1).
- (7) Pipette 100 μ L of Peroxidase-conjugated Antibody Solution to all wells, and shake as step (4).
- (8) Stick a plate seal (*③) on the plate and incubate the plate for 1 hour at 20°C - 25°C.
- (9) Discard the reaction mixture and rinse wells as step (1).
- (10) Pipette 100 μ L of TMB solution to wells, and shake as step (4).
- (11) Stick a Plate Seal (*③) on the plate and incubate the plate for 20 minutes at 20°C - 25°C.
- (12) Add 100 μ L of the Stop Solution to all wells and shake as step (4).
- (13) Measure the absorbance of each well at 450 nm (reference wavelength, 620 nm*) using a plate reader within 30 minutes.

*Refer to 12. Summary of Assay Procedure for notes on *①, *② and *③.

10. Calculations

- (1) Prepare a standard curve using two-way logarithmic section paper by plotting absorbance (Y-axis) against standard concentration (ng/mL) on X-axis.
- (2) Using the standard curve, read the Der f 2 concentration of a sample at its absorbance, and multiply the assay value by dilution factor if the sample has been diluted. Though the assay range is wide enough, in case the absorbance of some samples is higher than that of the highest standard, please repeat the assay after proper dilution of samples with the buffer solution. We recommend the use of 3rd order regression curve for log-log plot, or 4 parameters method for log-normal plot in computer calculation.

11. Trouble shooting

- Low absorbance in all wells

Possible explanations :

- 1) The standard or samples might not be added.
- 2) Reagents necessary for coloration such as Peroxidase-conjugated Antibody Solution, or TMB Solution might not be added.
- 3) Wrong reagents related to coloration might have been added. Wrong dilution of Biotin-conjugated Antibody Solution or Peroxidase-conjugated Antibody Solution.
- 4) Contamination of enzyme inhibitor(s).
- 5) Influence of the temperature under which the kits had been stored.
- 6) Excessive hard washing of the well plate.

7) Addition of TMB Solution soon after taking out from a refrigerator might cause poor coloration owing to low temperature.

- Blank OD is higher than that of the lowest standard concentration (0.78 ng/mL)

Possible explanations :

Improper or inadequate washing. (Change washing repetition from 3 times to 4 - 6 times after the reaction with Peroxidase-conjugated Antibody Solution.)

- High coefficient of variation (CV)

Possible explanation :

- 1) Improper or inadequate washing.
- 2) Improper mixing of standard or samples.
- 3) Pipetting at irregular intervals.

- Q-1 : Can I divide the plate to use it for the other testing?

A-1 : Yes, cut off the clear seal on the plate with cutter along strip. Put the residual plate, which is still the seal on, in a refrigerator soon.

- Q-2 : I found 96 well-plate is empty when I opened the box.

A-2 : As this kit is dried type.


Summary of Assay Procedure ☐ : Use as a check box

*First, read this instruction manual carefully and start your assay after confirmation of details.

- ☐ Bring the well-plate and all reagents back to 20°C - 25°C for 2 hours.

- ☐ Wash Solution (10×) concentrate must be diluted to 10 times by deionized water (or distilled water) that returned to 20°C - 25°C.

- ☐ Der f 2 Standard dilution example :


- ☐ Antibody-coated Plate
- ☐ ↓ Washing 3 times (*①) (*⑤)
- ☐ Diluted Samples or Standards 100 µL
- ☐ ↓ Shaking (*②), Incubation for 1 hour at 20°C - 25°C. (Standing (*③))
- ☐ Peroxidase-conjugated Antibody Solution (D) Dilute to 1 : 100 by using Buffer (C) and use. Dilute reagents during the first reaction.
- ☐ ↓ Washing 3 times (*②)
- ☐ Peroxidase-conjugated Antibody Solution 100 µL
- ☐ ↓ Shaking (*②), Incubation for 1 hour at 20°C - 25°C. (Standing (*③))
- ☐ ↓ Washing 3 times (*①) (*⑤)
- ☐ TMB Solution 100 µL
(After dispense, the color turns to blue depending on the concentration.)
- ☐ ↓ Shaking (*②), Incubation for 20 minutes at 20°C - 25°C. (Standing (*③))
- ☐ Stop Solution 100 µL
(After dispense, the color turns to yellow depending on the concentration.)
- ☐ ↓ Shaking (*②) (Immediately shake.)
- ☐ Measurement of absorbance (450 nm, Ref 620 nm (*④))
(Ref. wave cancels the dirt in the back of plate.)

- *① After dispensing wash buffer to wells, lightly shake the plate on your palm for 10 seconds and remove the buffer. Guideline of washing volume : 300 μ L/well for an automatic washer and for a pipette if the washing buffer is added by pipette. In case of washing by using 8 channel pipette, sometimes the back ground tends to be high. If so, change washing frequency from 3 times to 4 - 6 times at the constant stroke after the reaction with Peroxidase-conjugated Streptavidin Solution.
- Standard of plate-washing pressure: 5 mL/min - 25 mL/min. (Adjust it depending on the nozzle's diameter.)
- *② Guideline of shaking : 600 rpm - 1200 rpm for 3 times for 10 seconds.
- *③ Seal the plate during the reaction after shaking. Peel off the protective paper from the seal and stick the seal on the plate. Do not reuse the plate seal used once.
- *④ 600 nm - 650 nm can be used as reference wavelength.
- *⑤ After removal of wash buffer, immediately dispense the next reagent.

Worksheet example

	Strip 1&2	Strip 3&4	Strip 5&6	Strip 7&8	Strip 9&10	Strip 11&12
A	50 ng/mL	Sample 1	Sample 9	Sample 17	Sample 25	Sample 33
B	25 ng/mL	Sample 2	Sample 10	Sample 18	Sample 26	Sample 34
C	12.5 ng/mL	Sample 3	Sample 11	Sample 19	Sample 27	Sample 35
D	6.25 ng/mL	Sample 4	Sample 12	Sample 20	Sample 28	Sample 36
E	3.13 ng/mL	Sample 5	Sample 13	Sample 21	Sample 29	Sample 37
F	1.56 ng/mL	Sample 6	Sample 14	Sample 22	Sample 30	Sample 38
G	0.78 ng/mL	Sample 7	Sample 15	Sample 23	Sample 31	Sample 39
H	0 (Blank)	Sample 8	Sample 16	Sample 24	Sample 32	Sample 40

Assay worksheet

LBIS™ Der f 2 ELISA Kit

- [Storage] Store the kit at 2°C - 10°C (Do not freeze).
- [Expiration data] Indicated on the container.
- [Package] For 96 tests
- [Cat #] 290-86901

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