LBIS™ Mouse MCP-1 (CCL2) ELISA Kit

Please, read this instruction carefully before use.

1. Intended use

Mouse Monocyte chemoattractant protein-1 (MCP-1) consisting of 125 amino acid residues, is a member of the C-C chemokine family with chemoattractant properties for monocytes/macrophages. Like many chemokines, MCP-1 expression is induced by inflammatory stimuli such as TNF- α and IL-1. It is secreted as an inflammatory mediator by immunocompetent cells such as monocytes and macrophages, in addition to many cells such as vascular endothelial cells, fibroblasts, uriniferous tubule epithelial cells, and smooth muscle cells.

MCP-1 and its cognate receptor, chemokine motif receptor 2 (CCR2), specifically recruit monocyte/ macrophages into an inflammatory site and subsequently prolong inflammatory response, and MCP-1/CCR2 systems are thought to be involved in the pathogenesis of chronic inflammatory diseases. Furthermore, in obesity-associated adipose tissues, enhanced MCP-1 expression promotes the infiltration of monocytes/macrophages into adipose tissue, which is suggested to be involved in the pathophysiology of insulin resistance. From these findings, MCP-1 has gained increasing attention as a therapeutic target molecule for chronic inflammatory diseases such as metabolic syndrome.

This kit is a sandwich enzyme system for quantitative measurement of mouse MCP-1 (CCL2). This kit can measure MCP-1 (CCL2) in mouse serum (plasma) and urine, specifically and sensitively. This product is intended for research use only. Not for use in diagnostics procedures.

2. Assay Principle

In FUJIFILM Wako Pure Chemical Corporation's LBIS™ Mouse MCP-1 (CCL2) ELISA Kit, standards or samples are incubated in polyclonal antibody-coated wells to capture MCP-1 (CCL2). After 2 hours incubation and washing, biotinylated anti-MCP-1 (CCL2) antibody is added and incubated further for 1 hour to bind with captured MCP-1 (CCL2). After washing, HRP (horse radish peroxidase)-conjugated streptavidin is added and incubated for 30 minutes. After washing, bound HRP-conjugated streptavidin is reacted with a TMB Solution for 20 minutes, and reaction is stopped by addition of acidic solution, and absorbance of yellow product is measured spectrophotometrically at 450 nm. The absorbance is nearly proportional to MCP-1 (CCL2) concentration. The standard curve is prepared by plotting absorbance against standard MCP-1 (CCL2) concentrations. MCP-1 (CCL2) concentrations in unknown samples are determined using this standard curve.

3. Performance Characteristics

- Precision of assay (Within assay variation)
2 serum samples, 5 replicates assay

n/ID	Sample1 (pg/mL)	Sample2 (pg/mL)
1	170	18.1
2	171	18.9
3	172	18.7
4	173	18.4
5	169	18.8
mean	171	18.6
SD	1.8	0.32
CV (%)	1.1	1.7

- Reproducibility (Between assay variation)
3 serum samples, 4 days, 3 replicates assay

Day/ID	Sample1 (pg/mL)	Sample2 (pg/mL)	Sample3 (pg/mL)
0	210	50.8	8.67
1	211	48.0	8.33
2	215	48.2	8.25
3	226	51.9	8.84
mean	215	49.7	8.52
SD	7.2	1.9	0.28
CV (%)	3.3	3.8	3.3

- Recovery test (Standard mouse MCP-1 (CCL2) was added in 4 concentrations to serum or plasma sample and assayed.) The recoveries were 95.0% - 104%.



Serum sample

Added (pg/mL)	Found (pg/mL)	Recovered (pg/mL)	Recovery (%)
0	25.7	-	-
25	50.7	25.0	100
50	75.5	49.8	99.6
100	124	98.3	98.3
200	223	197	98.5

Plasma sample (EDTA)

Added (pg/mL)	Found (pg/mL)	Recovered (pg/mL)	Recovery (%)
0	14.7	-	-
25	40.5	25.8	103
50	66.5	51.8	104
100	115	100	100
200	205	190	95.0

- Dilution test (Serum sample and plasma sample were serially diluted by 3 steps.) The dilution curves showed linearity with $R^2=0.9992$ - 0.9998.

4. Technical Tips / Precautions

- ELISA can be affected by assay environment. Ensure that room temperature at places for assay operation and incubation is strictly controlled at 20°C - 25°C. Avoid performing assay under the airstream velocity over 0.4 m/sec. and the humidity less than 30%, seal the well plate with a plate seal and place the well plate in an incubator or a styrofoam box in each step of incubation.
- To avoid denaturation of the coated antibody, do not let the plate go dry.
- In order to avoid dryness of wells, contamination of foreign substances and evaporation of dispensed reagents, never forget to cover the well plate with a plate seal supplied, during incubation.
- Be careful to avoid any contamination of assay samples and reagents. We recommend the use of disposal pipette tips, and 1 tip for 1 well.
- As the Antibody-coated Plate is module type of 8 wells $\times$ 12 strips, each strip can be separated by cutting the cover sheet with a knife and used independently.
- The TMB Solution should be almost clear pale yellow before use. It turns blue during reaction and gives yellowish color after addition of Stop Solution. Greenish color means incomplete mixing.
- Time the reaction from the pipetting of the reagent to the first well.
- Prepare a standard curve for each assay.
- For professional use only, beginners are advised to use this kit under the guidance of experienced person. In manual operation, proficiency in pipetting technique is recommended.
- Wear protective gloves, clothing, eye, and face protection.
- Avoid contact of skin and mucous membranes with kit reagents or specimens. If any reagents come in contact with eyes, skin, or mucous membranes, wash with copious amounts of water and contact a physician.
- Do not drink, eat, or smoke in a place where this kit is used.
- The materials must not be pipetted by mouth.
- Do not use reagents with different lot numbers together.
- The reagents are prepared to give accurate results only when used in combination within the same box. Therefore, do not combine the reagents from kits with different lot numbers. Even if the lot number is the same, it is best not to mix the reagents with those that have been preserved for some period.
- Handle the sample with proper care, being aware that the sample may have an infection risk. This kit contains animal-derived ingredients.

- Residual samples and used tips should be sterilized before disposal.
- Dispose consumable materials and unused contents in accordance with applicable regional/national regulatory requirements.

5. Reagents supplied

Components	Use Status	Amount
(A) Antibody-coated Plate	Use after washing.	96 wells/1 plate
(B) Mouse MCP-1 (CCL2) Standard	Freeze-dried. Use after reconstitution.	1 bottle
(C) Buffer	Ready for use.	60 mL/1 bottle
(D) Biotin-conjugated Antibody	Freeze-dried. Use after reconstitution.	1 bottle
(E) Peroxidase-conjugated Streptavidin Solution	Concentrated. Use after dilution.	100 μ L/1 bottle
(F) TMB Solution	Ready for use.	12 mL/1 bottle
(H) Stop Solution	Ready for use.	12 mL/1 bottle
(I) Wash Solution (10 $\times$)	Concentrated. Use after dilution.	100 mL/1 bottle
(J) Plate Seal	–	4 sheets

[Storage and Stability]

[(A) Antibody-coated Plate]

If seal is not removed, put the strip back in a plastic bag with zip lock originally used for well-plate container and store at 2°C - 10°C.

[(B) Mouse MCP-1 (CCL2) Standard]

The reconstituted standard solution (original standard solution) should be used within 2 weeks if stored in a refrigerator. Dispose the remaining diluted standard solutions after use.

[(C) Buffer] & [(F) TMB Solution]

Use only volume you need for your assay. Remaining reagents should be stored at 2°C - 10°C closing the cap tightly. Once opened, we recommend using as soon as possible to avoid influence by environmental condition.

[(D) Biotin-conjugated Antibody]

The reconstituted (D) solution should be used within 2 weeks if stored in a refrigerator. Dispose the remaining diluted working solution.

[(E) Peroxidase-conjugated Streptavidin Solution]

Remaining working solution (already diluted) should be disposed. The rest of the undiluted solution (unused) : close the cap tightly and store at 2°C - 10°C. Once opened, we recommend using as soon as possible to avoid influence by environmental condition.

[(H) Stop Solution]

Close the cap tightly and store at 2°C - 10°C.

[(I) Wash Solution (10 $\times$)]

The rest of undiluted solution (unused) : close the cap tightly and store at 2°C - 10°C. Dispose any remaining diluted buffer.

6. Equipments or supplies required but not supplied ☐ Use as a check box

☐ Deionized water (or Distilled water) ☐ Test tubes for preparation of standard solution series. ☐ Glassware for dilution of Wash Solution (10 $\times$) (a graduated cylinder, a bottle) ☐ Pipettes (disposable tip type). One should be able to deliver 10 μ L - 100 μ L precisely, and another for 100 μ L - 1000 μ L. ☐ Syringe-type repeating dispenser like Eppendorf multipette plus which can dispense 100 μ L. ☐ Paper towel to remove washing buffer remaining in wells. ☐ A vortex-type mixer. ☐ A shaker for 96 well-plate (500 rpm - 1200 rpm) ☐ An automatic washer for 96 well-plate (if available), or a wash bottle with a jet nozzle. ☐ A 96 well-plate reader (450 nm $\pm$ 10 nm, 620 nm : 600 nm - 650 nm) ☐ Software for data analysis.

7. Preparation of Samples

- This kit is intended to measure MCP-1 (CCL2) in mouse serum or plasma or urine.
- We recommend using EDTA-2Na, EDTA-2K or heparin as anticoagulant.
- Samples should be immediately assayed or stored below -35°C until assay. Before starting assay, shake thawed samples sufficiently. Do not repeat freeze-and-thaw cycles.
- Centrifuge the sample to remove turbidity or insoluble matters where necessary before use for the assay.
- If presence of interfering substance is suspected, examine by dilution test at more than 2 points.
- Dilution of a sample should be made in a PP or PE test tube using (C) Buffer prior to adding them to wells.

8. Preparation of Reagents

- ◆ Bring all reagents of the kit to room temperature (20°C - 25°C) before use.
- ◆ Prepare reagent solutions in appropriate volume for your assay. Do not store the diluted reagents.

[Concentrated reagents]

[(B) Mouse MCP-1 (CCL2) Standard]

Reconstitute the (B) Mouse MCP-1 (CCL2) Standard with purified water in the amount specified in "Reconstitution of standard" *to prepare a original standard solution (5000 pg/mL).


Then prepare using the kit-attached (C) Buffer that has been allowed to warm up to room temperature.

*Please confirm "Reconstitution of standard" from this product page. Since the amount of purified water added varies by lot, please be sure to check each lot.

Below is an example of preparing each standard solution.

Volume of standard solution	(C) Buffer	Concentration
Original standard solution 50 μ L	450 μ L	500 pg/mL
500 pg/mL solution 200 μ L	250 μ L	222 pg/mL
222 pg/mL solution 200 μ L	250 μ L	98.8 pg/mL
98.8 pg/mL solution 200 μ L	250 μ L	43.9 pg/mL
43.9 pg/mL solution 200 μ L	250 μ L	19.5 pg/mL
19.5 pg/mL solution 200 μ L	250 μ L	8.67 pg/mL
8.67 pg/mL solution 200 μ L	250 μ L	3.85 pg/mL
Blank	350 μ L	0

Standard Curve

[(D) Biotin-conjugated Antibody]

Add 100 μ L of deionized water (or distilled water) of room temperature to (D) Biotin-conjugated Antibody (Freeze-dried) (Original (D) solution). Prepare working solution by dilution of original (D) solution with the (C) Buffer to **1 : 100**. 10 mL of the diluted solution is enough for 96 wells.

[(E) Peroxidase-conjugated Streptavidin Solution]

Prepare working solution by dilution of (E) with the (C) Buffer to **1 : 100**. 10 mL of the diluted solution is enough for 96 wells.

[(I) Wash Solution (10 $\times$)]

Dilute 1 volume of the concentrated (I) Wash Solution (10 $\times$) to **10 volumes** with deionized water (or distilled water) to prepare working solution. Example : 100 mL of concentrated (I) Wash Solution (10 $\times$) and 900 mL of deionized water (or distilled water).


9. Assay Procedure

Remove the cover sheet of the Antibody-coated Plate after bringing up to room temperature.

- (1) Wash the (A) Antibody-coated Plate by filling the wells with washing buffer and discard 4 times (*①), then strike the plate upside-down onto folded several sheets of paper towel to remove residual buffer in the wells (*⑤).
 - (2) Pipette **50 μ L** of standard solution to the wells designated for standards.
 - (3) Pipette **40 μ L** of (C) Buffer to the sample wells, and then add **10 μ L** of sample to the wells.
 - (4) Shake the plate gently on a plate shaker (*②).
 - (5) Stick a Plate Seal (*③) on the plate and incubate for 2 hours at 20°C - 25°C.
 - (6) Discard the reaction mixture and rinse wells as step (1).
 - (7) Pipette **50 μ L** of Biotin-conjugated Antibody solution to all wells and shake as step (4).
 - (8) Stick a Plate Seal (*③) on the plate and incubate the plate for 1 hour at 20°C - 25°C.
 - (9) Discard the reaction mixture and rinse wells as step (1).
 - (10) Pipette **50 μ L** of Peroxidase-conjugated Streptavidin Solution to all wells and shake as step (4).
 - (11) Stick a Plate Seal (*③) on the plate and incubate the plate for 30 minutes at 20°C - 25°C.
 - (12) Discard the reaction mixture and rinse wells as step (1).
 - (13) Pipette **100 μ L** of (F) TMB Solution to wells and shake as step (4).
 - (14) Stick a Plate Seal (*③) on the plate and incubate the plate for 20 minutes at 20°C - 25°C.
 - (15) Add **50 μ L** of the (H) Stop Solution to all wells and shake as step (4).
 - (16) Measure the absorbance of each well at 450 nm (reference wavelength, 620 nm (*④)) immediately using a plate reader.
- *Refer to 12. Summary of Assay Procedure for notes on *①, *②, *③, *④, and *⑤.

10. Calculations

- (1) Create a standard curve by reducing the data using computer software capable of generating 3rd order regression curve, or 4 or 5 parameters. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the mouse MCP-1 (CCL2) concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.
- (2) Using the standard curve, read the MCP-1 (CCL2) concentration of samples at its absorbance, and multiply the assay value by dilution factor if the sample has been diluted. The standard operation method is 5-fold dilution.

11. Trouble Shooting

- Low absorbance in all wells

Possible explanations :

- 1) The standard or samples might not be added.
- 2) Reagents necessary for coloration such as Biotin-conjugated Antibody, Peroxidase-conjugated Streptavidin Solution, or TMB Solution might not be added.
- 3) Wrong reagents related to coloration might have been added. Wrong dilution of Biotin-conjugated Antibody or Peroxidase-conjugated Streptavidin Solution.
- 4) Contamination of enzyme inhibitor(s).
- 5) Influence of the temperature under which the kits had been stored.
- 6) Excessive hard washing of the well plate.
- 7) Addition of chromogenic substrate reagent soon after taking out from a refrigerator might cause poor coloration owing to low temperature.

- Blank OD was higher than that of the lowest standard concentration (3.85 pg/mL).

Possible explanations :

Improper or inadequate washing. (Change washing frequency from 4 times to 5 - 8 times at the constant stroke after the reaction with Peroxidase-conjugated Streptavidin Solution.)

- High coefficient of variation (CV)

Possible explanation :

- 1) Improper or inadequate washing.
- 2) Improper mixing of standard or samples.
- 3) Pipetting at irregular intervals.

- Q-1 : Can I divide the plate to use it for the other testing?

A-1 : Yes, cut off the clear seal on the plate with cutter along strip. Put the residual plate, which is still the seal on, in a refrigerator soon.

- Q-2 : I found 96 well-plate is empty when I opened the box.

A-2 : As this kit is dried type, not preservation stabilizer is added.

12. Summary of Assay Procedure ☐ : Use as a check box

*First, read this instruction manual carefully and start your assay after confirmation of details.

☐ Bring the well-plate and all reagents back to **20°C - 25°C for 2 hours**.

☐ (I) Wash Solution (10×) must be diluted to **10 times** by deionized water (or distilled water) that returned to 20°C - 25°C.

☐ Standard solution dilution example

Reconstitute the (B) Mouse MCP-1 (CCL2) Standard with purified water in the amount specified in "Reconstitution of standard" * to prepare a original standard solution (5000 pg/mL).

Then prepare using the kit-attached (C) Buffer that has been allowed to warm up to room temperature.

* Please confirm "Reconstitution of standard" from this product page. Since the amount of purified water added varies by lot, please be sure to check each lot.

Preparation of diluted MCP-1 (CCL2) standard solutions :

Conc. (pg/mL)	500	222	98.8	43.9	19.5	8.67	3.85	0
Std. sol. (μL)	50	200*	200*	200*	200*	200*	200	—
Buffer (μL)	450	250	250	250	250	250	250	350

*One rank higher standard.

☐	Antibody-coated Plate	
☐	↓ Washing 4 times (*①), (*⑤)	
☐	Samples (diluted sample : Buffer 40 μL + Sample 10 μL), or Standards	50 μL
☐	↓ Shaking (*②), Incubation for 2 hours at 20°C - 25°C. (Standing (*③))	
☐	Add deionized water (or distilled water) of room temperature to (D) Biotin-conjugated Antibody (Freeze-dried) (Original (D) solution). Prepare working solution by dilution of original (D) solution to 1 : 100 with the (C) Buffer returned to 20°C - 25°C. (This should be prepared during incubation.)	
☐	↓ Washing 4 times (*①), (*⑤)	
☐	Biotin-conjugated Antibody solution (Diluted)	50 μL
☐	↓ Shaking (*②), Incubation for 1 hour at 20°C - 25°C. (Standing (*③))	
☐	Dilute Peroxidase-conjugated Streptavidin Solution to 1 : 100 with the (C) Buffer returned to 20°C - 25°C. (This should be prepared during incubation.)	
☐	↓ Washing 4 times (*①), (*⑤)	
☐	Peroxidase-conjugated Streptavidin Solution (Diluted)	50 μL
☐	↓ Shaking (*②), Incubation for 30 minutes at 20°C - 25°C. (Standing (*③))	
☐	↓ Washing 4 times (*①), (*⑤)	
☐	TMB Solution	100 μL
☐	After dispensing, the color turns to blue depending on the concentration.	
☐	↓ Shaking (*②), Incubation for 20 minutes at 20°C - 25°C. (Standing (*③))	
☐	Stop Solution	50 μL
☐	After dispensing, the color turns to yellow depending on the concentration.	
☐	↓ Shaking (*②), Immediately shake.	
☐	Measurement of absorbance (450 nm, Ref 620 nm (*④)) immediately. Ref. wave cancels the dirt in the back of plate.	

*① After dispensing wash buffer to wells, lightly shake the plate on your palm for 10 sec. and remove the buffer. Guideline of washing volume : 300 μL/well for an automatic washer and for a pipette if the washing buffer is added by pipette. In case of washing by using 8 channel pipettes, sometimes the background tends to be high. If so, change washing frequency from 4 times to 5 - 8 times at the constant stroke after the reaction with Peroxidase-conjugated Streptavidin Solution.

Standard of plate-washing pressure : 5 mL/min. - 25 mL/min. (Adjust it depending on the nozzle's diameter.)


*② Guideline of shaking : 600 rpm - 1200 rpm 3 times for 10 seconds.

*③ Seal the plate during the reaction after shaking. Peel off the protective paper from the seal and stick the seal on the plate. Do not reuse the plate seal used once.

*④ 600 nm - 650 nm can be used as reference wavelength.

*⑤ After removal of wash buffer, immediately dispense the next reagent.

Assay Worksheet

13. Storage and Expiration

The complete kit is stored at 2°C - 10°C (Do not freeze). Opened reagents should be used as soon as possible to avoid less than optimal assay performance caused by storage environment.

LBIS™ Mouse MCP-1 (CCL2) ELISA Kit

[Storage] Store the kit at 2°C - 10°C (Do not freeze)

[Expiration date] Indicated on the container

[Package] For 96 tests

[Cat #] 290-88101

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://flwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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