## DNA Extractor® WB Kit

### (Genomic DNA Extraction/NaI Method)

Obtaining super macromolecular genomic DNA from whole blood and culture cell is a primary requisite for genetic analyses of genome and molecular biology procedures. A variety of DNA purification procedures involve either hazardous processes such as the phenol/chloroform extraction or time-consuming and inefficient procedures such as lengthy dialysis and the transfer of DNA from one tube to another. The "DNA Extractor® WB Kit" employs a new extraction procedure for DNA purification after the lysis of cells. This procedure, using Sodium Iodide (NaI) as a chaotropic agent, realizes intact DNA isolation of both purity and high yield. Furthermore, extraction is done through several brief microfuge centrifugations from a single microfuge tube. This kit is for laboratory use only.

#### 【Reagents (50 Assays)】

One kit includes ;

(1)	296-50513	Lysis Solution	65 mL	2 bottles
(2)	297-50521	Enzyme Reaction Solution	10 mL	1 bottle
(3)	294-50531	Sodium Iodide Solution	15 mL	1 bottle
(4)	291-50541	Washing Solution (A)	50 mL	1 bottle
(5)	298-50551	Washing Solution (B)	50 mL	1 bottle
(6)	295-50561	Protease	10 mg	1 bottle

#### 【Materials and Apparatus】

Reagent	1) Isopropyl alcohol	25 mL
	2) Distilled deionized water	0.6 mL
Apparatus	1) Microcentrifuge (Max. 12,000 g)	
	2) Microcentrifuge tube with screwed cap (1.5 mL)	
	3) Microfuge tube mixer	

**Note :** During the storage, in such a case, Enzyme Reaction Solution and the NaI Solution may crystalize out, solubilize the precipitate completely by warming the solution at 50°C.

#### 【Principle of DNA Extraction】

The procedure includes three major steps : 1) Isolation of nucleus fraction. 2) Liberation of DNA. 3) Purification of DNA.

In the first step, cells are lysed by a non-ionic surfactant, polyoxyethylene (10) oxyphenyl ether, followed by a brief centrifugation. By repeating this operation, the majority of cellular macromolecules, including RNAs, are removed from the nucleus fraction. In the next step, the nucleus fraction is treated with protease in the presence of 1% sodium dodecyl sulfate for 1 hour. By this treatment DNA packed in the nucleus is expanded and liberated from nucleus protection. Simultaneously, protein in the solution is digested into polypeptide. In the final step, the resultant solution is subjected to NaI extraction of DNA as

described in literature. Polypeptide and other biological molecules remain soluble in a high concentration of NaI, a chaotropic salt, even after isopropanol is added to the solution to precipitate the nucleic acid.

#### 【Storage】

Store at 2 ~ 10°C.

#### 【Protocol】

##### 〈Preparation of Reagents〉

10 mg of Protease is reconstituted by adding 0.6 mL of distilled water. Protease solution reconstituted can be stored at -20°C.

##### 〈Scheme I : DNA Extraction Procedure from Whole Blood〉

- 0.5 mL of whole blood, prepared with EDTA-2Na as anti-coagulant, is dispensed into a 1.5 mL microcentrifuge tube with screwed cap and then placed standing on ice.
- After adding 0.5 mL of Lysis Solution and capping the tube, mix gently several times by inversion. (Follow this for the mixing operations below.)
- After a centrifugation at 10,000 × g, for 20 seconds at 4°C, carefully remove the supernatant, keeping the dark red precipitate intact.
- Add 1 mL of Lysis Solution to the tube and set the capped tube on a microfuge tube mixer for 30 seconds at a moderate speed.
- After a centrifugation of the mixture at 10,000 × g for 20 seconds at 4°C, remove the supernatant from the tube.
- Repeat operation step 4-5.
- Suspend the resultant pellet in 200 μL of Enzyme Reaction Solution. Add 10 μL of Protease Solution to the suspension and mix gently by inversion.
- Incubated at 37°C for 1 hour. (During the incubation, mix the solution several times by inversion.)
- After the 1 hour-incubation, add 0.3 mL of NaI Solution to the mixture and mix well by inversion.
- Add 0.5 mL of isopropyl alcohol to the mixture and mix well until DNA, a whitish material, appears.
- After centrifugation at 10,000 × g for 10 minutes at room temperature, gently removed the resultant supernatant. Put the tube upside down to remove the solution remaining on the surface of the tube.
- After adding 1 mL of Washing Solution (A), mix thoroughly so that the pellet is removed from the tube wall.
- After a centrifugation of the mixture at 10,000 × g for 5 minutes at room temperature, remove the supernatant from the tube.
- After adding 1 mL of Washing Solution (B), mix thoroughly so that the pellet is removed from the tube wall.
- After a centrifugation of the mixture at 10,000 × g for 5 minutes at room temperature, remove the supernatant from the tube.
- The resultant pellet is vacuum-dried for about 3 minutes. (Drying should be made within 3 minutes because DNA completely dried may be hard to be dissolved).

##### 〈Scheme II : DNA Extraction Procedure from Culture Cell (Example : HeLa, SQ-5)〉

- Culture cell pellet ( $10^3 \sim 10^6$  cells) such as HeLa and SQ-5, collected in a 1.5 mL microfuge tube is suspended in 1 mL Lysis Solution and then set the 1.5 mL microcentrifuge tube with screwed cap in a

microfuge tube mixer.

2. After a centrifugation of the mixture at  $3,000\times g$  for 5 minutes at  $4^{\circ}\text{C}$ , remove the supernatant from the tube. Follow with operations through step 4-16 as described in "Scheme I DNA Extraction Procedure from Whole Blood".

#### 【Application of "DNA Extractor WB Kit" for volume reduction procedure】

Genomic DNA was extracted from  $100\mu\text{L}$  of whole blood with the "DNA Extractor WB Kit" using smaller volumes of reagents than the standard procedure. In these procedures, 1/5 or 1/4 of the original reagent volumes were employed. Microfuge tubes were  $0.6\text{ mL}$ . As a reference, genomic DNA was extracted from  $500\mu\text{L}$  of whole blood using the standard procedure.

#### Recovery and purity of genomic DNA

Sample Volume	Reagents	Recovery ( $\mu\text{g DNA/mL}$ whole blood)	260/280
$100\mu\text{L}$ wholeblood	1/5	23.8	1.97
$100\mu\text{L}$ wholeblood	1/4	22.2	1.98
$500\mu\text{L}$ wholeblood	1/1	22.0	1.92

#### Note

1. DNA extracted from this kit is suitable for several applications, including restriction enzyme digestion and Southern blot analysis.
2. Whole blood anti-coagulated by EDTA is appropriate as a sample for this kit, but DNA isolated from heparinized whole blood with greater than 20 units/mL may influence the result of some experiments.
3. This kit is applicable to  $100\mu\text{L} \sim 500\mu\text{L}$  of whole blood by proportionally varying the volumes of solutions used in the standard method described below. (1 mL of whole blood can be applied by using a 2 mL microcentrifuge tube and doubling the amount of each reagent.)
4. Whole blood stored refrigerated for 3 weeks did not show a significant decrease in the recovery of DNA with this kit.
5. This kit is not compatible for DNA extraction from smear preparation samples.

## FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://fjwkw.fujifilm.co.jp>

#### FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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