Code No. 291-55203 (for 5,000 cm²)

ImmunoStar™ Reagents (for chemiluminescence detection)

ImmunoStar reagents from Wako are designed for a simple and highly sensitive immunoblotting utilizing detection by enhanced chemiluminescence. Detection levels comparable to those reached with radioactive labels are achieved by use of a unique enhancer. Sensitivity is further enhanced by use of ABC Solution (Streptavidin and Biotin-conjugated peroxidase complex).

【Features】

1. The immunoblotting assay using these chemiluminescent reagents can detect sub-pg amounts of the target protein within 1-3 minute(s) of exposure, which is 10 to 100 times more sensitive than those utilizing chromogenic detection.
2. The chemiluminescence emission continues for several hours, enabling the detection of additional small amounts of the target protein by prolonging the exposure.
3. The membrane is reusable after washing.

【Contents】

| Reagent | 5,000 cm ² |
|-------------------------|-----------------------|
| Luminescence Solution A | 330 mL |
| Luminescence Solution B | 330 mL |
| Luminescence Solution C | 120 mL |

【Materials and apparatuses to be prepared for immunoblotting】

- 1) Primary antibody
- 2) Blocking solution
Usually, 5% skim milk solution is used for blocking of the blotting membrane. For detection of phosphotyrosine, BSA or gelatin is recommended as the blocker.
- 3) Secondary antibody (Biotin-conjugated secondary antibody or peroxidase-conjugated secondary antibody)
- 4) ABC solution
- 5) Wash solution [0.05 w/v% Tween 20 and 50 mmol/L TBS, pH 7.2.]
- 6) Diluent [50 mmol/L TBS, pH 7.2] (as necessary)
- 7) Membrane
Nitrocellulose membrane or PVDF.
- 8) Detection apparatuses
X-ray film or luminescence image analyzer

Please prepare commercial or homemade products as needed. The dilution rate of the antibody depends on the titer of the antibody used. In addition, the performance of ABC solutions varies depending on the manufacturer and preparation method. It is recommended to consider appropriate conditions of use by performing dot-blot analysis, etc.

【Procedure】

1. Preparation of working solution

Just before use, chemiluminescent solutions are prepared by mixing chemiluminescent reagents A, B, and C at a ratio of 6 mL, 6 mL, and 2 mL per 100 cm² of membrane, respectively.

2. Chemiluminescent staining

- 1) Membrane infiltration according to the instructions for use of the membrane.
- 2) Transfer from gel to membrane by transfer equipment (semi-dry or tank type).
- 3) Shake the membrane in 50 mL of Wash solution for 5 minutes.
- 4) Incubate the membrane in 12 mL of blocking solution at 37°C for 1 hour.
- 5) Shake the membrane in 50 mL of Wash solution for 5 minutes. Change Wash solution and repeat the procedure.
- 6) Incubate the membrane in 10 mL of primary antibody solution at 37°C for 1 hour or under pre-determined antibody conditions.
- 7) Shake the membrane in 50 mL of Wash solution for 5 minutes. Change Wash solution and repeat the procedure twice.
- 8) Secondary antibody reaction:
 - 8-1) Incubate the membrane in 10 mL of biotin-labeled secondary antibody at room temperature or 37°C for 20 minutes.
 - 8-2) If a peroxidase-labeled secondary antibody is used, incubate the membrane in 10 mL of antibody solution for 20 minutes, usually at room temperature or 37°C. Since the following steps 9) and 10) are not necessary, proceed to 11) for washing step.
- 9) Shake the membrane in 50 mL of Wash solution for 5 minutes. Change the Wash solution and repeat the process.
- 10) Perform ABC reaction ; incubate the membrane in 10 mL of ABC solution at room temperature or 37°C for 20 minutes.
- 11) Shake the membrane in 50 mL of Wash solution for 5 minutes. Replace Wash solution and repeat the procedure twice.
- 12) Immerse the membrane in 14 mL of the chemiluminescence solution prepared immediately before for 1 minute. Place the membrane between plastic wrap, remove air, and expose to the X-ray film in a dark room. 30 seconds to 1 hour later, develop the X-ray film.

(Precaution)

1. The exposure time is determined by the strength of the detected image and background. Usually, good results can be obtained in about 1 minute.
2. When using a device such as LAS 1000, first expose for 5 minutes. Adjust the exposure time according to the results.

3. Reuse of used membranes

Once immunologically stained, the membrane can be reused after removing the antibodies used for detection. The procedure is described below.

- 1) Place the membrane in 50 mL of the solution containing 7 mol/L guanidine hydrochloride, 50 mmol/L glycine, 0.05 mol/L EDTA-3Na, 0.1 mol/L KCl, 20 mmol/L 2-mercaptoethanol. Place the solution and membrane in a rocker for 10 minutes.
- 2) Shake the membrane in 50 mL of Wash solution for 10 minutes. Change Wash solution and repeat the procedure.
- 3) The membrane is ready for use in the immunochemiluminescence method described above.

[Notes]

1. Before starting the procedure, calculate the volume of solution required in each step according to the size of the membrane (cm²) based on the instructions provided in the protocol.
2. To save the primary antibody, place a few drops of antibody solution on parafilm, cover with a membrane, and allow the membrane to react with the antibody.
3. Avoid the use of peroxidase inhibitors such as sodium azide.

[Trouble Shooting]

| Possible causes | Treatments |
|------------------------------------|--|
| 1. No signal detected | |
| • Failure blotting | Check the transfer process |
| • Proteolytic digestion of sample | Add protease inhibitor in sample preparation |
| • Property of the primary antibody | Determination the titer or optimal conditions for the immunoreaction |
| 2. Low signal | |
| • Insufficient sample | Increase the sample quantity used |
| • Insufficient exposure time | Prolong the exposure time |
| • The sample causes as in 1 | See the treatment in 1 |

3. Signal Too high

- | | |
|---|-----------------------------|
| • Excess sample amount | Decrease the sample amount |
| • Conditions for immunoreaction of the primary antibody | Optimize the immunoreaction |

4. High background

- | | |
|---|---|
| • Conditions for immunoreaction of the primary antibody | Optimize the immunoreaction |
| • Failure of blocking | Prolong the blocking time or try other blocker |
| • Contamination of apparatuses used | Use clean apparatuses |
| • Longer exposure | Diminish the exposure time or expose the membrane after 5 min on the reaction |
| • Insufficient washing in step 11) | Repeat the wash several times |

5. Non specific signal

- | | |
|--|--|
| • The property of the primary antibody | Optimize the immunoreaction |
| • Insufficient blocking | Prolong the blocking or try another blocking reagent |

[Storage]

2-10°C, in the dark.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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