mtDNA Extractor CT Kit

Mitochondrial DNA (mtDNA) is a replicon containing in eukaryotic cytoplasm and its variants are known to arise independently of genomic DNA and participate in appearance of some diseases as diabetes, Alzheimer disease, and muscle disorders. Distribution of mtDNA variants is thought to be quite different among each tissue, and even each cells, and there are cells containing mitochondria at different portions between the normal and its variant mtDNAs. Therefore, qualified mtDNA from various tissues and cells are required to study the relationships between the diseases and mtDNA variants. However, existing methods for purification of mtDNA are consisted of very tedious manipulations, inappropriate for simultaneous isolation of mtDNAs from many samples for screening analyses, and many cases, total DNA fraction is used to analyze mtDNA, although mtDNA occupies at a quite small portion of the fraction.

“mtDNA Extractor CT Kit” is designed to isolate mtDNA from a variety of cells and tissues in high yield and high quantity with about 1 hour.

[Features]

1. mtDNA isolation with easy handling and time-saving.
2. Isolated mtDNA is pure enough to apply to the subsequent experiments as PCR and restriction enzyme digestion.
3. mtDNA isolated from as much as 5 mg of tissue is enough to amplify given DNA fragment by PCR. mtDNA isolated from 50 mg of tissues besides muscle tissues (250 mg) can be detected by EtBr-staining after agarose gel electrophoresis.
4. The method is applicable to frozen tissues as well as fresh ones.
5. The preparation is performed without hazardous reagents such as chloroform and aqueous phenol.

[Contents]

1. Buffer for Homogenate	25 mL	1 bottle
2. DNA Extraction Solution I	1.3 mL	1 vial
3. DNA Extraction Solution II (A)	1.3 mL	1 vial
4. DNA Extraction Solution II (B)	1.3 mL	1 vial
5. DNA Extraction Solution III	1.9 mL	1 vial
6. Sodium Iodide Solution	7.5 mL	1 bottle
7. Washing Solution	50 mL	1 bottle

[Materials and apparatuses to be prepared]

Reagents

1. Isopropanol, analytical grade
2. TE buffer at pH 8.
3. PBS (–), for cultured cells

Apparatuses and Lobo wares

1. Motor-driven Glass homogenizer or the equivalent
2. Automixer (MAZELA S type Z-1100, Tokyorikakikai IK. K., Tokyo)
3. 1.5mL-microtubes (sterilized)
4. Vortex mixer
5. Microfuge (max. 12,000 × g)
6. A centrifuge under reduced pressure
7. Cell scribe for cultured cells

[Procedures]

A. Isolation of mtDNA from mammalian tissues

1. Homogenize 50 mg of minced tissue with scissors in 1 mL of the homogenizing buffer with a glass homogenizer by up-and down of the pestle at 1,000 r.p.m. several times on ice until the pestle reaches the bottom of the glass well.

〈Notes〉

- 1) The sample tissues should be as fresh as possible. Qualified mtDNA could be prepared from 50 mg of frozen mammalian tissues at – 80 degrees C. for 1 month, and that from those frozen at – 20 degrees C. for 1 week.
- 2) mtDNA obtained from 50 mg of tissues besides muscles can be detected as a doublet bands on EtBr-staining an agarose gel electrophoresis, while mtDNA isolated from 250 mg of skeletal muscle can be detected on the same agarose gel.
- 3) Homogenization of rigid tissues such as skeletal and heart muscles should be done using glass homogenizer but not with glass Teflon homogenizer due to inefficient fragmentation of the tissues.
2. After transfer of the homogenate to another tube, centrifuge it at 1,000 × g at 4 degrees C. for 1 minute to precipitate large materials like nuclei, and then transfer the supernatant to another microtube by papering.
3. Centrifuge the microtube at 10,000 × g at 4 degrees C. for 10 minute to precipitate the membrane fraction from the supernatant, and then remove the resultant supernatant from the tube.
4. Add 50 micro L of DNA Extractor Solution I to the tube to mix well by papering or vortexing.
5. Add 100 micro L of pre-mixed solution of DNA extractor solution II (A) and (B) at equal volumes to the tube, and vortex it a couple of times, leaving it on ice for 5 minutes.

- 〈Note〉 Pre-mixed solution of DNA extractor Solution II (A) and (B) can be used at most for 3 weeks on storage at 2~10 degrees C after mixing. Although DNA Extractor Solution II (B) may form precipitate during the storage. It can be used after the precipitate is dissolved by warming followed by cooling.
6. Add 75 micro L of chilled DNA Extractor Solution III to the tube, and briefly vortex the tube, leaving the tube on ice for 5 minutes.
 7. Centrifuge the tube at $12,000 \times g$ for 5 minutes at 4 degrees C, followed by transfer of the supernatant to another tube by papering.
 〈Note〉 Transfer of the supernatant from the tube should be done, avoiding precipitate material to transfer to the new tube. On taking brownish material, try again from the centrifugation.
 8. Add 300 micro L of Sodium Iodide Solution to the tube, and mixing.
 9. Add 500 micro L of Isopropanol to the tube, and mix well.
 〈Note〉 The solution in the tube may be colored after addition of some isopropanol reagent, probably, due to deterioration of the reagent or its insufficient quantity. Analytical grade of Isopropanol should be used.
 10. Centrifuge the tube at $12,000 \times g$ for 10 minutes at room temperature, followed by decantation to remove the supernatant.
 〈Note〉 The precipitate can be seen due to co-precipitant containing in the solution provided.
 11. Add 1 mL of the Washing Solution to the tube, and is followed by brief vortexing.
 12. Centrifuge the tube at $12,000 \times g$ for 5 minutes at room temperature, followed by decantation.
 13. Repeat the wash consisted of Step 11 and 12.
 14. After drying the precipitate under reduced pressure, dissolve it in an appropriate volume of TE, at pH 8.
 〈Note〉 The final product contains RNA, which can be removed by digestion of RNase treatment : incubation with 10 micro g/mL RNase for 1 hour at room temperature.

B. Isolation of mtDNA from cultured cell

1. Rinse culture cells twice with PBS on the culture dish, followed by collection of the cells using cell scribe into glass homogenizer.
2. Proceed the whole steps from 1-14 described in A.


[Storage]

Keep at 2~10 degrees C.

[References]

- 1) 安田和基 : 医学のあゆみ, **174**, 420 (1995).
- 2) Davis, R. E. *et al.* : *Proc. Natl. Acad. Sci. USA*, **94**, 4526 (1997).
- 3) Poulton, J. *et al.* : *Prenatal Diagnosis*, **16**, 1247 (1996).
- 4) Jazin, E. E. *et al.* : *Proc. Natl. Acad. Sci. USA*, **93**, 12382 (1996).
- 5) 田中雅嗣, 小澤高将 : 日本臨床, **47**, 304 (1989).

mtDNA Extractor CT Kit (catalog No. : 291-55301, 25 tests)
【Protocol】

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
 Telephone : + 81-6-6203-3741
 Facsimile : + 81-6-6201-5964
<http://ffwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
 Richmond, VA 23237
 U.S.A.
 Telephone : + 1-804-271-7677
 Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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