## LabAssay™ Creatinine

Please, read this instruction carefully before use.

### 1. Introduction

Creatinine is a breakdown product of creatine phosphate in muscle. It is mainly filtered by the kidneys, though a small amount is actively secreted. LabAssay™ Creatinine is based on an *in vitro* colorimetric Jaffé method for the quantitative determination of creatinine in mouse serum or urine. It is a simultaneous multi-sample assay format using a microplate, but measurements can also be made using a test tube.

### 2. Kit Contents

	Components	State	Amount
(1)	Deproteinizing Reagent	liquid	150 mL/1 bottle
(2)	Picric Acid Reagent	liquid	50 mL/1 bottle
(3)	0.75 mol/L Sodium Hydroxide Solution	liquid	50 mL/1 bottle
(4)	Standard Solution	liquid	15 mL/1 bottle

※ Assay in a microplate ; 500 tests. Assay in a test tube ; 50 tests.

### 3. Storage and Expiration

When the complete kit is stored at 2°C - 10°C, the kit is stable until the expiration date shown on the label on the container.

### 4. Assay Principle

The supernatant is separated, after adding deproteinizing reagent to a sample and centrifuging. Creatinine in the sample reacts with picric acid in alkaline solution to produce tangerine condensate (Jaffé reaction).

Quantitation of total creatinine in the sample can be made by measurement of the absorbance.

### 5. Materials and Apparatuses to be prepared

- 96 wells microplate (transparent type) • Micropipette • Incubator maintaining at 25°C - 30°C \*
- Plate mixer \*
- Microplate reader with 520 nm wavelength filter
- ( \* if the microplate reader is not equipped.)

#### (For Test Tube method)

- Test tube • Pipette • Spectrophotometer or colorimeter with 500 nm - 530 nm wavelength filter

### 6. Preparation of Reagents to be used

- ① Deproteinizing Reagent : ready-to-use
- ② Picric Acid Reagent : ready-to-use
- ③ 0.75 mol/L Sodium Hydroxide Solution : ready-to-use
- ④ Standard Solution

Standard solution is prepared by dilution of the provided Creatinine Standard

No.	Provided Std.	Distilled water for dilution	Sampling volume of diluted std.	Creatinine concentration
1	50 $\mu$ L	150 $\mu$ L	50 $\mu$ L	2.5 mg/dL
2	100 $\mu$ L	100 $\mu$ L	50 $\mu$ L	5.0 mg/dL
3	150 $\mu$ L	50 $\mu$ L	50 $\mu$ L	7.5 mg/dL
4	undiluted	—	50 $\mu$ L	10.0 mg/dL

### 7. Preparation of Samples

- Heparin, citrate, oxalate and EDTA as an anticoagulant and sodium fluoride as an antiglycolysis may not affect the assay.
- Toluene should be used as the preservative for urine samples.
- A substance with an active methylene group causes an increase of the value.
- Ascorbic acid, bilirubin and hemolysis may not significantly affect the assay.

## 8. Procedure

### (1) Assay in a Microplate

Perform the assay according to the following table scheme.

	Test	Standard	Blank
<b>First operations should be performed in the sample tubes.</b>			
Sample	Sample 50 $\mu$ L	Standard solution 50 $\mu$ L	Water 50 $\mu$ L
Deproteinizing Reagent	300 $\mu$ L	300 $\mu$ L	300 $\mu$ L
	Mix well, and stand at room temperature for 10 min. Centrifuge over 2500 rpm for 10 min. <b>Subsequent operations should be performed in the wells.</b>		
Supernatant	100 $\mu$ L	100 $\mu$ L	100 $\mu$ L
	Stand at 25°C - 30°C		
Picric Acid Reagent	50 $\mu$ L	50 $\mu$ L	50 $\mu$ L
0.75 mol/L Sodium Hydroxide Solution	50 $\mu$ L	50 $\mu$ L	50 $\mu$ L
	Stand at 25°C - 30°C for 20 min. Measure the absorbance of the test sample and standard solution at 520 nm with the blank solution as the control within the next 30 min.		


### (2) Assay in a Test Tube

Perform the assay in a test tube according to the following table scheme.

	Test	Standard	Blank
Sample	Sample 0.5 mL	Standard solution 0.5 mL	Water 0.5 mL
Deproteinizing Reagent	3.0 mL	3.0 mL	3.0 mL
	Add with shaking, and stand at room temperature for 10 min. Centrifuge over 2500 rpm for 10 min.		
Supernatant	2.0 mL	2.0 mL	2.0 mL
	Stand at 25°C - 30°C		
Picric Acid Reagent	1.0 mL	1.0 mL	1.0 mL
0.75 mol/L Sodium Hydroxide Solution	1.0 mL	1.0 mL	1.0 mL
	Stand at 25°C - 30°C for 20 min. Measure the absorbance of the test sample and standard solution with the blank solution as the control within the next 30 min. Colorimeter with 500 nm - 530 nm wavelength filter. Spectrophotometer : 520 nm		

## 9. Standard Curve

[assay in a microplate]

## 10. Performance

### (1) Sensitivity [assay in a Test Tube]

- The absorbance is 0.010 to 0.020, when measuring purified water as a sample.
- The absorbance is 0.400 to 0.500, when measuring 10 mg/dL creatinine as a sample

### (2) Specificity

The creatinine concentration is less than  $\pm 10\%$ , when measuring the known concentration of a control serum as a sample.

## 11. Usage Notes

### (1) Notes on the assay

- Do not use the reagents described above in any procedures other than those described herein. Performance cannot be guaranteed if the reagents are used in other procedures or for other purposes.
- Operate the measurement apparatuses according to operator's manuals under appropriate conditions. Consult the apparatus manufacturer for details.
- Keep the provided reagents under the indicated conditions before the expiration.
- The measurement should be done at an absorbance of 520 nm, although maximum absorption of the pigment is 480 nm, because picric acid may cause a high blank.
- Reaction temperature should be constant, because the temperature affects the assay value.
- Picric Acid Reagent and 0.75 mol/L Sodium Hydroxide Solution can be mixed beforehand. However, the mixed solution should be used within the day.
- 0.75 mol/L Sodium Hydroxide Solution should not be used if it is contaminated by Picric Acid Solution. Pay especially attentions not to contaminate with pipette.
- Urine creatinine samples below 10 mg/dL should be measured without the dilution.
- Acetone should not be used for drying of tubes, because it effects the absorbance.
- In the case of urine, if there are no opacity after addition of Deproteinizing Reagent, it is not necessary to centrifuge.
- This kit is for research use only, but not for diagnostic use.

### (2) Safety Precautions

- If a reagent comes into contact with the mouth, eyes, or skin, immediately wash with a lot of water. Consult a physician if necessary.
- A pipette should be used with a safety pipette filler.
- Vials should be opened carefully

### (3) Waste

- The waste should be processed appropriately according to the law.
- All the devices including reagents and vials that come into contact with specimens should be considered potentially infectious

## [References]

1. Bonsnes, R. W. and Taussky, H. H. : *J. Biol. Chem.*, 158, 581 (1945).
2. Henry, R. J. : *Clinical Chemistry*, 287 (Harper & Row), New York (1966).

LabAssay™ Creatinine

[Storage]	Store at 2°C - 10°C
[Expiration date]	Indicated on the label
[Package]	For 500 tests
[Cat #]	291-93901

---

## FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://fwk.fujifilm.co.jp>

### FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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