LabAssay™ Glucose

Please, read this instruction carefully before use.

1. Introduction

α -D-Glucose and β -D-glucose in solution maintain equilibrium in a constant ratio. Glucose oxidase reacts only with β -D-glucose and does not react with α -D-Glucose. Therefore, α -D-Glucose is converted into β -D-glucose using mutarotase.

LabAssay™ Glucose is the reagent kit for assay of glucose based on an enzymatic method with a combination of mutarotase and glucose oxidase.

The kit is used for the quantitative determination of glucose in mouse serum, plasma or urine. It is a simultaneous multi-sample assay format using a microplate, but measurements can also be made using a test tube.

2. Kit Contents

	Components	State	Amount
(1)	Buffer Solution	liquid	150 mL/1 bottle
(2)	Chromogen Substrate	freeze dry	For 150 mL/1 bottle
(3)	Standard Solution I	liquid	3 mL/1 bottle
(4)	Standard Solution II	liquid	3 mL/1 bottle


※ Assay in a microplate ; 500 tests. Assay in a test tube ; 50 tests.

3. Storage and Expiration

When the complete kit is stored at 2°C - 10°C, the kit is stable until the expiration date shown on the label on the container.

4. Assay Principle

When a sample is mixed with the Chromogen Reagent, the α -form of glucose in the sample is converted to β -form by mutarotase. β -D-Glucose is oxidized and yields hydrogen peroxide by glucose oxidase (GOD). In the presence of peroxidase (POD), the formed hydrogen peroxide yields a red pigment by quantitative oxidation condensation with phenol and 4-aminoantipyrine. The glucose concentration is obtained by measuring absorbance of the red pigment.

5. Materials and Apparatuses to be prepared

- 96 wells microplate (transparent type) • Micropipette • Incubator maintaining at 37°C *
- Plate mixer*
- Microplate reader with 505 nm wavelength filter (* if the microplate reader is not equipped.)

(For Test Tube method)

- Test tube • Pipette • Incubator maintaining at 37°C
- Spectrophotometer or colorimeter with 505 nm wavelength filter

6. Preparation of Reagents to be used

①Chromogen Reagent :

Dissolve 1 bottle of Chromogen Substrate with 150 mL of Buffer Solution. After reconstitution, the solution should be stored at 2°C - 10°C and used within 1 month.

②Standard Solution (Microplate method)

Standard solution is prepared by dilution of the provided Standard Solution.

No.	Standard Solution I	Standard Solution II	Distilled or deionized water	Sample volume	Concentration
1	50 μ L	–	150 μ L	2 μ L	50 mg/dL
2	100 μ L	–	100 μ L	2 μ L	100 mg/dL
3	undiluted	–	–	2 μ L	200 mg/dL
4	–	150 μ L	100 μ L	2 μ L	300 mg/dL
5	–	200 μ L	50 μ L	2 μ L	400 mg/dL
6	–	Undiluted	–	2 μ L	500 mg/dL

7. Preparation of Samples

Sample

- Samples should be used immediately after collecting.
- The glucose in whole blood decreases due to glycolysis in erythrocytes when a sample is kept after collecting blood. Therefore, separation of blood cells should be done as quickly as possible.
- Ascorbic acid and bilirubin cause a slightly negative effect on the assay.
- Hemolysis gives a slightly positive effect on the assay.
- Anticoagulants such as heparin, citrate, oxalate and EDTA and sodium fluoride as a glycolytic inhibitor, should not affect the assay when they are used in their usual amounts.

8. Procedure

(1) Assay in a Microplate

Perform the assay in the wells according to the following table scheme.

	Test	Standard	Blank
Chromogen reagent	300 μ L	300 μ L	300 μ L
Sample	2 μ L	2 μ L	– * ¹
	Mix well and incubate at 37°C for 5 min. Measure the absorbance at 505 nm ^{*2} of the test sample and standard solution with the blank solution as the control.		

*1 The addition of water is omitted because the absorbance difference in the presence of the addition of 2 μ L of water does not exhibit any practical use.

*2 Use 505 nm (Main) and 600 nm (Sub) for two wavelengths.

(2) Assay in a Test Tube


Perform the assay in a test tube according to the following table scheme.

	Test	Standard	Blank
Chromogen reagent	3.0 mL	3.0 mL	3.0 mL
Sample	20 μ L	20 μ L	– * ¹
	Mix well and incubate at 37°C for 5 min. Measure the absorbance of the test sample and standard solution with the blank solution as the control. Colorimeter with 505 nm. Spectrophotometer : 505 nm ^{*2}		

*1 The addition of water is omitted because the absorbance difference in the presence of the addition of 20 μ L of water does not exhibit any practical use.

*2 Use 505 nm (Main) and 600 nm (Sub) for two wavelengths.

9. Standard Curve

10. Performance

(1) Sensitivity [assay in a test tube]

- The absorbance is below 0.03, when measuring purified water as a sample.
- The absorbance is 0.40 - 0.55, when measuring 200 mg/dL glucose as a sample.

(2) Specificity

The glucose concentration is less than $\pm 12\%$, when measuring the known concentration of a control serum or a urine as a sample.

11. Usage Notes

(1) Notes on the Assay

- Do not use the reagents described above in any procedures other than those described herein. Performance cannot be guaranteed if the reagents are used in other procedures or for other purposes.
- Operate the measurement apparatuses according to operator's manuals under appropriate conditions. Consult the apparatus manufacturer for details.
- Keep the provided reagents under the indicated conditions before the expiration.
- Do not use reagents, which were frozen in error. Such reagents may give false results.
- After opening the reagents, it is recommended to use them immediately. The opened vials should be capped and kept under the specified conditions.
- Do not use the containers and other materials in the package for any purposes other than those described herein.
- Do not use the reagents at other reaction temperatures and reaction times than described herein.
- The reaction completes in about 4 min at 37°C. When the incubation is continued, the absorbance gradually decreases. Therefore, the incubation should be done within 20 min.
- Glucose can be also measured by a reaction at room temperature (15°C or more) for 15 min.
- Color level of the reaction solution should show very little change within 1 hour at room temperature.
- This kit is for research use only. Not for diagnostic use.

(2) Safety Precautions

- If a reagent comes into contact with the mouth, eyes, or skin, immediately wash with a lot of water. Consult a physician if necessary.
- A pipette with a safety pipette filler should be used.
- Be careful not to cut yourself with the aluminum cap when removing it from the vial.

(3) Waste

- The waste should be processed appropriately according to regulations.
Buffer contains 500 mg/L phenol.
- All the devices including reagents and vials that come into contact with specimens should be considered potentially infectious.

[Reference]

Miwa, I., Okuda, J., Maeda, K. and Okuda, G. : *Clin. Chim. Acta.*, 37, 538 (1972).

LabAssay™ Glucose

[Storage] Store at 2°C - 10°C

[Expiration date] Indicated on the container

[Cat #] 291-94001

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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