## LabAssay™ GLDH (AA-type)

### 1. Intended Use

Glutamate dehydrogenase (GLDH) is present in many microorganisms and in the mitochondria of eukaryotes. In mammals, it is located in the liver, kidneys, heart, and brain, and is particularly abundant in hepatocytes. Clinically, while scarcely present in normal serum, GLDH is released into the extracellular space when cells are damaged, acting as a leakage enzyme. It is also an essential enzyme for urea synthesis, facilitating the interconversion of glutamate and  $\alpha$ -ketoglutaric acid.


In non-clinical toxicity studies, elevated ALT levels in the blood have been used for a long time as an indicator of liver injury. However, blood ALT does not change in cases of mild liver injury, and there are reports that ALT levels can rise due to factors other than liver injury. Therefore, measuring GLDH which is abundant in the liver and highly specific as a liver injury marker in combination with ALT and other markers has drawn attention for its potential to improve the detection sensitivity of liver injury.

This kit is a research reagent for measuring GLDH in human, mouse, rat, dog, rabbit, and monkey serum and plasma. This kit contains research reagents for use in automated analyzers only.

### ◆Product Features

- Enables rapid measurement using an automated analyzer.
- This kit is used to measure GLDH in human, mouse, rat, dog, rabbit, and monkey serum and plasma.

### 2. Assay Principle

### 3. Kit Performance

- Sensitivity : 1.0 U/L
- Linearity : 150.0 U/L
- Precision

Within-run (intra-assay variation)

With three replicates, the average CV values were less than 20.0% for the low concentration samples and less than 15.0% for the high concentration samples.

### 4. Precautions

- This kit should be used by performers who have completed training in operating automated analyzers or are under the supervision of an instructor.
- Do not use for purposes or methods other than those described in the instruction manual.
- Be aware that assay procedures may vary depending on the type of automated analyzer used.
- Wear gloves, goggles, and protective clothing during preparation and operation.
- Do not allow reagents to come into contact with the skin. If reagents from this kit accidentally come into contact with the eyes, mouth, wounds, or skin, immediately perform first aid such as thoroughly rinsing with tap water, and seek medical attention if necessary.
- Do not eat, drink, or smoke in areas where this kit is being used.
- Handle samples with caution, treating them as infectious materials.
- Do not mouth-pipette reagents.
- Do not mix reagents from different lot numbers; performance and properties may change. Avoid adding reagents from the same lot, as it may affect assay results.
- Always store the kit refrigerated (2 - 10°C) and do not freeze.

## 5. Kit Component Reagents

Components	Use Status	Amount
(A) GLDH Standard	Freeze-dried. Use after reconstitution.	2 bottles
(B) Reaction Solution 1a (R1a)	Add the entire liquid volume to R1b and use it.	17.1 mL/2 bottles
(C) Reaction Reagent 1b (R1b)*	Freeze-dried. Use after reconstitution.	2 bottles
(D) Reaction Solution 2 (R2)	Ready for use	15 mL/1 bottle

\* (C) Reaction Reagent 1b (R1b) is in lyophilized form attached to the bottom of the bottle.

### [Storage]

(A) GLDH Standard

Store unopened GLDH standard at 2 - 10°C.

Reconstituted original standard solution (60.0 U/L) may be aliquoted into microtubes and stored frozen at -35°C or below.

Use within 45 days of freezing.

Repeated freeze-thaw cycles of the original standard solution (60.0 U/L) are strictly prohibited.

(B) Reaction Solution 1a (R1a)

Store at 2 - 10°C.

(C) Reaction Reagent 1b (R1b)

Store Reaction Reagent 1b (R1b) at 2 - 10°C.

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Store at 2 - 10°C.

## 6. Equipment Required But Not Included in the Kit

- Automated analyzer
- Saline solution
- Purified water (distilled water)
- Pipettes with replaceable tips
- Vortex mixer

## 7. Sample Preparation

### Serum/Plasma

Measure immediately after blood collection or freeze at -35°C or below for long-term storage. If high values are expected, dilute appropriately before the assay. Do not use severely hemolyzed or highly lipemic samples.

### [Precautions Regarding Samples]

- When diluting samples, please dilute them with purified water or saline. Multiply the obtained value by the dilution factor to calculate the measured value.
- Thaw frozen samples just before the assay and mix thoroughly. Avoid repeated freeze-thaw cycles, as they may compromise the accuracy of results.
- Ensure that the pH of the specimen is between 6.0 and 7.8.
- When the presence of interfering materials is suspected, confirm dilution linearity by testing at two or more different dilutions of the sample.

## 8. Reagent Preparation

- Do not use reagents that have passed the expiration date indicated on the kit's outer label.

### Reaction Solution 1 (R1)

#### Preparation Method

Add the entire contents of “(B) Reaction Solution 1a (R1a) (17.1 mL)” to “(C) Reaction Reagent 1b (R1b)” and dissolve. This is “Reaction Solution 1 (R1)”.

\* (C) Reaction Reagent 1b (R1b) is in lyophilized form attached to the bottom of the bottle.

\* Use within 3 weeks after preparation. Refer to “10. Technical Tips”.

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(D) Use Reaction Solution 2 (R2) as is.

\* Use within 3 weeks after opening. Refer to “10. Technical Tips”.

## (A) GLDH Standard

### Preparation Method

Reconstitute (A) GLDH Standard with purified water in the amount specified in "Reconstitution of standard"\* to prepare the original standard solution (60.0 U/L). Allow the reagent to dissolve at room temperature, then gently mix by inversion to avoid foaming (use of a rolling mixer or gentle pipetting is recommended).

- Avoid vigorous agitation, as it may reduce enzyme activity.
- Use the solution on the day of dissolution.
- The reconstituted original standard solution (60.0 U/L) may be aliquoted into microtubes and stored frozen at -35°C or below.
- Aliquot the original standard solution prior to measurement. Do not store any the original standard solution that has been used.
- Use within 45 days of freezing. Repeated freeze-thaw cycles of the original standard solution (60.0 U/L) are strictly prohibited.
- If the measured value of LabAssay™ GLDH Control (128-07131) falls outside  $\pm 15\%$  of the assigned value, perform calibration.

\*Find and check "Reconstitution of standard" on this product page.

**As the amount of purified water to be added varies by lot, be sure to check it for every lot.**


### Create a calibration curve

Create a calibration curve using two standard concentrations (0.0 and 60.0 U/L). Use purified water for the 0.0 U/L standard concentration.

## 9. Assay Procedure

### Operation Method with Automated Analyzer (e.g., Hitachi Automatic Analyzer 3500)

**\*Be aware that the operation method varies depending on the automated analyzer used.**

## 10. Technical Tips

- Avoid contaminating samples and reagents with impurities. Using disposable tips is recommended.
- Handle samples and reagents carefully to reduce foaming. Before placing samples in the automated analyzer, ensure there are no bubbles. If there are any bubbles, please remove them before starting the measurement.
- Always perform calibration before measuring samples.
- Soak used samples and consumables in 1% formalin, 2% glutaraldehyde, or  $\geq 0.1\%$  sodium hypochlorite solution for at least 1 hour. Alternatively, sterilize by autoclaving before disposal. Dispose of used consumables and unused chemicals in accordance with your institution's regulations and local laws.
- Although Reaction Solutions 1 and 2 can be stored after opening for up to 3 weeks when refrigerated (2 - 10°C), storage duration may be shorter depending on conditions. For example, if stored in the automated analyzer after measurement, the period may be shortened due to the analyzer's operating conditions. If not in use, remove from the analyzer, close the cap, and store.

## 11. Storage and Expiration Date of Kit

Store at 2 - 10°C and do not freeze. Once opened, use each reagent promptly, as storage conditions may affect reagent performance. Do not use reagents that have passed their expiration date.

## 12. Contact Information

For analytical parameters, please contact :

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku,

Osaka 540-8605, Japan

Email : ffwk-labchem-tec@fujifilm.com

LabAssay™ GLDH (AA-type)

[Storage] Store at 2 - 10°C (Do not freeze)

[Term of validity] Indicated on the label.

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