QCdetect™ dsRNA ELISA Kit

Please read this instruction carefully before use.

1. Introduction

In vitro transcription (IVT) is a key step in the development of mRNA-based therapeutics, including mRNA vaccines. During this process, double-stranded RNA (dsRNA) is generated as an unintended byproduct. Because dsRNA is highly immunogenic, accurate measurement and quantification of residual dsRNA are essential in the manufacturing of mRNA-based pharmaceuticals and vaccines.

This product is an ELISA kit specifically designed to measure the concentration of dsRNA with high sensitivity. This kit is intended for research use only.

2. Assay Principle

Each well of the assay plate is pre-coated with anti-dsRNA antibody clone J2. Standards or samples are added to the wells and incubated. Next, biotin-conjugated anti-dsRNA antibody clone J2 is added and allowed to react. Peroxidase-conjugated streptavidin is then added and incubated. The peroxidase remaining in the wells is allowed to react with the TMB solution. The reaction is stopped by adding the Stop Solution, and the resulting yellow product is measured colorimetrically at 450 nm (reference wavelength : 620 nm). A calibration curve is generated by plotting absorbance against standard concentrations and used to determine the dsRNA concentration in each sample.

3. Performance Characteristics

Assay range	UTP : 7.7 pg/mL - 1,000 pg/mL m1ψ : 15.4 pg/mL - 2,000 pg/mL
Target Analyte	dsRNA
Analysis sample	Purified mRNA solution
Measurement duration	Approx. 4.5 hours
Detection Method	Colorimetric assay
Measurement Wavelengths	Primary wavelength : 450 nm Secondary wavelength : 620 nm

4. Precautions

1. Bring all reagents of the kit to room temperature (20°C - 25°C) for approximately 2 hours before use.
2. Before beginning each step, prepare all reagents required for the next step.
3. For professional use only, beginners are advised to use this kit under the guidance of experienced person.
4. Avoid contact with the acidic Stop Solution and TMB Solution, which contains hydrogen peroxide and tetramethylbenzidine. Wear gloves and eye and clothing protection when handling these reagents.
5. Be careful not to allow the reagent solutions of the kit to touch the skin, eyes and mucus membranes. Especially be careful for the Stop Solution because it is sulfuric acid. The Stop Solution and the substrate solution may cause skin/eyes irritation. In case of contact with these wash skin/eyes thoroughly with water and seek medical attention, when necessary.
6. Do not drink, eat or smoke in the areas where assays are carried out.
7. Unused samples and used tips should be rinsed in more than 0.1% sodium hypochlorite solution for more than 1 hour, or be treated by an autoclave before disposal.
8. Dispose consumable materials and unused contents in accordance with applicable regional/national regulatory requirements.
9. Do not combine reagents from different lots.
10. In order to avoid dryness of wells, contamination of foreign substances and evaporation of dispensed reagents, never forget to cover the well plate with a Plate Seal supplied, during incubation.
11. ELISA can be easily affected by your laboratory environment. Room temperature should be at 20°C - 25°C strictly. Avoid airstream velocity over 0.4 m/sec. (including wind from air conditioner), and humidity less than 30%.
12. Use the same type of tubes or deep well plates when preparing standard solutions and mRNA samples. Differences in material or surface treatment can affect RNA adsorption and influence measurement results. Always use the same containers and addition methods for each measurement.

5. Kit Contents

5-1. Reagents supplied

Components	Use Status	Amount
(A) Antibody-coated Plate	Use after washing.	96 wells (8×12)/1 plate
(B) UTP dsRNA Standard	Concentrated. Use after dilution.	50 μ L/1 bottle
(C) m1 ψ dsRNA Standard	Concentrated. Use after dilution.	50 μ L/1 bottle
(D) Sample Reaction Buffer	Ready for use.	25 mL/1 bottle
(E) Antibody Reaction Buffer	Ready for use.	25 mL/1 bottle
(F) Biotin-conjugated Antibody	Lyophilized. Use after reconstitution.	1 bottle
(G) Peroxidase-conjugated Streptavidin Solution	Concentrated. Use after dilution.	120 μ L/1 bottle
(H) TMB Solution	Ready for use.	12 mL/1 bottle
(I) Stop Solution	Ready for use.	12 mL/1 bottle
(J) Wash Solution (10×)	Concentrated. Use after dilution.	100 mL/1 bottle
Plate Seal	Ready for use.	4 sheets

5-2. Storage and Stability

When using only a part of the kit, refer to the storage and stability information below for unused reagents.

(1) (A) Antibody-coated Plate

The plate can be used in individual strips.

If seal is not removed, put the strip back in a plastic bag with zip-seal originally used for well-plate container and store at 2°C - 10°C.

Stable until the expiration date.

(2) (B) UTP dsRNA Standard / (C) m1 ψ dsRNA Standard

If unused, store at 2°C - 10°C. Stable until the expiration date.

Note : Do not store standard solutions after dilution.

(3) (D) Sample Reaction Buffer / (E) Antibody Reaction Buffer

When using a portion of the solution, transfer slightly more than the required volume into a separate container.

Immediately cap the remaining solution tightly and store at 2°C - 10°C ; do not allow it to reach room temperature.

Stable until the expiration date.

(4) (F) Biotin-conjugated Antibody

Immediately cap the remaining solution tightly and store at 2°C - 10°C ; do not allow it to reach room temperature. Stable until the expiration date.

Store the reconstituted biotin-conjugated antibody stock solution at 2°C - 10°C and use within 2 weeks.

For longer-term storage, freeze and store at -20°C or below. Do not subject to repeated freeze-thaw cycles.

Note : Discard any remaining diluted solution.

(5) (G) Peroxidase-conjugated Streptavidin Solution

Immediately cap the remaining solution tightly and store at 2°C - 10°C ; do not allow it to reach room temperature. Stable until the expiration date.

When using only part of the solution, remove only the volume required for dilution from the refrigerator. Immediately cap the remaining undiluted stock tightly and store at 2°C - 10°C ; do not allow it to reach room temperature. Stable until the expiration date.

Discard any unused working solution (already diluted). The remaining undiluted solution may be stored tightly closed at 2°C - 10°C.

(6) (H) TMB Solution

When using only a part of the solution, transfer slightly more than the required volume into a separate container.

The remaining undiluted solution may be stored tightly closed at 2°C - 10°C.

Stable until the expiration date.

(7) (I) Stop Solution

When using only a part of the solution, transfer slightly more than the required volume into a separate container.

The remaining undiluted solution may be stored tightly closed at 2°C - 10°C.

Stable until the expiration date.

(8) (J) Wash Solution (10×)

Cap tightly and store at 2°C - 10°C. Stable until the expiration date.

The remaining undiluted buffer may be stored tightly closed at 2°C - 10°C. Dispose any unused diluted buffer.

6. Equipment or supplies required but not supplied ☐ : Use as a check box

- ☐ Deionized water (ultrapure water, RNase-free recommended)
- ☐ RNase-free tubes and/or deepwell plates for standards and antibody dilution
- ☐ Glassware for dilution of Wash Solution (10×) (a graduated cylinder, a bottle)
- ☐ Micropipettes (disposable tip type) (One capable of accurately pipetting 10 μ L, and one capable of accurately pipetting 200 μ L - 500 μ L)
- ☐ Repeater pipette capable of continuously dispensing 100 μ L (recommended)
- ☐ Paper towel to remove washing buffer remaining in wells.
- ☐ A vortex-type mixer
- ☐ A shaker for 96-well plate (approximately 600 rpm to 800 rpm)
- ☐ An automatic washer for 96-well plate (recommended), or a wash bottle with a jet nozzle
- ☐ A 96-well plate reader (450 nm $\pm$ 10 nm, 620 nm : 600 nm - 650 nm)
- ☐ Software for data analysis

7. Preparation of Samples

Dilute the sample as necessary to ensure that the concentration falls within the assay range. Use the same containers and addition methods for each measurement.

(1) Method for Measuring mRNA Sample Concentration

For mRNA samples after *in vitro transcription* (IVT) reactions, we recommend quantification using RNA-binding fluorescent dyes.

Our company has successfully performed quantification using "RNA Quantification Kit for Purified EV" (Product Code : 297-97901).

For purified samples (after impurity removal), in addition to the above method, quantification by absorbance measurement is also possible.

(2) Concentration Adjustment Method for UTP mRNA Samples and m1 ψ mRNA Samples

We recommend adjusting the concentration of UTP mRNA samples to 100 ng/mL or less, and m1 ψ mRNA samples to 500 ng/mL or less.

8. Preparation of Reagents

Bring all reagents of the kit to room temperature (20°C - 25°C) for approximately 2 hours before use.

8-1. [(B) UTP dsRNA Standard / (C) m1 ψ dsRNA Standard]

Prepare the standard solutions by mixing each concentration with the (D) Sample Reaction Buffer included in the kit, equilibrated to room temperature. Use the same containers and addition methods for each measurement.

Use (B) UTP dsRNA Standard when measuring UTP samples, and (C) m1 ψ dsRNA Standard when measuring m1 ψ samples. The following is provided as an example.

[(B) UTP dsRNA Standard]

Concentration (pg/mL)	Volume of UTP dsRNA standard solution	(D) Sample Reaction Buffer
1,000,000	1 μ g/mL original solution : 50 μ L	—
100,000	1 μ g/mL original solution : 10 μ L	90 μ L
1,000	100,000 pg/mL solution : 6 μ L	594 μ L
444	1,000 pg/mL solution : 200 μ L	250 μ L
198	444 pg/mL solution : 200 μ L	250 μ L
87.8	198 pg/mL solution : 200 μ L	250 μ L
39.0	87.8 pg/mL solution : 200 μ L	250 μ L
17.3	39.0 pg/mL solution : 200 μ L	250 μ L
7.71	17.3 pg/mL solution : 200 μ L	250 μ L
0	0 (Blank)	250 μ L

[(C) m1ψ dsRNA Standard]

Concentration (pg/mL)	Volume of m1ψ dsRNA standard solution	(D) Sample Reaction Buffer
1,000,000	1 μg/mL original solution : 50 μL	–
200,000	1 μg/mL original solution : 10 μL	40 μL
2,000	200,000 pg/mL solution : 6 μL	594 μL
889	2,000 pg/mL solution : 200 μL	250 μL
395	889 pg/mL solution : 200 μL	250 μL
176	395 pg/mL solution : 200 μL	250 μL
78.0	176 pg/mL solution : 200 μL	250 μL
34.7	78.0 pg/mL solution : 200 μL	250 μL
15.4	34.7 pg/mL solution : 200 μL	250 μL
0	0 (Blank)	250 μL

8.2. [(F) Biotin-conjugated Antibody]

Before beginning the procedure, reconstitute (F) Biotin-conjugated Antibody (lyophilized) in 120 μL of purified water equilibrated to room temperature (RNase-free recommended) and allow to stand for **at least 10 minutes**. Then dilute the required volume **100-fold** with (E) Antibody Reaction Buffer to prepare the biotin-conjugated antibody solution.

The reconstituted stock solution can be stored for up to 2 weeks at 2°C - 10°C. For longer-term storage, keep at -20°C or below.

Do not subject to repeated freeze-thaw cycles. Do not store the diluted solution; discard after use.

8.3. [(G) Peroxidase-conjugated Streptavidin Solution]

During the 1-hour incubation with biotin-conjugated antibody, dilute (G) Peroxidase-conjugated Streptavidin Solution **100-fold** with (E) Antibody Reaction Buffer equilibrated to room temperature.

8.4. [(J) Wash Solution (10×)]

Dilute (J) Wash Solution **10-fold** with purified water (RNase-free recommended) equilibrated to room temperature.

* All other reagents are ready for use as supplied.

9. Assay Procedure

- Before beginning each step, prepare all reagents to be dispensed in the following step.
 - Use of a multichannel pipette for continuous dispensing is recommended for the dispensing steps (2), (3), (6), (9), (12), and (15).
- (1) Wash the (A) Antibody-coated Plate by filling the wells with 300 μL of Wash Solution (1×) and discard 3 times, then strike the plate upside-down onto several sheets of paper towel to remove residual buffer in the wells.
If bubbles remain during manual washing using a wash bottle or other method, try allowing the wash solution to overflow only during the final wash cycle to remove bubbles remaining in the wells.
 - (2) Pipette 100 μL each of the standard dilution solutions into the respective wells.
 - (3) Pipette 100 μL of each sample, diluted with (E) Sample Reaction Buffer, into the sample wells.
 - (4) Cover the plate with a plate seal and incubate at room temperature (20°C - 25°C) for 2 hours with gentle shaking at approximately 500 rpm using a microplate shaker.*
 - (5) Discard the reaction mixture and rinse the wells as described in step (1).
 - (6) Pipette 100 μL of biotin-conjugated antibody solution into each well.
 - (7) Cover the plate with a plate seal and incubate at room temperature (20°C - 25°C) for 1 hour with gentle shaking at approximately 500 rpm using a microplate shaker.*
 - (8) Discard the reaction mixture and rinse the wells as described in step (1).
 - (9) Pipette 100 μL of diluted (G) Peroxidase-conjugated Streptavidin Solution into each well.
 - (10) Cover the plate with a plate seal and incubate at room temperature (20°C - 25°C) for 1 hour with gentle shaking at approximately 500 rpm using a microplate shaker.*
 - (11) Discard the reaction mixture and rinse the wells as described in step (1).
 - (12) Pipette 100 μL of (H) TMB solution (equilibrated to room temperature) into each well.
 - (13) Shake using a microplate shaker. (Recommended shaking conditions: approximately 600 rpm for 10 seconds, repeated 3 times.)
 - (14) Cover the plate with a plate seal and incubate at room temperature (20°C - 25°C) for 30 minutes under light-protected conditions.*
 - (15) Add 100 μL of (I) Stop Solution (equilibrated to room temperature) to each well, and shake using a microplate shaker. (Recommended shaking conditions: approximately 600 rpm for 10 seconds, repeated 3 times.)

(16) Promptly measure the absorbance of each well at 450 nm and at the reference wavelength of 620 nm (within the range of 600 nm - 650 nm).

[Supplementary notes]

Peel off the protective backing from the plate seal and apply the adhesive side to the plate.


Do not reuse plate seals.

10. Calculation

- (1) Create a calibration curve for each measurement, with the standard concentration (pg/mL) on the X-axis and absorbance on the Y-axis.
- (2) From the calibration curve, read the dsRNA (pg/mL) corresponding to the absorbance of the sample. If a diluted sample is used, the concentration reading is multiplied by the sample dilution ratio to obtain the measurement value.

*In computer software calculations, we recommend using 4 parameters or 5 parameters.

11. Standard Curve (Example)

12. Summary of Assay Procedure ☐ Use as a check box

- *First, read this instruction manual carefully and start your assay after confirmation of details.
- ☐ Bring the well-plate and all reagents back to 20°C - 25°C for 2 hours.
- ☐ (J) Wash Solution (10×) : Dilute **10-fold** with purified water (RNase-free recommended) equilibrated to room temperature.
- ☐ (B) UTP dsRNA Standard and (C) m1ψ dsRNA Standard dilution example :
Prepare dilutions using (D) Sample Reaction Buffer supplied with the kit, equilibrated to room temperature.
Use (B) UTP dsRNA Standard when measuring UTP samples, and (C) m1ψ dsRNA Standard when measuring m1ψ samples. The following is provided as an example.

Concentration (pg/mL)	100000	1000	444	198	87.8	39.0	17.3	7.71	0
(B)UTP dsRNA Standard (μL)	Orig.sol. : 10	6*	200*	200*	200*	200*	200*	200*	—
(D)Sample Reaction Buffer (μL)	90	594	250	250	250	250	250	250	250

Concentration (pg/mL)	200000	2000	889	395	176	78.0	34.7	15.4	0
(C)m1ψ dsRNA Standard (μL)	Orig.sol. : 10	6*	200*	200*	200*	200*	200*	200*	—
(D)Sample Reaction Buffer (μL)	40	594	250	250	250	250	250	250	250

*One rank higher standard.

- ☐ **Antibody-coated Plate**
- ☐ ↓ Washing 3 times (*①)
- ☐ **Samples/Standards** **100 μ L**
- ☐ ↓ Incubate at 20°C - 25°C for 2 hours with gentle shaking at approximately 500 rpm using a microplate shaker.
(F) Biotin-conjugated Antibody :
- ☐ Reconstitute in 120 μ L of purified water, then dilute **100-fold** with (E) Antibody Reaction Buffer equilibrated to room temperature.
- ☐ ↓ Washing 3 times (*①)
- ☐ **Biotin-conjugated Antibody Solution** **100 μ L**
- ☐ ↓ Incubate at 20°C - 25°C for 1 hour with gentle shaking at approximately 500 rpm using a microplate shaker.
(G) Peroxidase-conjugated Streptavidin Solution :
- ☐ Dilute (G) Peroxidase-conjugated Streptavidin Solution **100-fold** with (E) Antibody Reaction Buffer equilibrated to room temperature.
- ☐ ↓ Washing 3 times (*①)
- ☐ **Peroxidase-conjugated Streptavidin Solution** **100 μ L**
- ☐ ↓ Incubate at 20°C - 25°C for 1 hour with gentle shaking at approximately 500 rpm using a microplate shaker.
- ☐ ↓ Washing 3 times (*①)
- ☐ **TMB Solution** **100 μ L**
- ☐ ↓ Shaking (*②). Incubation for 30 minutes at 20°C - 25°C. (Standing *③)
- ☐ **Stop Solution** **100 μ L**
- ☐ ↓ Shaking (*②)
- ☐ **Measurement of absorbance (450 nm, Ref 620 nm) immediately.**
(Ref. wave cancels the dirt in the back of plate.)

(*①) Wash the (A) Antibody-coated Plate by filling the wells with 300 μ L of washing buffer and discard 3 times, then strike the plate upside-down onto several sheets of paper towel to remove residual buffer in the wells. If bubbles remain during manual washing using a wash bottle or other method, try allowing the wash solution to overflow only during the final wash cycle to remove bubbles remaining in the wells.

(*②) Recommended shaking conditions: approximately 600 rpm for 10 seconds, repeated 3 times.

(*③) Seal the plate during the reaction after shaking. Peel off the protective paper from the seal and stick the seal on the plate. Do not reuse the Plate Seal used once.

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