Presep[®] PFAS (150 mg/6 mL)

This product is a solid-phase extraction (SPE) column packed with a divinylbenzene-polymethacrylate-based polymer modified with anion-exchange functional groups. It is suitable for the extraction and recovery of per- and polyfluoroalkyl substances (PFAS), including perfluorooctanesulfonic acid (PFOS) and perfluorooctanoic acid (PFOA). The column offers high recovery rates for a wide range of PFASs, from short-chain to long-chain PFAS, which have been difficult to recover using conventional SPE columns.

【Example of Solid-Phase Extraction Procedure*】

Column Conditioning

- 1) CH₃OH with 0.1% NH₄OH 5 mL
- 2) CH₃OH 5 mL
- 3) H₂O 5 mL

Sample Loading/Concentration

Water sample (sample volume : -1 L)

↓

Load (flow rate: 10-20 mL/min.)

↓

Remove residual water (Air purge or N₂ purge)

↓

Elute (CH₃OH with 0.1% NH₄OH : 5 mL)

↓

Concentration (N₂ gas purge : → 0.5 mL)

↓

LC/MS(/MS) measurement (injection volume : 5 μL)

* This procedure is provided as an example only. Optimization of conditions may be required depending on the sample matrix and target analytes.

【Instruction for Use and Handling】

1. Prior to use, be sure to condition the column according to the recommended procedure. Drying of the packing agent during column conditioning may adversely affect recovery and reproducibility.
2. To ensure efficient recovery of target analytes, use appropriate solvents for column conditioning and elution. As an example, CH₃OH with 0.1% NH₄OH can be used.
3. Load the sample solution under vacuum or positive pressure so that a constant flow rate is maintained.
4. When processing samples containing large amounts of suspended solids, filter the sample with an appropriate filter prior to SPE to prevent column clogging.
5. Particular care should be taken to avoid PFAS contamination from laboratory equipment, containers, reagents, and airborne dust. Avoid the use of PTFE and other fluoropolymer-based laboratory ware, tubing, and sealing tape.

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